



cITm Instrument Compatibility Matrix

Document no.	Document type	Part document	Version	Status
10452469	SWC	000	16	FR

Electronically signed by

Phase	Role	Name	UserID	Date	Time
Author(s)	Signee	J . Martinez Gomar	MARTJ144	11.08.2015	18:23:26
Reviewer	Signee	M . Nabel	NABELM	12.08.2015	07:13:03
Reviewer	Signee	P . Tan	TANP9	12.08.2015	09:59:55
Release	Signee	J . Jaeger	JAEGERJ2	12.08.2015	11:28:05

Document validity A printed copy of this document is valid at the time of printing, afterwards the validity of the document has to be checked in the document management system or in the process manual.

Document classification

Characteristic	Characteristic value
Department	SW
Paper documents	No
filing	DHF (R&D)
process	CONFIG & CHANGE MGMT
document type	VARIOUS

cobas IT middleware

Instrument Driver Compatibility Matrix

Document	ID-Number	10452469 SWC	Version	16
	Filename	cobas_IT_middleware_Instrument_Driver_Compatibility_Matrix v16.xls		

	Project Role	Name	Signature	Date
Creation	SPL Instrument Drivers	see front page	see front page	see front page
Review	TPM cobas IT middleware	see front page	see front page	see front page
Review	IPM cobas IT middleware	see front page	see front page	see front page
Approval	Quality Management	see front page	see front page	see front page

	Version	Name	Reason For Change	Date
Changes	16	Jose Martínez	Added new driver releases: • ARCHITCTv2 v3.1.0 • bep2000 v3.0.0 • BIOFLASH v3.0.0 • cobas6kv2 v3.2.0 • cobas8100 v3.1.0 • cobasIT1k v3.0.0 • cobasLIAT v3.0.0 • COMBI13 v3.0.1 • Desaga2k v2.0.1 • IRMA v3.0.1 • LIAISON v3.0.0 • MAXMATPL v3.0.0 • SEBIACAP v3.2.0 • SYSSUITV3 v3.0.0 Release Notes updated: • u411 v2.1.2 Driver releases retired: • cobas8800 v3.0.0 • uri1100 v2.2.0 (replaced by u411 v2.1.2) • SYSSUIT v3.0.1 (replaced by SYSSUITV3 v3.0.0) • ilactop v2.0.1 (replaced by BIOFLASH v3.0.0) • liaison v2.4.0 (replaced by LIAISON v3.0.0) • phoresis v2.2.2 (replaced by SEBIACAP v3.2.0)	See approval date
	15	Jose Martínez	Added new driver releases: • Act5Diff v3.0.0 • AELEVEN v3.0.0 • amplilink v3.0.0 • au680 v2.2.1 • bc5800 v2.1.0 • cobas8k v3.2.0 • cobas4800 v3.1.0 • cobas6500 v3.1.0 • COMBI13 v3.0.0 • CP2000 v3.0.0 • ELECSYSv2 v3.1.0 • ImmuLite v2.3.0 • IQ200 v2.2.0 • IRMA v3.0.0 • LABUMAT v3.0.0 • liaison v2.4.0 • mek7222 v2.1.0 • p501-p701 v3.1.0 • PIANO v3.0.0 • Platinum v3.0.0 • POCT1A v3.1.0 • quanta v2.2.0 • RSADYN v3.3.0 • SYSSUIT v3.0.1 • TCA v3.1.0 • TEST1THL v2.0.2	19.05.2015
	14	Jose Martínez	Added new driver releases: • ABL700 v2.5.3 • AU680 v2.2.0 • BMSTAv2 v3.1.0 • cobas4800 v3.0.0 (SR v2.1) • miditron v2.0.2 Release Notes updated: • coasys v2.2.0	26.01.2015
	13	Jose Martínez	Added new driver releases: • cobas8800 v3.0.0	19.11.2014

12	Jose Martínez	Added new driver releases: <ul style="list-style-type: none"> • architctv2 v3.0.3 • biodiast v2.2.0 • ILAcITop v2.0.1 • MPA v3.2.0 • quanta v2.1.0 • RSA v3.3.0 • RSADYN v3.2.0 • spaplus v2.1.0 • SysSuit v2.3.1 • u411 v2.1.2 • URISED2 v2.1.0 • vidas v2.0.1 Release Notes updated: <ul style="list-style-type: none"> • c4800 v2.0.0 • CA6000 v2.2.2 	11.11.2014
11	Jose Martínez	Added new driver releases: <ul style="list-style-type: none"> • CA400 v2.0.1 • chorus v3.0.0 • mek6400 v2.0.1 • OCDiana v2.0.1 • AU680 v2.1.0 • JUNIOR v3.0.0 Release Notes updated: <ul style="list-style-type: none"> • CA6000 v2.2.2 	05.08.2014
10	Jose Martínez	Added new driver releases: <ul style="list-style-type: none"> • cobas8100 v3.0.1 • CobasInt v2.5.2 • EVOLIS v2.1.0 • hit912 v2.3.0 • ImmuLite v2.2.3 • IQ200 v2.1.0 • liaison v2.2.1 • RSA v3.2.0 • RSADYN v3.1.0 • PHORESIS v3.0.1 • SYSSUIT v2.3.0 • Uri2400 v2.2.0 • architctv2 v3.0.2 • ELECSYSV2 v3.0.2 	06.06.2014
09	Jose Martinez	Added new driver releases: <ul style="list-style-type: none"> • COBAS8K v3.0.5 • FreeEVO v2.0.0 • EMERALD v2.0.0 • AU400 v2.1.1 • MPA v3.1.1 • PHORESIS v2.2.2 • ABL700 v2.5.2 Release Notes updated: <ul style="list-style-type: none"> • OMNIS v2.3.2 • K1000 v2.1.1 	24.03.2014
08	Jose Martínez	cobasp612d 3.0.2 renamed to RSADYN 3.0.3	29.01.2014
07	Carlos Santin	Added new driver releases: <ul style="list-style-type: none"> • EVOLIS 2.0.0 • QUANTA 2.0.0 • SYSSUIT 2.2.1 • COBAS6500 3.0.0 • TEST1THL 2.0.1 • URI2400 2.1.1 • cobasp612d 3.0.2 • CA6000 2.2.2 	14.01.2014
06	Jose Martínez	Added new driver releases: <ul style="list-style-type: none"> • OSMOMTR 2.0.0 • C3112 3.0.0 • PENTASTM 2.0.2 • XT2000i 2.0.2 • KONELAB 2.1.0 • SPAPLUS 2.0.0 • TEST1THL 2.0.0 • REFLOTRN 2.0.1 • PHORESIS 3.0.0 • SRM 3.0.2 • ABL700 2.5.1 • OMNIS 2.3.2 • HIT912 2.2.0 • K1000 2.1.1 • COBAS8K 3.0.4 • COBAS8100 3.0.0 	30.10.2013

05	Carlos Santín	Added new driver releases: • COBAS6K2 3.0.1 • SYSXASTM 3.0.0 • IMMUCAP 2.2.0 • IMMULITE 2.2.2 • U500 2.0.0 • CENTAUR 2.0.2	01.08.2013
04	Carlos Santín	Added new driver releases: • RSA 3.1.2 • RSADYN 3.0.1 • ARCHITCT2 3.0.1 • ELECSYS2 3.0.1 • COBASINT 2.5.1 • COBAS8K 3.0.3 • IMMULITE 2.2.1 • LIAISON 2.2.0 • BMSTA2 3.0.0 • OMNIS 2.3.1 • AU680 2.0.0 • STKS2A 2.0.3 • URI1100 2.2.0 • COASYS 2.2.0 • AMPLILINK 2.2.3	29.05.2013
03	Carlos Santín	Added new driver releases: • C4800 2.0.0 • COBASINT 2.5.0 • HA8180 2.0.0 • MASSLYNX 2.0.0 • PHORESIS 2.2.1 • SRM 3.0.1 • URISYS 3.0.0	16.03.2013
02	Carlos Santín	Added new driver releases: • ARCHITCT2 3.0.0 • COBAS8K 3.0.1 • COBAS6K2 3.0.0 • ELECSYS2 3.0.0 • RSA 3.1.1 • RSADYN 3.0.0 • SRM 3.0.0	20.12.2012
01	David Barba	Added new driver releases: • MPA 3.1.0 • RSA 3.1.0 • IMMUCAP 2.1.0 • C111 2.1.0 • OMNIS 2.3.1 • URI2400 2.1.0	05.11.2012
00	David Barba	First released version	13.02.2012

Instrument			Driver							CM Version
Vendor	Instrument(s)	Instrument FW version (if applic.)	Driver Name	Driver Version	Version of cobas IT that driver is validated against		Version of cobas IT required for this driver			
					cIT 3000/5000	cITm	cIT 3000/5000	cITm		
77 Elektronika Kft	Urised2	2.0.0	urised2	2.1.0	2.03.11	1.01	All versions	>=1.0	12	
77 Elektronika Kft	LabUMat		LABUMAT	3.0.0	2.06.00	1.04	See Rel. Notes	>=1.04	15	
Abbott	Cell-Dyn Sapphire		sapphire	2.1.0	-	-	See Rel. Notes	>=1.0	0	
Abbott	Cell-Dyn Emerald		emerald	2.0.0	2.03.11	1.01.01	All versions	>=1.0	9	
Abbott	Architect system		ARCHITCTv2	3.1.0	2.07.00	1.04.00	See Rel. Notes	>=1.04	16	
Adams	A1c HA-8180		ha8180	2.0.0	-	-	All versions	>=1.0	3	
Alifax	Test1 THL		TEST1THL	2.0.2	-	1.04.01	All versions	>=1.0	15	
Amelung	Amax 200		amax200	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Amelung	Amax Destiny		amaxdest	2.1.0	-	-	See Rel. Notes	>=1.0	0	
Amelung	KC10		kc10	2.1.1	-	-	See Rel. Notes	>=1.0	0	
Amelung	KC4		kc10	2.1.1	-	-	See Rel. Notes	>=1.0	0	
Arkray Inc.	Aution Eleven AE-4020		AELEVEN	3.0.0	2.06.00	1.04.01	See Rel. Notes	>=1.03.03	15	
Aution	MAX AX-4280		ax4280	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Baumann Medical	Mediff 9000 LC		mediff9k	2.1.0	-	-	See Rel. Notes	>=1.0	0	
Bayer	Clinitek 500		bayclini	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Bayer	Clinitek 200+		bayclini	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Bayer	Clinitek Advantus		bayclini	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Bayer	Chiron 850		chiro850	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	ACT5 Diff Hematology Analyzers	1.03	Act5Diff	3.0.0	2.06.01	1.04.01	See Rel. Notes	>=1.03.03	15	
Beckman Coulter	Act5 diff with ASTM protocol		act5astm	2.0.4	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	Act5 with ASTM protocol		act5astm	2.0.4	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	AU480		au680	2.2.1	2.03.11	1.01.01	All versions	>=1.0	15	
Beckman Coulter	AU680		au680	2.2.1	2.03.11	1.01.01	All versions	>=1.0	15	
Beckman Coulter	AU5800		au680	2.2.1	2.03.11	1.01.01	All versions	>=1.0	15	
Beckman Coulter	DXC600 System		codxc600	2.1.0	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	Immage Immunochemistry System		couimage	2.0.2	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	ACT8		coum8	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	MD2		coum8	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	MD8		coum8	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	HMX		stks2a	2.0.3	2.03.09	-	All versions	>=1.0	4	
Beckman Coulter	LX 20 PRO		syncix	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	Access 2 System		caccess2	2.0.2	-	-	See Rel. Notes	>=1.0	0	
Beckman Coulter	ONYX Analyzer		onyx	2.0.1	-	-	See Rel. Notes	>=1.0	0	
Becton Dickinson	SediSystem		sedisys	2.0.1	-	-	See Rel. Notes	>=1.0	0	
Becton Dickinson	SediSys 15		sedisys	2.0.1	-	-	See Rel. Notes	>=1.0	0	
Biomerieux	MiniVidas		minvidas	2.0.3	-	-	See Rel. Notes	>=1.0	0	
Biomerieux	Thrombolyzer		coasys	2.2.0	2.03.07	1.04	All versions	>=1.0	14	
Biomerieux	Vidas	v4.7.0	vidas	2.0.1	2.03.06	1.04.02	See Rel. Notes	>=1.0	12	
Biorad	Variant II		biovari	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Biorad	Diastat D10		biodiast	2.2.0	2.03.04	1.01	All versions	>=1.0	12	
Biorad	Evolis		evolis	2.1.0	2.03.11	1.01	All versions	>=1.01	10	
Care Diagnostica	Economic		economic	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Caretium	XI-921B		xi921	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Dade Behring	BCT		behbct	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Dade Behring	BEP 3		behbep3	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Dade Behring	Fibrintimer		behfra	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Dade Behring	Dimension Clinical Chemistry System	7.3.1	dimencc	2.3.0	-	-	See Rel. Notes	>=1.0	0	
Dade Behring	BEP 2000		BEP2000	3.0.0	-	1.04.00	See Rel. Notes	>=1.03.03	16	
Desaga GmbH	Desaga Sediplus S2000		desaga2k	2.0.1	2.07.00	1.04.01	See Rel. Notes	>=1.0	16	
Diagnova RFCL Limited	Best 2000		best2000	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Diamond Diagnostics	Furuno CA-400 Chemistry Analyzer		ca400	2.0.1	2.03.11	1.01	All versions	>=1.0	11	
Diamond Diagnostics	Furuno CA-90 Chemistry Analyzer		ca400	2.0.1	2.03.11	1.01	All versions	>=1.0	11	
DiaSorin	Liaison System		LIAISON	3.0.0	2.07.01	1.04.01	See Rel. Notes	>=1.04	16	
DiaSorin	Liaison XL		LIAISON	3.0.0	2.07.01	1.04.01	See Rel. Notes	>=1.04	16	
Diesse	Chorus TRIO		chorus	3.0.0	-	1.03.03	See Rel. Notes	>=1.03.03	11	
DPC Biermannn / Siemens	Immolute one		immulte	2.3.0	2.03.11	1.02	All versions	>=1.0	15	
DPC Biermannn / Siemens	Immolute 1000		immulte	2.3.0	2.03.11	1.02	All versions	>=1.0	15	
DPC Biermannn / Siemens	Immolute 2000		immulte	2.3.0	2.03.11	1.02	All versions	>=1.0	15	
DPC Biermannn / Siemens	Immolute 2500		immulte	2.3.0	2.03.11	1.02	All versions	>=1.0	15	
Eiken Chemical	OC-Diana		ocdiana	2.0.1	2.03.04	1.01.01	See Rel. Notes	>=1.0	11	
EKF Diagnostic	BIOSSEN 5040		bios5040	2.0.0	-	-	See Rel. Notes	>=1.0	0	
EKF Diagnostic	BIOSSEN 5140	1.00	bios5040	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Eppendorf	Ebio		ebioplus	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Erba Lachema	Maxmat PL		MAXMATPL	3.0.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	16	
Erba Lachema	Maxmat PL II		MAXMATPL	3.0.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	16	
Erba Lachema	Multiflex		MAXMATPL	3.0.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	16	
General Electric	Celldyn 3200		ab3200	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Greiner	SRS 20/II		srs20ii	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Helena Biosciences Europe	Platinum		PLATINUM	3.0.0	-	1.03.03	See Rel. Notes	>=1.03.03	15	
Hitado Diagnostic Systems	LabURReader	02000243	laburead	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Horiba	MICROS		micros	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Horiba	Penta 60C+		pentastm	2.0.2	2.03.11	-	All versions	>=1.0	6	
Horiba	Pentra 80 / Pentra XL80		pentastm	2.0.2	2.03.11	-	All versions	>=1.0	6	
Horiba	Pentra DX 120		pentastm	2.0.2	2.03.11	-	All versions	>=1.0	6	
HUMAN GmbH	Combilyzer 13		COMBI13	3.0.1	2.07.00	1.04	See Rel. Notes	>=1.04	16	
HUMAN GmbH	Elisys Duo		evolis	2.1.0	2.03.11	1.01	All versions	>=1.01	10	
IL	ACL-1000		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-300		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-300/R		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-3000		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-3000+		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-100		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-200		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-8000		acl8000	2.1.0	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-9000		acl8000	2.1.0	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-10000		acl8000	2.1.0	-	-	See Rel. Notes	>=1.0	0	
IL	ACL-2000		acl1000	2.0.1	-	-	See Rel. Notes	>=1.0	0	
IL	ACL Futura, Elite, Advanced		acl8000	2.1.0	-	-	See Rel. Notes	>=1.0	0	
IL	ACL TOP		BIOFLASH	3.0.0	2.07.00	1.04.01	See Rel. Notes	>=1.03.03	16	
IL	BIOFLASH		BIOFLASH	3.0.0	2.07.00	1.04.01	See Rel. Notes	>=1.03.03	16	
INOVA Diagnostics	QUANTA Link	v3.1 / v3.2	quanta	2.2.0	2.03.11	1.01	All versions	>=1.01	15	
Instrumentation Laboratory SpA	Leucodiff 700		leuco700	2.0.1	-	-	See Rel. Notes	>=1.0	0	
Instrumentation Laboratory SpA	Leucodiff 800plus		leuco700	2.0.1	-	-	See Rel. Notes	>=1.0	0	
Iris International	iQ200		iq200	2.2.0	2.03.11	1.01	All versions	>=1.0	15	
ITC	IRMA SL		IRMA	3.0.1	2.07.00	1.04	See Rel. Notes	>=1.04	16	
ITC	IRMA SL 2000		IRMA	3.0.1	2.07.00	1.04	See Rel. Notes	>=1.04	16	
Kabe	ESR 48 System		kaesr48	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Kabe	Glucoseanlysatör GA1		kabega1	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Menarini	MAX AX-4280		ax4280	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Menarini	HA-8160		ha8160	2.0.0	-	-	See Rel. Notes	>=1.0	0	
Menarini	AM-6050		osmomtr	2.0.0	-	1.01.01	-	>=1.0	6	
Mindray	BC-5800		bc5800	2.1.0	2.05.00	1.01.01	All versions	>=1.0	15	

Mission	U500		u500	2.0.0	2.03.11	1.01.01	All versions	>=1.0	5
MLA	Electra 1000C		el1600c	2.0.0	-	-	See Rel. Notes	>=1.0	0
MLA	Electra 1400C		el1600c	2.0.0	-	-	See Rel. Notes	>=1.0	0
MLA	Electra 1600C		el1600c	2.0.0	-	-	See Rel. Notes	>=1.0	0
MLA	Electra 1800C		el1600c	2.0.0	-	-	See Rel. Notes	>=1.0	0
Nihon Kohden	MEK6400	02-06	mek6400	2.0.1	2.03.01	1.01.01	See Rel. Notes	>=1.0	11
Nihon Kohden	MEK7222	05-02	mek7222	2.1.0	2.03.05	1.01	All versions	>=1.0	15
Nihon Kohden	MEK8222		mek7222	2.1.0	2.03.05	1.01	All versions	>=1.0	15
Nova	Electrolyte Analyzer		nov4plus	2.0.0	-	-	See Rel. Notes	>=1.0	0
Olympus	AU400		au400	2.1.1	2.01.04	1.01.01	All versions	>=1.0	9
Olympus	AU640		au400	2.1.1	2.01.04	1.01.01	All versions	>=1.0	9
Olympus	AU2700		au2700	2.0.0	-	-	See Rel. Notes	>=1.0	0
Orgentec	Alegria		alegria	2.0.0	-	-	See Rel. Notes	>=1.0	0
Ortho Clinical Diagnostics	Vitros ECiQ Immunodiagnostic System		vitros	2.2.0	-	-	See Rel. Notes	>=1.0	0
Ortho Clinical Diagnostics	Vitros Chemistry 950/700/700C/700R/700S/500/350/250		vitrosc	2.0.2	-	-	See Rel. Notes	>=1.0	0
Ortho Clinical Diagnostics	Vitros 5.1 FS		vitros51	2.0.4	-	-	See Rel. Notes	>=1.0	0
Ortho Clinical Diagnostics	Vitros 3600		vitros36	2.0.0	-	-	See Rel. Notes	>=1.0	0
PersonalLab	Mini-BOS		minibos	2.0.5	-	-	See Rel. Notes	>=1.0	0
Phadia	Immucapac Information Data Manager		immucap	2.2.0	2.03.09	1.01.01	>=2.03.07	>=1.0	5
Pharmacia	MasterCap		unicap	2.0.1	-	-	See Rel. Notes	>=1.0	0
Pharmacia	UniCAP 100 RM		unicap	2.0.1	-	-	See Rel. Notes	>=1.0	0
Radim	Alisei		alisei	2.0.0	-	-	See Rel. Notes	>=1.0	0
Radiometer	ABL555		ABL700	2.5.3	2.04.00	1.01	All versions	>=1.0	14
Radiometer	ABL700		ABL700	2.5.3	2.04.00	1.01	All versions	>=1.0	14
Radiometer	ABL800 FLEX		ABL700	2.5.3	2.04.00	1.01	All versions	>=1.0	14
Radiometer	ABL80 COOX		ABL700	2.5.3	2.04.00	1.01	All versions	>=1.0	14
Roche Diagnostics	MPA	06-11	MPA	3.2.0	2.07.00	1.01	>=2.07.00	>=1.01	12
Roche Diagnostics	cobas p 512 (Standard Interface)	1.4.1	RSA	3.3.0	-	1.01.01	-	>=1.01	12
Roche Diagnostics	cobas p 612 (Standard Interface)	1.4.1	RSA	3.3.0	-	1.01.01	-	>=1.01	12
Roche Diagnostics	cobas p 312 (Dynamic Interface)	1.1.1	RSADYN	3.3.0	2.07.00	1.01.01	See Rel. Notes	>=1.01	15
Roche Diagnostics	cobas p 512 (Dynamic Interface)	1.4.1	RSADYN	3.3.0	2.07.00	1.01.01	See Rel. Notes	>=1.01	15
Roche Diagnostics	cobas p 612 (Dynamic Interface)	1.4.1	RSADYN	3.3.0	2.07.00	1.01.01	See Rel. Notes	>=1.01	15
Roche Diagnostics	cobas p 471 (Dynamic Interface)	1.3.0	RSADYN	3.3.0	2.07.00	1.01.01	See Rel. Notes	>=1.01	15
Roche Diagnostics	cobas p 671 (Dynamic Interface)	1.3.0	RSADYN	3.3.0	2.07.00	1.01.01	See Rel. Notes	>=1.01	15
Roche Diagnostics	STA		BMSTAv2	3.1.0	2.07.00	1.04	>=2.07.00	>=1.04	14
Roche Diagnostics	STA compact		BMSTAv2	3.1.0	2.07.00	1.04	>=2.07.00	>=1.04	14
Roche Diagnostics	STA-R		BMSTAv2	3.1.0	2.07.00	1.04	>=2.07.00	>=1.04	14
Roche Diagnostics	cobas 1411 (ASTM)		BMSTAv2	3.1.0	2.07.00	1.04	>=2.07.00	>=1.04	14
Roche Diagnostics	Coasys 110		coasys	2.2.0	2.03.07	1.04	All versions	>=1.0	14
Roche Diagnostics	cobas 1411		coasys	2.2.0	2.03.07	1.04	All versions	>=1.0	14
Roche Diagnostics	Elecsys 1010		ELECSYSv2	3.1.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	Elecsys 2010		ELECSYSv2	3.1.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	cobas e 411		ELECSYSv2	3.1.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	Hitachi 912		hit912	2.3.0	2.01.04	1.01.01	All versions	>=1.0	10
Roche Diagnostics	Hitachi 917		hit917	2.0.0	-	-	See Rel. Notes	>=1.0	0
Boehringer	Miditron		miditron	2.0.2	2.03.11	1.01.01	>=2.03.11	>=1.01	14
Boehringer	Miditron M		miditron	2.0.2	2.03.11	1.01.01	>=2.03.11	>=1.01	14
Boehringer	Miditron Junior		miditron	2.0.2	2.03.11	1.01.01	>=2.03.11	>=1.01	14
Roche Diagnostics	Urisys 1800 - miditron mode		miditron	2.0.2	2.03.11	1.01.01	>=2.03.11	>=1.01	14
Roche Diagnostics	cobas u 411 (ASTM)		u411	2.1.2	2.03.11	1.01.01	All versions	>=1.0	16
Roche Diagnostics	Urisys 1100 (ASTM)		u411	2.1.2	2.03.11	1.01.01	All versions	>=1.0	16
Roche Diagnostics	Urisys 1800 (ASTM)		u411	2.1.2	2.03.11	1.01.01	All versions	>=1.0	16
Roche Diagnostics	Integra 400		cobasint	2.5.2	-	1.01.01	All versions	>=1.0	10
Roche Diagnostics	Integra 700		cobasint	2.5.2	-	1.01.01	All versions	>=1.0	10
Roche Diagnostics	Integra 800		cobasint	2.5.2	-	1.01.01	All versions	>=1.0	10
Roche Diagnostics	Amplink		amplink	3.0.0	2.06.00	1.01.01	See Rel. Notes	>=1.0	15
Roche Diagnostics	902 System		hi902	2.0.0	-	-	See Rel. Notes	>=1.0	0
Roche Diagnostics	cobas c 111		c111	2.1.0	-	-	See Rel. Notes	>=1.0	1
Roche Diagnostics	cobas c 311		C311v2	3.0.0	-	1.01.00	-	>=1.0	6
Roche Diagnostics	UriSys 2400		uri2400	2.2.0	2.03.11	1.01	All versions	>=1.0	10
Roche Diagnostics	cobas 8000		COBAS8K	3.2.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	cobas 6000		COBAS6Kv2	3.2.0	2.07.00	1.04.00	See Rel. Notes	>=1.04	16
Roche Diagnostics	cobas p 501		p501-p701	3.1.0	-	1.01	-	>=1.01	15
Roche Diagnostics	cobas p 701		p501-p701	3.1.0	-	1.01	-	>=1.01	15
Roche Diagnostics	cobas 4800 System	SR v1.1	c4800	2.0.0	2.03.07	1.04.02	>=2.03.07	>=1.0	12
Roche Diagnostics	cobas 4800 System	SR v2.1	cobas4800	3.1.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	cobas 6500	2.0.1	cobas6500	3.1.0	2.06.01	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	u601	2.0.1	cobas6500	3.1.0	2.06.01	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	u701	2.0.1	cobas6500	3.1.0	2.06.01	1.04.01	See Rel. Notes	>=1.04.01	15
Roche Diagnostics	Urisys1800	4.0.0.0704	URISYS1800	3.0.0	-	1.0	>=2.05	>=1.0	3
Roche Diagnostics	cobas b 121 (Omni C)		omnis	2.3.2	2.03.11	1.01.01	All versions	>=1.01	9
Roche Diagnostics	cobas b 221 (Omni S)		omnis	2.3.2	2.03.11	1.01.01	All versions	>=1.01	9
Roche Diagnostics	cobas b 123		omnis	2.3.2	2.03.11	1.01.01	All versions	>=1.01	9
Roche Diagnostics	cobas 8100		cobas8100	3.1.0	-	1.02	-	>=1.02	16
Roche Diagnostics	Reflotrn Sprint		reflotrn	2.0.1	-	1.01.01	All versions	>=1.0	6
Roche Diagnostics	ACCU-Check Inform II		POCT1A	3.1.0	2.07.00	1.04	See Rel. Notes	>=1.04	15
Roche Diagnostics	cobas h 232		POCT1A	3.1.0	2.07.00	1.04	See Rel. Notes	>=1.04	15
Roche Diagnostics	cobas b 123		POCT1A	3.1.0	2.07.00	1.04	See Rel. Notes	>=1.04	15
Roche Diagnostics	cobas b 101		POCT1A	3.1.0	2.07.00	1.04	See Rel. Notes	>=1.04	15
Roche Diagnostics	cobas IT 1000 application		cobasIT1k	3.0.0	-	1.04.00	See Rel. Notes	>=1.03.03	16
Roche Diagnostics	cobas LIAT System (HL7)		cobasLIAT	3.0.0	2.07.00	1.04.01	See Rel. Notes	>=1.04.01	16
Sebia	CAPILLARYS / PHORESIS		SEBIACAP	3.2.0	2.06.01	1.04.02	See Rel. Notes	>=1.03.03	16
Sekisui Medical Co	Coapresta 2000		CP2000	3.0.0	2.06.01	-	See Rel. Notes	>=1.03.03	15
Siemens	Centaur System		centaur	2.0.2	2.03.06	1.01.01	All versions	>=1.0	5
Siemens	BEP III		behbep3	2.0.0	-	-	See Rel. Notes	>=1.0	0
Siemens	BCS XP		behra	2.0.0	-	-	See Rel. Notes	>=1.0	0
Siemens	Clinitek Atlas		uri2400	2.2.0	2.03.11	1.01	All versions	>=1.0	10
Sigma	Amax 200		amax200	2.0.0	-	-	See Rel. Notes	>=1.0	0
Sigma	Amax Destiny		amaxdest	2.1.0	-	-	See Rel. Notes	>=1.0	0
SpinReactor	Spinlab 100		junior	3.0.0	-	1.03.03	See Rel. Notes	>=1.03.03	11
Sysmex	KX21		k1000	2.1.1	-	1.01.01	All versions	>=1.0	9
Sysmex	SE 9000		se9000	2.0.0	-	-	See Rel. Notes	>=1.0	0
Sysmex	SE 9500		se9000	2.0.0	-	-	See Rel. Notes	>=1.0	0
Sysmex	UF-50		sysuf50	2.0.1	-	-	See Rel. Notes	>=1.0	0
Sysmex	XS-800i		xt2000i	2.0.2	2.03.11	1.01.01	All versions	>=1.0	6
Sysmex	XT-1800i		xt2000i	2.0.2	2.03.11	1.01.01	All versions	>=1.0	6
Sysmex	XT2000i		xt2000i	2.0.2	2.03.11	1.01.01	All versions	>=1.0	6
Sysmex	CA500 (Family)		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CA540		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CA560		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CA600		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CA1500		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CA6000		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CA7000		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CS2000i		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CS2100i		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	CS5100		ca6000	2.2.2	2.03.11	1.01.01	All versions	>=1.0	12
Sysmex	SF3000		sf3000	2.0.0	-	-	See Rel. Notes	>=1.0	0
Sysmex	K-4500		k4500	2.0.0	-	-	See Rel. Notes	>=1.0	0
Sysmex	K1000		k1000	2.1.1	-	1.01.01	All versions	>=1.0	9

Sysmex	HC LAB		syshclab	2.1.0	-	-	See Rel. Notes	>=1.0	0
Sysmex	XS-1000i		SYSXASTM	3.0.0	-	1.02	-	>=1.02	5
Sysmex	UF-1000i (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XE-2100 (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XE-2100D (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XE-2100L (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XE-5000 (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XN-Series (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XS-800i (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XS-1000i (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XT-1800i (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Sysmex	XT-2100i (through eiPU)		SYSSUITV3	3.0.0	-	1.04.01	See Rel. Notes	See Rel. Notes	16
Tecan	Freedom EVO		freeevo	2.0.0	2.03.11	1.01.01	All versions	>=1.0	9
Thermo Scientific	Konelab series		konelab	2.1.0	2.04.022	1.01	All versions	>=1.0	6
Thermo Scientific	TCAutomation		TCA	3.1.0	-	1.05	-	>=1.05	15
Tokio Boeki Machinery	SPAPLUS		spaplus	2.1.0	2.03.11	1.01	>=2.03.09	>=1.0	12
TOSOH	PIANO TURBO PROGRAM		PIANO	3.0.0	-	1.04.01	See Rel. Notes	>=1.03.03	15
Trinity Biotech (formerly Amelung)	Amax 200		amax200	2.0.0	-	-	See Rel. Notes	>=1.0	0
Trinity Biotech (formerly Amelung)	Amax Destiny		amaxdest	2.1.0	-	-	See Rel. Notes	>=1.0	0
Waters	MassLynx		masslynx	2.0.0	2.03.07	-	All versions	>=1.0	3
YD Diagnostic	Uriscan Pro+		uriscan	2.0.1	-	-	See Rel. Notes	>=1.0	0
YD Diagnostic	Uriscan Super		uriscans	2.0.0	-	-	See Rel. Notes	>=1.0	0