

DECLARATION OF CONFORMITY

370 West 1700 South
Logan, Utah 84321 USA
435 752 6011
800 453 2725
fax 435 752 4127
web www.elitechgroup.com

Manufacturer's Name: ELITechGroup Inc.
Manufacturer's Address: 370 West 1700 South
Logan, Utah 84321-5294
USA

European Authorized Representative: **MT Promedt Consulting, GmbH**
Altenhofstrasse 80
D-66386 St. Ingbert / Germany
Tel.: +49 6894 - 58 10 20
Fax: +49 6894 - 58 10 21
www.mt-procons.com

Declares, the product

Product Name: Aerospray® TB Slide Stainer
Model Number: 7722
Type of Equipment: AC line-powered, motorized slide stainer


Conforms to the requirements of the following Standards and Specifications:**Safety**

CAN/CSA-C22.2 No. 61010-1/R:2009
CAN/CSA-C22.2 No. 61010-2-020:2009
CAN/CSA-C22.2 No. 61010-2-081/R:2009
CAN/CSA-C22.2 No. 61010-2-101:2004
EN 61010-1:2010
EN 61010-2-101:2002
EN 61010-2-081+A1:2003
EN 61010-2-020:2006
UL 61010-2-020:2002
UL 61010-1:2012

Thereby complying with the In-vitro Diagnostic Directive 98/79/EC and therefore bearing the IVD / CE Marking.



Dawn T. Perdue,
Director of Quality Assurance and Regulatory Affairs



Date

This declaration applies to the instrument itself and to the accessories and supplies listed hereafter.

Accessories and Supplies:

| | |
|------------------|--|
| AC-028 | 12-slide vented carousel for Aerospray |
| AC-034 | Nozzle tool for Aerospray |
| AC-035 | Hex wrench for Aerospray nozzles |
| AC-038 | Space-Saver Container with Cap, 5 L, Aerospray |
| AC-039 | Spigot for AC-038 container, Aerospray |
| AC-041 | Aerospray Drain Tube, 1.8 m |
| AC-043-05 | Reagent bottle with cap, 500 mL, 5/PK, Aerospray |
| AC-043-160 | Reagent bottle with cap, 500 mL, 160/PK, Aerospray |
| AC-057 | 30-slide vented carousel for Aerospray |
| AC-059 | Aerospray Nozzle Orifice Cleaning Wire, 3 m |
| AC-069 | Pump primer for Aerospray |
| AC-072 | 5L reagent bottle assembly w/barb fitting Aerospray |
| AC-110 | Safety Shield for Aerospray/Cytopro |
| AC-160 | Cytopro Rotor |
| AC-170 | Waste container for Aerospray |
| AC-181 | 1D Barcode Reader with Stand |
| AC-182 | Waste Container with level detection |
| AC-183 | 5L reagent bottle assembly w/level detection |
| AC-184 | Aerospray Nozzle Maintenance Kit |
| RP-464 | Manual, application, for Aerospray 7722 (TB) |
| SS-029 | Aerospray Nozzle Cleaning Solution, 355mL |
| SS-029C | Aerospray Nozzle Cleaning Solution, 250 mL Concentrate |
| SS-029C-EU | Aerospray Nozzle Cleaning Solution, 250 mL Concentrate |
| SS-029CG | Nozzle Cleaning Solution, 1.89 liter Concentrate |
| SS-029CG-EU | Nozzle Cleaning Solution, 2.5 liter Concentrate |
| SS-061A | Aerospray TB Reagent A: Carbol Fuchsin Decolorizer, 500 mL |
| SS-061A-EU | Aerospray TB Reagent A: Carbol Fuchsin Decolorizer, 500 mL |
| SS-061AF | Aerospray TB Reagent A: Fluorescence Decolorizer, 500 mL |
| SS-061AF-EU | Aerospray TB Reagent A: Fluorescence Decolorizer, 500 mL |
| SS-061BBG | Aerospray TB Reagent B: Brilliant Green, 500 mL |
| SS-061BBG-EU | Aerospray TB Reagent B: Brilliant Green, 500 mL |
| SS-061BCS-EU | Aerospray TB Reagent B: Carbol Solution, 500 mL |
| SS-061BMB | Aerospray TB Reagent B: Methylene Blue, 500 mL |
| SS-061BMB-EU | Aerospray TB Reagent B: Methylene Blue, 500 mL |
| SS-061BP | Aerospray TB Reagent B: Potassium Permanganate, 500 mL |
| SS-061BP-EU | Aerospray TB Reagent B: Potassium Permanganate |
| SS-061BRT-EU | Aerospray TB Reagent B: Thiazin Red, 500 mL |
| SS-061CA | Aerospray TB Reagent C: Auramine O, 500 mL |
| SS-061CA-EU | Aerospray TB Reagent C: Auramine, 500 mL |
| SS-061CAR | Aerospray TB Reagent C: Auramine Rhodamine (mod), 500 mL |
| SS-061CAR-EU | Aerospray TB Reagent C: Auramine Rhodamine, 500 mL |
| SS-061CEZ | Aerospray TB Reagent C: ELITe ZN™, 500mL |
| SS-061CEZ-EU | Aerospray TB Reagent C: ELITe ZN™, 500mL |
| SS-061M | Aerospray TB Mycohold Cell Adhesive, 500 mL |
| SS-061M-EU | Aerospray TB Mycohold Cell Adhesive, 500 mL |
| SS-105 | Aerospray Nozzle Caps, package of 6 |
| SS-133 | Decontamination Solution Concentrate, 15 mL |
| SS-161A | Aerospray TB Carbol Fuchsin Decolorizer Conc., 30 mL |
| SS-161AF | Aerospray TB Fluorescence Decolorizer Conc., 30 mL |
| SS-161BBG | Aerospray TB Brilliant Green Concentrate, 130 mL |
| SS-161BMB | Aerospray TB Methylene Blue Concentrate, 130 mL |
| SS-161BP | Aerospray TB Potassium Permanganate Conc., 200 mL |
| SS-161M | Aerospray TB Mycohold Cell Adhesive Conc., 100 mL |
| SS-266 | PM Cleaning Solution, 500 mL |
| SS-266-EU | PM Cleaning Solution, 500 mL |
| SS-267 | PM Chart for Aerospray TB (7722), pad of 24 sheets |
| SS-MeOH | Aerospray Reagent-Grade Methanol, 500 mL |