

## Olympus EndoTherapy - Your Picks

### HEMOSTASIS

#### SINGLE-USE INJECTION NEEDLES



- Extra-stiff sheath to enhance penetration and prevent kinking
- Large range for wide variety of applications
- Ergonomic handle with positive “click” action
- Different needle bevels (upper GI: regular; lower GI: short) ensure optimal results in both applications

Article Number	Article Name	Quantity	Min. Working Channel ø	Working Length	Needle Bevel	Needle Length (A)	Needle ø (B)	Area of Compatibility
N6138750	NM-400L-0421	5	2.8 mm	1650 mm	Regular	4 mm	21G = 0.8 mm	
N6139250	NM-400L-0621	5	2.8 mm	1650 mm	Regular	6 mm	21G = 0.8 mm	
N6138850	NM-400L-0423	5	2.8 mm	1650 mm	Regular	4 mm	23G = 0.6 mm	
N6139050	NM-400L-0523	5	2.8 mm	1650 mm	Regular	5 mm	23G = 0.6 mm	
N6139350	NM-400L-0623	5	2.8 mm	1650 mm	Regular	6 mm	23G = 0.6 mm	
N6138950	NM-400L-0425	5	2.8 mm	1650 mm	Regular	4 mm	25G = 0.5 mm	
N6139150	NM-400L-0525	5	2.8 mm	1650 mm	Regular	5 mm	25G = 0.5 mm	
N6139450	NM-400L-0625	5	2.8 mm	1650 mm	Regular	6 mm	25G = 0.5 mm	
N6139650	NM-400U-0423	5	2.8 mm	2300 mm	Short	4 mm	23G = 0.6 mm	
N6139850	NM-400U-0523	5	2.8 mm	2300 mm	Short	5 mm	23G = 0.6 mm	
N6140050	NM-400U-0623	5	2.8 mm	2300 mm	Short	6 mm	23G = 0.6 mm	
N6139750	NM-400U-0425	5	2.8 mm	2300 mm	Short	4 mm	25G = 0.5 mm	

**Olympus EndoTherapy - Your Picks**

N6139950	NM-400U-0525	5	2.8 mm	2300 mm	Short	5 mm	25G = 0.5 mm	
N6140150	NM-400U-0625	5	2.8 mm	2300 mm	Short	6 mm	25G = 0.5 mm	