

RB / WB

WB

∅ 3.5 mm  
★ RB (Regular Base)  
BLX



6 mm	–
8 mm	061.3308
10 mm	061.3310
12 mm	061.3312
14 mm	061.3314
16 mm	061.3316
18 mm	061.3318

∅ 3.75 mm  
★ RB (Regular Base)  
BLX



6 mm	061.4306
8 mm	061.4308
10 mm	061.4310
12 mm	061.4312
14 mm	061.4314
16 mm	061.4316
18 mm	061.4318

∅ 4.0 mm  
★ RB (Regular Base)  
BLX



6 mm	061.5306
8 mm	061.5308
10 mm	061.5310
12 mm	061.5312
14 mm	061.5314
16 mm	061.5316
18 mm	061.5318

∅ 4.5 mm  
★ RB (Regular Base)  
BLX



6 mm	061.6306
8 mm	061.6308
10 mm	061.6310
12 mm	061.6312
14 mm	061.6314
16 mm	061.6316
18 mm	061.6318

∅ 5.0 mm  
★ WB (Wide Base)  
BLX



6 mm	061.7306
8 mm	061.7308
10 mm	061.7310
12 mm	061.7312
14 mm	061.7314
16 mm	061.7316
18 mm	061.7318

∅ 5.5 mm  
★ WB (Wide Base)  
BLX



6 mm	061.8306
8 mm	061.8308
10 mm	061.8310
12 mm	061.8312
14 mm	–
16 mm	–
18 mm	–

∅ 6.5 mm  
★ WB (Wide Base)  
BLX



6 mm	061.9306
8 mm	061.9308
10 mm	061.9310
12 mm	061.9312
14 mm	–
16 mm	–
18 mm	–

**Roxolid® SLActive®**

**Roxolid® SLA®**

**Closure Cap**

RB (Regular Base)

WB (Wide Base)



H 0.4 mm  
064.41005



H 0.5 mm  
064.81025

**Healing Abutment**



AH = abutment height  
GH = gingiva height  
☒ = diameter

Healing Abutments for Crown

Healing Abutments for Crown

Healing Abutments for Bridge

Healing Abutments for Crown - Wide Profile

Healing Abutments for Crown - Wide Profile

GH	AH 2 mm		
	∅ 4.0 mm	∅ 5.0 mm	∅ 6.5 mm
1.5 mm	064.4202S	064.4212S	064.4222S
2.5 mm	064.4204S	064.4214S	064.4224S
3.5 mm	064.4206S	064.4216S	064.4226S

GH	AH 4 mm		
	∅ 4.0 mm	∅ 5.0 mm	∅ 6.5 mm
1.5 mm	064.4203S	064.4213S	064.4223S
2.5 mm	064.4205S	064.4215S	064.4225S
3.5 mm	064.4207S	064.4217S	064.4227S

GH	AH 2 mm	AH 4 mm
	∅ 5.0 mm	∅ 5.0 mm
1.5 mm	064.4232S	064.4233S

GH	AH 2 mm	
	∅ 6.0 mm	∅ 7.5 mm
0.75 mm	064.8201S	064.8203S
1.5 mm	064.8212S	064.8214S

GH	AH 4 mm	
	∅ 6.0 mm	∅ 7.5 mm
0.75 mm	064.8202S	064.8204S
1.5 mm	064.8213S	064.8215S

**Auxiliary Instruments**



## RB/WB prosthetics

Fit on RB and WB implants  
(Ø 3.5 mm - Ø 6.5 mm)



### Impression Posts; Implant Analogs



Open-tray 065.0031  
Closed-tray 065.0033  
Analog 065.4310  
Analog 065.0021

### Scanbodies and Repositionable Analogs



Scanbody 065.0035  
Analog 065.0023

### Temporary Abutments



For Crown, AH 11 mm (incl. screw)  
Ø 3.8 mm Ø 4.5 mm Ø 6.0 mm

	1.5 mm	064.4362	064.4372	—
GH	2.5 mm	064.4363	064.4373	064.4382
	3.5 mm	064.4364	064.4374	064.4383

For Bridge, AH 5.5 mm (incl. screw)  
Ø 4.5 mm

GH	1.5 mm	064.4352
----	--------	----------

## WB prosthetics

Fit only on WB implants  
(Ø 5.0 mm - Ø 6.5 mm)



### Impression Posts; Implant Analogs



Open-tray 065.0032  
Closed-tray 065.0034  
Analog 065.4810  
Analog 065.0022

### Scanbodies and Repositionable Analogs



Scanbody 065.0035  
Analog 065.0024


### Temporary Abutments



Temporary Abutments for Crown, AH 11 mm (incl. screw)  
Ø 5.5 mm


	0.75 mm	064.4391
GH	1.5 mm	064.4392

## Variobase®, Burn-out copings



Variobase® for Crown, AH 5.5 mm (incl. screw)

	Ø 3.8 mm	Ø 4.5 mm
	062.4934	062.4944
GH	2.5 mm	062.4935
	3.5 mm	062.4936
	062.4945	062.4946



Burn-out Copings for Crown

	Ø 3.8 mm	Ø 4.5 mm
	065.0014	065.0015

Variobase® Basal Screw

Basal Screw	Basal Screw AS
065.0036	065.0037



Variobase® for Crown AS, AH 5.5 mm (incl. screw)

	Ø 4.5 mm
GH	1.5 mm
	062.4972

Burn-out Copings for Crown AS

	Ø 4.5 mm
	065.0018



Variobase® for Bridge/Bar, AH 3.5 mm (incl. screw)

	Ø 4.5 mm
GH	1.5 mm
	062.4961

Burn-out Copings for Bridge/Bar

	Ø 4.5 mm
	065.0017

## Variobase®, Burn-out copings




Variobase® for Crown, AH 5.5 mm (incl. screw)

	Ø 5.5 mm
GH	0.75 mm
	062.4953
	1.5 mm
	062.4954

Variobase® for Crown AS, AH 5.5 mm (incl. screw)

	Ø 5.5 mm
GH	1.5 mm
	062.4971




Burn-out Copings for Crown

	Ø 5.5 mm
	065.0016

Burn-out Copings for Crown AS


	Ø 5.5 mm
	065.0019



Variobase® Basal Screw

Basal Screw	Basal Screw AS
065.0036	065.0037

## Multi-unit Abutments



Titanium Abutments (incl. screw)

	Ø 4.6 mm		
	0°	17°	30°
	062.4722S	—	—
	062.4723S	—	—
GH	1.5 mm	062.4724S	062.4733S
	2.5 mm	062.4724S	062.4743S
	3.5 mm	062.4725S	062.4744S
	4.5 mm	062.4735S	062.4745S
	5.5 mm	—	—

Titanium Copings (incl. screw)

	Ø 4.6 mm
	023.4747
Crown	023.4751
Bridge	023.4751




Protective cap

	H 4.5 mm	H 5.1 mm
	024.0020-04	024.4323-04

Temporary TAN Copings (incl. screw)

	Ø 4.6 mm
	024.0023
Crown	024.0024
Bridge	024.0024



Impression posts

	Ø 4.6 mm	
	Closed-tray	Open-tray
	025.0014	025.0012

Variobase® for Bridge/Bar (incl. screw)

	Ø 4.6 mm
	023.0028



Analog for Screw-retained Abutments

	Ø 4.6 mm
	023.4756


Scanbody

025.0001

Repositionable Analog

025.0008

## Novaloc



Novaloc® Abutments

	Ø 3.8 mm		
	Angle	0°	15°
	1.5 mm	062.4501	—
	2.5 mm	062.4502	062.4507
	3.5 mm	062.4503	062.4508
GH	4.5 mm	062.4504	062.4509
	5.5 mm	062.4505	062.4510
	6.5 mm	062.4506	062.4511
	7.5 mm	—	062.4512

Processing Package titanium 2010.601-STM

## Auxiliary Instruments



Screwdriver 15mm 046.400  
Screwdriver 21mm 046.401  
Screwdriver 27mm 046.402  
AS Screwdriver 15mm 046.786  
AS Screwdriver 21mm 046.787  
AS Screwdriver 27mm 046.788

Straumann®Dental Implant System  
**BONE LEVEL IMPLANTS**

**SC**

Ø 2.9 mm  
SC (Small CrossFit)

Bone Level  
Tapered (BLT)



**NC**

Ø 3.3 mm  
NC (Narrow CrossFit)

Bone Level (BL)      Bone Level Tapered (BLT)



**RC**

Ø 4.1 mm  
RC (Regular CrossFit)

Bone Level (BL)      Bone Level Tapered (BLT)



**RC**

Ø 4.8 mm  
RC (Regular CrossFit)

Bone Level (BL)      Bone Level Tapered (BLT)



**Roxolid® SLActive®**

8 mm	—
10 mm	021.0010
12 mm	021.0012
14 mm	021.0014
16 mm	—
18 mm	—

021.2308	8 mm	021.3308
021.2310	10 mm	021.3310
021.2312	12 mm	021.3312
021.2314	14 mm	021.3314
—	16 mm	021.3316
—	18 mm	021.3318

021.4308	8 mm	021.5308
021.4310	10 mm	021.5310
021.4312	12 mm	021.5312
021.4314	14 mm	021.5314
—	16 mm	021.5316
—	18 mm	021.5318

021.6308	8 mm	021.7308
021.6310	10 mm	021.7310
021.6312	12 mm	021.7312
021.6314	14 mm	021.7314
—	16 mm	021.7316
—	18 mm	021.7318

**Roxolid® SLA®**

8 mm	—
10 mm	021.0110
12 mm	021.0112
14 mm	021.0114
16 mm	—
18 mm	—

021.2508	8 mm	021.3508
021.2510	10 mm	021.3510
021.2512	12 mm	021.3512
021.2514	14 mm	021.3514
—	16 mm	021.3516
—	18 mm	021.3518

021.4508	8 mm	021.5508
021.4510	10 mm	021.5510
021.4512	12 mm	021.5512
021.4514	14 mm	021.5514
—	16 mm	021.5516
—	18 mm	021.5518

021.6508	8 mm	021.7508
021.6510	10 mm	021.7510
021.6512	12 mm	021.7512
021.6514	14 mm	021.7514
—	16 mm	021.7516
—	18 mm	021.7518

**Cover Screw**

**SC**



H 0.5 mm  
024.0006S

**NC**



H 0.0 mm      H 0.5 mm  
024.2100S      024.2105S

**RC**



H 0.0 mm      H 0.5 mm  
024.4100S      024.4105S

**Healing Abutments**

H	Conical, oval
2.0 mm	024.0007S
3.5 mm	024.0008S
5.0 mm	024.0009S
6.5 mm	024.0010S

H	Ø 3.6 mm	Ø 4.8 mm	Ø 3.3 mm
2.0 mm	024.2222S	024.2242S	—
3.5 mm	024.2224S	024.2244S	024.2234S
5.0 mm	024.2226S	024.2246S	024.2236S

H	Ø 5.0 mm
7.0 mm	024.2270S

H	Ø 5.0 mm	Ø 6.5 mm	Ø 4.4 mm	Ø 4.7 mm
2.0 mm	024.4222S	024.4242S	—	—
4.0 mm	024.4224S	024.4244S	024.4234S	—
6.0 mm	024.4226S	024.4246S	—	024.4236S

H	Ø 7.0 mm
7.0 mm	024.4270S

**Auxiliary Instruments**



Torque Control Device 046.049      Ratchet 046.119      Screwdriver 15mm 046.400      Screwdriver 21mm 046.401      Screwdriver 27mm 046.402      AS Screwdriver 15mm 046.786      AS Screwdriver 21mm 046.787      AS Screwdriver 27mm 046.788

BL Bone Profiler 1 026.0022      BL Bone Profiler 2 026.0023      BL Bone Profiler 3 026.0024      BL Bone Profiler SC 026.0070      NC Guiding Pin for Bone Profiler 026.0025S      RC Guiding Pin for Bone Profiler 026.0026S      SC Guiding Pin for Bone Profiler 026.0071S

# Straumann® Dental Implant System

## BONE LEVEL PROSTHETICS

### Straumann® Bone Level Implant CrossFit® Connection

#### Implant-level Impression

Open-tray Impression

	SC	NC	RC
Crown short	025.0021	025.2202	025.4202
Crown long	025.0022	025.2205	025.4205
Bridge short	—	025.0057	025.0060

Closed-tray Impression

	SC	NC	RC
	025.0020	025.2201	025.4201

Analog

	SC	NC	RC
	025.0023	025.2101	025.4101

Temporary Abutments

	SC	
1.0 mm	024.0011	
2.0 mm	024.0015	
3.0 mm	024.0016	
Incl. screw	025.0031	

Temporary Abutments

	NC	RC
Crown TAN	024.2371	024.4371
Bridge TAN	024.2375	024.4375
TAN/PMMA	024.2372	024.4372
TAN Incl. screw	025.2900	025.4900
PMMA Incl. screw	025.2908	025.4908

Repositionable Analog

	SC	NC	RC
	025.0024	025.2102	025.4102

CARES® Mono-Scanbody

	SC	NC	RC
	025.0025	025.2915	025.4915

#### Implant level Cementable / Variobase® Abutments

Variobase® for Crown

	SC		NC		RC	
AH	3.5 mm	3.5 mm	5.5 mm	3.5 mm	5.5 mm	
GH 1.0 mm	022.0038	025.2921	022.0027	025.4921	022.0026	
GH 2.0 mm	022.0039	022.0102	022.0106	022.0103	022.0107	
GH 3.0 mm	022.0040	022.0104	022.0108	022.0105	022.0109	
Incl. screw	025.0031	025.2900		025.4900		

Variobase® for bridge/bar Cylindrical

	NC	RC
AH	3.5 mm	
GH 1.0 mm	022.0110	022.0111
Incl. screw	025.2926	

Variobase® for Crown AS

	NC		RC	
AH	3.5 mm	5.5 mm	3.5 mm	5.5 mm
GH 1.0 mm	022.0084	022.0093	022.0087	022.0096
Incl. screw	025.0055			

Variobase® for Cerec

	NC	RC
AH	4.7 mm	
GH 1.0 mm	022.0043	022.0044
Incl. screw	025.2900	025.4900

Burn-out Copings for Variobase®

	SC		NC		RC	
AH	3.5 mm	3.5 mm	5.5 mm	3.5 mm	5.5 mm	
Crown	023.0011V4	023.2756V4	023.0018V4	023.4759V4	023.0017V4	

Burn-out Copings for Variobase® Cylindrical

	NC / RC
Bridge/bar	023.4759

Burn-out Copings for Variobase® for Crown AS

	NC	RC
	023.0025	023.0026

#### Hybrid

##### Novaloc®

Processing Package	Impression Cap	Analog
2010.601-NOV	2010.722-NOV	2010.721-NOV

#### Abutment level Cementable / Multi-unit Abutments

Titanium Abutments (incl. screw)

	NC		RC		NC		RC	
	0°				17°			
∅	4.6 mm				4.6 mm			
GH	1.5 mm	022.01285	022.01325	—	—	—	—	—
	2.5 mm	022.01295	022.01335	Type A	022.01365	022.01505	—	—
				Type B	022.01375	022.01515	—	—
	3.5 mm	022.01305	022.01345	Type A	022.01385	022.01525	022.01445	022.01585
				Type B	022.01395	022.01535	022.01455	022.01595
	4.5 mm	—	—	Type A	022.01405	022.01545	022.01465	022.01605
				Type B	022.01415	022.01555	022.01475	022.01615
	5.5 mm	—	—	Type A	022.01425	022.01565	022.01485	022.01625
				Type B	022.01435	022.01575	022.01495	022.01635

Variobase® Coping (incl. screw 023.4763)

	NC / RC
Cylindrical	023.0028

Analog

	NC / RC
	023.4756

Open-tray Impression

	NC / RC
Crown	025.2244
Bridge/bar	025.0012

Titanium Copings (incl. screw 023.4763)

	NC / RC
Bridge	023.4751
Crown	023.4747

Protective caps

	NC / RC
H 4.5 mm	024.0020-04
H 5.1 mm	024.4323-04

CARES® Scanbody

	NC / RC
	025.0001

Temporary TAN Copings (incl. screw 023.4763)

	NC / RC
Crown	024.0023
Bridge	024.0024

Burn-out Copings (incl. screw 023.4763)

	NC / RC
Crown	023.4748
Bridge	023.4758

Repositionable Analog

	NC / RC
	025.0008

Novaloc® Abutments

		NC		RC			
		0°					
0°	—	1.0 mm	022.0046	022.0052			
		2.0 mm	022.0047	022.0053			
		3.0 mm	022.0048	022.0054			
		4.0 mm	022.0049	022.0055			
		5.0 mm	022.0050	022.0056			
		6.0 mm	022.0051	022.0057			
		15°	Type A	2.0 mm	022.0062	022.0067	
				3.0 mm	022.0063	022.0068	
				4.0 mm	022.0064	022.0069	
				5.0 mm	022.0065	022.0070	
6.0 mm	022.0066			022.0071			
Type B	—			2.0 mm	022.0072	022.0077	
		3.0 mm	022.0073	022.0078			
		4.0 mm	022.0074	022.0079			
		5.0 mm	022.0075	022.0080			
		6.0 mm	022.0076	022.0081			

System Overview  
**TLX SURGERY**

Ø 3.75 mm 

Ø 4.5 mm 

Ø 5.5 mm 

Ø 6.5 mm 

Standard Plus (SP)  
RT Ø4.8 mm



Standard Plus (SP)  
RT Ø4.8 mm



Standard Plus (SP)  
WT Ø6.5 mm



Standard Plus (SP)  
WT Ø6.5 mm



**Roxolid® SLActive®**

6 mm	035.3106S
8 mm	035.3108S
10 mm	035.3110S
12 mm	035.3112S
14 mm	035.3114S
16 mm	035.3116S
18 mm	035.3118S

6 mm	035.3506S
8 mm	035.3508S
10 mm	035.3510S
12 mm	035.3512S
14 mm	035.3514S
16 mm	035.3516S
18 mm	035.3518S

6 mm	035.3706S
8 mm	035.3708S
10 mm	035.3710S
12 mm	035.3712S
14 mm	-
16 mm	-
18 mm	-

6 mm	035.3806S
8 mm	035.3808S
10 mm	035.3810S
12 mm	035.3812S
14 mm	-
16 mm	-
18 mm	-

**Roxolid® SLA®**

6 mm	035.2106S
8 mm	035.2108S
10 mm	035.2110S
12 mm	035.2112S

6 mm	035.2506S
8 mm	035.2508S
10 mm	035.2510S
12 mm	035.2512S

6 mm	035.2706S
8 mm	035.2708S
10 mm	035.2710S
12 mm	035.2712S

6 mm	035.2806S
8 mm	035.2808S
10 mm	035.2810S
12 mm	035.2812S

**Closure Caps**

RT

RT

WT

WT

RT Closure Cap		Ø 2.7 mm	Ø 5.5 mm
	H 0 mm	036.3200S	-
	H 1.5 mm	-	036.1201S

WT Closure Cap		Ø 2.7 mm
	H 0 mm	036.3200S

**Healing Caps**

RT Healing Cap		Ø 5.5 mm
	H 2.0 mm	036.1202S
	H 3.0 mm	036.1203S
	H 4.5 mm	036.1204S

WT Healing Cap		Ø 7.2 mm
	H 2.0 mm	036.2202S
	H 3.0 mm	036.2203S
	H 4.5 mm	036.2204S

**Important Information**

	Prosthetic platform	Implant Diameters	Implant platform
RT ••	Two dots  	Ø 3.75 mm Ø 4.5 mm	RT 

Components laser-marked RT (Regular, TorcFit™) and with two dots can be used on all TLX implants with the RT neck.

	Prosthetic platform	Implant Diameters	Implant platform
WT •••	Three dots  	Ø 5.5 mm Ø 6.5 mm	WT 

Components laser-marked WT (Wide, TorcFit™) and with three dots can be used on all TLX implants with the WT neck.

**Auxiliary Instruments**



Torque Control Device  
066.1100



Ratchet  
046.119



Surgical Handle  
066.4000



RT/RN Bone Profiler 28 mm  
044.304



WT/WN Bone Profiler 25 mm  
044.305



Guiding Pin for TLX Implant  
049.406S

# TLX PROSTHETICS


## Impression Posts Digital Impression Posts Implant Analogs Scanbodies

### RT - Regular TorcFit™ Ø 4.8 mm



Impression Post Open Tray	Impression Post Closed Tray	Implant Analog	Repositionable Implant Analog	CARES® Mono Scanbody
036.1000 036.1001	036.1002	036.0102	036.2100	036.3220

### WT - Wide TorcFit™ Ø 6.5 mm



Impression Post Open Tray	Impression Post Closed Tray	Implant Analog	Repositionable Implant Analog	CARES® Mono Scanbody
036.2000	036.2002	036.1100	036.2102	036.3220

## Temporary Restorations

Temporary Abutment, crown	Temporary Abutment, bridge	Replacement Screw
037.1000	037.1001	036.3110

Temporary Abutment, crown	Temporary Abutment, bridge	Replacement Screw
037.2000	037.2001	036.3110

## Variobase® Burn-out Copings

Variobase® for Crown	Variobase® for Crown AS	Variobase® for Bridge/Bar Cylindrical	Variobase® C
037.1201 037.1211	037.1203 037.1212	037.1204 037.1213	037.1205
Replacement Screw		Replacement Screw for Variobase® for Crown AS	
036.3110		036.3111	

Variobase® for Crown	Variobase® for Crown AS	Variobase® for Bridge/Bar Cylindrical	Variobase® C
037.2201 037.2211	037.2203 037.2212	037.2204 037.2213	037.2205
Replacement Screw		Replacement Screw for Variobase® for Crown AS	
036.3110		036.3111	

## Novaloc®

Novaloc Abutment	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm
0°	037.1500	037.1501	037.1502	037.1503	037.1504	037.1505
15°	-	037.1510	037.1511	037.1512	037.1513	037.1514

Processing Package titanium 2010.601-STM

Novaloc Abutment	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm
0°	037.2500	037.2501	037.2502	037.2503	037.2504	037.2505
15°	-	037.2510	037.2511	037.2512	037.2513	037.2514

Processing Package titanium 2010.601-STM


## Auxiliary Instruments

Screwdriver 15mm 046.400	Screwdriver 21mm 046.401	Screwdriver 27mm 046.402

AS Screwdriver 15mm 046.786	AS Screwdriver 21mm 046.787	AS Screwdriver 27mm 046.788

**Straumann®**  
Dental Implant System  
**TISSUE LEVEL**  
**PROSTHETICS**

RN  
Ø 4.8 mm



synOcta®

WN  
Ø 6.5 mm



synOcta®

**Impression Posts**  
**Impression Caps**  
**Implant Analogs**

Open-tray Impression      Closed-tray Impression      Analog

Implant Level

048.090 / 048.010      048.862      048.124

Open-tray Impression      Closed-tray Impression      Analog

Implant Level

048.091      048.870      048.171

**Digital Impression Posts**  
**Scanbodies and**  
**Repositionable Analogs**

048.129      048.168

048.172      048.169

**Temporary Restorations**

Crown 048.715      Bridge 048.716      Crown 048.664

Crown 048.717      Bridge 048.718      Crown 048.242

**Prosthetic Restoration**

**Cementable / Screw retained Prosthetics Parts**

**Cementable / Screw retained Prosthetics Parts**

**Titanium Abutments**

**Variobase® for Crown**  
AH 4.0 mm | AH 6.0 mm  
048.710 | 022.0022

**Variobase® for bridge/bar**  
**Cylindrical**  
AH 4.0 mm  
048.378

**Variobase® for Crown AS**  
AH 4.0 mm | AH 6.0 mm  
048.877 | 048.880

**Variobase® for Cerec**  
AH 7.7 mm  
022.0019

**Variobase® for Crown**  
AH 4.5 mm | AH 6.5 mm  
048.711 | 022.0023

**Variobase® for bridge/bar**  
**Cylindrical**  
AH 4.5 mm  
048.379

**Variobase® for Crown AS**  
AH 4.5 mm | AH 6.5 mm  
048.878 | 048.881

**Variobase® for Cerec**  
AH 8.2 mm  
022.0020

**Plastic Copings**

**Burn-out Coping for Variobase® for Crown**  
AH 3.5 mm | AH 5.5 mm  
048.268V4 | 023.0015V4

**Burn-out Coping for Variobase® for bridge/bar**  
**Cylindrical**  
AH 3.5 mm  
048.268V4

**Burn-out Coping for Variobase® for Crown AS**  
048.897

**Burn-out Coping for Variobase® for Crown**  
AH 3.5 mm | AH 5.5 mm  
048.269V4 | 023.0016V4

**Burn-out coping for Variobase® for bridge/bar**  
**Cylindrical**  
AH 3.5 mm  
048.269  
048.269V4

**Burn-out coping for Variobase® for Crown AS**  
048.898

**Hybrid Abutments**

Novaloc®

Novaloc®

Novaloc Abutment	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm
0°	048.812	048.813	048.814	048.815	048.816	048.817
15°	—	048.832	048.833	048.834	048.835	048.836

Processing Package titanium  
2010.601-STM



Novaloc Abutment	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm
0°	048.818	048.819	048.820	048.821	048.822	048.823
15°	—	048.837	048.838	048.839	048.840	048.841

Processing Package titanium  
2010.601-STM



**Straumann® Dental Implant System**

**System Overview**  
**TISSUE LEVEL IMPLANTS**

**Ø 3.3 mm**

**Ø 4.1 mm**

**Ø 4.8 mm**



**Roxolid® SLActive®**

4 mm	–
6 mm	–
8 mm	033.511S
10 mm	033.512S
12 mm	033.513S
14 mm	033.514S

4 mm	033.043S
6 mm	033.560S
8 mm	033.561S
10 mm	033.562S
12 mm	033.563S
14 mm	033.564S

4 mm	033.044S
6 mm	033.590S
8 mm	033.591S
10 mm	033.592S
12 mm	033.593S
14 mm	033.594S

4 mm	033.045S
6 mm	033.610S
8 mm	033.611S
10 mm	033.612S
12 mm	033.613S
14 mm	–

**Roxolid® SLA®**

6 mm	–
8 mm	043.260S
10 mm	043.261S
12 mm	043.262S
14 mm	043.263S

6 mm	043.161S
8 mm	043.162S
10 mm	043.163S
12 mm	043.164S
14 mm	043.165S

6 mm	043.066S
8 mm	043.067S
10 mm	043.068S
12 mm	043.069S
14 mm	043.070S

6 mm	043.655S
8 mm	043.656S
10 mm	043.657S
12 mm	043.658S
14 mm	–

**Ti SLA® Loxim™**

6 mm	–
8 mm	043.911S
10 mm	043.912S
12 mm	043.913S
14 mm	043.914S

6 mm	043.925S
8 mm	043.926S
10 mm	043.927S
12 mm	043.928S
14 mm	043.929S

6 mm	043.939S
8 mm	043.940S
10 mm	043.941S
12 mm	043.942S
14 mm	043.943S

6 mm	043.948S
8 mm	043.949S
10 mm	043.950S
12 mm	043.951S
14 mm	–

**Cover Screw**

**RN**

**RN**

**RN**

**WN**



**Healing Caps**

