


PremiCron®

Non-absorbable sutures

DESCRIPTION

PremiCron® is intended for use in soft tissue approximation of the wound edges to render possible an undisturbed wound healing.

PremiCron® is indicated for use in general surgery, cardiovascular and vascular indications and reconstructive and plastic surgery, when surgical practice requires the use of a non-absorbable suture material. PremiCron® is intended for adult patients regardless of gender who do not meet any contraindication.



PremiCron®

PremiCron®

Product overview

Polyester, braided, coated



Colour	Green or white
Origin	Synthetic
Sizes	USP 6/0 (metric 0.7) to USP 5 (metric 7)
Type of absorption	Non-absorbable
Sterilisation	Ethylene oxide / Gamma-irradiation

PremiCron®

Needle length mm	Thread length cm colour	USP (metric)						
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle								
	 75 90 green			C2026094		C2026091		
	 Pledgets 6x3mm Firm 75 green					C0026907		
	 Pledgets 6x3mm Firm 75 white					C0027010		
	 8 x 75 green/white					M0027707		
	  Pledgets 3x3mm Firm 8 x 75 green/white					M0027953		
	  Pledgets 6x3mm Firm 4 x 75 8 x 75 green/white					M0027905 M0027706		
	  Pledgets 6x3mm Soft 8 x 75 green/white					M0027525		
	 	45 75 green			4.2 p.d. C0026024	4.1 p.d. C0026025	C0026026	C0026007

C Box of 36 units	B Box of 24 units	G Box of 12 units	M Box of 6 units
-------------------	-------------------	-------------------	------------------