



CONTINENCE CARE & UROLOGY PRODUCTS

WITH PATIENTS & SPECIALISTS WE BRING SAFETY TO UROLOGY

All together, for the future



Safe solutions for patients with chronic and acute urinary disorders

Urinary disorders affect millions of people all over the World. There are no real age limitations, no specific genders. It affects children, teenagers, adults, seniors, men, women, disabled people, and otherwise healthy people that might be treated in hospital facilities or at home.

Continence care products are meant to deal with chronic disorders for home care patients looking for a better quality of life on a daily basis for an active and independent life. Urology products are a reliable and Safe partners for health care professionals to deal with acute disorders in hospitals.

For adults or chronic care patients, B. Braun provides a complete solution for safety products in terms of protocol and materials in compliance with urinary management international guidelines for health care and patient safety.



For a better understanding...

This catalogue has been designed to help identify our users.


They could be home care patients represented by the "big family" members according to the product, or the health care professionals, represented by the smaller characters here under.

In addition, a distinction between Home care and a Hospital solutions is made, choosing the small house logo for Home care & "H" logo for Hospital products.


We hope you will find our brochure instructive and enjoyable to read...

Home care solutions

The big family



Hospital solutions



Summary



URINARY CATHETERIZATION

04



URINARY COLLECTION & MEASUREMENT

37



BOWEL MANAGEMENT & URINARY INCONTINENCE

43



WINNING COMBINATIONS

48



ORDERING INFORMATION

50





URINARY CATHETERIZATION

Intermittent Catheterization

8 ACTREEN TECHNOLOGY

More than 15 years of experience

Intermittent Catheters

10 ACTREEN® MINI RANGE (CATH & SET)

Small things can make a big difference!

12 ACTREEN® HI LITE RANGE (CATH & SET)

To use Anywhere, Anytime!

14 ACTREEN SAFESET®

The innovative solution for Urinary Tract Infection prevention - (UTI)

Suprapubic Catheters

16 CYSTOFIX®

Geared toward safety solutions

18 CYSTOFIX® SG

Allows a supra pubic catheter placement with minimal risk

Percutaneous Nephrostomy Set

19 NEPHROFIX®

Percutaneous - Nephrostomy sets

Ureteral Stents and Catheters

20 URECATH®

Ureteral stents and catheters

Urinary Dry Catheters

21 NELATON AND TIEMANN CATHETERS

Single use dry catheters for intermittent urinary catheterization

Indwelling Transurethral Catheters

22 URIMED® CATH

Silicone Foley Catheters

Indwelling Catheter Maintenance

24 URO-TAINER®

Improving Patient Quality of Life

ACTREEN® TECHNOLOGY

More than 15 years of experience


1

A unique technology based on hydrophilic lubricant

Components: Water and glycerine based formula maintained throughout catheterization.


Patient's Benefits

- **No mess, no splash:** gel lubricant properties
- **Ready to use product:** no additional water required for preparation
- **No risk of drying out during insertion:** the level of lubrication stays the same during the insertion⁽¹⁾



With time, most catheters dry out within the first 10min. Actreen® remains at the same level of lubrication for up to 60 minutes⁽¹⁾.

Increased stickiness



Report N° 7.5AC059.2

No risk of residue:

When the catheter is immersed in synthetic urine, Actreen® does not present any visually detectable particules in suspension⁽²⁾.

No effect on sperm mobility

An in vitro comparative study performed by an external expert demonstrates that Actreen® has negligibly toxicological effect on sperm motility and motion⁽³⁾.



⁽¹⁾ Report N°7.5AC059.2

⁽²⁾ Report N°7.5.AA.034.3

⁽³⁾ Effect of short-term exposure to two hydrophilic-coated and one gel prelubricated urinary. Dr Auger et al. Paris, Hôpital Cochin Minerva Urologica e nefrologica, 2007;59(2):115-124

⁽⁴⁾ «Enquiry about user acceptability of Actreen® Hi-Lite Set». Data on file.

⁽⁵⁾ Report n° 7.5.AA.100.1


⁽⁶⁾ For safety and hygienic reasons, we advise you to use the product in the shortest possible time after opening

2 The catheter features

- PVC free or DEHP free
- New material for PVC free catheter: TPO (Thermo-Plastic Polyolefin)
- Non traumatic eyes > Smooth insertion⁽⁵⁾



Smooth
Drainage eyes⁽⁵⁾



3 B.Braun innovative packaging No Touch technique




The catheter has no direct contact
with hands and fingers

- > Minimises risks of catheter contamination during handling



A specific connector design

- > Providing an easy grip without touching the catheter⁽⁴⁾



Environment oriented packaging

Packaging of Actreen® is made from recyclable materials. The polyethylene plastic material turns into water vapour and carbon dioxide when incinerated.

Discard the catheter according to the current local regulation.



ACTREEN® MINI RANGE

Small things can make a big difference!

Mini catheters and sets for female urinary catheterization

Common benefits...


The Mini catheters are designed with the same concept as the other actreen products.

- Pre-lubricated with Hydrophilic lubricant, ready to use catheter, no water needed
- Safe catheter material, PVC free, DEHP free
- Non traumatic eyes⁽²⁾
- Straight forward use
- Easy to carry thanks to their low profile and lightweight



The optimal discretion for woman

The Actreen® Mini range is your perfect indoor and outdoor companion, thanks to its feminine low profile packaging and famous discreet pouch that comes with each box.



For who? For what?

Used by female patients for intermittent self urinary catheterization.

For more information

www.actreen.com



Scan to watch instructional videos.



Brochure

⁽¹⁾ Report n° 7.5.AC.067.2
⁽²⁾ Report n° 7.5.AC.020.1

⁽³⁾ Report n° 7.5.AB.043.2
⁽⁴⁾ Report n° 7.5.AB.042.2




Actreen® Mini Cath to stay




Easy to open

- Clearly identifiable
- Easy open packaging⁽¹⁾
- Direct access to the catheter



Easy to use

- Universal connector: easy to grip and handle during use
- Hygienic function: the catheter's size avoids contact with the toilets



A feminine design

- Small and discreet packaging
- A feminine design for more discretion

Actreen® Mini Set to go



All in one lightweight product

- Catheter + bag sealed together and folded in a low profile shape



Anti-reflux valve

- Integrated in the bag to avoid leakage during storage⁽³⁾




Graduated bag of 700 mL

- Approximate measurement of urine
- Easy to use and empty thanks to its tear off corner⁽⁴⁾

Ordering information

Product description	CH	Ø (mm)	Units/box	Reference
ACTREEN® MINI CATH - 9 cm (catheter)				
	■ 10	3,3	30 + 1 discreet pouch	228010*
	□ 12	4,0	30 + 1 discreet pouch	228012*
	■ 14	4,7	30 + 1 discreet pouch	228014*
ACTREEN® MINI SET - 9 cm (catheter)				
	■ 10	3,3	30 + 1 discreet pouch	239010*
	□ 12	4,0	30 + 1 discreet pouch	239012*
	■ 14	4,7	30 + 1 discreet pouch	239014*

*For each country, a specific letter is added at the end of the reference.



ACTREEN® HI-LITE RANGE



To use Anywhere, Anytime!

Catheter and Sets for intermittent urinary catheterization

Common benefits...

The Hi Lite range is designed with the common features below:

- Safe catheter material, PVC free, DEHP free
- Non traumatic eyes for smooth bladder drainage⁽¹⁾
- Pre-lubricated with hydrophilic lubricant, ready to use catheter, no need for water
- Opening design for good and reduced dexterity
- Safe handling with two «No touch» options
- Foldable product for more discretion**



2 options to avoid touching the catheter during use.



Movable no touch



Accordion no touch



The invisible solution all day long

Its foldable shape allows to carry the product anywhere in the discreet pouch provided with each box**.



For who? For what?

Used by male, female and pediatric patients for intermittent self urinary catheterization.

For more information

www.actreen.com



Scan to watch instructional video.



Actreen® Hi-Lite cath



Actreen® Hi-Lite set

** The CH18 and the 20 cm version can not be folded and no discreet pouch is included inside these boxes.

⁽¹⁾ Report n° 7.5.AA.100.1

⁽²⁾ Report n°7.5.AB.054.1

⁽³⁾ CE MARK in progress - Product available soon



Actreen® Hi-Lite Cath to stay



Designed to be open by everyone

- 2 options for opening adapted to the level of patient's hand dexterity
- Direct access to the catheter



Unique connector shape

- Allows to easily handle the catheter during use

Actreen® Hi-Lite Set to go



All in one lightweight

- Catheter + bag sealed together and folded in a compact size**
- Fits in a pocket



Anti-reflux valve

- Integrated into the bag to avoid leakages during storage⁽²⁾



Graduated bag of 1000 mL

- Approximate measurement of urine
- Easy to read
- Bag with a friendly design
- Easy to empty thanks to its tear off corner⁽²⁾

Ordering information

Actreen® Hi-Lite Cath

	CH	ø (mm)	Units/box	Ref.
TIEMANN - 41 cm				
	08	2,7	30 + 1 discreet pouch	238108*
	10	3,3	30 + 1 discreet pouch	238110*
	12	4,0	30 + 1 discreet pouch	238112*
	14	4,7	30 + 1 discreet pouch	238114*
	16	5,3	30 + 1 discreet pouch	238116*
NELATON - 41 cm				
	08	2,7	30 + 1 discreet pouch	238208*
	10	3,3	30 + 1 discreet pouch	238210*
	12	4,0	30 + 1 discreet pouch	238212*
	14	4,7	30 + 1 discreet pouch	238214*
	16	5,3	30 + 1 discreet pouch	238216*
NELATON - 20 cm				
	06	2,0	30 + 1 discreet pouch	238306*
	08	2,7	30 + 1 discreet pouch	238308*
	10	3,3	30 + 1 discreet pouch	238310*
	12	4,0	30 + 1 discreet pouch	238312*
	14	4,7	30 + 1 discreet pouch	238314*
	16	5,3	30 + 1 discreet pouch	238316*

Actreen® Hi-Lite Set

	CH	ø (mm)	Units/box	Ref.
TIEMANN - 37 cm				
	08	2,7	30 + 1 discreet pouch	242108 ⁽³⁾
	10	3,3	30 + 1 discreet pouch	242110*
	12	4,0	30 + 1 discreet pouch	242112*
	14	4,7	30 + 1 discreet pouch	242114*
	16	5,3	30 + 1 discreet pouch	242116*
	18	6,0	30	242118*
NELATON - 37 cm				
	08	2,7	30 + 1 discreet pouch	242208 ⁽³⁾
	10	3,3	30 + 1 discreet pouch	242210*
	12	4,0	30 + 1 discreet pouch	242212*
	14	4,7	30 + 1 discreet pouch	242214*
	16	5,3	30 + 1 discreet pouch	242216*
	18	6,0	30	242218*

* For each country, a specific letter is added at the end of the reference radical.
⁽³⁾ CE MARK in progress - Product available soon

ACTREEN SAFESet®

The innovative solution for Urinary Tract Infection prevention - (UTI)

Closed system concept⁽¹⁾⁽²⁾⁽³⁾ – all in one with integrated needle free sample port for single catheterization at hospital

Benefits for:

- **Patient:**
 - Slipperiness maintained during catheter insertion⁽⁴⁾
 - Prevents risk of contracting UTI^(1,2,3)
- **Personnel care:**
 - Reduced risk of injury thanks to the needle free urine sample port
 - Time saving thanks to the all in one concept
- **Hygienists:**
 - Reliable partner for UTI reduction^(1,2,3)
- **Pharmacist/Purchaser:**
 - Number of rationalized products and orders
 - Simplified urinary catheterization procedure and reduced waste ⁽¹⁾



Focus on main handling steps



Product Opening: guaranteed catheter sterility thanks to the «No touch» technique.



Desinfect the sample port before taking urine sampling, Sterility is guaranteed inside the pouch.



The urine can be taken with different urine sampling devices. Ready for urinalysis.



For who? For what?

Used by nurses in hospitals for single intermittent urinary catheterisation.

Winning combination



Examination and protective gloves. Vasco® Nitril long gloves.

For more information



Scan to watch instructional videos.



Brochure



Poster

⁽¹⁾ Infections urinaires nosocomiales de l'adulte. Conférence de Consensus SPILF & AFU, Novembre 2002, page 11

⁽²⁾ ANAES (National Agency for health evaluation & accreditation) « Qualité de la pose et de la surveillance des sondes urinaires, décembre 1999»

⁽³⁾ SFHH Recommendations (French Hospital Hygienists society) - «Surveiller et prévenir les infections associées aux soins, septembre 2010»

⁽⁴⁾ Report N°7.5AC059.2

⁽⁵⁾ Report n°7.5.AB.045.1



Urinary catheterization



All in one system

- Catheter + bag integrating an urine sample port, all sealed together: no possible detachment⁽⁵⁾



Ready to use

- Pre-lubricated catheter, no need of extra lubricant



Guaranteed sterility

- «No touch» system to reduce the risk of catheter contamination
- Sterile double packaging for safe handling in operating room



No touch

Urine collection and sampling



Graduated bag of 1500 mL

- Approximate measurement of urine
- Tear off to empty the bag



Anti-reflux valve

- To avoid urine backflow from the bag⁽⁵⁾
- Prevents the urine splash that could contaminate the environment



B. Braun patented

Integrated needle free sample port


- Suitable with all luer sampling devices
- Designed to maintain a closed system until urine sampling⁽⁵⁾

Ordering information


	CH	Ø (mm)	Units/box	Reference
ACTREEN SAFESET® NETATON - 47 cm				
	10	3,3	50	234210
	12	4,0	50	234212
	14	4,7	50	234214
	16	5,3	50	234216
	18	6,0	50	234218

Needle free sample port

Suitable with different luer devices.



Syringe - Vacuum systems



CYSTOFIX®

Geared toward safety solutions

Suprapubic puncture & exchange sets

User benefits

Puncture sets:

- **Ready to use** setting: pre assembled catheter placed into the cannula, to reduce risks of sharp injury and save time
- **Thin and sharp stainless steel cannula**
 - 3-wings handling grip for easy prehension
- **Wide range** of catheter designs and sizes
 - Silicone catheters with integrated balloon for an easy replacement
 - Polyurethane pigtail catheters
- **Open tip** catheters to facilitate replacement with a guide wire

Replacement sets:

- **Integrated balloon silicone catheters** to facilitate replacement
- **Wide choice** of sizes and sets to adapt to the patient's situation



Puncture sets: how to use



The catheter is already in the cannula. Keep the protection tube on the cannula.



Pull the catheter until the tip is right below the cannula bevel.

Remove the protection tube to begin the puncture.

Winning combinations



To be used with **Urimed® SP** or **Urimed® SP+**

Replacement sets and dilation sets



Integrated balloon for easy replacement.



Short tip for patient's comfort.



Range of catheter sizes.



Range of dilators' size.



Ordering information

Puncture sets



PUNCTION SETS WITH BALLOON FOR ALL USE. NEW RANGE CYSTOFIX® READY-TO-USE

Each set contains: one split sheath cannula, protection tube, catheter pre inserted into the cannula, stopper.

Catheter design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
Integral balloon straight tip Ready to Use	Silicone	12 cm	5,4 mm	■ 12	44 cm	10	4450200
		12 cm	5,4 mm	■ 14	44 cm	10	4450220
Balloon & J tip (classical balloon) Ready to Use	Silicone	12 cm	5,4 mm	■ 12	53 cm	10	4450210

PUNCTION SETS WITH PIGTAIL CATHETER FOR SHORT TERM USE. 7 EYES. NEW RANGE CYSTOFIX® READY-TO-USE

Each set contains: one split sheath cannula, protection tube, catheter pre inserted into the cannula, sliding clamp, urine bag according to model.

Catheter design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
Complete set with 2 L urine bag Ready to Use	Polyurethane	8 cm	3,4 mm	■ 10	65 cm	10	4450150
		12 cm	3,4 mm	■ 10	65 cm	10	4450160
		12 cm	5,4 mm	■ 15	65 cm	10	4450170
Eco set (no bag) Ready to Use	Polyurethane	8 cm	3,4 mm	■ 10	65 cm	10	4450100
		12 cm	3,4 mm	■ 10	65 cm	10	4450120
		12 cm	5,4 mm	■ 15	65 cm	10	4450130

PUNCTION SET WITH PIGTAIL CATHETER FOR PEDIATRIC AND NEONATES

Each set contains: one split sheath cannula, protection tube, catheter. Paediatric set: + urine bag & sliding clamp. Minipaed set: + soft mandril & connectors.


Catheter design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
Paediatric complete set with 2 L urine bag Ready to Use	Polyurethane	5 cm	3,4 mm	■ 10	50 cm	10	4450180
Minipaed set	Polyurethane	5 cm	2,2 mm	■ 5	50 cm	10	4441036

PUNCTION SET WITH PIGTAIL CATHETER FOR PER AND POST OPERATIVE USE. 28 EYELTS

Each set contains: one split sheath cannula, protection tube, catheter sliding clamp.


Catheter design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
TUR set	Polyurethane	12 cm	5,4 mm	■ 15	65 cm	10	4441354


Exchange sets and dilation sets



EXCHANGE SETS WITH INTEGRAL BALLOON CATHETER

Set design	Catheter design	Catheter material	CH	Catheter length	Units/box	Reference
Dilation and replacement with guidewire and dilator	Balloon straight tip	Silicone	■ 12	44 cm	5	4450512
			■ 14			4450514
			■ 16			4450516
Replacement with guidewire	Balloon straight tip	Silicone	■ 10	44 cm	5	4450010
			■ 12			4450012
			■ 14			4450014
			■ 16			4450016
			■ 18			4450018
Replacement without guidewire	Balloon straight tip	Silicone	■ 12	44 cm	5	4450712
			■ 14			4450714
			■ 16			4450716
			■ 18			4450718
			■ 20			4450720
Replacement single pack	Balloon straight tip	Silicone	■ 12	44 cm	1	4450812
			■ 14			4450814
			■ 16			4450816
			■ 18			4450818
			■ 20			4450820






CYSTOFIX® SG

Allows a suprapubic catheter placement with minimal risk
Suprapubic sets according to the seldinger technique

User benefits

- Ready to use under local anesthesia
- Minimal pressure required to insert the 18G puncture needle in comparison with a trocar
- Gentle dilation of the tissues thanks to the use of a needle followed by a dilator
- Withdrawal and replacement of the catheter facilitated by the integrated silicone balloon
- Seldinger technique is well known for the insertion of Central Venous Catheter and Percutaneous Nephrostomy Catheters



18G 12 cm needle

- Local anesthesia facilitated



Dilator with split sheath

- Designed for a gentle dilation of the tissues




Silicone catheter with integrated balloon

- Facilitated withdrawal and replacement thanks to the flat balloon

Ordering information

PUNCTION SETS WITH INTEGRAL BALLOON SILICONE CATHETERS FOR ALL USE. NEW RANGE CYSTOFIX® SG
Each set contains: one scalpel, 18G needle, guidewire, dilator, catheter, 2 syringes, spigot.

Catheter design	Catheter material	Dilator length	Dilator diameter	CH	Catheter length	Units/box	Reference
Integrated balloon straight tip	Silicone	15 cm	5.1 mm	■ 10	44 cm	5	4450410
		15 cm	5.9 mm	□ 12	44 cm	5	4450412
		15 cm	6.6 mm	■ 14	44 cm	5	4450414
		15 cm	6.6 mm	■ 16	44 cm	5	4450416




NEPHROFIX®

Percutaneous Nephrostomy sets

User benefits

- Complete sets with dilators and working sheath for short and long term
- 22 cm long needle
- 1, 2 or 3 dilators for a gentle dilation of the tissues
- Choice of pigtail and balloon catheters



Pigtail tip catheters

- Certon® (polyurethane)
- Radiopaque and flexible
- Low risk of incrustation
- Choice of 3 sizes
- Large drainage eyes



Balloon catheter

- Silicone
- Short tip and short balloon
- 2 eyes above balloon and 1 eye below balloon

Ordering information

Complete sets

Product	GUIDE					PIGTAIL CATHETER			BALLOON CATHETER	Units/box
	Style	Tip	Ø (inch)	Ø (mm)	Length (cm)	CH8	CH11	CH14	CH12	
Nephrofix®	Seldinger	Straight	0,035"	0,9	80	4435087	4435117	4435141	-	5
Nephrofix® Certo	Schüller	J3	0,032"	0,8	90	4435583	4435613	4435648	4435745	5

All sets contain: one catheter, 3 pieces needle, guide wire, dilators, scalpel, syringe, ruler.

Exchange sets

Product	GUIDE					PIGTAIL CATHETER		BALLOON CATHETER	Units/box
	Style	Tip	Ø (inch)	Ø (mm)	Length (cm)	CH8	CH11	CH12	
Nephrofix® Nephrofix® Certo	Seldinger	Straight	0,035"	0,9	80	4436083	4436113	4436741	5

All exchange sets contain: one catheter, guide wire & dilator of the appropriate size.

Puncture needles: • Nephrofix®: 3 pieces 17G x 220 mm • Nephrofix® Certo: 3 pieces 18G x 220 mm – **Schüller guide wire:** rigid part 300 mm




URECATH®

Ureteral stents and catheters

User benefits

- Made of polyurethane (Certon®)
- X-ray visible
- Stents with smooth pigtail on bladder side designed to improve comfort



Ordering information

Stents

Stent Tip	CH	Length (mm)				*
		200	250	275	300	
OPEN / CLOSE STENT DESIGN						
Tiemann	5	4430921	4430930	4430948	4430956	⁽¹⁾
	6	4431057	4431073	4431081	4431090	⁽²⁾
	7,5	4431103	4431111	4431120	4431138	⁽²⁾
	9	-	4431316	-	-	⁽²⁾
Round	7,5	4431200	4431219	4431227	-	⁽²⁾

CH	Length (mm)			*
	250	275		
OPEN / OPEN STENT DESIGN				
6	4437497	4437500		⁽³⁾
7,5	4437535	4437543		⁽³⁾

CH	Length (mm)			
	200	250	275	*
CLOSE / CLOSE STENT DESIGN				
7,5	4431600	4431618	4431626	⁽³⁾

Catheters

Tip shape	CH	Length (cm)	Number of eyes	Eyes location			Reference	*
				Tip	Loop	Shaft		
URETERIC CATHETERS								
Straight	4	100	4	x			4432690	⁽⁴⁾
	6	100	4	x			4432703	⁽⁴⁾
	8	100	4	x			4432711	⁽⁴⁾
Loop	6	70	8		x	x	4432002	⁽⁵⁾
	7,5	70	8		x	x	4432010	⁽⁵⁾
	9	70	8		x	x	4432029	⁽⁵⁾
Loop	7	70	4		x		4432800	⁽⁵⁾

Packaging: Boxes of 10.

⁽¹⁾ Delivered with guide wire 900 mm x 0,6 mm (0,024"), pusher, cap, clamp and patient's card.

⁽²⁾ Delivered with guide wire 900 mm x 1,1 mm (0,043"), pusher, cap, clamp and patient's card.

⁽³⁾ Delivered with guide wire 1500 mm x 0,9 mm (0,035"), 1 side curled, pusher, cap, clamp and patient's card.

⁽⁴⁾ Delivered with introducer.

⁽⁵⁾ Delivered with introducer, cap and clamp.



UROLOGICAL CATHETERS

Single use dry catheters for intermittent urinary catheterization

User benefits

- Medical PVC material - Latex Free
- Transparent PVC – allowing urine visibility through it
- Wide range composed of Nelaton and Tiemann tips and 2 sizes for male (40 cm) and female (18 cm)
- Individually packaged – sterile
- Colour coded Charrière enabling immediate size identification




2 lateral eyelets



The universal colour-coded connectors provide an immediate identification of the sizes.

Ordering information

	Length (cm)	CH	Ø (mm)	Units/box	Reference
TIEMANN CATHETER					
	40	08	2,7	100	4500200
		10	3,3	100	4500201
		12	4,0	100	4500202
		14	4,7	100	4500203
		16	5,3	100	4500204
		18	6,0	100	4500205
		20	6,6	100	4500206
NELATON CATHETER					
	40	08	2,7	100	4500207
		10	3,3	100	4500208
		12	4,0	100	4500209
		14	4,7	100	4500210
		16	5,3	100	4500211
		18	6,0	100	4500212
		20	6,6	100	4500213
NELATON CATHETER					
	18	08	2,7	100	4500214
		10	3,3	100	4500215
		12	4,0	100	4500216
		14	4,7	100	4500217
		16	5,3	100	4500218
		18	6,0	100	4500219



URIMED® CATH

Silicone Foley Catheters

Indwelling balloon catheters


User benefits

Silicone catheters are hypoallergenic⁽¹⁾ and therefore preferable for long-term settings:

- Less tissue irritation
- Reduced risk of encrustation in long-term catheterized patients⁽¹⁾
- Less traumatic catheterization and less risk of blockage compared to latex thanks to thinner walls and larger lumen, for the same charrière size⁽¹⁾⁽²⁾
- Visible on X-Rays



Double sterile packaging 2 openings: lateral and longitudinal



Wrinkle the primary peel-pack.



Tear it off with the side pre-cuts...



or with the pre-cuts down the catheter.

Winning combination

- **Urine Measurement**
Ureofix® 500 Classic Et Compact



- **Urine Collection**
2 Way: Urimed® SP bags



- **3 Way Prostatic catheters:**
4 L bag



- **Bladder washout, irrigation and instillation**
Uro-Tainer®




For more information



Brochure



Adapted to every need

**2 Way Nelaton**

- 1st intention product
- Suitable for both Male & Female
- For routine bladder drainage
- Closed tip
- 2 drainage eyes

**2 Way Tiemann**


- Designed to ease catheterization in Male
- For routine bladder drainage
- Closed tip
- 1 drainage eye

**3 Way Dufour Prostatic Catheters**

- A 3rd channel enables irrigation after prostate surgery
- For routine bladder drainage and irrigation and/or instillation
- Open tip
- 3 drainage eyes + 1 bevel

Ordering information

	CH	Balloon capacity (mL)	Units/box	Reference
URIMED® CATH FOLEY NELATON				
	■ 12	10	10	4563250
	■ 14	10	10	4563251
	■ 16	10	10	4563252
	■ 18	10	10	4563253
	■ 20	10	10	4563254
	■ 22	10	10	4563255
	■ 24	10	10	4563256
URIMED® CATH FOLEY TIEMANN				
	■ 12	10	10	4563257
	■ 14	10	10	4563258
	■ 16	10	10	4563259
	■ 18	10	10	4563260
	■ 20	10	10	4563261
URIMED® CATH FOLEY DUFOUR				
	■ 16	40	10	4563262
	■ 18	50	10	4563263
	■ 20	60	10	4563264
	■ 22	70	10	4563265
	■ 24	80	10	4563266



⁽¹⁾ Geng V et al. Evidence-based Guidelines for Best Practice in Urological Health Care : Catheterisation - Indwelling catheters in adults - urethral and suprapubic. EAUN, 2012

⁽²⁾ Morris NS et al. Which indwelling urethral catheters resist encrustation by *Proteus mirabilis* biofilms? BJU 1997; 80:58-63

URO-TAINER®

Improving Patient Quality of Life

User benefits

Uro-Tainer® technology

- Patients with indwelling urinary catheters are at high risk of developing complications such as infection and encrustation
- There is a strong rationale for on-going catheter maintenance to prevent complications and maximize the life of the catheter together with patient comfort and quality of life
- For over 25 years the Uro-Tainer® technology and solutions have been used in clinical practice to manage encrustation and debris formation

Uro-Tainer® solutions

- **TWIN:**
Prevention & removal of blockage
- **Uro-Tainer® Polihexanide:**
Bacterial decolonisation
- **Uro-Tainer® NaCl:**
Bladder instillations



For who?

50% of catheterized patients affected due to encrustation and blockage of indwelling urinary catheters⁽¹⁾.

Winning combination

Urimed® Cath
Silicone Foley Catheters



For more information



Brochure

⁽¹⁾ Getliffe KA. The characteristics and management of patients with recurrent blockage of long-term urinary catheters. J Adv Nurs 1994; 20: 140-149



Uro-Tainer® Solutions

PVC
free



Prevention & removal of blockage



Bacterial decolonisation



Bladder instillations

Twin

• Suby G

Indication: A mildly hypotonic fluid that is less irritating as a result of the addition of magnesium. This fluid is specially designed to prevent phosphate crystallization and dissolve existing calcification in indwelling catheters.

• Solutio R

Indication: A mildly hypotonic fluid that is specifically designed for catheters with stubborn calcification where Suby G does not provide a sufficient result. This is due to its higher concentration of citric acid and gluconolactone.

In addition, this fluid minimises trauma when removing an indwelling catheter.

Uro-Tainer® Polihexanide

Indication: Uro-Tainer® Polihexanide 0.02% is used for routine decolonisation (bacteria removal) of the catheter.

Uro-Tainer® NaCl

• NaCl 0.9%

Indication: This isotonic fluid is used primarily for mechanically cleaning the bladder and catheter. eg. in the case of debris forming in the bladder.

• M NaCl 0.9% with injection port

Indication: The Uro-Tainer® M has an injection port that is specifically for administering drugs.

Products manufactured by:
B.Braun Medical AG
Seesatz 17
6204 Sempach
Switzerland

Ordering information

URO-TAINER® Twin


	Size	Units/box	Reference
URO-TAINER® TWIN SUBY G			
Standard	100 mL	10	FB99839
Twin	60 mL	10	9746609
URO-TAINER® TWIN SOLUTION R			
Standard	100mL	10	FB99841
Twin	60 mL	10	9746625

URO-TAINER® Polihexanide

	Size	Units/box	Reference
URO-TAINER® POLIHEXANIDE			
	100 mL	10	FB99965

URO-TAINER® NaCl

	Size	Units/box	Reference
URO-TAINER® NAACL 0.9%			
	100 mL	10	FB99833
URO-TAINER® M NAACL 0.9% with injection port			
	100 mL	10	FB99853





URINARY COLLECTION & MEASUREMENT

Urine Measurement Systems

- 28 **UREOFIX® 500 CLASSIC**
Accurate and safe urine output management
- 30 **UREOFIX® COMPACT**
Accurate and safe urine output management

Closed Urine Drainage Bags

- 32 **UREOFIX® 112 PLUS**
Closed urine drainage system with triple safety concept
- 33 **POST OPERATIVE 4 L BAG**
Closed urine drainage system with large capacity for post TUR use
- 34 **URIMED® SP**
Quality within everyone's reach

Adult Urine Bags

- 36 **URIMED® BAGS**
A sophisticated range of urine bags for daily use
- 37 **URIMED® ACCESSORIES**
A complementary range of accessories
- 38 **URIMED® B'BAGS**
Adult urine collection bag

Pediatric Urine Bags

- 40 **URINOCOL®**
Looking for adapted solutions
- 41 **PEDIATRIC URINE BAGS®**
Cost effective range for urine pediatric collection


UREOFIX® 500 CLASSIC

Accurate and safe urine output management

3 chambers Urine Measurement System

User benefits


- Precise measurement with 3 chambers
- For reduced risk of retrograde migration of bacteria thanks to the Pasteur chamber and 2 non return valves
- Safety for the user thanks to the needle free sample port
- User convenience with a wide choice of bags sizes and tube lengths
- Direct transfer of urine from chamber to bag with the drainage tab



User convenience



Needle Free sample port



Clear marking drainage tab



Cross drainage outlet with dog house



Technical aspects



Closed System

- Dry Pasteur chamber
- With non return valve



Accurate measurement

- Precise measuring chamber up to 50 mL
- 3 chambers for precision in all situation



Choice of fixation

- 2 straps for permanent fixation
- 1 built-in hook
- 1 pin for use of optional pin strap



Non return valve

- Between chamber and bag



Choice of bag


- 2 L bottom outlet bag
- 1.5 L & 3.5 L closed bag



Winning combination

- Urimed® Cath Silicone Foley Catheters

Ordering information



	Tubing length (cm)	Units/box	Reference
UREOFIX® 500 CLASSIC SETS WITH CLOSED BAGS			
Ureofix® 500 Classic with exchange bag 1,5 L	120	10	4417910
	170	10	4417920
Ureofix® 500 Classic with exchange bag 3,5 L	170	10	4417960
UREOFIX® 500 CLASSIC SETS WITH BOTTOM OUTLET BAGS			
Ureofix® 500 Classic with premounted b.o bag 2 L	120	10	4417930
	170	10	4417940
	200	10	4417950
UREOFIX® 500 CLASSIC SPARE BAGS			
	Volume		
Ureofix® 500 Classic closed bag	1.5 L	100	4417535
	3.5 L	100	4417543
Ureofix® 500 Classic Bottom Outlet Bag	2 L	100	4417551
UREOFIX® 500 CLASSIC ACCESSORIES			
Fixation strap		10x10	4417755


UREOFIX® COMPACT

Accurate and safe urine output management

2 chambers compact Urine Measurement System

User benefits

- Precise measurement with 2 separated chambers
- For reduced risk of urinary tract infections thanks to the Pasteur chamber and 2 non return valves
- Safety of use thanks to the needle free sample port
- Convenience of use with 3 practical fixation means and the compact size
- Quick transfer of urine from chamber to bag. (5" for 80 mL – Test report 7.10W019.1)



User convenience



Double hooks for quick and stable fixation



2 straps for permanent fixation



Pin strap for permanent fixation and quick repositioning



Needle free sample port to avoid sharp injury



Technical aspects



- Pasteur chamber with non return valve



- Non return valve between chamber and bag



- Large diameter tube between chamber and bag for quick urine transfer



- Accurate urine measurement. Precise and easy to read measuring scale in the critical range 0-50 mL




- Compact design to facilitate the placement on the bedframe and avoid contamination by the floor



Winning combination

- Urimed® Cath Silicone Foley Catheters



Ordering information

Production description	Tubing/Lenght (cm)	Units/Box	Reference
UREOFIX® COMPACT			
UREOFIX® Compact - 2 Chambers - 2 L bag	150	10	4417901
	120	10	4417902
Fixation strap	-	10x10	4417755



UREOFIX® 112 PLUS

Closed urine drainage system with triple safety concept

User benefits


Triple safety concept :

- **Infection prevention** thanks to the drip chamber and the two non return valves
- **Safety for the users** with the needle free sample port
- **Bottom outlet stored** in a vented protective cap to avoid contamination with the floor

- **Convenient fixation** means with 2 hooks and 1 cord



User convenience



Safe drip chamber

Designed to reduce the risk of retrograde bacteria migration.

Air-Vent for a pressure balance during drainage.

Double Safety System:

Anti-Reflux valves within the Drip Chamber and within the bag.

Needle free sample port

Latex-free with universal connector.

Eliminates the risk of needle stick injury during urine sampling.

Vented protective cap on the bag

For cross outlet fixation to avoid contamination.

Ordering information

	Units/Box	Reference
UREOFIX® 112 PLUS		
	25	28431S



+ Features:

- 2000 mL bag capacity – Sterile
- Graduated scale from 25 to 2000 mL
- Tube of 110 cm
- 2 hanging systems (rope hanger + universal double-hook for bed fixation)
- Cross outlet with stopper
- Latex free gloves included



POST OPERATIVE 4 L BAG

Closed urine drainage system with large capacity for post TUR use

User benefits

- 4 L capacity for use after TURP and TURB
- **Prevention of infection** thanks to the drip chamber and the 2 non return valves
- **Safety for the users** with the needle free sample port
- **Bottom outlet stored** in a vented protective cap to avoid contamination with the floor
- **Convenient fixation** means with 2 hooks and 1 cord



User convenience



Needle free sample port on the staged connector.
Eliminates the risk of needle stick injury during urine sampling.



Tube without and with pump

- Standard model without pump for use with high irrigation flow
- Model with pump to allow washouts without disconnecting the catheter



2 non return valves:
Anti-Reflux valves in the Drip Chamber and within the bag.



Winning combination

To be used with the 3 way silicone catheter Urimed® Cath Dufour.

Ordering information

Production description	Units/Box	Reference
4 L POST OPERATIVE URINE BAG		
4 L bag TUR*	20	4801020
4 L bag TUR with pump*	20	4801030

*Delivered with clamp, bed sheet clip in double pouch.



URIMED® SP

Quality within everyone's reach

Urine bags with needle free sample port.

User's benefits

- **Reduced risk of infections** thanks to the integrated anti-reflux valve
- **Safety of use** with the needle free sample port
- **Cross contamination** with the floor is limited thanks to the outlet storage on the bag
- **Non kinking tube** with large diameter for a better flow of urine
- **Universal connection**




Urimed® SP




Urimed® SP+

Focus on main handling steps



Clamp the tube in order to stop the urine flow



Desinfect the needle free sampling port located on the connector



Take the urine from the needle free sample port with a luer outlet device such as a syringe needle free

For whom? For what?

Used by nurses for urine collection and sampling during a short term period of time.

Winning combination



CYSTOFIX® Ready to use Suprapubic puncture sets.

For more information



Brochure



Technical characteristics




Clamp

- Used to stop urine flow during sampling



Needle Free Sample Port

- Eliminates the risk of needle sharp injury during urine sampling



Air filter

- Allows a complete filling of the bag
- Provides additional safety when the pressure is applied on the bag




Cross outlet

- Easy to manage, with only one hand
- Free flow, avoids urine stagnation and limits risk of infection



Outlet Storage

- The outlet can be stored to avoid contamination with the floor risk



Additional feature on URIMED SP+

Double Hanging System:

- Universal double-hook for bed fixation
- Rope hanger facilitating patient movement

Ordering information

	Length tube (cm)	Ø Tube (mm)	Units/box	Reference
URIMED® SP 2 L				
	100	9	10	28610
	100	9	50	28650
URIMED® SP+ 2 L				
	130	9	30	28630



URIMED® BAG

A sophisticated range of urine bags for daily use

Adult urine collection bag

Patient benefits

- A highly sophisticated product range for a large range of application
- Designed to be easy to use
- Safe usability thanks to the anti-reflux valve
- Secure drainage: outlets designed to be easy to handle
- Non kinking tube with large diameter for a better flow of urine
- Universal connection



Anti-reflux valve on all bags

Comfort

Main features

Night and Day bags (2 L)



Large capacity of urine collection (2 L)



Free flow cross outlet with stopper



A soft back foil on the bags to avoid skin irritation



Free flow lever tap designed to be easy and secure to drain with one hand

Leg bags (0,5 L - 0,8 L - 1,5 L)

Urimed® Tribag Plus 800 mL



Specific 3 chamber pouch for a comfortable leg fit.

Ordering information

Night and Day bags (2 L)

	Length tube	Ø Tube (mm)	Units/box	Reference
URIMED® BAG 2 L				
Non sterile	120 cm connected connector	9	10	28300A
MEDICARE® 2 L				
Sterile	90 cm	6.5	10	28115
Non sterile	90 cm	6.5	10	28114

Leg bags

	Length tube	Ø Tube (mm)	Units/box	Reference
URIMED® BAG PLUS 0,5 L				
Non sterile	60 cm cut to fit	9	10	28501A
URIMED® TRIBAG PLUS 0,8 L				
3 Chamber - Sterile	60 cm connected connector	9	10	28305
Non sterile	60 cm cut to fit	9	10	28306
URIMED® BAG PLUS 1.5 L				
Non sterile	90 cm cut to fit	9	10	28150A



URIMED® ACCESSORIES

A complementary range of accessories


Patient benefits

- A useful range of accessories to use in complement with the urine bags
- Practical and adapted to different applications



Urimed® Fix

- Can be worn on the thigh or lower leg



Urimed® Straps

- Adapted for leg bag



Urimed® Klett

- Fixation for tube and catheter



Urimed® Handle

- To hang at the bedside or wheelchair

Ordering information

Accessories

	Tube	Units/box	Reference
URIMED® FIX*			
	Small	2	68520R
	Medium	2	68530R
	Large	2	68540R
URIMED® STRAPS			
		1 pair	68550R

	Tube	Units/box	Reference
URIMED® KLETT			
	Small	10	28260
	Large	10	28261
URIMED® HANDLE			
		1	68548R

* Distributed by B. Braun

URIMED® B'BAGS

Adult urine collection bag

User benefits

- **Wide range** of urine bags:
 - Closed and Drainable bags
 - Large capacity: 2 L and 1,5 L
 - Sterile – Non Sterile
 - With or without Anti-reflux valve
- Latex free material
- **Graduated scale** on all bags
- **Patient data** can be filled in the multilingual table
- Universal connection



Types of bags

Closed bags



- For collection of urine and secretion

Drainable bags



- For urine collection

For who? For what?

Used by the nurses or by the patients for urine collection and secretion.

Winning combination



Urimed® Cath
Indwelling balloon catheters,
Silicone Foley Catheters

For more information



Brochure



Closed bags



Tear off
for emptying the bag.

Drainable bags



Free flow cross outlet
designed to be managed
with only one hand.



Outlet can be stored
to limit risk of
contamination.



x2



x4

Holes for universal hanging system



Anti-reflux system

Closed bags exist with and without anti-reflux valve and all drainable bags include an anti-reflux valve to avoid urine backflow during use.

Ordering information

Closed bags

		Capacity (L)	Units/box	Reference
URIMED® B'BAG - Tube length: 90 cm - Tube diameter intR : 5.05 Et extR 6.5				
Sterile bags	No ARV*	2	100	29400
		1,5	100	29420
Non sterile bags	No ARV*	2	250	29415
		1,5	250	29435
	With ARV*	2	250	29410
		1,5	250	29430
Non sterile bags	No ARV*	2	30 + 1 hanging accessory**	29416
Non sterile bags	With ARV*	2	30 + 1 hanging accessory**	29411

*ARV= Anti-reflux valve



**Hanging accessory

Some references exist with an hanging accessory included in the box.

Drainable bags

		Capacity (L)	Units/box	Reference
URIMED® B'BAG - Tube length: 90 cm - Tube diameter intR : 5.05 Et extR 6.5				
Sterile bags	With ARV*	2	100	29440
Non sterile bags	With ARV*	2	200	29450
		1,5	200	29455
Non sterile bags	With ARV*	2	30 + 1 hanging accessory**	29451




URINOCOL®

Looking for adapted solutions

Pediatric urine bags

Patient benefits

- **Wide sterile range** designed to be adapted to all babies including premature ones
- **Shape of adhesive adapted** to the boy and girl anatomy
- High skin **tolerance adhesive**: gentle on application and removal
- **Safe material**: PVC Free / Latex free
- **Permanent visual control of urine** thanks to the transparent material
- **Bags available** with or without drain and anti-reflux valve



Main features

Latex free

PVC free


For baby > 2500g

Designed to fit the baby anatomy



One **specific shape** for boys, one for girls

Suitable for different indications



1. **No drainage Closed bag**
Urine sampling for biochemical, cytological or bacteriological examination
2. **With drainage Opened bag**
collecting bag for diuresis management



Sampling with syringe

Connection to urine bag

For baby from 1200g to 2500g Premature

Depending on the baby weight – **exists in 2 sizes**: Standard et Mini



Parachute shape
Designed to maintain the bag in place more effectively.

Available with extension line: adapted for incubator.

Ordering information

URINOCOL PAEDIATRIC				Units/box	Reference	URINOCOL PREMATURE				
		Weight (g)	Units/box	Reference			Reference			
	Sterile open paediatric bags	Boy	50	227546*		Sterile premature bags	Standard	1200-2500	50	227571*
		Girl	50	227566*			Mini	<1200	50	227570*
	Sterile closed paediatric bags	Boy	100	227550*		Sterile premature bags (with extension line**)	Standard	1200-2500	50	227573*
		Girl	100	227560*			Mini	<1200	50	227572*

*For each country, a specific letter is added at the end of the reference.

**made of PVC: DEHP free




STANDARD PEDIATRIC URINE BAGS

Cost effective range for urine pediatric collection

User benefits

- Self adhesive **unisex shape**
- Size of the stick surface is **variable**
- **Graduated scale** at 10 mL up to 100 mL
- **Bags available** with or without drain and anti-reflux valve
- **Transparent material film**



- A unisex shape
- Self adhesive, the size of the stick surface is variable for boys and girls



- A graduated scale on the bag up to 100 mL to easily measure the urine through the film (approximate measurement)

Ordering information - To be confirmed

	Sterile / Non sterile	ARV* / Without ARV*	Closed / Drainable	Units/box	Reference
STANDARD PEDIATRIC URINE BAGS					
For Children	Non sterile	Without ARV*	Drainable - With drainage tube of 90 cm - With slide clamp	20 - Individually packed	4423046
	Non sterile	With ARV*	Drainable - With bottom outflow nozzle	10 - Bulk packed	4423135
	Sterile – Individually wrapped	Without ARV*	Closed	10 - Individually packed	4425030
	Non sterile	Without ARV*	Closed	10 - Bulk packed	4423038
For Neonates	Sterile – Individually wrapped	Without ARV*	Closed	10 - Individually packed	4425049

* ARV= Anti-reflux valve





BOWEL MANAGEMENT URINARY & INCONTINENCE

Colonic Lavage - Colonic Irrigation

44 **RECTOBAG®**

Enema and Diagnosis Closed Systems

Accessories

45 **RECTAL CATHETERS**

For enema purposes and intestine evacuation
or for the injection of contrast media

Male External Catheters

46 **URIMED® VISION**

Connect to freedom

Female Intravaginal Tampon

47 **DIVEEN®**

Feel free again

RECTOBAG®

Enema and Diagnosis Closed Systems

User benefits

- Single use closed enema kit always at hand
- 3 different kits available
- Hygienic handling: Rectal balloon catheter designed to be lodged in the rectum
- 2 L bag capacity with eyelets
- 100% silicone rectal balloon catheter
- Latex-Free



- 2 L lavage bag to be closed with liquid-tight cap



- Silicone soft balloon rectal catheter



- Special set for X-Ray examination procedures

Ordering information

	Product description	Units/box	Reference
RECTOBAG® FLUID COLLECTING SYSTEMS			
	Closed System for intestinal lavage <ul style="list-style-type: none"> • 2 L lavage bag • 150 cm large bore tubing with side clamp • Rectal tube with balloon • Omnifix® 50/60 mL Luer Lock syringe Sterile / single-packed	20	4423364
	Closed System for intestinal lavage and irrigation <ul style="list-style-type: none"> • 2 L lavage bag • 150 cm large bore tubing with side clamp • Y-piece with connector and side clamp • Rectal tube with balloon • Omnifix® 50/60 mL Luer Lock syringe Sterile / single-packed	20	4423348
	Closed System for lavage bag <ul style="list-style-type: none"> • 2 L lavage bag • Connector • 150 cm large bore tubing with side clamp Non sterile / bulk packed	40	4423356



RECTAL CATHETERS

For enema purposes and intestine evacuation
or for the injection of contrast media

User benefits

- Latex-free
- Individually packaged in box of 100
- Nelaton - 40 cm length
- Immediate size identification
- Colour coded Charrière enabling immediate size identification




- 2 large lateral eyelets for even instillation of liquids



- The universal colour-coded connectors provide an immediate size identification.

Ordering information



RECTAL CATHETERS	Length (cm)	CH	Units/box	Reference
	40	22	100	4500300
		25	100	4500301
		28	100	4500302
		30	100	4500303



URIMED® VISION

Connect to freedom


100% Silicone Male External Catheter

User benefits

Male External Catheters offer a discreet and reliable solution to male urinary incontinence.

The real benefit for patients is to continue having an active life with minimal inconvenience in a secured way.

- **Security**
 - Skin-friendly: excellent skin tolerance thanks to the silicone material
 - Silicone material enables high performance in water vapour permeability to avoid maceration
- **Autonomy**
 - Connection to all major types of collecting bags
 - Connected to urine bag, Urimed® leg & bed bags
- **Convenience**
 - One piece system



Kink proof silicone funnel

- Helps to ensure uninterrupted urine flow

Perfect Fit!

Choose the right connection!



24h
autonomy

Ordering information

	Ø Tube (mm)	Units/box	Reference
URIMED® VISION STANDARD - 7,5 cm length with 4,5 cm adhesive area			
	25	30	IH2525
	29	30	IH2529
	32	30	IH2532
	36	30	IH2536
	41	30	IH2541

	Ø Tube (mm)	Units/box	Reference
URIMED® VISION ULTRA - 7,5 cm length with full adhesive area			
	25	30	IH3525
	29	30	IH3529
	32	30	IH3532
	36	30	IH3536
	41	30	IH3541

	Ø Tube (mm)	Units/box	Reference
URIMED® VISION SHORT - 4,5 cm length with full adhesive area			
	25	30	IH4525
	29	30	IH4529
	32	30	IH4532
	36	30	IH4536
	41	30	IH4541

Description	Units/box	Reference
URINE BAGS		
Urimed® Bag Plus 500 mL	10	28501A
Urimed® Bag Plus 800 mL	10	28306A
Urimed® Bag Plus 1500 mL	10	28150A
Urimed® Bag 2000 mL	10	28300A





What is Diveen®?

Diveen® is a intravaginal device indicated for women suffering from stress or mixed urinary incontinence.


As simple as a tampon to apply!⁽¹⁾

B.Braun patent⁽¹⁾

DIVEEN®

Feel free again

For stress or mixed urinary incontinence⁽¹⁾



User benefits

- A unique mechanism enabling significant reduction of the risk of urine loss.⁽¹⁾
- **100% discreet** protection.
- Biocompatible material⁽¹⁾.
- **12 hours of incontinence control** to enjoy physical activity with total peace of mind.⁽²⁾
- **Easy to use.** No need to remove Diveen® during the day, even when going to the toilet.
- Can be inserted without applicator.

Innovative and patented concept⁽¹⁾

Diveen® has been designed to support the anterior vaginal wall in the area below the bladder neck. During an effort, the weight of the bladder puts pressure on the flexible ring of the device. This causes the rigid part to move and ensures closure of the urethra, thus significantly reducing the risk of urine loss.⁽²⁾

Clinical study

Clinically proven efficacy in the management of stress urinary incontinence in women.^(2a)

- Diveen® reduces 69% of urine loss on average.^(2a)

Ordering information

	Size	Units/box	Reference
	Small	5 Devices + 1 Applicator	3132
	Medium	5 Devices + 1 Applicator	3133
	Small	15 Devices + 1 Applicator	3022
	Medium	15 Devices + 1 Applicator	3023

⁽¹⁾ Patent EP 1 924 219 B1

⁽²⁾ Cornu J.N. et al. 75nc007 device for noninvasive stress urinary incontinence management in women : a randomized controlled trial (55 patients were controlled and analyzed - 26 controls and 29 treated). Total of 41 analyzed (24 in the control group and 17 in the treatment group). International Urogynecology Journal 2012, 23(12):1727-34.

^(2a) After 14 days, reduction of incontinence episodes frequency (primary parameter) of -68,8% ± 22% versus -8,25% ± 25,4% in the control group (p<0,0001) in per-protocol analysis.



WINNING COMBINATIONS



Suprapubic Catheters

50 **CYSTOFIX®**

- URIMED® SP 2 L
- or
- URIMED® SP+ 2 L
WITH BED
FIXATION

Vasco® Nitril long gloves - Sterile Hand protection
Softaskin®

Male External Catheters

50 **URIMED® VISION**

- URIMED® TRIBAG
- or
- URIMED® BAG
1.5 L AND 2 L

Urimed® Handle - Urimed® Straps

Indwelling and Prostatic Catheters

50 **URIMED® CATH**

- UREOFIX® 500 CLASSIC
- or
- POST OPERATIVE 4 L BAG
- or
- UREOFIX® COMPACT
- or
- UREOFIX® 112
- or
- URO-TAINER®

Softaskin®

CREATE YOUR SOLUTION

To make it concrete, B. Braun is proposing solutions called «The winning combinations». «Winning» in terms of time saving, efficiency in the results, safety for patients & health care professionals.

Suprapubic Catheters

CYSTOFIX®



Vasco® Nitril long gloves - Sterile
Hand protection

Softaskin®
Hands cleansing

URIMED® SP 2 L
Urine Collecting & Sampling

**URIMED® SP+ 2L
WITH BED FIXATION**
Urine Collection & Sampling



Male External Catheters

URIMED® VISION



Urimed® Handle

Urimed® Straps

URIMED® TRIBAG
Leg bag urine collection

**URIMED® BAG
1.5 L AND 2 L**
Urine collection



Indwelling Transurethral Catheters

URIMED® CATH



UREOFIX® 500 CLASSIC

Urine measurement in Intensive Care Unit



UREOFIX® COMPACT

Urine measurement in Intensive Care Unit



POST OPERATIVE 4 L BAG

TUR with Dufour type



UREOFIX® 112

Long-term catheterization



URO-TAINER®

Catheter maintenance & Bladder instillation



ORDERING INFORMATION





ACTREEN® MINI RANGE

Product description	CH	Ø (mm)	Units/box	Reference
ACTREEN® MINI CATH - 9 cm (catheter)				
	■ 10	3,3	30 + 1 discreet pouch	228010*
	□ 12	4,0	30 + 1 discreet pouch	228012*
	■ 14	4,7	30 + 1 discreet pouch	228014*
ACTREEN® MINI SET - 9 cm (catheter)				
	■ 10	3,3	30 + 1 discreet pouch	239010*
	□ 12	4,0	30 + 1 discreet pouch	239012*
	■ 14	4,7	30 + 1 discreet pouch	239014*

*For each country, a specific letter is added at the end of the reference.



ACTREEN® HI-LITE RANGE

Actreen® Hi-Lite Cath

	CH	Ø (mm)	Units/box	Ref.
TIEMANN - 41 cm				
	■ 08	2,7	30 + 1 discreet pouch	238108*
	■ 10	3,3	30 + 1 discreet pouch	238110*
	□ 12	4,0	30 + 1 discreet pouch	238112*
	■ 14	4,7	30 + 1 discreet pouch	238114*
	■ 16	5,3	30 + 1 discreet pouch	238116*
NELATON - 41 cm				
	■ 08	2,7	30 + 1 discreet pouch	238208*
	■ 10	3,3	30 + 1 discreet pouch	238210*
	□ 12	4,0	30 + 1 discreet pouch	238212*
	■ 14	4,7	30 + 1 discreet pouch	238214*
	■ 16	5,3	30 + 1 discreet pouch	238216*
NELATON - 20 cm				
	■ 06	2,0	30 + 1 discreet pouch	238306*
	■ 08	2,7	30 + 1 discreet pouch	238308*
	■ 10	3,3	30 + 1 discreet pouch	238310*
	□ 12	4,0	30 + 1 discreet pouch	238312*
	■ 14	4,7	30 + 1 discreet pouch	238314*
	■ 16	5,3	30 + 1 discreet pouch	238316*

Actreen® Hi-Lite Set

	CH	Ø (mm)	Units/box	Ref.
TIEMANN - 37 cm				
	■ 08	2,7	30 + 1 discreet pouch	242108 ⁽³⁾
	■ 10	3,3	30 + 1 discreet pouch	242110*
	□ 12	4,0	30 + 1 discreet pouch	242112*
	■ 14	4,7	30 + 1 discreet pouch	242114*
	■ 16	5,3	30 + 1 discreet pouch	242116*
	■ 18	6,0	30	242118*
NELATON - 37 cm				
	■ 08	2,7	30 + 1 discreet pouch	242208 ⁽³⁾
	■ 10	3,3	30 + 1 discreet pouch	242210*
	□ 12	4,0	30 + 1 discreet pouch	242212*
	■ 14	4,7	30 + 1 discreet pouch	242214*
	■ 16	5,3	30 + 1 discreet pouch	242216*
	■ 18	6,0	30	242218*

* For each country, a specific letter is added at the end of the reference radical.
⁽³⁾ CE MARK in progress - Product available soon



ACTREEN SAFESET®

	CH	Ø (mm)	Units/box	Reference
ACTREEN SAFESET® NETATON - 47 cm				
	■ 10	3,3	50	234210
	□ 12	4,0	50	234212
	■ 14	4,7	50	234214
	■ 16	5,3	50	234216
	■ 18	6,0	50	234218



CYSTOFIX®

Puncture sets



PUNCTION SETS WITH BALLOON FOR ALL USE. NEW RANGE CYSTOFIX® READY-TO-USE

Each set contains: one split sheath cannula, protection tube, catheter pre inserted into the cannula, stopper.

Catheter design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
Integral balloon straight tip Ready to Use	Silicone	12 cm	5,4 mm	■ 12	44 cm	10	4450200
		12 cm	5,4 mm	■ 14	44 cm	10	4450220
Balloon Et J tip (classical balloon) Ready to Use	Silicone	12 cm	5,4 mm	■ 12	53 cm	10	4450210

PUNCTION SETS WITH PIGTAIL CATHETER FOR SHORT TERM USE. 7 EYES. NEW RANGE CYSTOFIX® READY-TO-USE

Each set contains: split sheath cannula, protection tube, catheter pre inserted into the cannula, sliding clamp, urine bag according to model.

Catheter design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
Complete set with 2 L urine bag Ready to Use	Polyurethane	8 cm	3,4 mm	■ 10	65 cm	10	4450150
		12 cm	3,4 mm	■ 10	65 cm	10	4450160
		12 cm	5,4 mm	■ 15	65 cm	10	4450170
Eco set (no bag) Ready to Use	Polyurethane	8 cm	3,4 mm	■ 10	65 cm	10	4450100
		12 cm	3,4 mm	■ 10	65 cm	10	4450120
		12 cm	5,4 mm	■ 15	65 cm	10	4450130

PUNCTION SET WITH PIGTAIL CATHETER FOR PEDIATRIC AND NEONATES

Each set contains: split sheath cannula, protection tube, catheter. Paediatric set: + urine bag & sliding clamp. Minipaed set: + soft mandril & connectors.

Set design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
Paediatric complete set with 2 L urine bag Ready to Use	Polyurethane	5 cm	3,4 mm	■ 10	50 cm	10	4450180
Minipaed set	Polyurethane	5 cm	2,2 mm	■ 5	50 cm	10	4441036

PUNCTION SET WITH PIGTAIL CATHETER FOR PER AND POST OPERATIVE USE. 28 EYELETS

Each set contains: split sheath cannula, protection tube, catheter sliding clamp.

Set design	Catheter material	Cannula length	Cannula diameter	CH	Catheter length	Units/box	Reference
TUR set	Polyurethane	12 cm	5,4 mm	■ 15	65 cm	10	4441354

Exchange sets and dilation sets 

EXCHANGE SETS WITH INTEGRAL BALLOON CATHETER

Set design	Catheter design	Catheter material	CH	Catheter length	Units/box	Reference
Dilation and replacement with guidewire and dilator	Balloon straight tip	Silicone	■ 12	44 cm	5	4450512
			■ 14			4450514
			■ 16			4450516
Replacement with guidewire	Balloon straight tip	Silicone	■ 10	44 cm	5	4450010
			■ 12			4450012
			■ 14			4450014
			■ 16			4450016
Replacement without guidewire	Balloon straight tip	Silicone	■ 12	44 cm	5	4450712
			■ 14			4450714
			■ 16			4450716
			■ 18			4450718
			■ 20			4450720
Replacement single pack	Balloon straight tip	Silicone	■ 12	44 cm	1	4450812
			■ 14			4450814
			■ 16			4450816
			■ 18			4450818
			■ 20			4450820



CYSTOFIX® SG

PUNCTION SETS WITH INTEGRAL BALLOON SILICONE CATHETERS FOR ALL USE. NEW RANGE CYSTOFIX® SG

Each set contains: scalpel, 18G needle, guidewire, dilator, catheter, 2 syringes, spigot.

Catheter design	Catheter material	Dilator length	Dilator diameter	CH	Catheter length	Units/box	Reference
Integrated balloon straight tip	Silicone	15 cm	5.1 mm	■ 10	44 cm	5	4450410
		15 cm	5.9 mm	■ 12	44 cm	5	4450412
		15 cm	6.6 mm	■ 14	44 cm	5	4450414
		15 cm	6.6 mm	■ 16	44 cm	5	4450416

URINARY CATHETERIZATION

Indwelling and Prostatic Catheters



NEPHROFIX®

Complete sets

Product	GUIDE					PIGTAIL CATHETER			BALLOON CATHETER	Units/box
	Style	Tip	Ø (inch)	Ø (mm)	Length (cm)	CH8	CH11	CH14	CH12	
Nephrofix®	Seldinger	Straight	0,035"	0,9	80	4435087	4435117	4435141	-	5
Nephrofix® Certo	Schüller	J3	0,032"	0,8	90	4435583	4435613	4435648	4435745	5

All sets contain: one catheter, 3 pieces needle, guide wire, dilators, scalpel, syringe, ruler.

Exchange sets

Product	GUIDE					PIGTAIL CATHETER		BALLOON CATHETER	Units/box
	Style	Tip	Ø (inch)	Ø (mm)	Length (cm)	CH8	CH11	CH12	
Nephrofix®	Seldinger	Straight	0,035"	0,9	80	4436083	4436113	4436741	5
Nephrofix® Certo	Schüller	J3	0,032"	0,8	90				

All exchange sets contain: one catheter, guide wire & dilator of the appropriate size.

Puncture needles: • Nephrofix®: 3 pieces 17G x 220 mm • Nephrofix® Certo: 3 pieces 18G x 220 mm - **Schüller guide wire:** rigid part 300 mm



URECATH®

Stents

Stent Tip	CH	Length (mm)				
		200	250	275	300	*
OPEN / CLOSE STENT DESIGN						
Tieman	5	4430921	4430930	4430948	4430956	(1)
	6	4431057	4431073	4431081	4431090	(2)
	7,5	4431103	4431111	4431120	4431138	(2)
	9	-	4431316	-	-	(2)
Round	7,5	4431200	4431219	4431227	-	(2)

CH	Length (mm)		
	250	275	*
OPEN / OPEN STENT DESIGN			
6	4437497	4437500	(3)
7,5	4437535	4437543	(3)

CH	Length (mm)			
	200	250	275	*
CLOSE / CLOSE STENT DESIGN				
7,5	4431600	4431618	4431626	(3)

Catheters

Tip shape	CH	Length (cm)	Number of eyes	Eyes location			Reference	*
				Tip	Loop	Shaft		
URETERIC CATHETERS								
Straight	4	100	4	x			4432690	(4)
	6	100	4	x			4432703	(4)
	8	100	4	x			4432711	(4)
Loop	6	70	8		x	x	4432002	(5)
	7,5	70	8		x	x	4432010	(5)
	9	70	8		x	x	4432029	(5)
Loop	7	70	4		x		4432800	(5)

Packaging: Boxes of 10.

(1) Delivered with guide wire 900 mm x 0,6 mm (0,024"), pusher, cap, clamp and patient's card.

(2) Delivered with guide wire 900 mm x 1,1 mm (0,043"), pusher, cap, clamp and patient's card.

(3) Delivered with guide wire 1500 mm x 0,9 mm (0,035"), 1 side curled, pusher, cap, clamp and patient's card.

(4) Delivered with introducer.

(5) Delivered with introducer, cap and clamp.



UROLOGICAL CATHETERS

	Length (cm)	CH	Ø (mm)	Units/box	Reference
TIEMANN CATHETER					
	40	08	2,7	100	4500200
		10	3,3	100	4500201
		12	4,0	100	4500202
		14	4,7	100	4500203
		16	5,3	100	4500204
		18	6,0	100	4500205
		20	6,6	100	4500206
NELATON CATHETER					
	40	08	2,7	100	4500207
		10	3,3	100	4500208
		12	4,0	100	4500209
		14	4,7	100	4500210
		16	5,3	100	4500211
		18	6,0	100	4500212
		20	6,6	100	4500213
NELATON CATHETER					
	18	08	2,7	100	4500214
		10	3,3	100	4500215
		12	4,0	100	4500216
		14	4,7	100	4500217
		16	5,3	100	4500218
		18	6,0	100	4500219



URIMED® CATH

	CH	Balloon capacity (mL)	Units/box	Reference
URIMED® CATH FOLEY NELATON				
	12	10	10	4563250
	14	10	10	4563251
	16	10	10	4563252
	18	10	10	4563253
	20	10	10	4563254
	22	10	10	4563255
	24	10	10	4563256
URIMED® CATH FOLEY TIEMANN				
	12	10	10	4563257
	14	10	10	4563258
	16	10	10	4563259
	18	10	10	4563260
	20	10	10	4563261
URIMED® CATH FOLEY DUFOUR				
	16	40	10	4563262
	18	50	10	4563263
	20	60	10	4563264
	22	70	10	4563265
	24	80	10	4563266



URO-TAINER®

URO-TAINER® Twin

	Size	Units/box	Reference
URO-TAINER® TWIN SUBY G			
Standard	100 mL	10	FB99839
Twin	60 mL	10	9746609
URO-TAINER® TWIN SOLUTION R			
Standard	100 mL	10	FB99841
Twin	60 mL (Twin)	10	9746625

URO-TAINER® NaCl

	Size	Units/box	Reference
URO-TAINER® NaCl 0.9%			
	100 mL	10	FB99833
URO-TAINER® M NaCl 0.9% with injection port			
	100 mL	10	FB99853

URO-TAINER® Polihexanide

	Size	Units/box	Reference
URO-TAINER® POLIHEXANIDE			
	100 mL	10	FB99965



UREOFIX® 500 CLASSIC

	Tubing length (cm)	Units/box	Reference
UREOFIX® 500 CLASSIC SETS WITH CLOSED BAGS			
Ureofix® 500 Classic with exchange bag 1,5 L	120	10	4417910
	170	10	4417920
Ureofix® 500 Classic with exchange bag 3,5 L	170	10	4417960
UREOFIX® 500 CLASSIC SETS WITH BOTTOM OUTLET BAGS			
Ureofix® 500 Classic with premounted b.o bag 2 L	120	10	4417930
	170	10	4417940
	200	10	4417950
UREOFIX® 500 CLASSIC SPARE BAGS			
Ureofix® 500 Classic closed bag	Volume		
	1.5 L	100	4417535
	3.5 L	100	4417543
Ureofix® 500 Classic Bottom Outlet Bag	2 L	100	4417551
UREOFIX® 500 CLASSIC ACCESSORIES			
Fixation strap		10x10	4417755

URINARY COLLECTION & MEASUREMENT

Urine Measurement Systems



UREOFIX® COMPACT

Production description	Tubing/Lenght (cm)	Units/Box	Reference
UREOFIX® COMPACT			
Ureofix® Compact - 2 Chambers 2 L bag	150	10	4417901
	120	10	4417902
Fixation strap	-	10x10	4417755



URINARY COLLECTION & MEASUREMENT

Closed Urine Drainage Bags



UREOFIX® 112 PLUS

	Units/Box	Reference
UREOFIX® 112 PLUS		
	25	284315



POST OPERATIVE 4 L BAG

Production description	Units/Box	Reference
4 L POST OPERATIVE URINE BAG		
4 L bag TUR*	20	4801020
4 L bag TUR with pump*	20	4801030

* Delivered with clamp, bed sheet clip in double pouch



URIMED® SP

	Length tube (cm)	Ø Tube (mm)	Units/box	Reference
URIMED® SP 2 L				
	100	9	10	28610
	100	9	50	28650
URIMED® SP+ 2 L				
	130	9	30	28630

URINARY COLLECTION & MEASUREMENT

Adult Urine Bags



URIMED® BAG

 Night and Day bags (2 L)

	Length tube	Ø Tube (mm)	Units/box	Reference
URIMED® BAG 2 L				
Non sterile	120 cm connected connector	9	10	28300A
MEDICARE® 2 L				
Sterile	90 cm	6.5	10	28115
Non sterile	90 cm	6.5	10	28114

Leg bags

	Length tube	Ø Tube (mm)	Units/box	Reference
URIMED® BAG PLUS 0,5 L				
Non sterile	60 cm cut to fit	9	10	28501A
URIMED® TRIBAG PLUS 0,8 L				
3 Chamber - Sterile	60 cm connected connector	9	10	28305
Non sterile	60 cm cut to fit	9	10	28306
URIMED® BAG PLUS 1.5 L				
Non sterile	90 cm cut to fit	9	10	28150A



URIMED® ACCESSORIES

Accessories

Tube	Units/box	Reference
URIMED® FIX*		
Small	2	68520R
Medium	2	68530R
Large	2	68540R
URIMED® STRAPS		
	1 pair	68550R

* Distributed by B. Braun

Tube	Units/box	Reference
URIMED® KLETT		
Small	10	28260
Large	10	28261
URIMED® HANDLE		
	1	68548R



URIMED® B'BAGS

Closed bags

Capacity (L)	Units/box	Reference
URIMED® B'BAG - Tube length: 90 cm - Tube diameter intR : 5.05 & extR 6.5		
Sterile bags	No ARV*	2, 100, 29400
		1,5, 100, 29420
Non sterile bags	No ARV*	2, 250, 29415
		1,5, 250, 29435
	With ARV*	2, 250, 29410
		1,5, 250, 29430
Sterile bags	No ARV*	2, 30 + 1 hanging accessory**, 29416
Non sterile bags	With ARV*	2, 30 + 1 hanging accessory**, 29411

Drainable bags

Capacity (L)	Units/box	Reference
URIMED® B'BAG - Tube length: 90 cm - Tube diameter intR : 5.05 & extR 6.5		
Sterile bags	With ARV*	2, 100, 29440
		2, 200, 29450
Non sterile bags	With ARV*	1,5, 200, 29455
		2, 30 + 1 hanging accessory**, 29451

* ARV= Anti- Reflux valve
**made of PVC: DEHP free

URINARY COLLECTION & MEASUREMENT

Pediatric Urine Bags



URINOCOL®

Units/box	Reference
URINOCOL PAEDIATRIC	
	Boy, 50, 227546*
	Girl, 50, 227566*
	Boy, 100, 227550*
	Girl, 100, 227560*

*For each country, a specific letter is added at the end of the reference.

Weight (g)	Units/box	Reference
URINOCOL PREMATURE		
	Standard, 1200-2500, 50, 227571*	
	Mini, < 1200, 50, 227570*	
	Standard, 1200-2500, 50, 227573*	
	Mini, < 1200, 50, 227572*	

**made of PVC: DEHP free



STANDARD PEDIATRIC URINE BAGS To be confirmed

Sterile/ Non sterile	ARV*/ Without ARV*	Closed / Drainable	Units/box	Reference
STANDARD PEDIATRIC URINE BAGS				
For Children	Non sterile	Without ARV*	Drainable - With drainage tube of 90 cm - With slide clamp	20 - Individually packed, 4423046
	Non sterile	With ARV*	Drainable - With bottom outflow nozzle	10 - Bulk packed, 4423135
	Sterile - Individually wrapped	Without ARV*	Closed	10 - Individually packed, 4425030
	Non sterile	Without ARV*	Closed	10 - Bulk packed, 4423038
For Neonates	Sterile - Individually wrapped	Without ARV*	Closed	10 - Individually packed, 4425049

* ARV= Anti-reflux valve



RECTOBAG®

Product description		Units/box	Reference
RECTOBAG® FLUID COLLECTING SYSTEMS			
	Closed System for intestinal lavage <ul style="list-style-type: none"> • 2 L lavage bag • 150 cm large bore tubing with side clamp • Y-piece with connector and side clamp Sterile / single-packed	20	4423364
	Closed System for intestinal lavage and irrigation <ul style="list-style-type: none"> • 2 L lavage bag • Connector Sterile / single-packed	20	4423348
	Closed System for lavage bag Non sterile / bulk packed	40	4423356



RECTAL CATHETERS

Length (cm)	CH	Units/box	Reference
RECTAL CATHETERS			
40		100	4500300
		100	4500301
		100	4500302
		100	4500303



URIMED® VISION

Ø Tube (mm)	Units/box	Reference	Ø Tube (mm)	Units/box	Reference
URIMED® VISION STANDARD - 7,5 cm length with 4,5 cm adhesive area			URIMED® VISION SHORT - 4,5 cm length with full adhesive area		
25	30	IH2525	25	30	IH4525
29	30	IH2529	29	30	IH4529
32	30	IH2532	32	30	IH4532
36	30	IH2536	36	30	IH4536
41	30	IH2541	41	30	IH4541
URIMED® VISION ULTRA - 7,5 cm length with full adhesive area			URINE BAGS		
25	30	IH3525	Description Units/box Reference		
29	30	IH3529	Urimed® Bag Plus 500 mL 30 28501A		
32	30	IH3532	Urimed® Bag Plus 800 mL 10 28306A		
36	30	IH3536	Urimed® Bag Plus 1500 mL 30 28150A		
41	30	IH3541	Urimed® Bag 2000 mL 30 28300A		



DIVEEN®

Size	Units/box	Reference
DIVEEN®		
	5 Devices + 1 Applicator	3132
	5 Devices + 1 Applicator	3133
Small	15 Devices + 1 Applicator	3022
Medium	15 Devices + 1 Applicator	3023

This document, its contents, including institutional data, information, trademarks and logos mentioned here are the exclusive property of B. Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun.

Commercialized by:

B. Braun Medical | 204 avenue du Maréchal Juin | 92100 Boulogne Billancourt - France

Tél. 01 41 10 53 00 | Fax 01 41 10 53 99 | www.bbraun.fr

B. Braun Medical SAS with a share capital of 31 000 000 € | RCS Nanterre 562050856



For most countries a specific letter is added at the end of the radical reference.
All volumes are measured according to ISO 8670-2 standard.
This medical device is a regulated health care product and on the basis of this regulation is certified CE.
Non-binding documents and photographs.
The products commercialized by B. BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.