

Killik

Embedding medium for cryostat sectioning.

| COLOUR | PACKAGING | CODE |
|---------|-----------|----------|
| Neutral | 4x100 ml | 05-9801 |
| Pink | 4x100 ml | 05-9801R |
| Blue | 4x100 ml | 05-9801B |
| Green | 4x100 ml | 05-9801V |



Crio Clor 0,3

Chlorexydine-based (0.3%) spray disinfectant for cryostat, effective against bacteria, fungi and other viruses (hepatitis B, poliomyelitis, herpes simplex).

| PACKAGING | CODE |
|-----------|---------|
| 1x125 ml | 05-9802 |



94 dalis

Cryo spray

Freezing spray for histology purposes: quickly freezes the tissue to be sectioned by cryostat, as well as cooling down paraffin embedded specimens before microtome cutting.

| PACKAGING | CODE |
|-----------|----------|
| 12x150 ml | 08-SPRAY |



Microtome oil

Microtome lubricant.

| PACKAGING | CODE |
|-----------|---------|
| 1x100 ml | 08-1721 |
| 2x100 ml | 08-1720 |

