

## Pipette stands

When the pipette is not in use, it should be stored in an upright position in order to avoid contamination from desks and benches. Sartorius Biohit provides stands for all of its pipettes. It is recommended that electronic pipettes be stored and charged on a charging stand whenever they are not in use. In this way, the batteries always remain charged when work begins. Linear Stand is designed for all Sartorius Biohit mechanical and electronic pipettes, particularly for mLINE, Proline Plus and Proline mechanical devices. This stand is also compatible with a wide range of other pipette makes.

Compact carousel stands are ideal for saving bench space and are available for mechanical pipettes, and for electronic pipettes for charging.

### Pipette Stands - Ordering details

Cat. No.	Item
730981	Charging Stand for one pipette*
730991	Charging Carousel for 4 pipettes*
725620	Linear Stand for all Sartorius Biohit pipettes
725600	Carousel Stand for 6 mechanical pipettes
725610	Holder for one mLINE/Proline Plus pipette
721259	Holder for Proline pipette

\* Includes Universal AC-adaptor (with Euro, US/Jpn, UK and China plugs)

113.7 pirkimo dalis →



Charging Stand



Charging Carousel



Linear Stand



Carousel Stand



mLINE/Proline Plus Holder



Proline Holder

← 113.7 pirkimo dalis

