

Calibration & Maintenance Solutions



Storage set for electrodes

ROTH SELECTION.

For safe storage of electrodes with shaft Ø 12 mm so that they are ready for use at all times.

Set 1

Delivery incl. GL 45 screw cap and 5x seals with hole. Without screw top bottle and storage solution, please order separately.

Version	Art. No.	Pack Qty.
Set 1	KH88.1	1 unit(s)

Set 2

Delivery incl. GL 45 screw cap, 5x seals with hole, 100 ml screw top bottle and 250 ml storage solution.

Version	Art. No.	Pack Qty.
Set 2	KCN1.1	1 unit(s)



Regeneration solution ammonium bifluoride

Mettler-Toledo.

To regenerate the pH-sensitive glass membrane.

Warning H315-H319

Suitable for	Art. No.	Pack Qty.
pH single-rod measuring cells	HTE6.1	25 ml



Cleaning solution

Mettler-Toledo.

Pepsin/HCl

To remove protein contaminations.

Danger H314-EUH208

Suitable for	Art. No.	Pack Qty.
pH single-rod measuring cells	HTE4.1	250 ml

Thiourea/HCl

For cleaning diaphragms with silver-sulphide contaminations.

Warning H351-H361d-H412

Suitable for	Art. No.	Pack Qty.
pH single-rod measuring cells	HTE5.1	250 ml



Electrode holder with cable duct

ROTH SELECTION.

For three electrodes, opening Ø 16 mm. One additional opening with Ø 11 mm. Three joints enable a high degree of adjustment flexibility. Base plate: 110 x 130 mm. Base plate weight: 0.6 kg. Arm weight: 0.1 kg. Max. arm height: 32 cm.

Version	Art. No.	Pack Qty.
Electrode holder	EHY1.1	1 unit(s)

Adapter for pH single-rod measuring cells

Mettler-Toledo.

Suitable for	Art. No.	Pack Qty.
pH single-rod measuring cells with BNC connector and pH meter with DIN 19262 socket	1EP5.1	1 unit(s)
pH single-rod measuring cells with BNC + Cinch connector and pH meter with DIN 19262 + 4 mm banana plug socket	1EP6.1	1 unit(s)