

**PAK<sup>®</sup>**

Paclitaxel coated peripheral  
angioplasty balloon catheter

Combination of confidence and effect

To keep your patients walking



Drug  
Treatment

POLIGRADE<sup>®</sup>  
Coating technology

Introducing **PAK<sup>®</sup>** microcrystalline paclitaxel coated balloon



CE<sub>2409</sub>




## Technology:


- **POLIGRADE<sup>®</sup>** special coating technology for superior drug adherence and improved transfer to the tissue.
- Microcrystalline drug particles with decreased embolization risk.
- **POLIGRADE<sup>®</sup>** 3 µg/mm<sup>2</sup> paclitaxel dosage.
- Special insertion device for protection of the drug coating during delivery through the introducer.
- Excellent pushability and trackability of the balloon catheter.

mcPCB

Control




### Late Lumen Loss at 6 months




### Evidence based confidence from comprehensive pre-clinical evaluation and clinical, controlled randomized trial:

- Superior efficacy with long-term, low-dose tissue drug retention.
- Optimal healing profile.
- Efficacy in neointimal inhibition as represented by nearly 3-fold reduction in late lumen loss.
- More than two fold significant reduction of restenosis, repeated revascularization rates and serious adverse events.

### Serious Adverse Events



### Serious Adverse Events and Repeated Revascularization at 12 months




Catalogue number	Nominal balloon diameter [mm]	Nominal balloon length [mm]	Use with guide wire	Use with introducer	Catheter length [cm]
		C			
PAKOTW3xC25110	3	20, 25, 30, 35, 40, 45, 50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100	.025"	5F	110
PAKOTW4xC25110	4		.025"	5F	
PAKOTW4xC35110			.035"	6F	
PAKOTW5xC25110	5		.025"	6F	
PAKOTW5xC35110			.035"	6F	
PAKOTW6xC25110	6		.025"	6F	
PAKOTW6xC35110			.035"	6F	
PAKOTW7xC35110	7		.035"	7F	
PAKOTW8xC35110	8		.035"	7F	

NOTE: On request we provide catheters with different balloon lengths and diameters.

The standard catheter length is 110 cm, however 80 cm catheters can be also ordered.

To order, please replace 110 by 80 in catalogue number – for example: PAKOTW8x203580.

Paclitaxel coated peripheral angioplasty balloon catheters with rapid exchange system **PAK<sup>®</sup>rx** are also available – please contact us for more information.



**Caution:**

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information on herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time. This publication is only a commercial information and cannot be handled the commercial offer in the meaning of the art.66 §1 of Polish Civil Code.