



COMPATIBLE WITH ACCURUS®, CONSTELLATION®, ASSOCIATE® 2500, EVA

EASY CONNECTION TUBING

VITRECTOMY SYSTEM CONNECTION

	ECT
Ingredients	<ul style="list-style-type: none"> ▶ medical grade PVC ▶ polycarbonate + acetal
Form	sterile solid
Packaging	pouch
Sterilization	ethylene oxide according to EN ISO 11135-1
Storage temperature	room temperature
Expiry date	60 months
Intended use	vitrectomy system connection
Class	medical device
GMDN code	P 47240
CND code	Q 0299
MD registration number	10222/R
Medical device class	Is
CE mark	CE 0546

ECT is connection tubing compatible with the following vitrectomy systems: Accurus®, Constellation®, Associate® 2500, EVA.

Ref.: ECT 003-00

Contents:

- ▶ 1 ECT compatible with Accurus®, Constellation®, Associate® 2500, EVA
- ▶ 1 package insert

GS1-128 barcode

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

VAT IT 00063370282

Phone +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com