


ENDOLINE® LOOP

30

POLYPECTOMY SNARE



**INNER CABLE
AVAILABLE IN 0.67 MM AND
0.9 MM DIAMETER**

**The use of a larger diameter
inner cable
(0.9 mm) enables the distal part
to be heated to a higher temper-
ature*

The **ENDOLINE® LOOP** polypectomy snare enables the resection of sessile/flat or pedicled polyps (≥ 6 mm) in the digestive tract via an upper or lower endoscopy. This snare is connected to the diathermy energy generator via the insert plug on the handle's movable slider.

ADVANTAGES

- Three fingers grip for easy manipulation.
- Radiopaque distal part for perfect visibility.
- Product compatible for use with the main endoscopes on the market.
- Connector for connection to a diathermic energy generator.

AVAILABLE FORMS

LARGE OVAL

Snare with **flexible** or rigid cable
Ø 30 x 50 mm

SMALL OVAL

Snare with **flexible** or rigid cable
Ø 15 x 25 mm

MONOWIRE OVAL

Snare with **flexible**
cable
Ø 25 x 50 mm
or 17 x 45 mm

HEXAGONAL

Snare with **flexible** or rigid cable
Ø 30 x 50 mm

ASYMMETRIC

Snare with **flexible** cable
Width, 15 mm or 25 mm

ABOUT THIS MEDICAL DEVICE

1 st release date	2001
Single use	Yes
Class	IIb
Sterile	Yes
Sterilisation method	Ethylene oxide sterilisation
Contains latex	No
Lifetime	5 years
Certification	CE
Packaging	Box of 5 units in individual sterile wrapping
Storage	In its original box in a cool and dry place

PRINCE MEDICAL SAS
ZA La Sente Du Moulin
64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

PRINCE MEDICAL
DEVICES

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

* Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

ENDOLINE® LOOP

31

LARGE OVAL

Reference	Type of snare	Diameter of the inner cable (mm)	Sheath diameter (mm)	Sheath length (cm)	Opening (mm)
AS2865023	Flexible	0.67	2.35	230	30 x 50
AS3755023	Rigid	0.9	2.35	230	30 x 50

SMALL OVAL

AS2922523	Flexible	0.67	2.35	230	15 x 25
AS3762523	Rigid	0.9	2.35	230	15 x 25

MONOWIRE OVAL

AS1961723	Flexible	0.67	2.35	230	17 x 45
AS1972523	Flexible	0.67	2.35	230	25 x 50

HEXAGONAL

AS2985023	Flexible	0.67	2.35	230	30 x 50
AS3775023	Rigid	0.9	2.35	230	30 x 50

ASYMMETRIC

AS3265023	Flexible	0.67	2.35	230	25
AS3961523	Flexible	0.67	2.35	230	15

PRINCE MEDICAL SAS
 ZA La Sente Du Moulin
 64 rue Gutenberg - 60530 Ercuis - France
www.prince-medical.com
 Tel: +33 (0)3 44 26 54 76
prince-medical@omerin.com

The information provided on this data sheet is indicative and subject to change without notice. PRINCE MEDICAL SAS shall not be held responsible for any incidents resulting from inappropriate use, failure to comply with the rules and regulations in force, or failure to follow the instructions for use provided. Medical devices sold by PRINCE MEDICAL SAS should only be used by a duly qualified practitioner or authorised person.

® Registered trademark of PRINCE MEDICAL SAS. Drawings and photos are not contractually binding. Reproduction is prohibited without the prior consent of PRINCE MEDICAL SAS.

PRINCE MEDICAL
 DEVICES