


PRODUCT DATA SHEET

Novosyn® CHD

Mid-term absorbable braided and coated suture made of polyglactin 910 coated with Chlorhexidine Diacetate



Novosyn® CHD is a mid-term absorbable braided synthetic suture made of polyglactin 910 coated with Chlorhexidine Diacetate. Novosyn® CHD degrades by hydrolysis and provides a predictable and reliable absorption.

Advantages

Novosyn® CHD has the same physical properties than Novosyn but with added antimicrobial protection.

It has a zone of inhibition against the most common pathogens that causes SSI.

- It helps on the prevention of microbial contamination of the implanted thread.

Features

- 40-50 % Tensile Strength Retention: 21 days.
- Complete Mass Absorption: 56 up to 70 days.

Indication

- Gastrointestinal surgery

- Gynecology / Obstetrics
- Orthopaedics
- Urology
- Skin closure (intracutaneous, subcutaneous)
- Ligatures

Article Description	Article Number
NOVOSYN CHD VIO 5/0 (1) 45CM DR12 (M)DDP	C1068134
NOVOSYN CHD VIOL 4/0(1,5)70CM DR12(M)DDP	C1068436
NOVOSYN CHD VIOL 4/0(1,5)70CM DR18(M)DDP	C1068150
NOVOSYN CHD VIOL 3/0 (2) 70CM DR18(M)DDP	C1068151
NOVOSYN CHD VIOL 3/0 (2) 70CM DR22(M)DDP	C1068160
NOVOSYN CHD VIOL 2/0 (3) 70CM DR26(M)DDP	C1068165
NOVOSYN CHD VIOL 3/0 (2) 70CM DR32(M)DDP	C1068170
NOVOSYN CHD VIOL 2/0 (3) 70CM DR32(M)DDP	C1068171
NOVOSYN CHD VIOL 0 (3,5) 70CM DR32(M)DDP	C1068172
NOVOSYN CHD VIOL 1 (4) 70CM DR32 (M) DDP	C1068173
NOVOSYN CHD VIO 5/0(1)45CM HR13F(M)DDP	C1068019
NOVOSYN CHD VIOLET 5/0(1)70CM HR13(M)DDP	C1068007
NOVOSYN CHD VIOL 4/0(1,5)70CM HR13(M)DDP	C1068008
NOVOSYN CHD UNDY 5/0 (1) 70CM HR13(M)DDP	C1069007
NOVOSYN CHD VIOL 0(3,5)70CM HR13SS(M)DDP	C1068111
NOVOSYN CHD VIOLET 5/0(1)70CM HR17(M)DDP	C1068012
NOVOSYN CHD VIOL 4/0(1,5)70CM HR17(M)DDP	C1068013
NOVOSYN CHD VIOL 3/0 (2) 70CM HR17(M)DDP	C1068014
NOVOSYN CHD VIOL 2/0 (3) 70CM HR17(M)DDP	C1068015
NOVOSYN CHD UNDY 5/0 (1) 70CM HR17(M)DDP	C1069012
NOVOSYN CHD UNDY 4/0(1,5)70CM HR17(M)DDP	C1069013
NOVOSYN CHD UNDY 3/0 (2) 70CM HR17(M)DDP	C1069014
NOVOSYN CHD VIO4/0(1,5)4X45CM HR17TOMDDP	C1088713
NOVOSYN CHD VIOL 3/0(2)4X45CM HR17TOMDDP	C1088714
NOVOSYN CHD VIO 4/0(1,5)8X45CM HR17TOMTD	M1088185

Article Description	Article Number
NOVOSYN CHD VIO 4/0(1,5)8X70CM HR17TOMTD	M1088255
NOVOSYN CHD VIOL 0 (3,5)70CM HR17S(M)DDP	C1068116
NOVOSYN CHD VIOL 1 (4) 70CM HR17S (M)DDP	C1068118
NOVOSYN CHD VIOL 3/0 (2) 45CM HR22(M)DDP	C1088325
NOVOSYN CHD VIOL 4/0(1,5)70CM HR22(M)DDP	C1068029
NOVOSYN CHD VIOL 3/0 (2) 70CM HR22(M)DDP	C1068030
NOVOSYN CHD VIOL 2/0 (3) 70CM HR22(M)DDP	C1068031
NOVOSYN CHD UNDY 4/0(1,5)70CM HR22(M)DDP	C1069029
NOVOSYN CHD VIO4/0(1,5)4X45CM HR22TOMDDP	C1088729
NOVOSYN CHD VIOL 3/0(2)4X45CM HR22TOMDDP	C1088730
NOVOSYN CHD VIOL 2/0(3)4X45CM HR22TOMDDP	C1088731
NOVOSYN CHD VIO 4/0(1,5)8X45CM HR22TOMTD	M1088177
NOVOSYN CHD VIOL 3/0(2)8X45CM HR22TO(MTD	M1088178
NOVOSYN CHD VIO4/0(1,5)4X70CM HR22TOMDDP	C1088629
NOVOSYN CHD VIOL 3/0(2)4X70CM HR22TOMDDP	C1088630
NOVOSYN CHD VIOL 2/0(3)70CM HR22S(M)DDP	C1068036
NOVOSYN CHD VIOL 0 (3,5)70CM HR22S(M)DDP	C1068037
NOVOSYN CHD VIOLET 1(4)70CM HR22S(M) DDP	C1068038
NOVOSYN CHD VIOLET 5/0(1)70CM HR26(M)DDP	C1068039
NOVOSYN CHD VIOL 4/0(1,5)70CM HR26(M)DDP	C1068040
NOVOSYN CHD VIOL 3/0 (2) 70CM HR26(M)DDP	C1068041
NOVOSYN CHD VIOLET 2/0(3)70CM HR26(M)DDP	C1068042
NOVOSYN CHD VIOLET 0(3,5)70CM HR26(M)DDP	C1068043
NOVOSYN CHD VIOLET 1(4)70CM HR26 (M) DDP	C1068044
NOVOSYN CHD VIOL 4/0(1,5)90CM HR26(M)DDP	C1068540
NOVOSYN CHD VIOL 3/0 (2) 90CM HR26(M)DDP	C1068541
NOVOSYN CHD VIOL 2/0 (3) 90CM HR26(M)DDP	C1068542
NOVOSYN CHD VIOL 0 (3,5) 90CM HR26(M)DDP	C1068543
NOVOSYN CHD VIOL 1 (4) 90CM HR26 (M) DDP	C1068544
NOVOSYN CHD UNDY 3/0 (2) 70CM HR26(M)DDP	C1069041
NOVOSYN CHD UNDY 2/0 (3) 70CM HR26(M)DDP	C1069042
NOVOSYN CHD VIOL 3/0(2)4X45CM HR26TOMDDP	C1088741

Article Description	Article Number
NOVOSYN CHD VIOL 2/0(3)4X45CM HR26TOMDDP	C1088742
NOVOSYN CHD VIO 4/0(1,5)8X45CM HR26TOMTD	M1088170
NOVOSYN CHD VIOL 3/0(2)8X45CM HR26TO(MTD	M1088171
NOVOSYN CHD VIOL 2/0(3)8X45CM HR26TO(MTD	M1088172
NOVOSYN CHD VIOL 3/0(2)4X70CM HR26TOMDDP	B1088641
NOVOSYN CHD VIOL 2/0(3)4X70CM HR26TOMDDP	B1088642
NOVOSYN CHD VIOL 3/0(2)8X70CM HR26TO(MTD	M1088239
NOVOSYN CHD VIOL 2/0(3)8X70CM HR26TO(MTD	M1088241
NOVOSYN CHD UNDY 2/0(3)8X45CM HR26TO(MTD	M1089117
NOVOSYN CHD VIOL 3/0(2)90CM 2XHR26(M)DDP	C1068901
NOVOSYN CHD VIOL 2/0(3)90CM 2XHR26(M)DDP	C1068902
NOVOSYN CHD VIOL 3/0 (2)70CM HR26S(M)DDP	C1068089
NOVOSYN CHD VIOL 2/0 (3)70CM HR26S(M)DDP	C1068090
NOVOSYN CHD VIOL 0 (3,5)70CM HR26S(M)DDP	C1068091
NOVOSYN CHD VIOL 1 (4) 70CM HR26S (M)DDP	C1068092
NOVOSYN CHD VIOL 2 (5) 70CM HR26S (M)DDP	C1068093
NOVOSYN CHD VIOL 2/0 (3)90CM HR26S(M)DDP	C1068590
NOVOSYN CHD VIOL 0 (3,5)90CM HR26S(M)DDP	C1068591
NOVOSYN CHD VIOL 1 (4) 90CM HR26S(M)DDP	C1068592
NOVOSYN CHD VIOL 2 (5) 90CM HR26S (M)DDP	C1068613
NOVOSYN CHD UNDY 2/0 (3)70CM HR26S(M)DDP	C1069090
NOVOSYN CHD UNDY 0 (3,5)70CM HR26S(M)DDP	C1069091
NOVOSYN CHD VIO 2/0(3)4X45CM HR26STOMDDP	C1088790
NOVOSYN CHD VIO 0(3,5)4X45CM HR26STOMDDP	B1088791
NOVOSYN CHD VIOL 1(4)4X45CM HR26STO(MDDP	B1088792
NOVOSYN CHD VIOL 2(5)4X45CM HR26STO(MDDP	B1088793
NOVOSYN CHD VIO 0(3,5)4X70CM HR26STOMDDP	B1088691
NOVOSYN CHD VIOL 1(4)4X70CM HR26STO(MDDP	B1088692
NOVOSYN CHD VIOL 2(5)4X70CM HR26STO(MDDP	B1088693
NOVOSYN CHD VIOLET 1(4)70CM HR26SS(M)DDP	C1068100
NOVOSYN CHD VIOL 1 (4) 90CM HR26SS(M)DDP	C1068580
NOVOSYN CHD VIOLET 3/0(2)70CM HR30(M)DDP	C1068046

Article Description	Article Number
NOVOSYN CHD VIOLET 2/0(3)70CM HR30(M)DDP	C1068047
NOVOSYN CHD VIOLET 0(3,5)70CM HR30(M)DDP	C1068048
NOVOSYN CHD VIOLET 1(4)70CM HR30 (M) DDP	C1068049
NOVOSYN CHD VIOLET 2(5)70CM HR30 (M) DDP	C1068045
NOVOSYN CHD VIOL 2/0(3) 90CM HR30 (M)DDP	C1068447
NOVOSYN CHD VIOL 2/0(3)4X70CM HR30TOMDDP	B1088661
NOVOSYN CHD VIOL 0 (3,5)70CM HR30S(M)DDP	C1068548
NOVOSYN CHD VIOL 1 (4) 70CM HR30S (M)DDP	C1068549
NOVOSYN CHD VIOLET 1(4)90CM HR30SS(M)DDP	B1068646
NOVOSYN CHD VIOL 1(4)4X70CM HR30SSTOMDDP	B1088803
NOVOSYN CHD VIOLET 3/0(2)70CM HR37(M)DDP	C1068050
NOVOSYN CHD VIOLET 2/0(3)70CM HR37(M)DDP	C1068051
NOVOSYN CHD VIOLET 0(3,5)70CM HR37(M)DDP	C1068052
NOVOSYN CHD VIOLET 1(4)70CM HR37 (M) DDP	C1068053
NOVOSYN CHD VIOLET 2(5)70CM HR37 (M) DDP	C1068054
NOVOSYN CHD VIOL 3/0 (2) 90CM HR37(M)DDP	C1068550
NOVOSYN CHD VIOL 2/0 (3) 90CM HR37(M)DDP	C1068551
NOVOSYN CHD VIOLET 0(3,5)90CM HR37(M)DDP	B1068552
NOVOSYN CHD VIOLET 1 (4) 90CM HR37(M)DDP	B1068553
NOVOSYN CHD UNDY 2/0 (3) 70CM HR37(M)DDP	C1069051
NOVOSYN CHD VIOL 2/0(3)4X45CM HR37TOMDDP	C1088770
NOVOSYN CHD VIOL 2/0(3)8X70CM HR37TO(MTD	M1088219
NOVOSYN CHD VIOL 1 (4) 8X45CM HR37S(M)TD	M1088138
NOVOSYN CHD VIOL 3/0 (2)70CM HR37S(M)DDP	C1068094
NOVOSYN CHD VIOL 2/0 (3)70CM HR37S(M)DDP	C1068095
NOVOSYN CHD VIOL 0 (3,5)70CM HR37S(M)DDP	C1068096
NOVOSYN CHD VIOL 1 (4)70CM HR37S (M) DDP	C1068097
NOVOSYN CHD VIOL 2 (5) 70CM HR37S (M)DDP	C1068098
NOVOSYN CHD VIOL 3/0 (2)90CM HR37S(M)DDP	C1068955
NOVOSYN CHD VIOL 2/0 (3)90CM HR37S(M)DDP	C1068595
NOVOSYN CHD VIOL 0(3,5)90CM HR37S(M)DDP	B1068956
NOVOSYN CHD VIOLET 1(4)90CM HR37S (M)DDP	B1068597

Article Description	Article Number
NOVOSYN CHD VIOLET 2(5)90CM HR37S (M)DDP	B1068598
NOVOSYN CHD UNDY 3/0 (2)70CM HR37S(M)DDP	C1069099
NOVOSYN CHD UNDY 2/0 (3)70CM HR37S(M)DDP	C1069095
NOVOSYN CHD UNDY 0 (3,5)70CM HR37S(M)DDP	C1069096
NOVOSYN CHD UNDY 2/0(3)90CM HR37S (M)DDP	C1069595
NOVOSYN CHD UNDY 0 (3,5)90CM HR37S(M)DDP	B1069596
NOVOSYN CHD UNDYED 1(4)90CM HR37S (M)DDP	B1069597
NOVOSYN CHD VIO 3/0(2)4X45CM HR37STOMDDP	C1088794
NOVOSYN CHD VIO 2/0(3)4X45CM HR37STOMDDP	C1088795
NOVOSYN CHD VIO 0(3,5)4X45CM HR37STOMDDP	B1088796
NOVOSYN CHD VIOL 1(4)4X45CM HR37STO(MDDP	B1088797
NOVOSYN CHD VIOL 2(5)4X45CM HR37STO(MDDP	B1088798
NOVOSYN CHD VIOL 1(4)4X70CM HR37STOMDDP	B1088697
NOVOSYN CHD VIOL 2(5)4X70CM HR37STO(MDDP	B1088698
NOVOSYN CHD VIOL 1 (4) 70CM HR37SS(M)DDP	C1068123
NOVOSYN CHD VIOLET 1(4)90CM HR37SS(M)DDP	B1068577
NOVOSYN CHD VIOL 2/0 (3)70CM HR40S(M)DDP	C1068055
NOVOSYN CHD VIOL 0 (3,5)70CM HR40S(M)DDP	C1068056
NOVOSYN CHD VIOL 1 (4) 70CM HR40S (M)DDP	C1068057
NOVOSYN CHD VIOL 2 (5) 70CM HR40S (M)DDP	C1068058
NOVOSYN CHD VIOL 3/0(2)90CM HR40S(M)DDP	C1068466
NOVOSYN CHD VIOL 2/0(3)90CM HR40S(M)DDP	C1068467
NOVOSYN CHD VIO 0 (3,5)90CM HR40S(M)DDP	B1068556
NOVOSYN CHD VIOL 1 (4) 90CM HR40S (M)DDP	B1068557
NOVOSYN CHD VIOLET 2(5)90CM HR40S (M)DDP	B1068558
NOVOSYN CHD VIOLET 1(4)120CM HR40S(M)DDP	B1068555
NOVOSYN CHD UNDYED 1 (4)90CM HR40S(M)DDP	B1069557
NOVOSYN CHD UNDYED 2(5)90CM HR40S (M)DDP	B1069558
NOVOSYN CHD VIOL 1(4)8X45CM HR40STO(M)TD	M1088117
NOVOSYN CHD VIOL 2(5)8X45CM HR40STO(M)TD	M1088118
NOVOSYN CHD VIOL 1(4)4X70CM HR40STO(MDDP	B1088667
NOVOSYN CHD VIOL 2(5)4X70CM HR40STO(MDDP	B1088668

Article Description	Article Number
NOVOSYN CHD VIO 1(4)150CM HR40S LOOPMDDP	B1068350
NOVOSYN CHD VIOL 0 (3,5) 70CM HR43(M)DDP	C1068946
NOVOSYN CHD VIOL 1 (4) 70CM HR43 (M) DDP	C1068947
NOVOSYN CHD VIOL 2 (5) 70CM HR43 (M) DDP	C1068948
NOVOSYN CHD VIOLET 1 (4) 90CM HR43(M)DDP	B1068847
NOVOSYN CHD VIOLET 2 (5) 90CM HR43(M)DDP	B1068848
NOVOSYN CHD VIOL 2/0 (3) 70CM HR48(M)DDP	C1068060
NOVOSYN CHD VIOL 0 (3,5) 70CM HR48(M)DDP	C1068061
NOVOSYN CHD VIOL 1 (4) 70CM HR48 (M) DDP	C1068062
NOVOSYN CHD VIOL 2 (5) 70CM HR48 (M) DDP	C1068063
NOVOSYN CHD VIOL 2/0 (3) 90CM HR48(M)DDP	C1068560
NOVOSYN CHD VIOLET 0(3,5)90CM HR48(M)DDP	B1068561
NOVOSYN CHD VIOLET 1 (4) 90CM HR48(M)DDP	B1068562
NOVOSYN CHD VIOLET 2 (5) 90CM HR48(M)DDP	B1068563
NOVOSYN CHD UNDYED 0(3,5)90CM HR48(M)DDP	B1069561
NOVOSYN CHD UNDYED 1(4)90CM HR48 (M) DDP	B1069562
NOVOSYN CHD UNDYED 2(5)90CM HR48 (M) DDP	B1069563
NOVOSYN CHD VIOL 1(4)8X45CM HR48TO(M)TD	M1088103
NOVOSYN CHD VIOL 2(5)8X45CM HR48TO(M)TD	M1088104
NOVOSYN CHD VIO 1(4)150CM HR48 LOP(M)DDP	B1068020
NOVOSYN CHD VIO 0 (3,5) 70CM HR60 (M)DDP	C1068065
NOVOSYN CHD VIOL 1 (4) 70CM HR60 (M) DDP	C1068066
NOVOSYN CHD VIOL 2 (5) 70CM HR60 (M) DDP	C1068068
NOVOSYN CHD VIO 1(4)150CM HR60 LOOPMDDP	B1068380
NOVOSYN CHD VIO 1 (4) 90CM HR65 (M) DDP	B1068639
NOVOSYN CHD VIOLET 2(5) 90CM HR65(M) DDP	B1068640
NOVOSYN CHD VIOLET 1(4)70CM HR70 (M) DDP	B1068183
NOVOSYN CHD VIOLET 2/0(3)70CM HR76(M)DDP	B1068074
NOVOSYN CHD VIOLET 1(4)90CM HR76 (M)DDP	B1068579
NOVOSYN CHD VIO 3/0 (2) 70CM FR26 (M)DDP	C1068414
NOVOSYN CHD VIOL 2/0 (3) 70CM FR26(M)DDP	C1068415
NOVOSYN CHD VIOL 0 (3,5) 70CM FR26(M)DDP	C1068416

Article Description	Article Number
NOVOSYN CHD VIOL 1 (4) 70CM FR26 (M) DDP	C1068417
NOVOSYN CHD VIOL 2/0(3)90CM 2XFR26(M)DDP	C1068820
NOVOSYN CHD VIOL 0(3,5)90CM 2XFR26(M)DDP	C1068828
NOVOSYN CHD VIOL 2/0 (3) 70CM FR36(M)DDP	C1068425
NOVOSYN CHD VIOL 0 (3,5) 70CM FR36(M)DDP	C1068426
NOVOSYN CHD VIO 1 (4) 70CM FR36 (M) DDP	C1068519
NOVOSYN CHD VIOL 0 (3,5) 70CM FR40(M)DDP	C1068490
NOVOSYN CHD VIO 2/0 (3) 90CM GR51 (M)DDP	C1068185
NOVOSYN CHD VIOL 2/0(3)70CM JR30(M) DDP	C1068475
NOVOSYN CHD VIOL 0(3,5)70CM JR36SS(M)DDP	C1068408
NOVOSYN CHD VIOL 1 (4) 70CM JR40S (M)DDP	C1068478
NOVOSYN CHD VIOL 4/0(1,5)70CM DRC18(MDDP	C1068124
NOVOSYN CHD VIOL 3/0(2)70CM DRC18(M)DDP	C1068128
NOVOSYN CHD VIOL 3/0(2)70CM DRC22(M)DDP	C1068827
NOVOSYN CHD VIOL 2/0(3)70CM DRC26(M)DDP	C1068131
NOVOSYN CHD VIOL 0(3,5)70CM DRC32(M)DDP	C1068304
NOVOSYN CHD VIOL 4/0(1,5)70CM HRC17(MDDP	C1068313
NOVOSYN CHD VIOL 3/0 (2)70CM HRC17(M)DDP	C1068314
NOVOSYN CHD VIOL 0 (3,5)70CM HRC17(M)DDP	C1068316
NOVOSYN CHD VIOL 1 (4) 70CM HRC17 (M)DDP	C1068081
NOVOSYN CHD VIOL 4/0(1,5)70CM HRC22(MDDP	C1068325
NOVOSYN CHD VIOL 3/0 (2)70CM HRC22(M)DDP	C1068326
NOVOSYN CHD VIOL 2/0 (3)70CM HRC22(M)DDP	C1068331
NOVOSYN CHD VIOL 3/0 (2)70CM HRC26(M)DDP	C1068335
NOVOSYN CHD VIOL 2/0 (3)70CM HRC26(M)DDP	C1068336
NOVOSYN CHD VIOL 0 (3,5)70CM HRC26(M)DDP	C1068343
NOVOSYN CHD VIO 2/0(3)4X45CM HRC26TOMDDP	C1088036
NOVOSYN CHD VIOL 2/0(3)70CM HRC26S(M)DDP	C1068441
NOVOSYN CHD UNDY 1 (4) 70CM HRC26S(M)DDP	C1069444
NOVOSYN CHD VIO 2/0 (3) 70CM HRC30(M)DDP	C1068372
NOVOSYN CHD VIOL 0 (3,5)70CM HRC30(M)DDP	C1068371
NOVOSYN CHD VIO 1 (4) 70CM HRC30 (M) DDP	C1068374

Article Description	Article Number
NOVOSYN CHD VIOL 0 (3,5)70CM HRC37(M)DDP	C1068362
NOVOSYN CHD VIOL 3/0 (2)90CM HRC37(M)DDP	C1068139
NOVOSYN CHD VIOL 2/0 (3)90CM HRC37(M)DDP	C1068140
NOVOSYN CHD VIO 0(3,5)90CM HRC37(M)DDP	B1068141
NOVOSYN CHD VIOLET 1(4)90CM HRC37(M)DDP	B1068142
NOVOSYN CHD VIOLET 2(5)90CM HRC37 (M)DDP	B1068143
NOVOSYN CHD VIOL 2/0(3)70CM HRC37S(M)DDP	C1068842
NOVOSYN CHD VIOL 0(3,5)70CM HRC37S(M)DDP	C1068843
NOVOSYN CHD VIOL 1 (4) 70CM HRC37S(M)DDP	C1068844
NOVOSYN CHD VIO 2/0(3) 90CM HRC37S(M)DDP	C1068017
NOVOSYN CHD VIOL 0(3,5)90CM HRC37S(M)DDP	B1068895
NOVOSYN CHD VIOL 1 (4) 90CM HRC37S(M)DDP	C1068018
NOVOSYN CHD VIOL 2 (5) 90CM HRC37S(M)DDP	B1068897
NOVOSYN CHD VIOL 0 (3,5)70CM HRC40(M)DDP	C1068909
NOVOSYN CHD VIOL 1 (4) 70CM HRC40S(M)DDP	C1068960
NOVOSYN CHD VIOL 1 (4) 70CM HRC43S(M)DDP	C1068627
NOVOSYN CHD VIOL 2 (5) 70CM HRC43S(M)DDP	C1068628
NOVOSYN CHD VIOL 1 (4) 90CM HRC43S(M)DDP	B1068982
NOVOSYN CHD VIOLET 2(5)90CM HRC43S(M)DDP	B1068983
NOVOSYN CHD VIOL 0 (3,5)70CM HRC48(M)DDP	C1068963
NOVOSYN CHD VIOL 1 (4) 70CM HRC48 (M)DDP	C1068964
NOVOSYN CHD VIO 0 (3,5)90CM HRC48(M)DDP	B1068573
NOVOSYN CHD VIOLET 1(4) 90CM HRC48(M)DDP	B1068574
NOVOSYN CHD VIOLET 2(5)90CM HRC48 (M)DDP	B1068575
NOVOSYN CHD VIO 1(4)150CM HRC48 LOOPMDDP	B1068987
NOVOSYN CHD VIOL 1 (4) 70CM JRC30S(M)DDP	C1068476
NOVOSYN CHD VIOL 0 (3,5)70CM HRN35(M)DDP	C1068518
NOVOSYN CHD VIOL 0 (3,5)90CM HRN35(M)DDP	B1068618
NOVOSYN CHD VIOLET 1 (4)90CM HRN35(M)DDP	B1068619
NOVOSYN CHD VIOLET 1(4) 90CM HRN40(M)DDP	B1068637
NOVOSYN CHD VIOLET 1(4)90CM HRN50 (M)DDP	B1068664
NOVOSYN CHD VIOL 4/0(1,5)70CM DS12(M)DDP	C1068508

Article Description	Article Number
NOVOSYN CHD UNDY 5/0 (1) 45CM DS12(M)DDP	C1069207
NOVOSYN CHD UNDY 4/0(1,5)45CM DS12(M)DDP	C1069208
NOVOSYN CHD VIOL 5/0 (1) 45CM DS16(M)DDP	C1068212
NOVOSYN CHD VIOL 4/0(1,5)45CM DS16(M)DDP	C1068213
NOVOSYN CHD VIOL 5/0 (1) 70CM DS16(M)DDP	C1068512
NOVOSYN CHD VIOL 4/0(1,5)70CM DS16(M)DDP	C1068513
NOVOSYN CHD VIOL 3/0 (2) 70CM DS16(M)DDP	C1068514
NOVOSYN CHD UNDY 4/0(1,5)45CM DS16(M)DDP	C1069213
NOVOSYN CHD VIOL 5/0 (1)45CM DS19 (M)DDP	C1068219
NOVOSYN CHD VIOL 4/0(1,5)45CM DS19(M)DDP	C1068220
NOVOSYN CHD VIOL 3/0 (2) 45CM DS19(M)DDP	C1068221
NOVOSYN CHD VIOL 5/0 (1) 70CM DS19(M)DDP	C1068419
NOVOSYN CHD VIOL 4/0(1,5)70CM DS19(M)DDP	C1068420
NOVOSYN CHD VIOL 3/0 (2) 70CM DS19(M)DDP	C1068421
NOVOSYN CHD VIOL 2/0 (3) 70CM DS19(M)DDP	C1068422
NOVOSYN CHD UNDY 5/0 (1) 45CM DS19(M)DDP	C1069219
NOVOSYN CHD UNDY 4/0(1,5)45CM DS19(M)DDP	C1069220
NOVOSYN CHD UNDY 3/0 (2) 45CM DS19(M)DDP	C1069221
NOVOSYN CHD UNDY 5/0 (1) 70CM DS19(M)DDP	C1069419
NOVOSYN CHD UNDY 4/0(1,5)70CM DS19(M)DDP	C1069420
NOVOSYN CHD UNDY 3/0 (2) 70CM DS19(M)DDP	C1069421
NOVOSYN CHD VIOL 4/0(1,5)70CM DS24(M)DDP	C1068234
NOVOSYN CHD VIOL 3/0 (2) 70CM DS24(M)DDP	C1068235
NOVOSYN CHD VIOL 2/0 (3) 70CM DS24(M)DDP	C1068236
NOVOSYN CHD VIOL 0 (3,5) 70CM DS24(M)DDP	C1068237
NOVOSYN CHD UNDY 4/0(1,5)70CM DS24(M)DDP	C1069234
NOVOSYN CHD UNDY 3/0 (2) 70CM DS24(M)DDP	C1069235
NOVOSYN CHD UNDY 2/0 (3) 70CM DS24(M)DDP	C1069236
NOVOSYN CHD VIOL 3/0 (2) 70CM DS30(M)DDP	C1068246
NOVOSYN CHD VIOL 2/0 (3) 70CM DS30(M)DDP	C1068247
NOVOSYN CHD VIOL 0 (3,5) 70CM DS30(M)DDP	C1068248
NOVOSYN CHD VIOL 1 (4) 70CM DS30 (M) DDP	C1068249

Article Description	Article Number
NOVOSYN CHD UNDY 3/0 (2) 70CM DS30(M)DDP	C1069246
NOVOSYN CHD UNDY 2/0 (3) 70CM DS30(M)DDP	C1069247
NOVOSYN CHD UNDY 0 (3,5) 70CM DS30(M)DDP	C1069248
NOVOSYN CHD VIOL 4/0(1,5)45CM HS21(M)DDP	C1068225
NOVOSYN CHD VIOL 3/0 (2) 45CM HS23(M)DDP	C1068230
NOVOSYN CHD VIOL 3/0 (2) 70CM HS23(M)DDP	C1068430
NOVOSYN CHD VIOL 2/0 (3) 70CM HS23(M)DDP	C1068431
NOVOSYN CHD UNDY 3/0 (2) 70CM HS23(M)DDP	C1069430
NOVOSYN CHD UN 2/0(3)4X45CM HS23SSTOMDDP	C1089680
NOVOSYN CHD VIOL 3/0 (2) 70CM HS26(M)DDP	C1068241
NOVOSYN CHD VIOL 2/0 (3) 70CM HS26(M)DDP	C1068242
NOVOSYN CHD VIOL 0 (3,5) 70CM HS26(M)DDP	C1068243
NOVOSYN CHD VIOL 2/0 (3)70CM HS26S(M)DDP	C1068290
NOVOSYN CHD VIOL 0 (3,5)70CM HS26S(M)DDP	C1068291
NOVOSYN CHD VIO 1 (4) 70CM HS26S (M) DDP	C1068292
NOVOSYN CHD UNDY 2/0 (3)70CM HS26S(M)DDP	C1069290
NOVOSYN CHD VIOL 2/0 (3) 70CM HS30(M)DDP	C1068251
NOVOSYN CHD VIOL 0 (3,5) 70CM HS30(M)DDP	C1068252
NOVOSYN CHD VIOL 1 (4) 70CM HS30 (M) DDP	C1068253
NOVOSYN CHD VIOL 3/0 (2) 90CM HS37(M)DDP	C1068494
NOVOSYN CHD VIOL 2/0 (3) 90CM HS37(M)DDP	C1068495
NOVOSYN CHD VIOLET 0(3,5)90CM HS37(M)DDP	B1068496
NOVOSYN CHD VIOLET 1 (4) 90CM HS37(M)DDP	B1068497
NOVOSYN CHD VIOLET 2 (5) 90CM HS37(M)DDP	B1068498
NOVOSYN CHD VIOL 2/0 (3)70CM HS37S(M)DDP	C1068294
NOVOSYN CHD VIOL 0 (3,5)70CM HS37S(M)DDP	C1068296
NOVOSYN CHD VIOL 1 (4) 70CM HS37S (M)DDP	C1068297
NOVOSYN CHD VIOL 2 (5) 70CM HS37S (M)DDP	C1068298
NOVOSYN CHD VIOL 2/0 (3)90CM HS37S(M)DDP	C1068596
NOVOSYN CHD VIOL 0 (3,5)90CM HS37S(M)DDP	B1068582
NOVOSYN CHD VIOLET 1(4)90CM HS37S (M)DDP	B1068583
NOVOSYN CHD VIOLET 2(5)90CM HS37S (M)DDP	B1068599

Article Description	Article Number
NOVOSYN CHD VIO 2/0 150CMHS37S/HR37SMDDP	B1068278
NOVOSYN CHD UN 0(3,5)4X45CM HS37SSTOMDDP	B1089718
NOVOSYN CHD VIO 0 (3,5) 90CM HS40(M)DDP	B1068486
NOVOSYN CHD VIOLET 1(4)90CM HS40 (M) DDP	B1068487
NOVOSYN CHD VIOLET 2 (5) 90CM HS40(M)DDP	B1068488
NOVOSYN CHD VIOL 2/0 (3)70CM HS40S(M)DDP	C1068267
NOVOSYN CHD VIOL 0 (3,5)70CM HS40S(M)DDP	C1068268
NOVOSYN CHD VIOL 1 (4)70CM HS40S (M) DDP	C1068269
NOVOSYN CHD VIOL 2/0 (3)90CM HS40S(M)DDP	C1068585
NOVOSYN CHD VIOL 0(3,5)90CM HS40S(M)DDP	B1068586
NOVOSYN CHD VIOLET 1(4)90CM HS40S (M)DDP	B1068587
NOVOSYN CHD VIOLET 2(5)90CM HS40S (M)DDP	B1068588
NOVOSYN CHD VIOL 0 (3,5) 70CM HS48(M)DDP	C1068261
NOVOSYN CHD VIOL 1 (4) 70CM HS48 (M) DDP	C1068262
NOVOSYN CHD VIOLET 0(3,5)90CM HS48(M)DDP	B1068461
NOVOSYN CHD VIOLET 1(4)90CM HS48 (M) DDP	B1068462
NOVOSYN CHD VIOLET 2 (5) 90CM HS48(M)DDP	B1068463
NOVOSYN CHD UNDYED 1 (4) 90CM HS48(M)DDP	B1069363
NOVOSYN CHD UNDYED 2 (5) 90CM HS48(M)DDP	B1069364
NOVOSYN CHD VIOL 3/0 (2) 70CM GS60(M)DDP	C1068265
NOVOSYN CHD VIOL 2/0 (3) 70CM GS60(M)DDP	C1068266
NOVOSYN CHDUNDY 4/0 (1,5)70CM GS60(M)DDP	C1069280
NOVOSYN CHD UNDY 3/0 (2) 70CM GS60(M)DDP	C1069265
NOVOSYN CHD UNDY 2/0 (3) 70CM GS60(M)DDP	C1069267
NOVOSYN CHD UNDY 5/0(1)45CM DSMP11(M)DDP	C1069608
NOVOSYN CHD VIOL 4/0(1,5)45CM DSMP13MDDP	C1068624
NOVOSYN CHD UNDY 5/0(1)45CM DSMP13(M)DDP	C1069619
NOVOSYN CHD UNDY 4/0(1,5)45CM DSMP13MDDP	C1069620
NOVOSYN CHD UNDY 5/0(1)45CM DSMP16(M)DDP	C1069525
NOVOSYN CHD UNDY 4/0(1,5)45CM DSMP16MDDP	C1069524
NOVOSYN CHD UNDY 4/0(1,5)45CM DSMP19MDDP	C1069629
NOVOSYN CHD UNDY 3/0(2)45CM DSMP19(M)DDP	C1069630

Article Description	Article Number
NOVOSYN CHD UNDY 4/0(1,5)70CM DSMP19MDDP	C1069529
NOVOSYN CHD UNDY 3/0(2)70CM DSMP19(M)DDP	C1069530
NOVOSYN CHD UNDY 2/0(3)70CM DSMP19(M)DDP	C1069531
NOVOSYN CHD UNDY 3/0(2)45CM DSMP24(M)DDP	C1069656
NOVOSYN CHD UNDY 2/0(3)45CM DSMP24(M)DDP	C1069657
NOVOSYN CHD UNDY 3/0(2)70CM DSMP24(M)DDP	C1069556
NOVOSYN CHD UNDY 2/0(3)70CM DSMP24(M)DDP	C1069554
NOVOSYN CHD VIOLET 4/0 (1,5) 3X45C DDP	C1058624
NOVOSYN CHD VIOLET 3/0 (2) 3X45CM DDP	C1058625
NOVOSYN CHD VIOLET 2/0 (3) 3X45CM DDP	C1058626
NOVOSYN CHD VIOLET 0 (3,5) 3X45CM DDP	C1058627
NOVOSYN CHD VIOLET 1 (4) 3X45CM DDP	C1058629
NOVOSYN CHD VIOLET 2 (5) 3X45CM DDP	C1058630
NOVOSYN CHD VIOLET 4/0 (1,5) 6X45CM DDP	C1058644
NOVOSYN CHD VIOLET 3/0 (2) 6X45CM DDP	C1058645
NOVOSYN CHD VIOLET 2/0 (3) 6X45CM DDP	C1058646
NOVOSYN CHD VIOLET 0 (3,5) 6X45CM DDP	B1058647
NOVOSYN CHD VIOLET 1 (4) 6X45CM DDP	B1058649
NOVOSYN CHD VIOLET 2 (5) 6X45CM DDP	B1058650
NOVOSYN CHD VIOLET 4/0 (1,5) 12X45CM DDP	C1058664
NOVOSYN CHD VIOLET 3/0 (2) 12X45CM DDP	C1058665
NOVOSYN CHD VIOLET 2/0 (3) 12X45CM DDP	C1058666
NOVOSYN CHD VIOLET 0 (3,5) 12X45CM DDP	B1058667
NOVOSYN CHD VIOLET 1 (4) 12X45CM DDP	B1058669
NOVOSYN CHD VIOLET 2 (5) 12X45CM DDP	B1058670
NOVOSYN CHD VIOLET 3/0 (2) 2X70CM DDP	C1058405
NOVOSYN CHD VIOLET 2/0 (3) 2X70CM DDP	C1058406
NOVOSYN CHD VIOLET 0 (3,5) 2X70CM DDP	C1058407
NOVOSYN CHD VIOLET 1 (4) 2X70CM DDP	C1058409
NOVOSYN CHD VIOLET 2 (5) 2X70CM DDP	C1058410
NOVOSYN CHD VIOLET 4/0 (1,5) 5X70CM DDP	C1058444
NOVOSYN CHD VIOLET 3/0 (2) 5X70CM DDP	C1058445

Article Description	Article Number
NOVOSYN CHD VIOLET 2/0 (3) 5X70CM DDP	C1058446
NOVOSYN CHD VIOLET 0 (3,5) 5X70CM DDP	C1058447
NOVOSYN CHD VIOLET 1 (4) 5X70CM DDP	B1058449
NOVOSYN CHD VIOLET 2 (5) 5X70CM DDP	B1058450
NOVOSYN CHD VIOLET 4/0 (1,5) 140CM DDP	C1058204
NOVOSYN CHD VIOLET 3/0 (2) 140CM DDP	C1058205
NOVOSYN CHD VIOLET 2/0 (3) 140CM DDP	C1058207
NOVOSYN CHD VIOLET 0 (3,5) 140CM DDP	C1058208
NOVOSYN CHD VIOLET 1 (4) 140CM DDP	C1058209
NOVOSYN CHD VIOLET 2 (5) 140CM DDP	C1058210
NOVOSYN CHD UNDYED 3/0 (2) 3X45CM DDP	C1059625
NOVOSYN CHD UNDYED 2/0 (3) 3X45CM DDP	C1059626
NOVOSYN CHD UNDYED 3/0 (2) 6X45CM DDP	C1059645
NOVOSYN CHD UNDYED 2/0 (3) 6X45CM DDP	C1059646
NOVOSYN CHD UNDYED 3/0 (2) 12X45CM DDP	C1059665
NOVOSYN CHD UNDYED 2/0 (3) 12X45Cm DDP	C1059666
NOVOSYN CHD UNDYED 2/0 (3) 5X70CM DDP	C1059446
NOVOSYN CHD UNDYED 0 (3,5) 5X70CM DDP	C1059447
NOVOSYN CHD UNDYED 1 (4) 5X70CM DDP	B1059449
NOVOSYN CHD UNDYED 3/0 (2) 140CM DDP	C1059205
NOVOSYN CHD UNDYED 2/0 (3) 140CM DDP	C1059207
NOVOSYN CHD UNDYED 0 (3,5) 140CM DDP	C1059208