

1 p.d.

For urinary tract decompression in fetuses 18-32 weeks.*





Illustration by Lisa Clark

*For the treatment of fetal postvesicular obstructive uropathy.

HARRISON FETAL
BLADDER STENT SET


sterili, vienkartinė, dvigubas stentas, ant pakuotės yra galiojimo laikas mėn, su numatyta pakuotės atidarymo vieta



HARRISON FETAL
BLADDER STENT SET

Used in fetal urinary tract decompression following the diagnosis of fetal postvesicular obstructive uropathy in fetuses of 18 to 32 weeks gestational age.

- The EchoTip® technology enhances the visualization of the needle tip under ultrasound.
- The usable length between pigtails is 1.5-3.5 cm.



Order Number	Reference Part Number	Guide Needle gage	Guide Needle Length cm	Wire Guide Diameter inch/mm	Wire Guide Length cm	Multilength Stent Fr
G16766	J-HFBS-503540	13	18	.038/.97	40	5.0

Humanitarian device: Authorized by federal (USA) law for use in the treatment of fetal postvesicular obstructive uropathy. The effectiveness of this device for this use has not been demonstrated.

CAUTION: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician with appropriate training and experience.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

















EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, IS.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
EDI: cookmedical.com/edi.do
Americas:
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com
Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4



Reference	Description	Features	Packing
TFB000AA160 available variants TFB000AA180 TFB000AA210	Standard blood giving set  PVC DEHP-free I.V. Blood giving set 160 cm long composed by non-vented spike with standard dropper, drip-chamber with 200 µm blood filter, tube Ø 3,0 x 4,1 mm, roller clamp, male luer lock connector with cap. 180 cm long 210 cm long	   20  200 µm    0425	Carton (29x29x49 cm) 125 pcs. Full pallet (80x120x220 cm) 42 cartons
TFB001AA150 available variants TFB000AA180 TFB000AA210	Blood giving set with vented spike Transfuzinė sistema su oro anga Be latekso  Hidrofilinis filtras PVC DEHP-free I.V. Blood giving set 150 cm long composed by vented spike with air-inlet (with hydrophobic filter and ON/OFF cap) with standard dropper, drip-chamber with 200 µm blood filter, tube Ø 3,0 x 4,1 mm, roller clamp, male luer lock connector with cap. Lašų kamera 200 µm dalelių filtras 180 cm long 180 cm ilgis 210 cm long	   20  200 µm    0425 Sterili	Carton (29x29x49 cm) 125 pcs. Full pallet (80x120x220 cm) 42 cartons

Caption: