

Sentrant™

Introducer Sheath  
with Hydrophilic Coating


# THE PROVEN CHOICE FOR ACCESS

The choice for superior hemostasis†



**Medtronic**  
Further. Together

# SUPERIOR SEAL IN STRAIGHTFORWARD AND COMPLEX ANATOMIES†




## Engineered to deliver procedural confidence

- EnsureSeal Technology delivers superior leak resistance versus competitors†
- Coil-reinforced tubing for added stability and kink resistance
- Maintains lubricity after multiple insertions
- Radiopaque dilator shaft and sheath tip for accurate visualization and guidance
- 64 cm configuration launched to service a broader range of anatomies and procedures
- Compatible with a wide range of endovascular portfolios

Vožtuvas užtikrina visišką hermetiškumą po daugkartinio maksimalaus dydžio priemonių įvedimo ir ištraukimo

Trumpiausias introduiserio ilgis - 28 cm, ilgiausias - 64 cm

## Superior leak resistance versus Cook Check-Flo Performer™ introducer sheath and Gore DrySeal™ Flex introducer sheath†



## PRODUCT ORDERING INFORMATION

Catalog Number	Inner Diameter (Fr)	Usable Length (cm)
SENSH1228W	12	28
SENSH1428W	14	28
SENSH1628W	16	28
SENSH1828W	18	28
SENSH2028W	20	28
SENSH2228W	22	28
SENSH2428W	24	28
SENSH2628W	26	28
SENSH1264W	12	64
SENSH1464W	14	64
SENSH1664W	16	64
SENSH1864W	18	64
SENSH2064W	20	64
SENSH2264W	22	64
SENSH2464W	24	64
SENSH2664W	26	64

†Leak Resistance Bench Test Data on file at Medtronic. Test data not indicative of clinical performance. Bench Test compared Cook Check-Flo Performer™ 16 Fr and Gore DrySeal™ Flex 16 Fr to Sentrant™ 16 Fr. Graph shows average leakage over an extrapolated 2.5 hour procedure.

**Important Information:** Prior to use, refer to the Instructions for Use for indications, contraindications, suggested procedure, warnings, and precautions.

**Indications for Use:** The Medtronic Sentrant™ Introducer Sheaths with Hydrophilic Coating are intended to provide a conduit for the insertion of diagnostic or endovascular devices into the vasculature and to minimize blood loss associated with such insertions.

**CAUTION:** Federal (USA) law restricts these devices to sale by or on the order of a licensed healthcare practitioner.

**Aortic | Peripheral | Venous**  
3576 Unocal Place  
Santa Rosa, CA 95403  
USA

**24-hour Technical Support**  
Toll free: +1.800.328.2518

**Orders**  
Toll free: +1.877.526.7890  
Fax: +1.800.838.3103

**CardioVascular LifeLine  
Customer Support**  
Tel: +1.763.526.7890  
Toll free: +1.877.526.7890

Medtronic and the Medtronic logo are trademarks of Medtronic.™ Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

[medtronic.com/aortic](http://medtronic.com/aortic)