

**HARMS TRABECULOTOME SET**

Product number: E0421

**BAUSCH + LOMB**  
**Instruments**

9 mm long pointed tips with 3 mm spread. Pair consists of left and right. Round knurled handle. Polished finish. Overall length 46 mm.

↓ [General reprocessing instructions](#)

UNITS PER BOX

› 1

STERILE


› No

© 2016 BAUSCH &amp; LOMB INCORPORATED. |

This site is intended for Healthcare Professionals in EMEA only. All information and materials on this site pertain to the Bausch & Lomb EMEA Region, unless otherwise indicated.

**HAEFLIGER PHACO CLEAVER**

Product number: 55306



left. enlarged surface and cutting edge facilitate the manipulation of the nucleus with less risk of damage to the posterior capsule. 100 mm

[↕ General reprocessing instructions](#)

INSTRUMENTGROUP

[➤ Manipulators, Hooks](#)

UNITS PER BOX

[➤ 1](#)

STERILE

[➤ No](#)

VARIANT

[➤ 55305 Right](#)**RELATED PRODUCTS****LUND IRIS HOOK**14 mm tip. Angled 45°. Overall length 108 mm**RENTSCH REPOSITION DOUBLE HOOK**0.2 mm hooks. Overall length 137 mm

## TROCAR PORTS GRASPING FORCEPS

Product number: 65365



Forceps for grasping trocar ports in 23G and 25G PPV surgery. For easy removal of trocar port from sclera. Compatible with the Bausch + Lomb Entry Site Alignment (ESA) valved cannula system.

↓ [General reprocessing instructions](#)

INSTRUMENTGROUP

› Forceps

UNITS PER BOX

› 1

STERILE

› No

## RELATED PRODUCTS



## 20GA ENDGRIPPING ILM FORCEPS

20G single use end grasping vitreoretinal forceps. Specifically designed to remove delicate ILM and epiretinal membrane. Low glare instrument shaft. Sterile single use. 12/box.



**UTRATA CAPSULORHEXIS FORCEPS**

Product number: E2002 F



delicate, fine angled tips ideal for grasping delicate capsule during 1.8mm MICS procedures. And 12 mm shafts. Overall length 87 mm

⚡ [General reprocessing instructions](#)

INSTRUMENTGROUP

➤ Forceps

UNITS PER BOX

➤ 1

STERILE

➤ No

**RELATED PRODUCTS**