

PremiCron[®]

The polyester suture generation designed for Cardiovascular Surgery
Non-absorbable, braided and coated suture made of polyester



Sutures

PremiCron®

The polyester suture generation designed for Cardiovascular Surgery



Passion for Cardiovascular Surgery

An increasingly ageing population combined with more complex and diffuse cardiac disease means there are more challenging cardiac cases in daily practice. For that reason B. Braun is committed to providing you with the highest standards of quality and innovative surgical products sharing with you the same Passion for Cardiovascular Surgery.



Designed for the most demanding surgical indications

PremiCron® is a non-absorbable braided suture made of polyester and uniformly silicone coated with all required features for heart valve and annuloplastic ring implantation.

Main indications:


- Cardiovascular surgery
- Heart valve replacement
- Annuloplastic ring implantation

Also indicated in:

- General surgery




A braided thread structure designed for high resistance and reliable knotting



For high resistance and reliable knotting

PremiCron®'s braided thread structure is designed with a central core covered by several small compact braids of polyester filaments for **high resistance and reliable knotting, with an excellent handling.**

PremiCron® sutures are coated with silicone to decrease surface friction for an easy knot run down and atraumatic suturing of the cardiac tissue.



Central braided core¹

Central core covered by several small compact braids coated with silicone for a secure and atraumatic knotting.



Silicone coating

To provide a smooth surface for easy knot run down and atraumatic suturing.

¹ - only sutures from 5 USP (7 metric) to 2/0 USP (3 metric)



PremiCron[®]

OVAL Pledgets



Superior anatomical adaptability and conformity with less foreign material

- Superior conformity to valve contour tissue
- Avoidance of overlapping effect
- Less foreign material



Oval pledget design for cardiac valve surgery

The selected geometry **offers high anatomical adaptability and conformity.**

Reduce overlapping effect

Bio-prosthetic and mechanical valve suture with pledgets are placed **with little to no space** between the sequential horizontal mattress sutures. For that reason **oval pledgets**, without corners, **reduce the overlapping problem** and in consequence the over inflammatory response as a defense reaction to foreign material.



Firm and soft pledgets



Firm: Compact PTFE pledget.

Soft: Superior flexibility with secure resistance properties for superior conformity to valve tissue.

Firm and soft PTFE pledgets:

- Firm**
 - Compact PTFE pledget for valve and ring replacement
- Soft**
 - Soft pledgets are also made of 100 % PTFE and laser cut with **increased porosity**
 - Higher flexibility than firm pledgets with similar resistance
- Advantages**
 - Superior conformity to the tissue

| Pledgets | Firm | Soft | |
|----------------|--|-------------------|--|
| Material | 100 % PTFE | 100 % PTFE | |
| Thickness | 1.5 mm +/- 0.2 mm | 1.5 mm +/- 0.2 mm | |
| Flexibility | 1500 –1600 mN | 400 – 500 mN | |
| Bridge rupture | 8.06 kg | 7.52 kg | |
| Design | <div><div></div><div></div><div></div><div></div></div> <div><div>Pledgets 3 x 3 mm</div><div>Pledgets 6 x 3 mm</div><div>Pledgets 3 x 3 mm</div><div>Pledgets 6 x 3 mm</div></div> <div><div></div><div></div><div></div></div> <div><div>Pledgets 9 x 5 mm</div><div>Pledgets 9 x 5 mm</div><div>Pledgets 16 x 7 mm</div></div> | | |



PremiCron[®]

CV Pass, Cardiovascular needle design

Precise penetration

Exclusive 300 stainless steel series laser drilled silicone coated needles which have specific tip and body design for Cardiovascular Surgery indication.

Non-reflective cardiovascular black needles: CV Pass Easyblack

Black needles avoid glare effect for higher precision and superior visual comfort during your procedure.

Black needles are identified by a 'b'.



■ Specific tip design for cardiac surgery for precise and smooth penetrations

- Taper Point
- Short Cutting Tip

■ Smooth and atraumatic penetrations

Silicone coating to ease
penetrations


■ Excellent needle to thread ratio

■ Semi square body design for a better stability in the needle holder




Specific cardiac valve needle design

Intra-annular technique



Supra-annular technique



Special curvature for valve replacement

PremiCron® offers several needle-thread combinations for heart valve repair and replacement, covering specific requirements depending whether supra-annular or intra-annular positioning of the valve prosthesis.

■ Intra-annular prosthesis positioning

Needle curvature $< 180^\circ$

Needle range: ■ HR18, HR23, HR27
■ HRC18, HRC23, HRC27

■ Supra-annular prosthesis positioning

Needle curvature $> 180^\circ$

Needle range: ■ HR17, HR22, HR26
■ HRC17, HRC22, HRC26


Specific cardiovascular needle tip design for a precise and smooth penetration.

Taper point

Finest tip design for Cardiac Surgery to reduce tissue trauma.

Short cutting point

Superior penetration for **calcified tissues** without being atraumatic.



PremiCron®

Product overview

















| | |
|----------------------|--|
| Structure | Polyfilament, braided and coated |
| Colour | Green or white |
| Chemical composition | Polyester |
| Coating | Silicone |
| Origin | Synthetic |
| Sizes | USP 6/0 (0.7 metric) to USP 5 (7 metric) |
| Type of absorption | Non-absorbable |
| Indication | Cardiovascular surgery General surgery |
| Sterilization | DDP* without Pledgets: gamma irradiation DDP* with Pledgets: ethylene oxide MLP*: ethylene oxide |

* DDP: Direct Dispense Packaging

**MLP: Multipack Packaging



Needle-thread combination




















| Needle | | | Thread | USP metric | | | | | | |
|---|---|---|------------------|------------|----------|----------|----------|----------|-------|-----|
| length mm | angle | diameter microns | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 1/4 circle lancet needle with micro-point | | | | | | | | | | |
|  |  | 90° | 45 | | G0026155 | | | | | |
| | | | Green | | | | | | | |
| | | | 330µ | | | | | | | |
| | | 90° | 45 | G0027114 | G0027115 | | | | | |
| | | 90° | White | | | | | | | |
| | | 90° | 45 | | | G0027116 | | | | |
| | | 90° | White | | | | | | | |
| 1/2 circle round bodied needle | | | | | | | | | | |
|  |  | | 75 | | C0026213 | | | | | |
| | | | Green | | | | | | | |
|  |  | | 45 | C0026008 | C0026009 | C0026010 | | | | |
| | | | Green | | | | | | | |
|  |  | | 75 | | C0026695 | | | | | |
| | | | Green | | | | | | | |
| | |  | 90 | | C0026792 | | | | | |
| | | | Green | | | | | | | |
|  |  | | 45 | | C0026029 | C0026030 | C0026031 | | | |
| | | | 75 | | | C0026152 | C0026032 | | | |
| | | | Green | | | | | | | |
| |  | | 4 x 75 | | | M0026699 | | | | |
| | | | Green | | | | | | | |
|  |  | | 75 | | C0026843 | C0026844 | C0026845 | | | |
| | | | Green | | | | | | | |
|  |  | | 75 | | C0026003 | C0026004 | C0026005 | C0026006 | | |
| | | | Green | | | | | | | |

C Box of 36 units













B Box of 24 units

G Box of 12 units

M Box of 6 units

| Needle | | Thread | USP metric | | | | | | |
|--|--|---|-------------|----------------------|----------------------|----------------------|----------------------|----------|--------|
| length mm | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 1/2 circle round bodied needle | | | | | | | | | |
|   "Take off" | | 4 x 45 Green | | | C0026704 | C0026740 | | | |
| | | 75 90 Green | C0026802 | C0026803 C0026903 | C0026804 C0026904 | C0026805 C0026905 | C0026806 C0026906 | | |
| | | 75 90 White | | | C0027455 | | C0027457 C0027906 | | |
| | | | | | | | | | |
|   |  | 4 x 75 Green | | | M0026708 | | | | |
| |  | 8 x 75 G/W | | | | | M0027775 | | |
| |   Pledgets 3x3 mm Firm | 4 x 45 8 x 75 G/W | | | | | M0027764 M0027802 | | |
| |   Pledgets 6x3 mm Firm | 8 x 75 G/W | | | | | M0027776 | | |
|   | | 75 Green | | | | | C0026211 | | |
| | |  Pledgets 3x3 mm Firm | 75 Green | | | | C0026252 | | |
| | |  Pledgets 3x3 mm Soft | 75 Green | | | | C0026257 | | |
| |  | 8 x 75 G/W | | | | | M0027151 | | |
| |   Pledgets 3x3 mm Firm | 4 x 75 8 x 75 G/W | | | | | M0027090 M0027119 | | |
| |   Pledgets 3x3 mm Soft | 4 x 75 8 x 75 G/W | | | | | M0027105 M0027133 | | |
| | | | | | | | | | |
| | | | | | | | | | |

Needle-thread combination

| Needle | | | Thread | USP metric | | | | | | | |
|---|---------|---|--|--|----------|------------|----------|----------|----------|----------|----------|
| length mm | | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | |
| 1/2 circle round bodied needle | | | | | | | | | | | |
| CV Pass | 2xHR 18 |  |  Multi pack |  Pledgets 3x3 mm Firm | 8 x 75 | | | | M0027178 | | |
| | | | | | G/W | | | | | | |
| CV Pass | 2xHR 20 |  |  Multi pack |  Pledgets 3x3 mm Firm | 8 x 75 | | | | M0027953 | | |
| | | | | | G/W | | | | | | |
| | | |  Multi pack |  Pledgets 6x3 mm Firm | 4 x 75 | | | | | M0027905 | |
| | | | | | 8 x 75 | | | | | M0027706 | |
| | | | | | G/W | | | | | | |
| | | |  Multi pack |  Pledgets 3x3 mm Soft | 8 x 75 | | | | | M0027945 | |
| | G/W | | | | | | | | | | |
|  Pledgets 6x3 mm Soft | 8 x 75 | | | | | | M0027525 | | | | |
| | | | G/W | | | | | | | | |
| CV Pass | HR 22 |  | | | 45 | | | | C0026007 | | |
| | | | | | 75 | | | C0026024 | | C0026025 | C0026026 |
| | | | | | Green | | | | | | |

C Box of 36 units









B Box of 24 units

G Box of 12 units

M Box of 6 units















PremiCron®

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | | |
|---|---|---|-------------------------|----------|------------|----------|----------------------|----------------------|--------|--|
| length mm | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | |
| 1/2 circle round bodied needle | | | | | | | | | | |
| <div>CV Pass</div> <div>2xHR 22</div> | | 90 Green | | | C0026730 | C0026731 | C0026732 | | | |
| |  Pledgets 6x3 mm Firm | 75 Green | | | C0026735 | C0026307 | C0026733 | | | |
| | | 75 White | | | | | C0027034 | | | |
| |  Pledgets 6x3 mm Firm | 2 x 75 G/W | | | | | B0027718 | | | |
| | <div>Multi pack</div> | 8 x 75 G/W | | | | | M0027705 | | | |
| | <div>Multi pack</div>  Pledgets 6x3 mm Firm | 8 x 75 8 x 90 G/W | | | | | M0027605 M0027704 | | | |
| | | <div>Multi pack</div>  Pledgets 6x3 mm Soft | 8 x 75 G/W | | | | | M0027867 | | |
| | <div>CV Pass</div> <div>2xHR 22b</div> | | 75 Green | | | | | C0026265 | | |
| | |  Pledgets 6x3 mm Firm | 75 Green | | | | | C0026460 | | |
| | |  Pledgets 6x3 mm Soft | 75 Green | | | | | C0026493 | | |
| <div>Multi pack</div> | | 8 x 75 G/W | | | | | M0027181 | | | |
| <div>Multi pack</div>  Pledgets 6x3 mm Firm | | 4 x 75 8 x 75 G/W | | | | | M0027157 M0027169 | | | |
| | | <div>Multi pack</div>  Pledgets 6x3 mm Soft | 4 x 75 8 x 75 G/W | | | | | M0027163 M0027175 | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|---|--|-----------------------------------|------------|------------|----------------------|----------------------------------|----------------------|----------|--------|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 |
| 1/2 circle round bodied needle | | | | | | | | | |
|  | 2xHR 22ss  "Take off" | 4 x 45 Green | | | | C0026736 | B0026729 | | |
|  | 2xHR 23  | 75 Green | | | | C0026598 | | | |
| |   | 8 x 75 G/W | | | | M0027232 | | | |
|  | HR 26  | 75 120 Green | | C0026014 | C0026015 C0026002 | C0026016 C0026111 | C0026017 C0026111 | C0026018 | |
| | | 75 White | | | | G0027016 | | | |
| | | 75 90 Green | | C0026814 | C0026738 C0026815 | C0026739 C0026816 | C0026817 | | |
| | | 75 90 White | | | | C0027839 C0027816 | C0027817 | | |
|  | 2xHR 26  | 90 Green | | | C0026915 | C0026916 | | | |
| |  | 90 White | | | | C0027306 | | | |
| |  | 4 x 75 8 x 75 8 x 90 G/W | | | | M0027700 M0027715 M0027713 | | | |
| |   | 8 x 75 8 x 90 G/W | | | | M0027716 M0027717 | | | |

C Box of 36 units











B Box of 24 units

G Box of 12 units

M Box of 6 units



















PremiCron®

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | | |
|--|--|---|-------------------------|------------|----------|----------------------|----------|----------|--------|----------|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | |
| 1/2 circle round bodied needle | | | | | | | | | | |
| <div><div>2xHR 26</div><div></div><div><div>Multi pack</div><div><div>Pledgets 6x3 mm Firm</div></div></div></div> | | 4 x 75 4 x 90 8 x 75 8 x 90 G/W | | | M0027752 | M0027766 M0027712 | M0027759 | | | |
| | <div><div>2xHR 26b</div><div></div><div><div>Multi pack</div><div><div>Pledgets 6x3 mm Firm</div></div><div><div>Multi pack</div><div><div>Pledgets 6x3 mm Soft</div></div></div></div></div> | | 75 Green | | | C0026500 | | | | |
| | | | 75 Green | | | C0026522 | | | | |
| | | | 75 Green | | | C0026527 | | | | |
| | | | 8 x 75 G/W | | | M0027205 | | | | |
| | | | 4 x 75 8 x 75 G/W | | | M0027187 M0027199 | | | | |
| | | | 4 x 75 8 x 75 G/W | | | M0027193 M0027211 | | | | |
| | <div><div>HR 26s</div><div></div><div><div>EasySlide</div></div></div> | | 75 Green | | | C0026052 | C0026050 | C0026051 | | B0026053 |
| | <div><div>HR 26s</div><div></div><div><div>EasySlide</div><div>"Take off"</div></div></div> | | 4 x 45 Green | | | C0026724 | B0026019 | | | |
| | <div><div>HR 26ss</div><div></div><div><div>EasySlide</div><div>"Take off"</div></div></div> | | 4 x 45 Green | | | | B0026020 | | | |
| <div><div>2xHR 27</div><div></div><div><div>Multi pack</div><div><div>Pledgets 6x3 mm Firm</div></div></div></div> | | 8 x 75 G/W | | | M0027280 | | | | | |

G/W = Green / White

Needle-thread combination

| Needle length mm | Thread length cm colour | USP metric | | | | | | |
|---|-----------------------------------|------------|----------|----------------------|----------------------|----------------------|----------------------|----------------------|
| | | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 |
| 1/2 circle round bodied needle | | | | | | | | |
|   | 75 Green | C0026056 | C0026057 | C0026058 | | C0026059 | | |
|   | 75 Green | C0026035 | C0026036 | C0026037 | C0026038 | C0026039 | | |
|   | 75 90 Green | C0026750 | C0026836 | C0026837 C0026752 | | | | |
|   | 75 90 Green | | | | C0026066 | C0026067 C0026267 | C0026068 C0026277 | C0026069 C0026759 |
|   "Take off" | 4 x 45 Green | | | B0026727 | | | | |
|   | 60 Green | | | | | C0026840 | | |
|   | 75 90 Green | | | | C0026074 C0026972 | B0026075 | B0026425 | B0026083 |
|   | 75 Green | | | C0026098 | C0026080 | C0026182 | B0026082 | |
| | 75 White | | | | | | B0027082 | |
|   "Take off" | 4 x 45 Green | | B0026065 | B0026721 | | | | |

C Box of 36 units








B Box of 24 units

G Box of 12 units

M Box of 6 units







PremiCron®

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|--------------------------------|--|---|------------|------------|----------|------------|----------|----------|----------|
| length mm | | length cm colour | 7/0 0.5 | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 |
| 3/8 circle round bodied needle | | | | | | | | | |
| CV Pass | 2xDR 8  | 45 Green | | | C0026822 | | | | |
| | | 45 White | | | C0027822 | | | | |
| CV Pass | 2xDR 10  | 75 Green | | C0026871 | | | | | |
| CV Pass | DR 12  | 45 Green | | | C0026121 | | | | |
| CV Pass | 2xDR 12  | 75 Green | | C0026852 | C0026853 | C0026854 | C0026859 | | |
| | | 75 White | | | C0027853 | C0027854 | | | |
| | |  Pledgets 3x3 mm Firm | | | | C0026954 | | | |
| CV Pass | DR 15  | 75 Green | | | C0026130 | C0026131 | | | |
| CV Pass | DR 18  | 75 Green | | | C0026104 | C0026105 | | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|--------------------------------|---|---------------------|------------|------------|----------|------------|----------|----------|----------|
| length mm | | length cm colour | 7/0 0.5 | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 |
| 3/8 circle round bodied needle | | | | | | | | | |
| CV/Pass |  | 75 | | C0026863 | C0026864 | C0026865 | | | |
| | | 90 | | | | C0026858 | C0026883 | | |
| | | Green | | | | | | | |
| | | 8 x 90 | | | | | M0027769 | | |
| | | G/W | | | | | | | |
| CV/Pass |  | 75 | | | | C0026226 | C0026227 | C0026228 | |
| | | Green | | | | | | | |
| CV/Pass |  | 75 | | | | C0026136 | C0026137 | | |
| | | Green | | | | | | | |
| CV/Pass |  | 90 | | | | | C0026873 | | |
| | | Green | | | | | | | |
| Needle | | Thread | USP metric | | | | | | |
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 |
| 3/8 circle round bodied needle | | | | | | | | | |
| Easyslide |  | 75 | | | C0026204 | C0026141 | C0026142 | C0026238 | C0026239 |
| | | Green | | | | | | | |
| Easyslide |  | 75 | | | | C0026328 | C0026071 | | |
| | | Green | | | | | | | |

C Box of 36 units














B Box of 24 units

G Box of 12 units

M Box of 6 units












PremiCron®

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|--|---|---------------------|------------|------------|----------------------|----------------------|----------|--------|--------|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 |
| Round bodied ski needle | | | | | | | | | |
|  |  SKR 19 | 30 Green | | | C0026535 | C0026536 | | | |
|  |  SKR 26 | 30 120 Green | | | C0026540 C0026546 | C0026550 C0026046 | C0026549 | | |
| Straight round bodied needle | | | | | | | | | |
|  |  GR 19 | 75 Green | | C0026113 | | | | | |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | |
|   2xHRC 17 | | 90 Green | | | C0026305 | C0026306 | | | |
| | | 90 White | | | | C0027317 | | | |
| |  | 75 Green | | | | C0026249 | | | |
| |  | 75 90 Green | | | C0026240 | C0026248 C0026308 | C0026251 | | |
| |  | 2 x 90 G/W | | | | B0027875 | | | |
| |  | 75 Green | | | | C0026560 | | | |
| |  | 75 Green | | | | C0026563 | | | |
| | | | | | | | | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | | |
|---|---|---|-------------------------|------------|----------|----------------------|----------------------|--------|--------|--|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | | |
| <div>CV Pass</div> <div>2xHRC 17</div>  | | 4 x 75 8 x 75 G/W | | | | M0027921 M0027722 | M0027585 | | | |
| | <div>Multi pack</div> | <div> Pledgets 3x3 mm Firm</div> | 8 x 75 G/W | | | M0027738 | | | | |
| | <div>Multi pack</div> | <div> Pledgets 6x3 mm Firm</div> | 4 x 75 8 x 75 G/W | | | M0027628 M0027723 | M0027579 M0027675 | | | |
| | <div>Multi pack</div> | <div> Pledgets 3x3 mm Soft</div> | 4 x 75 8 x 75 G/W | | | M0027824 M0027819 | | | | |
| | <div>Multi pack</div> | <div> Pledgets 6x3 mm Soft</div> | 4 x 75 8 x 75 G/W | | | M0027813 M0027811 | | | | |
| | | | 75 Green | | | | C0026538 | | | |
| | | <div> Pledgets 3x3 mm Firm</div> | 75 Green | | | | C0026556 | | | |
| | | <div> Pledgets 6x3 mm Firm</div> | 75 Green | | | | C0026579 | | | |
| | | <div> Pledgets 3x3 mm Soft</div> | 75 Green | | | | C0026568 | | | |
| | | <div> Pledgets 6x3 mm Soft</div> | 75 Green | | | | C0026587 | | | |
| | <div>Multi pack</div> | | 8 x 75 G/W | | | | M0027253 | | | |
| | <div>Multi pack</div> | <div> Pledgets 3x3 mm Firm</div> | 4 x 75 8 x 75 G/W | | | | M0027217 M0027241 | | | |
| <div>Multi pack</div> | <div> Pledgets 6x3 mm Firm</div> | 4 x 75 8 x 75 G/W | | | | M0027229 M0027259 | | | | |

C Box of 36 units




















B Box of 24 units

G Box of 12 units

M Box of 6 units





PremiCron®

Needle-thread combination

| Needle | | | Thread | USP metric | | | | | | |
|---|---|---|-------------------------|------------|----------|------------|----------|----------------------|----------|--------|
| length mm | | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | | |
|  2xHRC 17b  |  |  Pledgets 3x3 mm Soft | 4 x 75 | | | | | M0027223 | | |
| | | | 8 x 75 G/W | | | | | M0027247 | | |
| |  |  Pledgets 6x3 mm Soft | 4 x 75 | | | | | M0027235 | | |
| | | | 8 x 75 G/W | | | | | M0027265 | | |
|  2xHRC 18  | |  Pledgets 3x3 mm Firm | 75 Green | | | | | C0026638 | | |
| | |  Pledgets 6x3 mm Firm | 75 Green | | | | | C0026646 | | |
| | |  Pledgets 3x3 mm Soft | 75 Green | | | | | C0026570 | | |
| | |  Pledgets 6x3 mm Soft | 75 Green | | | | | C0026577 | | |
| |  | | 8 x 75 G/W | | | | | M0027498 | | |
| | |  Pledgets 3x3 mm Firm | 8 x 75 G/W | | | | | M0027480 | | |
| |  |  Pledgets 6x3 mm Firm | 8 x 75 G/W | | | | | M0027486 | | |
| | |  Pledgets 3x3 mm Soft | 8 x 75 G/W | | | | | M0027828 | | |
| |  |  Pledgets 6x3 mm Soft | 4 x 75 8 x 75 G/W | | | | | M0027832 M0027834 | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | | |
|---|---|-------------------------|---------------|------------|----------|----------------------|----------|----------|--------|----------|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | | |
| <div><div>CV Pass</div><div>Easyblock</div><div>2xHRC 18b</div></div> | | 75 Green | | | | C0026592 | | | | |
| |  | 75 Green | | | | C0026596 | | | | |
| |  | 75 Green | | | | C0026608 | | | | |
| |  | 75 Green | | | | C0026600 | | | | |
| |  | 75 Green | | | | C0026612 | | | | |
| | <div>Multi pack</div> | 8 x 75 G/W | | | | M0027353 | | | | |
| | <div>Multi pack</div> | 4 x 75 8 x 75 G/W | | | | M0027271 M0027295 | | | | |
| | <div>Multi pack</div> | 4 x 75 8 x 75 G/W | | | | M0027283 M0027360 | | | | |
| | <div>Multi pack</div> | 4 x 75 8 x 75 G/W | | | | M0027277 M0027302 | | | | |
| | <div>Multi pack</div> | 4 x 75 8 x 75 G/W | | | | M0027289 M0027368 | | | | |
| | | | | | | | | | | |
| | <div><div>CV Pass</div><div>HRC 22</div></div> | | 100 Green | | | | C0026296 | C0026309 | | |
| | | | 4 x 90 G/W | | | | | | | B0027978 |










C Box of 36 units

B Box of 24 units

G Box of 12 units

















M Box of 6 units

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | | |
|--|---|---|---|-------------------------|----------|----------------------------------|----------------------|----------------------|--------|--|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | | |
| <div>CV Pass</div> <div>2xHRC 22</div>  | | 90 Green | | | | C0026323 | | | | |
| |  Pledgets 6x3 mm Firm | 75 90 Green | | | | C0026263 C0026321 | B0026899 | | | |
| | | 90 White | | | | | C0027537 | | | |
| | | | | | | | | | | |
| |  Pledgets 6x3 mm Soft | 75 Green | | | | C0026585 | | | | |
| | <div>Multi pack</div> | 4 x 75 4 x 90 8 x 75 G/W | | | | M0027549 M0027616 M0027102 | M0027107 M0027610 | | | |
| | |  Pledgets 6x3 mm Firm | 4 x 75 4 x 90 G/W | | | | M0027949 M0027685 | M0027108 M0027618 | | |
| | | |  Pledgets 6x3 mm Firm | 4 x 90 8 x 75 G/W | | | | M0027617 M0027100 | | |
| |  Pledgets 6x3 mm Soft | | | 4 x 75 8 x 75 G/W | | | | M0027836 M0027841 | | |
| | | | | | | | | | | |
| | | <div>CV Pass</div> <div>2xHRC 22b</div>  | | 75 Green | | | | C0026616 | | |
| |  Pledgets 6x3 mm Firm | | 75 Green | | | | C0026620 | | | |
|  Pledgets 6x3 mm Soft | | | 75 Green | | | | C0026625 | | | |
| <div>Multi pack</div> | 8 x 75 G/W | | | | | M0027414 | | | | |
| | | | | | | | | | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|--|---|---|------------|------------|----------|----------|----------|--------|--------|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | |
|   2xHRC 22b |  |  Pledgets 6x3 mm Firm | 4 x 75 | | | | M0027375 | | |
| | | | 8 x 75 | | | | M0027391 | | |
| |  |  Pledgets 6x3 mm Soft | 4 x 75 | | | | M0027383 | | |
| | | | 8 x 75 | | | | M0027402 | | |
|   HRC 26 |  | | 8 x 75 | | | | M0027915 | | |
| | | | G/W | | | | | | |
| | | 90 | | | | C0026313 | | | |
| | | Green | | | | | | | |
| | | 90 | | | | C0027326 | | | |
| | | White | | | | | | | |
| | | 75 | | | | C0026319 | | | |
| | | 90 | | | | C0026318 | | | |
| | | Green | | | | | | | |
| |  Pledgets 6x3 mm Firm | 90 | | | | G0027336 | | | |
| | | White | | | | | | | |
| |  Pledgets 6x3 mm Soft | 75 | | | | C0026594 | | | |
| 90 | | | | | C0026676 | | | | |
| | Green | | | | | | | | |
|  | | 4 x 75 | | | | M0027371 | M0027925 | | |
| | | 4 x 90 | | | | M0027555 | M0027686 | | |
| | G/W | | | | | | | | |
|  |  Pledgets 3x3 mm Firm | 4 x 75 | | | | M0027372 | | | |
| | | 8 x 75 | | | | M0027801 | | | |
| | G/W | | | | | | | | |
|  |  Pledgets 6x3 mm Firm | 4 x 75 | | | | | M0027859 | | |
| | | G/W | | | | | | | |

C Box of 36 units




































B Box of 24 units

G Box of 12 units

M Box of 6 units


















PremiCron®

Needle-thread combination

| Needle | | | Thread | USP metric | | | | | | |
|---|---|---|---|------------|------------|----------|----------|----------|--------|--------|
| length mm | | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | | |
| CV Pass |  |  |  Pledgets 6x3 mm Firm | 4 x 75 | | | M0027543 | | | |
| | | | | 4 x 90 | | | M0027687 | | | |
| | |  |  Pledgets 6x3 mm Firm | 8 x 90 | | | M0027660 | M0027661 | | |
| | | | | G/W | | | | | | |
| CV Pass easyblock |  |  |  Pledgets 6x3 mm Soft | 8 x 75 | | | M0027845 | | | |
| | | | | 8 x 90 | | | M0027531 | | | |
| | |  | | | | | | | | |
| | | | | | | | | | | |
| | |  |  Pledgets 6x3 mm Firm | 75 | | | C0026634 | | | |
| | | | | Green | | | | | | |
| | |  |  Pledgets 6x3 mm Firm | 75 | | | C0026643 | | | |
| | | | | Green | | | | | | |
| | |  |  Pledgets 6x3 mm Soft | 75 | | | C0026650 | | | |
| | | | | Green | | | | | | |
| CV Pass |  |  | | 8 x 75 | | | M0027465 | | | |
| | | | | G/W | | | | | | |
| | |  |  Pledgets 6x3 mm Firm | 4 x 75 | | | M0027424 | | | |
| | | | | 8 x 75 | | | M0027451 | | | |
| | |  |  Pledgets 6x3 mm Firm | 4 x 75 | | | M0027443 | | | |
| | | | | 8 x 75 | | | M0027471 | | | |
| | |  | | | | | | | | |
| | | | | | | | | | | |
| CV Pass |  |  | | 75 | | | C0026982 | | | |
| | | | | 90 | | | C0026618 | | | |
| | |  |  Pledgets 6x3 mm Firm | 75 | | | C0026968 | | | |
| | | | | Green | | | | | | |
| | |  |  Pledgets 6x3 mm Soft | 75 | | | C0026602 | | | |
| | | | | Green | | | | | | |
|  | | | | | | | | | | |
| | | | | | | | | | | |
| CV Pass |  |  |  Pledgets 6x3 mm Soft | 4 x 75 | | | M0027594 | | | |
| | | | | G/W | | | | | | |
| CV Pass |  |  |  Pledgets 6x3 mm Soft | 8 x 75 | | | M0027851 | | | |
| | | | | G/W | | | | | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|---|---|-------------------------|------------|------------|----------|----------------------|----------------------|--------|--------|
| length mm | | length cm colour | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 |
| 1/2 circle round bodied needle with short cutting point | | | | | | | | | |
|   | | 75 Green | | | | C0026658 | | | |
| |  | 75 Green | | | | C0026665 | | | |
| |  | 75 Green | | | | C0026674 | | | |
| |  | 8 x 75 G/W | | | | M0027501 | | | |
| |  | 4 x 75 8 x 75 G/W | | | | M0027477 M0027489 | | | |
| |  | 4 x 75 8 x 75 G/W | | | | M0027483 M0027495 | | | |
| |  | | | | | | | | |
| 3/8 circle round bodied needle with short cutting point | | | | | | | | | |
|   | | 75 Green | C0026352 | | | | | | |
| | | 75 Green | | | | C0026465 | C0026466 | | |
| |  | 75 Green | | | | | C0026476 | | |
|   |  | 75 Green | | | | | C0026636 | | |
| |  | 4 x 75 8 x 75 G/W | | | | | M0027601 M0027311 | | |
| |  | 4 x 75 G/W | | | | | M0027600 | | |
| |  | | | | | | | | |

C Box of 36 units
















B Box of 24 units

G Box of 12 units

M Box of 6 units






PremiCron®

Needle-thread combination

| Needle | | | Thread | USP metric | | | | | |
|---|--|--|------------|------------|------------|----------|----------|----------|----------|
| length mm | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 3/8 circle round bodied needle with short cutting point | | | | | | | | | |
| CV Pass |  |   Pledgets 3x3 mm Soft | 4 x 75 | | | | M0027394 | | |
| | | | 8 x 75 | | | | M0027507 | | |
| | | | G/W | | | | | | |
| |   Pledgets 6x3 mm Soft | 4 x 90 | | | | M0027513 | | | |
| | | 8 x 90 | | | | M0027519 | | | |
| | | G/W | | | | | | | |
| CV Pass EasyBlock |  |  Pledgets 6x3 mm Firm | 75 | | | | C0026161 | | |
| | | | Green | | | | | | |
| | |  Pledgets 6x3 mm Soft | 75 | | | | C0026184 | | |
| | | | Green | | | | | | |
| | |   Pledgets 6x3 mm Firm | 4 x 75 | | | | M0027002 | | |
| | | | 8 x 75 | | | | M0027046 | | |
| | | G/W | | | | | | | |
| |   Pledgets 6x3 mm Soft | 4 x 75 | | | | M0027013 | | | |
| | | 8 x 75 | | | | M0027079 | | | |
| | | G/W | | | | | | | |
| Hook curved round bodied needle with short cutting point | | | | | | | | | |
| EasySlide |  | | 75 | | | | | | C0026987 |
| | | | Green | | | | | | |
| Progressive curved round bodied needle with short cutting point | | | | | | | | | |
| EasySlide |  | | 120 | | | | C0026102 | C0026099 | |
| | | | Green | | | | | | |
| 1/2 circle round bodied needle with trocar point | | | | | | | | | |
| EasySlide |  | | 75 | | C0026379 | C0026380 | | | |
| | | | Green | | | | | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|--|---|---------------------|------------|----------|----------|--------|----------|--------|----------|
| length mm | | length cm colour | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 |
| 1/2 circle round bodied needle with trocar point | | | | | | | | | |
| Easyslide |  | 75 Green | | | C0026392 | | C0026394 | | |
| | | 4 x 75 Green | | | | | B0026717 | | |
| Easyslide |  | 75 Green | | | C0026401 | | | | |
| Easyslide |  | 4 x 75 Green | | | | | B0026718 | | |
| Easyslide |  | 4 x 75 Green | | | | | | | G0026081 |
| Easyslide |  | 4 x 75 Green | | | | | | | G0026723 |

C Box of 36 units













B Box of 24 units

G Box of 12 units

M Box of 6 units


PremiCron®









Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|--|--|---------------------|------------|----------|------------|----------|----------|----------|----------|
| length mm | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 3/8 circle round bodied with trocar point | | | | | | | | | |
|   2xDRT 26 |  Pledgets 6x3 mm Firm | 90 Green | | | | | C0026357 | | |
| |  Multi pack Pledgets 6x3 mm Firm | 8 x 90 G/W | | | | | M0027747 | | |
| 3/8 circle round bodied needle with trocar and micro-point | | | | | | | | | |
|   DRTM 8 | | 45 Green | C0026368 | | | | | | |
| Needle | Thread | USP metric | | | | | | | |
| length mm | length cm colour | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 | |
| 1/2 circle round bodied needle with sternum point | | | | | | | | | |
|   HRS 40 | 4 x 75 Green | | | | | B0026766 | | | |
|   HRS 48 | 4 x 75 Green | | | | | B0026767 | | | G0026702 |
|   HRS 60 | 4 x 75 Green | | | | | | | | G0026705 |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|-----------------------------------|--|---------------------|------------|----------|------------|----------|----------|----------|--------|
| length mm | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 1/2 circle reverse cutting needle | | | | | | | | | |
| Easyslide | HS 15  | 45 75 Green | | | C0026197 | | C0026342 | | |

| Needle | | Thread | USP metric | | | | | | |
|-----------------------------------|--|---------------------|--------------|----------|----------|----------|----------|----------|----------|
| length mm | | length cm colour | 4/0 1.5.7 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 |
| 1/2 circle reverse cutting needle | | | | | | | | | |
| Easyslide | HS 21  | 45 Green | C0026199 | C0026200 | | | | | |
| Easyslide | HS 23  | 45 Green | | C0026202 | | | | | |
| Easyslide | HS 23ss  | 75 Green | | | | C0026403 | C0026404 | C0026405 | |
| Easyslide | HS 26  | 75 Green | | C0026420 | C0026421 | C0026423 | | | |
| Easyslide | HS 30  | 75 Green | | | C0026541 | C0026542 | C0026543 | C0026544 | B0026545 |
| Easyslide | HS 37  | 75 Green | | | C0026471 | C0026472 | C0026473 | | |
| Easyslide | HS 37s  | 75 Green | | | C0026467 | C0026468 | C0026469 | | |
| Easyslide | HS 40  | 75 Green | | | | C0026479 | | | B0026482 |

C Box of 36 units















B Box of 24 units

G Box of 12 units

M Box of 6 units









PremiCron®

Needle-thread combination

| Needle | | Thread | USP metric | | | | | | |
|---|---|-------------------------|------------|----------|----------------------|----------------------|----------------------|----------------------|----------|
| length mm | | length cm colour | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 |
| 1/2 circle reverse cutting needle | | | | | | | | | |
|  |  | 75 Green | | | | | C0026390 | | |
|  |  | 75 Green | | | | C0026375 | C0026376 | B0026377 | B0026399 |
| | | 4 x 75 Green | | | | | | | G0026203 |
|  |  | 75 Green | | | | | | | B0026409 |
| Needle | | Thread | USP metric | | | | | | |
| length mm | | length cm colour | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 3/8 circle reverse cutting needle | | | | | | | | | |
|  |  | 45 Green | C0026171 | C0026172 | C0026173 | | | | |
|  |  | 45 75 Green | | C0026176 | C0026177 C0026245 | C0026246 | | | |
|  |  | 45 75 Green | | C0026551 | C0026552 C0026231 | C0026553 C0026232 | C0026554 C0026233 | | C0026235 |
| | | 45 White | | | | C0027974 | | | |
|  |  | 45 75 90 Green | | | C0026514 | C0026189 C0026515 | C0026190 C0026516 | C0026517 C0026340 | C0026518 |
| | | 75 White | | | | | C0027965 | | |

G/W = Green / White

Needle-thread combination

| Needle | | Thread | USP metric | | | | | |
|---|--|---------------------|------------|----------------------|----------|----------|----------|----------|
| length mm | | length cm colour | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 |
| 3/8 circle reverse cutting needle | | | | | | | | |
|  DS 30 | | 75 Green | | C0026531 | C0026532 | C0026533 | C0026534 | |
|  DS 35 | | 75 Green | | C0026281 | C0026282 | C0026283 | C0026284 | |
|  DS 39 | | 75 Green | | C0026329 | C0026330 | C0026331 | C0026332 | B0026333 |
|  DS 48 | | 75 Green | | | C0026353 | C0026354 | C0026355 | |
|  DS 60 | | 75 Green | | | | | C0026278 | B0026279 |
| Straight reverse cutting needle | | | | | | | | |
|  GS 51 | | 75 Green | C0026505 | C0026506 | C0026507 | C0026508 | C0026509 | |
|  GS 60 | | 75 90 Green | C0026314 | C0026315 C0026441 | | | | |
|  GS 65 | | 75 Green | | C0026287 | C0026288 | C0026289 | | |

C Box of 36 units









B Box of 24 units

G Box of 12 units

M Box of 6 units

PremiCron®

Needle-thread combination

| Needle length mm | Thread length cm colour | USP metric | | | | | | |
|---|--|-------------|----------|------------|----------|----------|----------|--------|
| | | 6/0 0.7 | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 |
| 3/8 reverse cutting needle with precision point | | | | | | | | |
|  | 2xDSMP 11  | 45 White | | C0027031 | | | | |
|  | DSMP 13  | 45 Green | | | C0026193 | | | |
| | | 45 White | | C0027086 | C0027087 | | | |
|  | DSMP 19  | 45 Green | | | G0026610 | | | |
| | | 45 White | | | C0027096 | C0027097 | | |
|  | DSMP 24  | 75 White | | | | C0027957 | | |

Needle-thread combination

| | Thread | USP metric | | | | | | | |
|---------------------|------------------|------------|------------|----------|----------|----------|----------|----------|----------|
| | length cm | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 1 | 0 3.5 | 1 4 | 2 5 | 3 6 |
| Pre-cut sutures | | | | | | | | | |
| 3 threads of 45 cm | 3 x 45 Green | C0120090 | C0120091 | C0120092 | C0120093 | C0120094 | C0120095 | C0120096 | C0120097 |
| 5 threads of 45 cm | 5 x 45 Green | C0120100 | C0120101 | C0120102 | C0120103 | C0120104 | B0120105 | B0120106 | |
| | 5 x 45 White | | | | | | | B0102030 | |
| 10 threads of 45 cm | 10 x 45 Green | C0120005 | C0120006 | C0120007 | C0120008 | B0120009 | B0120010 | B0120011 | |
| 15 threads of 45 cm | 15 x 45 Green | | | C0120040 | | | | | |
| 17 threads of 45 cm | 17 x 45 Green | | | C0120050 | C0120051 | | | | |
| 5 threads of 60 cm | 5 x 60 Green | | | | | C0120110 | B0120111 | | |
| 10 threads of 60 cm | 10 x 60 Green | | C0120015 | C0120016 | C0120017 | | | | |
| | Thread | USP metric | | | | | | | |
| | length cm | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 |
| Pre-cut sutures | | | | | | | | | |
| 2 threads of 75 cm | 2 x 75 Green | | C0120080 | C0120081 | C0120082 | C0120083 | C0120084 | B0120085 | B0120087 |
| | 2 x 75 White | | | C0100017 | | C0100008 | | | |
| 5 threads of 75 cm | 5 x 75 Green | | C0120119 | C0120120 | C0120121 | B0120122 | B0120123 | | |
| 10 threads of 75 cm | 10 x 75 Green | | C0120020 | C0120021 | B0120022 | | B0120024 | | |

C Box of 36 units

B Box of 24 units

G Box of 12 units




M Box of 6 units





PremiCron®

Needle-thread combination

| | Thread length cm | 4/0 1.5 | 3/0 2 | 2/0 3 | USP metric | | | | |
|--------------------|------------------------|------------|------------|----------|------------|----------|----------|----------|----------|
| | | | | | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 |
| Pre-cut sutures | | | | | | | | | |
| 1 thread of 150 cm | 1 x 150 Green | C0120025 | C0120026 | C0120027 | C0120028 | C0120029 | C0120030 | C0120031 | B0120033 |
| 1 thread of 150 cm | 1 x 150 White | | | | | C0102020 | C0102025 | C0102035 | |
| | Thread length cm | 5/0 1 | 4/0 1.5 | 3/0 2 | 2/0 1 | 0 3.5 | 1 4 | 2 5 | 3 6 |
| Ligature reels | | | | | | | | | |
| 1 thread of 250 cm | 250 Green | G0120055 | G0120056 | G0120060 | G0120061 | G0120062 | G0120063 | G0120064 | G0120065 |
| | Thread length cm | 4/0 1.5 | 3/0 2 | 2/0 3 | 0 3.5 | 1 4 | 2 5 | 3 6 | 5 7 |
| Ligature reels | | | | | | | | | |
| 1 thread of 250 cm | 250 Green | | | | | | | | G0120067 |

Needle-thread combination

| Pledgets FIRM | | Units / foil | Code |
|----------------|---|--------------|----------|
| 3 x 3 x 1.5 mm |  | 6 | C0027996 |
| 6 x 3 x 1.5 mm |  | 6 | C0027995 |
| 9 x 5 x 1.5 mm |  | 6 | C0027994 |

| Pledgets SOFT | | Units / foil | Code |
|-----------------|---|--------------|----------|
| 3 x 3 x 1.5 mm |  | 6 | C0027990 |
| 6 x 3 x 1.5 mm |  | 6 | C0027991 |
| 9 x 5 x 1.5 mm |  | 6 | C0027992 |
| 16 x 7 x 1.5 mm |  | 6 | C0027993 |



B. Braun Surgical SA | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Phone +34 93 5 86 62 00 | Fax +34 93 6 99 63 30 | www.bbraun.es

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

Subject to technical changes. All rights reserved.
This brochure may only be used for the exclusive
purpose of obtaining information about our
products. Reproduction in any form partial or
otherwise is not permitted.

Brochure No. B21602

0212/1/1