

Endo GIA™ reloads with Tri-Staple™ technology and Endo GIA™ Ultra universal staplers



Endo GIA™ Ultra universal handle

45° of articulation, approved for 25 firings. **Compatible with Tri-Staple™ technology and Endo GIA™ Ultra universal single use loading units (SULUs).**


3.1.8 Su tuo pačiu aparatu galima panaudoti įvairaus ilgio/kabučių ankščių pasisukančias kasetes, kuriose integruotas vienkartinis peilis.

Order Code	Description	Box Qty
EGIAUSHORT	Endo GIA™ Ultra universal short stapler 6 cm shaft length	3
EGIAUSTND	Endo GIA™ Ultra universal standard stapler 16 cm shaft length	3
EGIAUXL	Endo GIA™ Ultra universal XL stapler 26 cm shaft length	3

3.1.4 Aparato darbinės dalies ilgis pasirinktinai 6 cm, 16 cm, 26 cm

Endo GIA™ loading units

The reload design with Tri-Staple™ technology, intended to be used over a wider range of tissue thicknesses. They are compatible with the Endo GIA™ Ultra and the standard Endo GIA™ stapling system.



Order Code	Color Code	Description	Box Qty
EGIA30AV	Grey	Endo GIA™ 30 mm articulating vascular reload	6
EGIA30AVM	Yellow	Endo GIA™ 30 mm articulating vascular/medium reloads With Tri-Staple™ technology	6
EGIA30AMT	Purple	Endo GIA™ 30 mm articulating medium/thick reloads With Tri-Staple™ technology	6
EGIA45AV	Grey	Endo GIA™ 45 mm articulating vascular reloads*	6
EGIA45AVM	Yellow	Endo GIA™ 45 mm articulating vascular/medium reloads With Tri-Staple™ technology	6
EGIA45AMT	Purple	Endo GIA™ 45 mm articulating medium/thick reloads With Tri-Staple™ technology	6
EGIA45AXT	Black	Endo GIA™ 45 mm articulating x-tra thick reloads With Tri-Staple™ technology	6
EGIA60AVM	Yellow	Endo GIA™ 60 mm articulating vascular medium reloads With Tri-Staple™ technology	6
EGIA60AMT	Purple	Endo GIA™ 60 mm articulating medium/thick reloads With Tri-Staple™ technology	6

*Code EGIA30AV and EGIA45AV do not come with Tri-Staple™ technology, have a slimmer fixed anvil

ENDO GIA™ Universal

3.1.6 darbinė dalis sukasi apie savo ašį 360°.



0° / 22° / 45°

3.1.5 siuvimo aparatas pasuka kasetę iki 45° į abi puses

ENDO GIA™ Universal

Key Features	Key Benefits
60mm Articulating loading units	Allows the surgeon to maintain similar stapling techniques as open surgery by using 60mm staple lines
Universal Handle for all loading units	Helps reduce procedure costs; surgeons may not need a second instrument in a procedure
Six rows of staples with all sizes	Staple line security
New sharp and sterile knife with each reload <small>3.1.8 p.d su tuo pačiu aparatu galima panaudoti įvairaus ilgio/kabučių aukščio pasisukančias kasetes, kuriose integruotas vienkartinis peilis.</small>	Never dull from cutting through staple lines, and no risk of cross contamination
4.8 mm staples	Hemostasis in thick tissue
Tissue Gap Control	Ensures precise staple formation along the entire staple line

ENDO GIA™ Universal

Key Features	Key Benefits
25 firings per instrument <i>3.1.3 p.d aparatas gali būti užtaisytas/“iššautas“ 25 kartus.</i>	Reduces procedure costs, most surgeons will never need more than one instrument
22 and 45 degree SULU articulation	Can reduce the need for additional trocars to gain access to difficult surgical sites
Wide jaw opening available for 30mm and 45mm sizes as well as 60mm	Easy placement on tissue
Lock-out feature	The instrument won't fire, if it has already been used. It won't fire if tissue thickness is well beyond color standards

ENDO GIA™ Universal

3.1.7 p.d Kasetė uždaroma/atidaroma viena ranka, kilpinę rankeną spaudžiant/pastumiant nuo savęs, galimas kasetės daugkartinis atidarymas/uždarymas prieš šovimą

