

# Pleuracan®

Set for closed pleura and thorax drainage



Pleura Drainage Systems

**B | BRAUN**  
SHARING EXPERTISE


# Pleuracan®

## Pleuracan®


Accumulations of air and/or fluid within the pleura-, pericardial or abdominal space can cause life-threatening disorders of the respiration as well as the cardiac function.

Therefore, the objective of every emergency therapy is to guarantee remedy for the air or fluid. This used to be done either by needle puncture or by placing a sterile, surgical suction drainage.

Large lumen catheters (so called Bülau Drainages) are to be avoided whenever possible, from the costs' point of view and the advantages of a "minimal invasive surgery" for the patient.



## Pleura Drainage Systems



### Presentations

	Code No. (REF)	Sales unit/Pcs.
<b>Pleuracan® A, Complete Set</b>	446 2556	10
<ul style="list-style-type: none"> <li>■ Thin-walled puncture needle with short bevel 3.35 x 78 mm</li> <li>■ Certon® catheter 2.7 x 450 mm</li> <li>■ Attachment needle</li> <li>■ Double anti-reflux valve</li> <li>■ Secretion bag 2.0 l treated with EO</li> <li>■ Omnifix® lock, 60 ml</li> <li>■ Discofix® three-way stopcock</li> </ul>		
<b>Pleuracan® B, Basic Set</b>	446 2505	10
<ul style="list-style-type: none"> <li>■ Thin-walled puncture needle with short bevel 3.35 x 78 mm</li> <li>■ Certon® catheter 2.7 x 450 mm</li> <li>■ Attachment needle</li> <li>■ Double anti-reflux valve</li> <li>■ Discofix® three-way stopcock</li> </ul>		
<b>Secretion bag</b>		
■ 2.0 l treated with EO	442 5200	250
■ 2.0 l single packed, sterile	442 3208	100

### Available on request

	Code No. (REF)	Sales unit/Pcs.
<b>Pleuracan® double anti-reflux valve</b>	446 2564	20

## Pleura Drainage Systems

### User benefits

- Universally applicable for drainage of discharged fluids
- Atraumatic during puncture and the course of therapy
- Connectable to a secretion bag or a thorax drainage system



**B | BRAUN**  
SHARING EXPERTISE

**B. Braun Melsungen AG**

34209 Melsungen

Germany

Tel +49(0) 56 61 71-0

[www.bbraun.com](http://www.bbraun.com)