

1 p.d. Bendras plotis 210mm,tinklelio zonos 200mm.
Juodo stačiakampio dydis 6x7mm. Suderinamas su
Zoncare iMAC120

ZONCARE iMAC 120
BIOCARE 1200 / 1210 / 1215 / 1230
SONOSCAPE IE 12 / IE 15

210 x 140 x 142

REORDER CODE 13108

PROPER STORAGE OF THERMAL RECORDING PAPER

Store in a cool, dry place. The finished ECG record should not be exposed to direct sunlight, temperatures over 100°F (38°C), relative humidities over 80% or placed in contact with adhesives, adhesive tapes or plasticizers such as those found in all PVC page protectors.

CE
0051

EDAN SE-1201

ZONCARE iMAC 120

BIOCARE 1200 / 1210 / 1215 / 1230

SONOSCAPE IE 12 / IE 15

210 x 140 x 142

REORDER CODE 13108

PROPER STORAGE OF THERMAL RECORDING PAPER

cool, dry place. The finished ECG record should not be exposed to temperatures over 100°F (38°C), relative humidities over 80% or placed in contact with adhesive tapes or plasticizers such as those found in all PVC page protectors.

