

# SPECIFICATIONS

CapSureFix Novus MRI™ SureScan™  
4076 Lead

## MR Conditional



Initial Implant  
Existing Implant

## Physical Characteristics

Polarity	Bipolar <b>Bipolaris</b>
Shape	Straight
Chamber	Ventricle or Atrium
MR-Conditional Lengths	35, 45, 52, 58, 65, 85 cm
Connector	IS-1 BI

## Materials

Insulator	Polyurethane 55D (outer), Silicone with Siloxane™ treatment (inner)
Conductor	MP35N Nickel Alloy
Helix Electrode	Platinum Alloy with Titanium Nitride coating
Ring Electrode	Platinum Alloy with Titanium Nitride coating
Connector Ring	Stainless steel
Connector Pin	Stainless steel

## Diameter

Body	1.9 mm (5.7 Fr)
------	-----------------

## Recommended Introducer Size **Galima įvesti naudojant 7F įvediklį**

Without Guide Wire	7 Fr
With Guide Wire	9 Fr

## Electrode Configuration **Pasirenkamos fiksacijos**

Extendable/Retractable Helix Screw
------------------------------------

## Electrode Surface Area

Helix	4.2 mm <sup>2</sup>
Ring	22 mm <sup>2</sup>



## Other Characteristics **Iskiriantis gliukokortikoidus**

Steroid	Dexamethasone acetate
Tip-to-Ring Spacing	10 mm
Unipolar Conductor Resistance	34 $\Omega$ (58 cm)
Bipolar Conductor Resistance	62 $\Omega$ (58 cm)
Serial Prefix	BBL

## Stylets

Inserted	Gray Straight
Packaged	Gray Straight Blue Straight Gray J-shaped* Blue J-shaped* White J-shaped*

\*J-shaped stylets are not available for leads 65 or 85 cm.

Medtronic

### **Brief Statement**

See the device manual for detailed information regarding the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events.



[www.medtronic.com/manuals](http://www.medtronic.com/manuals)

Consult instructions for use on this website. Manuals can be viewed using a current version of any major Internet browser. For best results, use Adobe Acrobat® Reader with the browser.

#### **Europe**

Medtronic International Trading Sàrl.  
Route du Molliat 31  
Case postale  
CH-1131 Tolochenaz  
[www.medtronic.eu](http://www.medtronic.eu)  
Tel: +41 0 21 802 70 00  
Fax: +41 0 21 802 79 00

#### **United Kingdom/Ireland**

Medtronic Limited  
Building 9  
Croxley Green Business Park  
Hatters Lane  
Watford  
Herts WD18 8WW  
[www.medtronic.co.uk](http://www.medtronic.co.uk)  
Tel: +44 0 1923 212213  
Fax: +44 0 1923 241004