



23. P.D.



FORMOKREZOL

Product category: **Dentistry and Endodontics**
Producent: Chema-Elektromet

FOR PROFESSIONAL USE ONLY

Packaging available

Bottle containing 10 g of liquid and a dropper.

Pakuotė: 10 g buteliuke

About the product

Application

Storage conditions

Application

It helps in removing dead tissue from instruments. Į sudėtį įeina formaldehido ir krezolio, kurie turi dezinfekuojančių ir antibakterinių savybių.

Contains: formaldehyde and cresol, which have disinfecting and antibacterial properties.

Use: Apply a small amount of the liquid on a cotton pad with a dropper and clean dirty endodontic instruments. Then the tools should be subjected to standard cleaning and disinfection.

Dokumentą elektroniniu parašu pasirašė

Data: 2025-03-06 10:04:15

Paskirtis: Pirkimo nr. 925806

Vieta: Energetikų g. 8, Kaunas

Kontaktinė informacija: Viešųjų pirkimų specialistė