

**20. P.D.**



RESTORATIVE

## Nexcomp

Nano Hybrid Dental Composite Resin

### FEATURE & BENEFIT



**Low 1.6% shrinkage  
for superior marginal  
integrity**



**Various shade**  
19 shades available



**Filled with Various  
nano size filler for  
excellent  
polishability**



**Optimum handling**  
Excellent Sculptability for  
restoring natural shape  
and function



**Non-Slumping**



**Non-Sticky**

### APPLICATION



**Class I,II,III,IV,V  
restorations**



**Core-buildup to replace  
missing tooth structure**



**Diastema  
closures**



**Direct veneers**



**Composite and  
porcelain repairs**

## SPECIFICATION

### Principle organic component of polymer base

Bis-GMA, Bis-EMA, Urethane dimethacrylate, Triethylene glycol dimethacrylate

### Range of dimensions of inorganic filler particles

0.01 to 0.7  $\mu\text{m}$

### percentage by volume of total inorganic fillers

wt 75%

### Recommended external energy source

Šviesoje kietėjančios

#### Specification required of Light curing

**device:** usual dental halogen/LED lamps with light intensity of more than  $400\text{mW}/\text{cm}^2$  and wavelength range of 430 – 480nm

## Nexcomp Refill



### Galimos spalvos:

x 1

#### 4g x 1 Syringe

A2E, B2E, A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, D2, D3, TL, WT, OP, OA2, OA3

21. P.D.



RESTORATIVE

## Nexcomp Flow

Nano Hybrid Flowable Dental Composite Resin

Nano hibridinis takus kompozitas

### FEATURE & BENEFIT



Filled with Various nano size filler for excellent polishability



Excellent Radiopacity



Available in 14 shades including opaque and translucent



No slumping and running

### APPLICATION



Restoration of invasive cavity



Pits and fissure sealing



Base and Liner for direct restoration



Repair acrylic temporary materials



Repair of esthetic indirect restorations

## SPECIFICATION

### Principle organic component of polymer base

Bis-GMA, Bis-EMA, Urethane dimethacrylate, Triethylene glycol dimethacrylate

### percentage by volume of total inorganic fillers

wt 75%

### Range of dimensions of inorganic filler particles

0.01 to 0.7  $\mu\text{m}$

## Šviesoje kietėjančios

### Recommended external energy source

Specification required of Light curing

**device:** usual dental halogen/LED lamps with light intensity of more than  $400\text{mW}/\text{cm}^2$  and wavelength range of 430 – 480nm



## Nexcomp Flow Refill



×2

■ 2g × 2 Syringe

Galimos spalvos: A1, A2, A3, A3.5, A4, B1, B2, B3, OA2, OA3, C2, D2, WT, TL



■ 10EA disposable tip

Dokumenta elektroniniu parašu pasirašė

Data: 2025-03-06 10:03:11

Paskirtis: Pirkimo nr. 925806

Vieta: Energetikų g. 8, Kaunas

Kontaktinė informacija: Viešųjų pirkimų specialistė