

Tepalas įvairiems odontologiniams antgaliams


22. P.D.

## MAINTENANCE PRODUCTS FOR STRAIGHT/CONTRA-ANGLE HANDPIECES/TURBINES

### HS Maxima® Oil Spray

purškiamas

Universal cleaning and maintenance spray for straight and contra-angle handpieces, air, and ball-bearing turbines. CFC free. Does not include a spray attachment.

Content	Barcode	Art.-no.
500 ml (without adapter)		9001570

Pakuotė: 500 ml



### HS S-Spray





Universal spray for cleaning and maintaining turbines, straight and contra-angle handpieces. Ball-bearing turbines. CFC free. Does not include a spray attachment.

Content	Barcode	Art.-no.
250 ml		9003292



### HS Maxima® Oil Spray Adapter

Maxima Oil Spray Adapter is suitable for various couplings and turbines.

Content	Barcode	Art.-no.
For coupling with INTRAmatic (contra-angle handpiece)		9799531
For coupling with RotoQuick (W&H)		9799532
For coupling with Multiflex (KaVo).		9799533
For ISO couplings (universal)		9799534




## ULTRASONIC CLEANING

### HS Ultrasonic Cleaning Device sweep 2000 (4.0 liters)

The HS Ultrasonic Cleaning Device sweep 2000 is a piece of equipment with state-of-the-art technology in accordance with European Standards. Particularly noteworthy are the automatically locking lid, the larger drain, the hanging basket, which allows draining without manual intervention, and the digital control panel. The housing is easy to clean. The device is ideal for cleaning processes in dental practices and laboratories. The lid ensures a closed cleaning cycle and very low noise. Sweep technology, 4-liter capacity.

**Technical data:** Power consumption: 95 W, heater: 55 W, timer: 30 minutes digital, dimensions (W x H x D): 31.8 x 31 x 17.8 cm, net weight: approx. 4 kg


Dimensions	Barcode	Art.-no.
(W x H x D): 31.8 x 31 x 17.8 cm		9006375



### HS Plaster Solvent

Ready-to-use cleaner with high solvent power. HS Plaster Solvent removes plaster and investment containing gypsum and is suitable for use in ultrasonic baths.

**Caution:** Discoloration may occur when used on surfaces containing nickel.

Content	Barcode	Art.-no.
2 liters		9883777



Dokumentą elektroniniu parašu  
pasirašė

Data: 2025-03-06 10:03:44

Paskirtis: Pirkimo nr. 925806

Vieta: Energetikų g. 8, Kaunas

Kontaktinė informacija: Viešųjų  
pirkimų specialistė