


bfd156_C dup 300g 148*30*323 read



Antimicrobial wound dressings
with copper *with adhesive contour*

Physical Microbial Barrier


Copper impregnated antimicrobial wound dressings
designed to address a wide variety of wounds including:

- Diabetic wounds
- Leg & foot ulcers
- Pressure ulcers
- First & second degree burns
- Surgical wounds
- Venous and arterial ulcers

Ext: 10 cm x 25 cm

Int: 5 cm x 20 cm

10Pcs



7 290018 406175

LBL-MSB-32-E, Rev 8, 14.11.22

- Do Not Use
If Sterile Package
Is Compromised
- Consult Instruction
For Use
- Use Once
And Discard
- Sterilized Using
Ethylene Oxide
- Keep
Dry
- Do Not Re-sterilize



REF: 2C-1025-01a

Store at room temperature

The maximum use of a single wound dressing is up to 7 days

Do not use if hypersensitivity to device ingredients is known

This device is to be used by or on the order of a physician or a


licensed practitioner. Rx only

*Note: OTC for superficial and minor wounds.

Apply the orange side of the wound dressings to the wound site



0425



MedCu Technologies Ltd.
11 Galgalei Haplada St.,
Herzliya 4672200, Israel
+972-9-7429242
info@medcu.com
www.medcu.com

Patented under one or more as follows:

U.S. Pat. No. 9403041
E.U. Pat. No. 1809306
China Pat. No. 2587036
Japan Pat. No. 5275629
S.A. Pat. No. 2007/03628
Korea Pat. No. 1528918



MedNet EC-REP GmbH
Borkstrasse 10
48163 Muenster
Germany