

## Kanghua Medical Devices

### Single – use catheter syringe

- 3-part
  - Single-use, sterile
  - Non-toxic, non-pyrogenic
  - Made of polypropylene
  - Latex-free
  - Clear, legible and permanent black graduating facilitating dosage
  - Hermetic plunger
  - Transparent barrel
- Size: 50 ml, 60 ml, 100 ml, 120 ml, 140(150)ml, 200 ml



## Kanghua Medical Devices

### Kateterio tipo ir perfuziniai švirkštai

- 3-jų dalių
  - Vienkartiniai, sterilūs
  - Netoksiški ir nepirogeniški
  - Pagaminti iš polipropileno
  - Su arba be latekso
  - Su aiškia ir antiseptikams atparia gradacija
  - Su hermetišku stūmokliu
  - Skaidraus korpuso
- Dydžiai: 50 ml, 60 ml, 100 ml, 120 ml, 140(150)ml, 200 ml

įpakuoti po 1vnt



**江苏康华医疗器材有限公司**  
**JIANGSU KANGHUA**  
**MEDICAL EQUIPMENT CO. LTD.**



Switch board(总机): 0086 519 88671038

Tel (电话): 0086 519 88670071, 88676715, 88676716

Fax (传真): 0086 519 88678968

E-mail: kanghua@chinakang-hua.com

Http://www.kanghuamedical.com

Add: Sanhekou, Changzhou, Jiangsu, 213115, China

地址: 中国, 江苏, 常州, 三河口, 邮编:  
213115

To whom it may concern

**Statement**

We, Jiangsu Kanghua Medical Equipment Co. LTD. (Add: Sanhekou, Changzhou, Jiangsu, 213115, China) hereby declare that all kind of syringes manufactured by us are sterile, non-toxic, hermetic and without latex components . visi švirkštai be latekso komponentų

Date: March 19, 2021

/sign/

*Shen Zhen*  
江苏康华医疗器材有限公司

/stamp/

JIANGSU KANGHUA MEDICAL EQUIPMENT CO. LTD.



**江苏康华医疗器材有限公司**  
**JIANGSU KANGHUA**  
**MEDICAL EQUIPMENT CO. LTD.**



Switch board(总机): 0086 519 88671038

Tel (电话): 0086 519 88670071, 88676715, 88676716

Fax (传真): 0086 519 88678968

E-mail: kanghua@chinakang-hua.com

Http://www.kanghuamedical.com

Add: Sanhekou, Changzhou, Jiangsu, 213115, China

地址: 中国, 江苏, 常州, 三河口, 邮编:  
213115

To whom it may concern

**STATEMENT**

We, **Jiangsu Kanghua Medical Equipment Co. LTD.** (Adress: Sanhekou, Changzhou, Jiangsu, 213115, China), **declare that our manufactured 50ml(60ml), 100ml, 150ml CatheterTip syringes meet this standard ISO 7886-1.**

Date: 31<sup>st</sup> August 2022

/stamp/



/sign/

**江苏康华医疗器材有限公司**  
**JIANGSU KANGHUA**  
**MEDICAL EQUIPMENT CO. LTD.**



Switch board(总机): 0086 519 88671038

Tel (电话): 0086 519 88670071, 88676715, 88676716

Fax (传真): 0086 519 88678968

E-mail: kanghua@chinakang-hua.com

Http://www.kanghuamedical.com

Add: Sanhekou, Changzhou, Jiangsu, 213115, China

地址: 中国, 江苏, 常州, 三河口, 邮编:  
213115

Kam tai gali būti aktualu,

**PREIŠKIMAS**

Mes, **Jiangsu Kanghua Medical Equipment Co. LTD.** (Adresas: Sanhekou, Changzhou, Jiangsu, 213115, Kinija), **patvirtiname, kad mūsų pagaminti 50 ml (60 ml), 100ml, 150 ml kateterio antgalio švirkštai atitinka šį standartą ISO 7886-1.**

Data: 2022m. rugpjūčio 31d.

Antspaudas//

Parašas//