

# Medtronic

The Signia™  
small diameter  
reload

75%  
greater access  
than 12 mm  
reloads<sup>2,3</sup>



3.1 ilgis 30mm

30 mm  
& 45 mm

jaw length with curved tip

6.3 mm

anvil

45 mm

4

rows of staples  
Comparable hemostasis  
with 6 row staplers.<sup>4,§,Ω</sup>

3.2 Kasetėje 4 eilės kabučiu

# A narrower profile.<sup>1,†</sup> For broader possibilities.

3.5 Kasetės darbinės dalies skersmuo 8 mm (matuojama nuo kasetės sąnario link rankenos), kasetės skersmuo 6,3 mm

8 mm

working shaft

- † Compared to the Ethicon Echelon Flex™\*\* powered vascular stapler with advanced placement tip.  
‡ Compared to the Echelon Flex™\*\* powered vascular stapler 35 mm white 2.5 mm cartridge reloads.  
§ Preclinical results may not correlate with clinical performance in humans.  
Ω Tested using the Endo GIA™ handle with Signia™ small diameter 30 mm grey 2.0 mm reload versus Ethicon Echelon Flex™\*\* powered vascular stapler with 35 mm white 2.5 mm reload. Preclinical results may not correlate with clinical performance.  
1. Based on report #RE00142825, Image creation for Signia™ small diameter reload. March 26, 2019.  
2. Based on report #RE00142825, Image creation for Signia™ small diameter reload. March 26, 2019.  
3. Based on internal report #RE00014682, Small diameter reload summative usability. June 19-20, 22, 2017.)  
4. Based on internal report #RE00190239, Signia™ small diameter reload chronic study statistical analysis. March 26, 2019.



Important: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.  
© 2023 Medtronic. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.  
All other brands are trademarks of a Medtronic company. EMEA-ST-2200043-signia-sd-infographic-en-we-7342569

[medtronic.com/covidien/uk](https://www.medtronic.com/covidien/uk)