


## SCBI for STEAM, Plasma & Ethylene-Oxide

SAP / Monitoring



### Quickly monitoring ..... with SCBI

The monitoring of sterilization condition at preferred places with the help of Self Contained Biological Indicators (or SCBI) is no longer practice of every day in European Hospitals and mainly used during validation of the sterilizers, whereas outside Europe the SCBI is widely used for daily routine checks of sterilizers in Hospitals.

REF: 107.001.0100	SCBI for STEAM, pop. 1x10E6	Unit: 100 pieces
REF: 107.002.0100	SCBI for Plasma, pop. 1x10E6	Unit: 100 pieces
REF: 107.003.0100	SCBI for EO, pop. 1x10E6	Unit: 100 pieces
REF: 120.001.0001	Dual Incubator for SCBI	Unit: 001 piece

#### Characteristics:

- ◆ Biological Indicators conform to ISO 11138
- ◆ SCBI EO: Bacillus Atrophaeus ATCC 9372 1x10E6 spores per vial
- ◆ SCBI STEAM/Plasma: Geobacillus Stearothermophilus ATCC 7953 1x10E6 spores per vial
- ◆ Color change: Green to yellow (EO) and purple to yellow (STEAM)
- ◆ Digital Incubator (220V) for 37°C / 48 hrs (EO) and 56-59°C / 24 hrs (STEAM & Plasma)

