





# CATALOGUE 2022



Biosense Webster®

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

## **Biosense Webster, a division of Johnson & Johnson Medical NV/SA**

---

Biosense Webster (BW) for EMEA is recognized worldwide as a leader in the technology for the diagnosis and treatment of cardiac arrhythmias.

The first deflectable tip catheter was pioneered by Will Webster more than three decades ago. Ever since, Biosense Webster has defined the state of design and craftsmanship in diagnostic and therapeutic cardiac catheters. To date, we offer more electrophysiology catheter designs than anyone else in the industry. Beyond the numbers, our continuous pursuit of practical innovation has led to a steady flow of breakthrough designs.

At the same time, Biosense Webster has delivered numerous mapping and navigation innovations to advance the science and practice of electrophysiology. Ongoing research and development at Biosense Webster has given electrophysiologists access to 3D, real-time, color-coded views into the heart's electrical activity, allowing diagnosis and treatment with increasing precision. Strong collaboration with leaders in related industries is pushing the clinical boundaries still further, thanks to such technologies as the CARTO™ 3 System, the catheters and software modules of the QDOT MICRO™ Technology, the CARTO VIZIGO™ Bi-Directional Guiding Sheath, the HELIOSTAR™ Balloon Technology as well as the CARTO PRIME™ Module with PENTARAY™ Catheter and OCTARAY™ Mapping Catheter.

Guided by the needs of electrophysiologists and cardiologists, Biosense Webster has pioneered innovative technological advancements for more than thirty years. Some might call this quite a heritage. At Biosense Webster, we say "It's just the beginning".

We are proud to present you the catalogue for 2022.

# Table of Contents

## Curve Types / Tips and Connectors

<b>Fixed Curves</b>	
Diagnostic Catheters .....	10
<b>Deflectable Curves</b>	
Special Diagnostic Catheters .....	11
Uni-Directional Diagnostic and Ablation Catheters .....	12
Bi-Directional Diagnostic and Ablation Catheters .....	13
<b>Tip Electrodes</b> .....	15
<b>Connectors</b>	
Connectors for CELSIUS™ Ablation Catheter .....	16
Connectors for NAVISTAR™ Catheter .....	16
Connectors to Generators .....	17

## Therapeutic Catheters

<b>Focal - Navigation Therapeutic Catheters</b>	
QDOT MICRO™ Uni-directional Catheter .....	20
QDOT MICRO™ Bi-directional Catheter .....	20
THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter .....	22
THERMOCOOL™ SF NAV Uni-Directional Catheter .....	22
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	23
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization .....	24
THERMOCOOL™ SF NAV Bi-Directional Catheter with curve visualization .....	24
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	25
NAVISTAR™ 4 mm Catheter .....	26
NAVISTAR™ DS Catheter .....	26
NAVISTAR™ RMT THERMOCOOL™ Catheter .....	27
NAVISTAR™ RMT 4 mm Catheter .....	27
<b>Focal - Non-Navigation Therapeutic Catheters</b>	
THERMOCOOL™ SF Uni-Directional Catheter .....	29
THERMOCOOL™ SF Bi-Directional Catheter .....	29
EZ STEER™ THERMOCOOL™ Catheter .....	30
EZ STEER™ 4 mm Bi-Directional Catheter .....	30
EZ STEER™ DS Bi-Directional Catheter .....	31
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter .....	31
CELSIUS™ RMT THERMOCOOL™ Catheter .....	32
CELSIUS™ 4 mm Catheter Thermocouple .....	32
CELSIUS™ 4 mm Catheter Thermistor .....	33
CELSIUS™ 4 mm Braided Tip Catheter .....	33
CELSIUS FLTR™ 8 mm Uni-Directional Catheter .....	34
CELSIUS FLTR™ 8 mm Bi-Directional Catheter .....	34
CELSIUS™ DS Catheter .....	35
CELSIUS™ RMT Catheter .....	35
<b>Balloon</b>	
HELIOSTAR™ Balloon Ablation Catheter .....	36
LASSOSTAR™ Catheter .....	36
LASSOSTAR™ NAV Catheter .....	36
GUIDESTAR™ Sheath .....	36

## Diagnostic Catheters

<b>Navigation Diagnostic Catheters</b>	
LASSO™ 2515 NAV eco Circular Mapping Catheter .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	38
OCTARAY™ Mapping Catheter .....	39
PENTARAY™ NAV eco High Density Mapping Catheter .....	39
DECANAV™ Catheter .....	40
CARTO™ FAM DX .....	40
<b>Auto-ID Diagnostic Catheters</b>	
LASSO™ Decapolar Catheter .....	41
LASSO™ Duo-Decapolar Catheter .....	41
WEBSTER™ CS Uni-Directional Catheter .....	42
WEBSTER™ CS Bi-Directional Catheter .....	42
WEBSTER™ Quadrapolar Catheter Deflectable .....	43
WEBSTER™ Decapolar Catheter Deflectable .....	43
WEBSTER™ HIS Catheter .....	44
WEBSTER™ Quadrapolar Catheter Fixed .....	44
WEBSTER™ Decapolar Catheter Fixed .....	44
HALO™ XP Catheter .....	45
ISMUS™ Catheter .....	45
<b>Non Auto-ID Diagnostic Catheters</b>	
LASSO™ 2515 Variable Mapping Catheter .....	46
LASSO™ Mapping Catheter .....	46
WEBSTER™ CS Bi-Directional Catheter .....	47
WEBSTER™ Octopolar Catheter Deflectable .....	47
WEBSTER™ Hexapolar Catheter Deflectable .....	48
WEBSTER™ Quadrapolar Catheter Deflectable .....	48
WEBSTER™ Decapolar Catheter Deflectable .....	49
WEBSTER™ Decapolar Catheter Fixed .....	50
WEBSTER™ Hexapolar Catheter Fixed .....	50
WEBSTER™ Quadrapolar Catheter Fixed .....	51

WEBSTER™ Duo-Decapolar Catheter .....	52
HALO™ XP Tricuspid Mapping Catheter .....	52
ISMUS™ Catheter .....	52
AVAIL™ Bipolar Fixed Shape Catheter .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	53
CRISTACATH™ Catheter .....	54
PARAHISIAN™ Catheter .....	54
Pediatric Catheter .....	54

## Access

<b>Guiding Sheaths - Fixed Guiding Sheaths</b>	
PREFACE™ Fixed Guiding Sheaths for EP .....	56
PREFACE™ Fixed Guiding Sheaths .....	56
CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath .....	57
<b>Guiding Sheaths - Steerable Guiding Sheaths</b>	
CARTO VIZIGO™ Sheaths .....	58
MOBICATH™ Sheaths .....	58
<b>Transseptal Needles</b>	
HEARTSPAN™ Transseptal Needles .....	59

## Implantable Devices

<b>The COHEREX WAVECREST™ LAA Implant</b>	
Device Description .....	61
Implant Characteristics .....	61
COHEREX WAVECREST™ LAA Implant .....	62
Sheath ordering information .....	62

## UltraSound Catheters

<b>Ultrasound Catheters</b>	
SOUNDSTAR™ 3D eco Catheter .....	64
SOUNDSTAR™ eco 8F Catheter .....	64
SOUNDSTAR™ 3D Ultrasound Catheter .....	65
ACUNAV™ 10F Ultrasound Catheter .....	65
ACUNAV™ 8F Ultrasound Catheter .....	66
<b>System Compatibility</b>	
CARTOSOUND® Module System Compatibility* .....	67
ACUNAV™ Catheter System Compatibility .....	67

## RF Generators and Pump Systems

<b>nGEN™ GENERATOR PLATFORM</b>	
nGEN™ GENERATOR PLATFORM .....	69
<b>SMARTABLATE™ System</b>	
SMARTABLATE™ System .....	70
SMARTABLATE™ System Diagram .....	71
SMARTABLATE™ System Service Contracts .....	72
<b>Stockert EP Shuttle Generator</b>	
Indifferent Electrodes and Corresponding Connection Cables .....	73
Spare Parts and Others .....	74
Stockert EP Shuttle Accessories .....	75
Service Contracts .....	75
<b>COOLFLOW™ Pump</b>	
COOLFLOW™ Pump Tubing <sup>1</sup> .....	76
COOLFLOW™ Pump Accessories .....	76
Diagram - Fully Integrated COOLFLOW™ Pump System .....	77

## Systems

<b>CARTO™ 3 System</b>	
CARTO™ 3 System .....	79
CARTO™ 3 System Accessories .....	79
CARTO™ 3 System Software Module .....	81
CARTO™ 3 System Service Contracts .....	83-84
CARTO™ 3 System External Reference Patch .....	85

## Cables

<b>Charts</b>	
CONNECTION OF CHOICE™ Chart .....	87-88
<b>Cables</b>	
CARTO™ 3 System Cables .....	89
Interface Cable EP Shuttle RF Generator .....	90
Interface SMARTABLATE™ System Cables .....	90
<b>Non-Navigation</b>	
Non-Navigation Diagnostic Cables .....	91
Non-Navigation Ablation Cables .....	91
<b>Reverse Cables</b>	
Reverse Cables .....	92

## Index


Please contact your local BW representative for availability.

\* Please contact your Biosense Webster representatives for further details on exact compatible versions.


<sup>1</sup> Tubing to be used with COOLFLOW™ Pump.





## THERAPEUTICS



# DIAGNOSTICS







\* Trademark of Siemens Medical Solutions USA, Inc.

# Addresses

---

## Austria

Johnson & Johnson  
Medical Products GmbH  
Vorgartenstraße 206B  
1020 Vienna, Austria  
Tel: +43 1360250

## Belgium

European HQ  
Biosense Webster  
Leonardo Da Vincilaan 15  
Diegem 1831, Belgium  
Tel: +32 27463401

## Czech Republic

Johnson & Johnson, s.r.o.  
Karla Engliše 3201/6  
15000 Praha 5, Czech Republic  
Tel: +420 227012669

## Denmark

Filial af Johnson & Johnson AB  
Bregnerødvej 133  
DK 3460 Birkerød, Denmark  
Tel: +45 45948200

## Dubai

Johnson & Johnson (Middle East), Inc  
Dubai Health Care City  
Building 14, Level 6  
PO Box 505080  
Dubai, UAE  
Tel: +971 (0)4 4297 200

## Egypt

Johnson & Johnson (Medical Egypt), Inc  
Florida Mall, Sheraton Heliopolis  
Level 5  
1229 El Sheikh Aly Gad El Hak Square  
Cairo, Egypt  
Tel: +2 (0)222 685 026

## Finland

Johnson & Johnson Finland  
Vaisalantie 2  
FI - 02130 Espoo, Finland  
Tel: +358 207531500

## France

Biosense Webster France  
Rue Camille Desmoulins TSA 71001  
92787 Issy-les-Moulineaux, France  
Tel: +33 155003756

## Germany

Cordis Medizinische Apparate GmbH  
Biosense Webster  
Robert-Koch-Straße 1  
22851 Norderstedt, Germany  
Tel: +49 1802999112

## Greece

Johnson & Johnson Greece  
Aegialias & Epidavrou 4, Maroussi  
15125 - Athens, Greece  
Tel: +30 2106875500

## Hungary

Johnson & Johnson  
Krisztina Palace  
1123 Budapest  
Nagyenyed u. 8-14, Hungary  
Tel: +36 23513800

## Ireland

Biosense Webster, Johnson & Johnson  
(Ireland) Limited Airton Road,  
Tallaght, Dublin 24, Ireland  
Tel: +353 14665200

## Israel

Johnson & Johnson Israel  
Kibbutz Shefayim, Kibbutz Shefayim  
60990, Israel  
Tel: +972 99601722

## Italy

Johnson & Johnson Medical  
S.p.A., BIOSENSE WEBSTER  
Via del Mare 56  
00040 Pomezia, Roma, Italy  
Tel: +39 06911941

## Latvia

Johnson & Johnson AB Latvia branch  
Mukusalas street 101  
Riga, LV-1004, Latvia  
Tel: +371 67103407

## The Netherlands

Johnson & Johnson Medical BV  
Computerweg 14, 3821 AB  
Amersfoort, Postbus 188, 3800 AD  
Amersfoort, The Netherlands  
Tel: +31 334500739

## Norway

Johnson & Johnson AB  
PB 144, 1325 Lysaker  
Drammensveien 288 N-0283, Norway  
Tel: +47 66981030

## Poland

Johnson & Johnson Poland  
ul. Ilzecka 24  
02135 Warszawa, Poland  
Tel: +48 667607139

## Portugal

Johnson & Johnson Medical Iberia  
Biosense Webster  
Estrada Consiglieri Pedroso, 69 A  
Queluz de Baixo  
2730-055 Barcarena, Portugal  
Tel: +351 214368835

## Russia

Johnson & Johnson Russia  
Krylatskaya Str. 17, bld. 3, floor 2  
121614 Moscow, Russia  
Tel: +7 4955807777

## Saudi Arabia

Johnson & Johnson Medical KSA  
Al Rosais commercial center  
appt. No. 308, Olaya st.  
PO Box 65305 Riyadh 11556, KSA  
Tel: +966 38874790

## Slovenia

Johnson & Johnson SE  
Smartinska cesta 53  
1000 Ljubljana, Slovenia  
Tel: +386 140118800

## South-Africa

Biosense Webster  
2 Medical street, Halfway House  
Johannesburg, 1685, South Africa  
Tel: +27 112651000

## Spain

Johnson & Johnson Medical Iberia  
Biosense Webster  
Pº de las Doce Estrellas 5-7  
Campo de Las Naciones  
28042 Madrid, Spain  
Tel: +34 917228262

## Sweden

Biosense Webster, Johnson & Johnson  
AB, Staffansväg 2  
191 84 Sollentuna, Sweden  
Tel: +46 86262200

## Switzerland

Johnson & Johnson AG,  
Biosense Webster  
Gubelstrasse 34  
6300 Zug, Schweiz  
Tel: +41 582312333

## Turkey

Johnson & Johnson Medical Turkey  
Ertürk Sokak Keceli Plaza No:13  
Kavacık Beykoz Istanbul, Turkey  
Tel: +90 2165382100

## UK

Biosense Webster, Johnson & Johnson  
Medical Ltd, Pinewood Campus  
Nine Mile Ride, Wokingham  
Berkshire RG40 3EW, United Kingdom  
Tel: +44 1344871000



# Curve Types / Tips and Connectors ● ● ● ● ● ●

## Curve Types / Tips and Connectors

### Fixed Curves

Diagnostic Catheters.....	10
---------------------------	----


### Deflectable Curves

Special Diagnostic Catheters .....	11
Uni-Directional Diagnostic and Ablation Catheters .....	12
Bi-Directional Diagnostic and Ablation Catheters.....	13

Tip Electrodes.....	15
---------------------	----

### Connectors

Connectors for CELSIUS™ Ablation Catheter .....	16
Connectors for NAVISTAR™ Catheter .....	16
Connectors to Generators.....	17





## Diagnostic Catheters




**A Curve**  
(Josephson Type / Yellow)



**D Curve**  
(Damato Type / Blue)



**F Curve**  
(Cournand Type / Black)




**G Curve**  
(Coronary Sinus / Orange)




**P Curve**  
(Coronary Sinus – From Above)


## Special Diagnostic Catheters




**HALO™ XP Catheter Curve**




**Posteroseptal D Curve**  
(Blue with Red band)




**ISMUS™ Cath Catheter Curve**  
(160° Type)



**LASSO™ Catheter Curve**  
(D Curves)



**PENTARAY™ Catheter Curve**  
(Shaft Deflection – D Curve)



**OCTARAY™ Mapping Catheter**  
(Shaft Deflection – D Curve)

## Uni-Directional Diagnostic and Ablation Catheters



The Biosense Webster unidirectional catheters are steered utilizing a push-pull mechanism and is well-known for its ease of use, reliability and immediate response.

Standard Insertion Lengths	
5F.....	80 cm
6F.....	92 cm
7F.....	115 cm
8F.....	90 cm

### A Curve

(Yellow / 3.8 cm)




### B Curve

(Red / 5.1 cm)




### C Curve

(Green / 6.4 cm)




### D Curve

(Blue / 6.4 cm)




### E Curve

(White / 6.4 cm)




### F Curve

(Orange / 7.6 cm)




### J Curve

(Black / 10.2 cm)




### K Curve

(Black / 7.6 cm)




## Bi-Directional Diagnostic and Ablation Catheters


EZ STEER™ Bi-Directional Catheter put amazing responsiveness, predictability, and comfort close at hand. Utilizing a patented Biosense Webster micro-tensioning system that enables micromovement abilities, they provide physicians with a new level of continuous control, smooth deflections, and procedural confidence.




**D-D Curve**  
(Blue / Blue)




**F-F Curve**  
(Orange / Orange)



**J-J Curve**  
(Black / Black)




**D-F Curve**  
(Blue / Orange)




**F-J Curve**  
(Orange / Black)




## Bi-Directional Diagnostic and Ablation Catheters




**B-B Curve**  
(Red / Red)



**B-D Curve**  
(Red / Blue)



**B-F Curve**  
(Red / Orange)




**D-J Curve**  
(Blue / Black)



## TIP ELECTRODES

**1 mm**


(Diagnostic Catheter)



**3.5 mm Irrigated**

(Ablation Catheter)


THERMOCOOL™ SF Catheter



**3.5 mm Irrigated**


(Ablation Catheter)

THERMOCOOL™ Catheter tip




**4 mm**

(Diagnostic and Ablation Catheter)



**8 mm Dual Sensor**

(Ablation Catheter)





## Connectors for CELSIUS™ Ablation Catheter


Connector on therapeutic EP Catheter

Connector on Connection Cable



Integrated Redel-10  
Standard connector for CELSIUS™ Ablation Catheter

Redel-10



Integrated Redel-10 with luer fitting  
CELSIUS™ THERMOCOOL™ Catheter

Redel-10

## Connectors for NAVISTAR™ Catheter

Connector on Catheter

Connector on Connection Cable

Connector on NAVISTAR™ Catheter



Hypertronics 25-pin  
for NAVISTAR™, NAVISTAR™ DS and REFSTAR™ Catheters

Hypertronics 25-pin



Hypertronics 25-pin with luer fitting  
for NAVISTAR™ THERMOCOOL™ Catheter

Hypertronics 25-pin

## Connectors to Generators\*

Connector


Generator



Redel-10 (THERMOCOOL™ Catheter and Thermistor)

Stockert  
EP Shuttle  
RF Generator

\* Refers to Stockert RF Generator connectors only.



## Therapeutic Catheters



<b>Focal - Navigation Therapeutic Catheters.....</b>	<b>20</b>
QDOT MICRO™ Uni-directional Catheter .....	20
QDOT MICRO™ Bi-directional Catheter .....	20
THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter.....	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	22
THERMOCOOL™ SF NAV Uni-Directional Catheter .....	22
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	23
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization .....	24
THERMOCOOL™ SF NAV Bi-Directional Catheter with curve visualization.....	24
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	25
NAVISTAR™ 4 mm Catheter .....	26
NAVISTAR™ DS Catheter .....	26
NAVISTAR™ RMT THERMOCOOL™ Catheter .....	27
NAVISTAR™ RMT 4 mm Catheter .....	27
<b>Focal - Non-Navigation Therapeutic Catheters .....</b>	<b>33</b>
THERMOCOOL™ SF Uni-Directional Catheter.....	35
THERMOCOOL™ SF Bi-Directional Catheter .....	35
EZ STEER™ THERMOCOOL™ Catheter.....	36
EZ STEER™ 4 mm Bi-Directional Catheter.....	36
EZ STEER™ DS Bi-Directional Catheter .....	37
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter .....	37
CELSIUS™ RMT THERMOCOOL™ Catheter .....	38
CELSIUS™ 4 mm Catheter Thermocouple .....	38
CELSIUS™ 4 mm Catheter Thermistor.....	39
CELSIUS™ 4 mm Braided Tip Catheter.....	39
CELSIUS FLTR™ 8 mm Uni-Directional Catheter .....	40
CELSIUS FLTR™ 8 mm Bi-Directional Catheter.....	40
CELSIUS™ DS Catheter.....	41
CELSIUS™ RMT Catheter .....	41
<b>Balloons .....</b>	<b>42</b>
HELIOSTAR™ Balloon Ablation Catheter.....	42
LASSOSTAR™ Catheter.....	42
LASSOSTAR™ NAV Catheter.....	42
GUIDESTAR™ Sheath.....	42

# Therapeutic Catheters



## NAV

<b>Focal - Navigation Therapeutic Catheters*</b> .....	<b>20</b>
QDOT MICRO™ Uni-directional Catheter .....	20
QDOT MICRO™ Bi-directional Catheter.....	20
THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter .....	21
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	22
THERMOCOOL™ SF NAV Uni-Directional Catheter .....	22
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	23
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization .....	24
THERMOCOOL™ SF NAV Bi-Directional Catheter with curve visualization.....	24
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	25
NAVISTAR™ 4 mm Catheter .....	26
NAVISTAR™ DS Catheter .....	26
NAVISTAR™ RMT THERMOCOOL™ Catheter .....	27
NAVISTAR™ RMT 4 mm Catheter .....	27

\* The size of the catheter varies depending on whether the diameter is measured at the shaft or at the electrodes. In the following pages, please note that the Outer Diameter reflects the diameter of the catheter measured at the ring electrodes.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## QDOT MICRO™ Uni-directional Catheter\*



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Spacing	Insertion Length
D139401	x	x	D ●	3.5 mm	1-6-2 mm	115 cm
D139402	x	x	F ●	3.5 mm	1-6-2 mm	115 cm
D139403	x	x	J ●	3.5 mm	1-6-2 mm	115 cm

## QDOT MICRO™ Bi-directional Catheter\*



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Spacing	Insertion Length
D139501	x	x	D-D ● ●	3.5 mm	1-6-2 mm	115 cm
D139502	x	x	F-F ● ●	3.5 mm	1-6-2 mm	115 cm
D139503	x	x	J-J ● ●	3.5 mm	1-6-2 mm	115 cm
D139504	x	x	FJ ● ●	3.5 mm	1-6-2 mm	115 cm
D139505	x	x	D-F ● ●	3.5 mm	1-6-2 mm	115 cm

### Associated Cable:

Catalog No.:	Connects
D140102	TX eco Cable (Dongle for QDOT MICRO™ catheter)
D135703	TX eco Extension Cable (for QDOT MICRO™ Catheter)

### Associated Tubing:

Catalog No.:	Connects
SAT001	SMARTABLATE™ System Irrigation Tubing Set

\* QDOT MICRO™ requires an nGEN™ generator and a software upgrade CARTO™ 3 system V7 Phase 2 Kit (KT-5400-230).

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter\*

UNI



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Spacing	Insertion Length
D134701	8F	≥8.5F	D ●	3.5 mm	1-6-2 mm	115 cm
D134702	8F	≥8.5F	F ●	3.5 mm	1-6-2 mm	115 cm
D134703	8F	≥8.5F	J ●	3.5 mm	1-6-2 mm	115 cm

## THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter\*

BI



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Spacing	Insertion Length
D134801	8F	≥8.5F	D-D ●●	3.5 mm	1-6-2 mm	115 cm
D134802	8F	≥8.5F	F-F ●●	3.5 mm	1-6-2 mm	115 cm
D134803	8F	≥8.5F	J-J ●●	3.5 mm	1-6-2 mm	115 cm
D134804	8F	≥8.5F	F-J ●●	3.5 mm	1-6-2 mm	115 cm
D134805	8F	≥8.5F	D-F ●●	3.5 mm	1-6-2 mm	115 cm

## THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter

UNI



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D133601	8F	≥8.5F	D ●	3.5 mm	6	1-6-2 mm	115 cm
D133602	8F	≥8.5F	F ●	3.5 mm	6	1-6-2 mm	115 cm
D133603	8F	≥8.5F	J ●	3.5 mm	6	1-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3434CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

\* THERMOCOOL™ SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.



# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter

BI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D132701	8F	≥8.5F	D-D ●●	3.5 mm	6	1-6-2 mm	115 cm
D132702	8F	≥8.5F	F-F ●●	3.5 mm	6	1-6-2 mm	115 cm
D132703	8F	≥8.5F	J-J ●●	3.5 mm	6	1-6-2 mm	115 cm
D132704	8F	≥8.5F	F-J ●●	3.5 mm	6	1-6-2 mm	115 cm
D132705	8F	≥8.5F	D-F ●●	3.5 mm	6	1-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3434CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

## THERMOCOOL™ SF NAV Uni-Directional Catheter\*

UNI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D131802†	8F	≥8F	D ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131803	8F	≥8F	F ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131804†	8F	≥8F	J ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

\* THERMOCOOL™ SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## THERMOCOOL™ SF NAV Bi-Directional Catheter\*

BI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
BBNI35BDH†	8F	≥8F	B-D ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35BFH†	8F	≥8F	B-F ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35DFH	8F	≥8F	D-F ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35FFH†	8F	≥8F	F-F ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BNI35FJH†	8F	≥8F	F-J ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

### Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

### Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

\* THERMOCOOL™ SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## THERMOCOOL™ SF NAV Uni-Directional Catheter\* with curve visualization

UNI



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D131504†	8F	≥8F	J ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3434CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle
D129704	PIU to SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

## THERMOCOOL™ SF NAV Bi-Directional Catheter\* with curve visualization

BI



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
BNI35BDCT†	8F	≥8F	B-D ●●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3434CT	Catheter to CARTO™ 3 System PIU

\* THERMOCOOL™ SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter

UNI



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
34H17M	8F	≥8F	B ●	3.5 mm	4	2-5-2 mm	115 cm
34H27M	8F	≥8F	C ●	3.5 mm	4	2-5-2 mm	115 cm
34H37M	8F	≥8F	D ●	3.5 mm	4	2-5-2 mm	115 cm
34H57M	8F	≥8F	F ●	3.5 mm	4	2-5-2 mm	115 cm
34HJ7M	8F	≥8F	J ●	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

## NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter

BI



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36H33M	8F	≥8F	D-D ●●	3.5 mm	4	2-5-2 mm	115 cm
36H35M	8F	≥8F	D-F ●●	3.5 mm	4	2-5-2 mm	115 cm
36H55M†	8F	≥8F	F-F ●●	3.5 mm	4	2-5-2 mm	115 cm
36H5JM	8F	≥8F	F-J ●●	3.5 mm	4	2-5-2 mm	115 cm
36HJJM	8F	≥8F	J-J ●●	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## NAVISTAR™ 4 mm Catheter

UNI



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
34A15M	7F	B ●	4 mm	4	1-7-4 mm	115 cm
34A25M†	7F	C ●	4 mm	4	1-7-4 mm	115 cm
34A35M	7F	D ●	4 mm	4	1-7-4 mm	115 cm
34A45M†	7F	E ○	4 mm	4	1-7-4 mm	115 cm
34A55M	7F	F ●	4 mm	4	1-7-4 mm	115 cm
34AJ5M†	7F	J ●	4 mm	4	1-7-4 mm	115 cm
34B35M	7F	D-F ● ●	4 mm	4	1-7-4 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

## NAVISTAR™ DS Catheter

UNI



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
34J27M	7F	C ●	8 mm	4	1-7-4 mm	115 cm
34J37M	7F	D ●	8 mm	4	1-7-4 mm	115 cm
34J57M	7F	F ●	8 mm	4	1-7-4 mm	115 cm
34JJ7M	7F	J ●	8 mm	4	1-7-4 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E68R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NAVIGATION THERAPEUTIC CATHETERS

## NAVISTAR™ RMT THERMOCOOL™ Catheter



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
NR7TCSIY	8F	≥8.5F	NA	3.5 mm	4	2-5-2 mm	130 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
D129704	PIU to SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

## NAVISTAR™ RMT 4 mm Catheter



1-7-4 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit*	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
NR7TCS4Y	8F	≥8.5F	4 mm	4	1-7-4 mm	125 cm

Associated Cable:

Catalog No.:	Connects
CR3425CT	Catheter to CARTO™ 3 System PIU
39E43R	PIU to Stockert EP Shuttle Generator
39E58M	Catheter to CARTO™ XP System

\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.



# Therapeutic Catheters ○○○○○○

## NON-NAV

<b>Focal - Non-Navigation Therapeutic Catheters .....</b>	<b>29</b>
THERMOCOOL™ SF Uni-Directional Catheter.....	29
THERMOCOOL™ SF Bi-Directional Catheter.....	29
EZ STEER™ THERMOCOOL™ Catheter.....	30
EZ STEER™ 4 mm Bi-Directional Catheter.....	30
EZ STEER™ DS Bi-Directional Catheter .....	31
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter .....	31
CELSIUS™ RMT THERMOCOOL™ Catheter .....	32
CELSIUS™ 4 mm Catheter Thermocouple .....	32
CELSIUS™ 4 mm Catheter Thermistor.....	33
CELSIUS™ 4 mm Braided Tip Catheter.....	33
CELSIUS FLTR™ 8 mm Uni-Directional Catheter .....	34
CELSIUS FLTR™ 8 mm Bi-Directional Catheter.....	34
CELSIUS™ DS Catheter .....	35
CELSIUS™ RMT Catheter .....	35
<b>Balloons.....</b>	<b>36</b>
HELIOSTAR™ Balloon Ablation Catheter .....	36
LASSOSTAR™ Catheter.....	36
GUIDESTAR™ Sheath.....	36

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## THERMOCOOL™ SF Uni-Directional Catheter\*

UNI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D131602†	8F	≥8F	D ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131603	8F	≥8F	F ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
D131604†	8F	≥8F	J ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing

## THERMOCOOL™ SF Bi-Directional Catheter\*

BI



European Ordering No.	Outer Diameter	Sheath Fit**	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
BDI35BDRT†	8F	≥8F	B-D ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35DDRT†	8F	≥8F	D-D ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35DFRT	8F	≥8F	D-F ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm
BDI35FJRT†	8F	≥8F	F-J ● ●	3.5 mm	4 + 2	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing

\* THERMOCOOL™ SF Catheter requires an upgrade of the Stockert EP shuttle RF Generator to Software Version 1.037. Please contact your local BW representative.

\*\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## EZ STEER™ THERMOCOOL™ Catheter



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36Q33M†	8F	D-D ●●	3.5 mm	4	2-5-2 mm	115 cm
36Q35M	8F	D-F ●●	3.5 mm	4	2-5-2 mm	115 cm
36Q55M	8F	F-F ●●	3.5 mm	4	2-5-2 mm	115 cm
36Q5JM	8F	F-J ●●	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

## EZ STEER™ 4 mm Bi-Directional Catheter



1-7-4 mm



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35I35R	7F	D-F ●●	4 mm	4	1-7-4 mm	115 cm
35I5JR	7F	F-J ●●	4 mm	4	1-7-4 mm	115 cm
35I33R*	7F	D-D ●●	4 mm	4	1-7-4 mm	115 cm
35I55R*	7F	F-F ●●	4 mm	4	1-7-4 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator

\* For availability, please contact your local Biosense Webster representative.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## EZ STEER™ DS Bi-Directional Catheter

BI



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35U35R	7F	D-F ● ●	8 mm	4	1-7-4 mm	115 cm
35U5JR	7F	F-J ● ●	8 mm	4	1-7-4 mm	115 cm
35U33R*	7F	D-D ● ●	8 mm	4	1-7-4 mm	115 cm
35U55R*	7F	F-F ● ●	8 mm	4	1-7-4 mm	115 cm
35UJJR*	7F	J-J ● ●	8 mm	4	1-7-4 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE™ Generator

## CELSIUS™ THERMOCOOL™ Uni-Directional Catheter

UNI



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35Q13R*	7F	B ●	3.5 mm	4	2-5-2 mm	115 cm
35Q33R	7F	D ●	3.5 mm	4	2-5-2 mm	115 cm
35Q53R	7F	F ●	3.5 mm	4	2-5-2 mm	115 cm
35Q73R*	7F	J ●	3.5 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## CELSIUS™ RMT THERMOCOOL™ Catheter

RMT



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Sheath Fit*	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
CR7TCSIRT**	8F	≥8.5F	3.5 mm	4	2-5-2 mm	130 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

Associated Tubing:

Catalog No.:	Description
CFT003	COOLFLOW™ Pump Tubing
SAT001	SMARTABLATE™ System Irrigation Tubing Set

## CELSIUS™ 4 mm Catheter Thermocouple

UNI



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35G07R†	6F	A ●	4 mm	4	2-5-2 mm	92 cm
35G17R	6F	B ●	4 mm	4	2-5-2 mm	92 cm
35G37R†	6F	D ●	4 mm	4	2-5-2 mm	92 cm
36N07R†	7F	A ●	4 mm	4	2-5-2 mm	115 cm
36N17R	7F	B ●	4 mm	4	2-5-2 mm	115 cm
36N27R	7F	C ●	4 mm	4	2-5-2 mm	115 cm
36N37R	7F	D ●	4 mm	4	2-5-2 mm	115 cm
36N47R†	7F	E ○	4 mm	4	2-5-2 mm	115 cm
36N57R	7F	F ●	4 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

\*\* For availability, please contact your local Biosense Webster representative.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product. Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## CELSIUS™ 4 mm Catheter Thermistor

UNI



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36O07R	7F	A ●	4 mm	4	2-5-2 mm	115 cm
36O17R	7F	B ●	4 mm	4	2-5-2 mm	115 cm
36O27R	7F	C ●	4 mm	4	2-5-2 mm	115 cm
36O37R	7F	D ●	4 mm	4	2-5-2 mm	115 cm
36O47R	7F	E ○	4 mm	4	2-5-2 mm	115 cm
36O57R	7F	F ●	4 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

## CELSIUS™ 4 mm Braided Tip Catheter

UNI



2-5-2 mm



Catheter handle

European Ordering No.	Outer Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
34D37R	7F	D ●	4 mm	4	2-5-2 mm	115 cm
34D47R	7F	E ○	4 mm	4	2-5-2 mm	115 cm
34D57R	7F	F ●	4 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## CELSIUS FLTR™ 8 mm Uni-Directional Catheter

UNI



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D135501†	7.5F	≥8F	F	8 mm	4	0.5-5-2 mm	115 cm
D135502†	7.5F	≥8F	XR	8 mm	4	0.5-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE™ Generator

## CELSIUS FLTR™ 8 mm Bi-Directional Catheter

BI



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
D135901	7.5F	≥8F	D-XR	8 mm	4	0.5-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE™ Generator

\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

† Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# FOCAL - NON-NAVIGATION THERAPEUTIC CATHETERS

## CELSIUS™ DS Catheter

UNI



European Ordering No.	Outer Diameter	Sheath Fit*	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35S17R	7F	≥8F	B <span style="color: red;">●</span> <span style="color: yellow;">●</span>	8 mm	4	1-6-2 mm	115 cm
35S37R	7F	≥8F	D <span style="color: blue;">●</span> <span style="color: yellow;">●</span>	8 mm	4	1-6-2 mm	115 cm
35S57R	7F	≥8F	F <span style="color: orange;">●</span> <span style="color: yellow;">●</span>	8 mm	4	1-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39E68R	Stockert EP Shuttle Generator
D130303	SMARTABLATE™ Generator

## CELSIUS™ RMT Catheter

RMT



European Ordering No.	Outer Diameter	Sheath Fit*	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
CR7TCS4RT	7F	≥8F	4 mm	4	2-5-2 mm	125 cm

Associated Cable:

Catalog No.:	Connects
39E43R	Stockert EP Shuttle Generator
D130302	SMARTABLATE™ Generator

\* Compatible sheath size. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.



## HELIOSTAR™ Balloon Ablation Catheter



European Ordering No.	French Size	Size
D138904	10.5F	28 mm

Associated Cable:

Catalog No.:	Connects
D133701	Cable for HELIOSTAR™ Catheter to nGEN™ Generator

## LASSOSTAR™ Circular Catheter



European Ordering No.	Loop Diameter	Number of Electrodes	Size	Spacing
D139001	15 mm	10	3F	4.5 mm
D139002	20 mm	10	3F	6 mm
D139003	25 mm	10	3F	8 mm

Associated Cable:

Catalog No.:	Connects
CB3434CT	Cable for LASSOSTAR™ to CARTO™ 3 System

## LASSOSTAR™ NAV Circular Catheter



European Ordering No.	Loop Diameter	Number of Electrodes	Size	Spacing
D140401	15 mm	10	3F	4.5 mm
D140402	20 mm	10	3F	6 mm
D140403	25 mm	10	3F	8 mm

Associated Cable:

Catalog No.:	Connects
	Cable for LASSOSTAR™ NAV to CARTO™ 3 System

## GUIDESTAR™ Deflectable Sheath



European Ordering No.	French Size	Usable Length
D141101	13.5F	71 cm
D141102	13.5F	71 cm
D141103	14F	71 cm

# Diagnostic Catheters



## Navigation Diagnostic Catheters


LASSO™ 2515 NAV eco Circular Mapping Catheter .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	38
OCTARAY™ Mapping Catheter .....	39
PENTARAY™ NAV eco High Density Mapping Catheter .....	39
DECANAV™ Catheter .....	40
CARTO™ FAM DX .....	40

## Auto-ID Diagnostic Catheters

LASSO™ Decapolar Catheter .....	41
LASSO™ Duo-Decapolar Catheter .....	41
WEBSTER™ CS Uni-Directional Catheter .....	42
WEBSTER™ CS Bi-Directional Catheter .....	42
WEBSTER™ Quadrapolar Catheter Deflectable .....	43
WEBSTER™ Decapolar Catheter Deflectable .....	43
WEBSTER™ HIS Catheter .....	44
WEBSTER™ Quadrapolar Catheter Fixed .....	44
WEBSTER™ Decapolar Catheter Fixed .....	44
HALO™ XP Catheter .....	45
ISMUS™ Catheter .....	45

## Non Auto-ID Diagnostic Catheters

LASSO™ 2515 NAV eco Circular Mapping Catheter .....	46
LASSO™ NAV eco Circular Mapping Catheter .....	46
WEBSTER™ CS Bi-Directional Catheter .....	47
WEBSTER™ Octopolar Catheter Deflectable .....	47
WEBSTER™ Hexapolar Catheter Deflectable .....	48
WEBSTER™ Quadrapolar Catheter Deflectable .....	48
WEBSTER™ Decapolar Catheter Deflectable .....	49
WEBSTER™ Decapolar Catheter Fixed .....	50
WEBSTER™ Hexapolar Catheter Fixed .....	50
WEBSTER™ Quadrapolar Catheter Fixed .....	51
WEBSTER™ Duo-Decapolar Catheter .....	52
HALO™ XP Tricuspid Mapping Catheter .....	52
ISMUS™ Catheter .....	52
AVAIL™ Bipolar Fixed Shape Catheter .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	53
CRISTACATH™ Catheter .....	54
PARAHISIAN™ Catheter .....	54
Pediatric Catheter .....	54




## LASSO™ 2515 NAV eco Circular Mapping Catheter



European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Diameter	Insertion Length
D134302	7F	D	10	8 mm	25-15 mm	115 cm
D134301	7F	D	20	2-6-2 mm	25-15 mm	115 cm

## LASSO™ NAV eco Circular Mapping Catheter




European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Diameter	Insertion Length	Interface Cable	Compatibility
D134901	7F	D	10	4.5 mm	15 mm	115 cm	D134401	CARTO™ 3 System
D134903	7F	D	10	6 mm	20 mm	115 cm	D134401	CARTO™ 3 System
D134905	7F	D	10	8 mm	25 mm	115 cm	D134401	CARTO™ 3 System

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Diameter	Insertion Length	Interface Cable	Compatibility
D134902	7F	D	20	4.5 mm	15 mm	115 cm	D134401	CARTO™ 3 System
D134904	7F	D	20	6 mm	20 mm	115 cm	D134401	CARTO™ 3 System
D134906	7F	D	20	8 mm	25 mm	115 cm	D134401	CARTO™ 3 System

Associated Cable:

Catalog No.:	Connects
D134401	Catheter to CARTO™ 3 System PIU
EM5050060	20 pole eco Cable

## OCTARAY™ Mapping Catheter



European Ordering No.	Diameter	Sheath Fit*	Curve	Number of electrodes	Spacing	Spline length
D160901	8F	≥8.5F Soft tip	D ●	48	2-2-2 mm	1.5 cm
D160902	8F	≥8.5F Soft tip	D ●	48	2-5-2 mm	2.0 cm
D160903	8F	≥8.5F Soft tip	D ●	48	3-3-3 mm	2.0 cm
D160904	8F	≥8.5F Soft tip	F ●	48	2-2-2 mm	1.5 cm
D160905	8F	≥8.5F Soft tip	F ●	48	2-5-2 mm	2.0 cm
D160906	8F	≥8.5F Soft tip	F ●	48	3-3-3 mm	2.0 cm

Associated Cable:

Catalog No.:	Connects
D134402	Catheter to Signal Processing Unit (SPU)

Associated Hardware:


Catalog No.:	Connects
KT-5070-00	Signal Processing Unit (SPU) to CARTO™ 3 System PIU

## PENTARAY™ NAV eco High Density Mapping Catheter

4-4-4 mm



Catheter handle



End View

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
D128207	7F	F ●	20	4-4-4 mm	115 cm
D128208	7F	F ●	20	2-6-2 mm	115 cm
D128210	7F	D ●	20	4-4-4 mm	115 cm
D128211	7F	D ●	20	2-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
D134401	Catheter to CARTO™ 3 System PIU
EM5050060	20 pole eco Cable

\* Approved sheaths: CARTO VIZIGO™ Bi-Birectional Guiding Sheath, HeartSpan, CardiaGuide® ML0, Zulpaz®, Swartz SL1, Agilis NXT



## DECANAV™ Catheter



Pairs at 10 mm (2-8-2 mm)



Catheter handle

European Ordering No.	Diameter	Curve	Tip Size	Number of Electrodes	Spacing	Insertion Length
R7D282CT	7F	D 	2 mm	10	2-8-2 mm	115 cm
R7F282CT	7F	F 	2 mm	10	2-8-2 mm	115 cm


Associated Cable:

Catalog No.:	Connects
CB3434CT	Catheter to CARTO™ 3 System PIU

## CARTO™ FAM DX

European Ordering No.	Description
KT5400165	FAM DX Module
KT5400165U	FAM DX Module OBL

### LASSO™ Decapolar Catheter




European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Loop Diameter	Insertion Length
35026RCT	7F	D ●	3F	10	4.5 mm	15 mm	115 cm
35036RCT	7F	D ●	3F	10	6 mm	20 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO™ 3 System PIU
D128622	Catheter to CARTO™ 3 System PIU

### LASSO™ Duo-Decapolar Catheter



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Loop Diameter	Insertion Length
D7L2015CT*	7F	D ●	3F	20	4.5 mm	15 mm	115 cm
D7L2020CT*	7F	D ●	3F	20	6 mm	20 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3434CT	Catheter to CARTO™ 3 System PIU
D128622	Catheter to CARTO™ 3 System PIU

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

## WEBSTER™ CS Uni-Directional Catheter



2-8-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrodes	Number of Electrodes	Spacing	Insertion Length	Auto ID	Compression Coil
D135303	6F	CS-D	2 mm	10	2-8-2 mm	115 cm	Yes	Yes
D135304	6F	CS-F	2 mm	10	2-8-2 mm	115 cm	Yes	Yes

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO™ 3 System PIU
CY1212CT	Catheter to CARTO™ 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox



## WEBSTER™ CS Bi-Directional Catheter



2-8-2 mm




Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36G35QCT	7F	D-F 	10	2-8-2 mm	115 cm
36G5JMCT	7F	F-J 	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO™ 3 System PIU
CY1212CT	Catheter to Carto™ 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER™ Quadrapolar Catheter Deflectable



European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36B38RCT	6F	D ●	4	2-5-2 mm	92 cm
36B33RCT	6F	D ●	4	5 mm	92 cm
36C38RCT	7F	D ●	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to Carto™ 3 System PIU
CY1212CT	Catheter to CARTO™ 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER™ Decapolar Catheter Deflectable



European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
35X33QCT*	6F	D ●	10	2-8-2 mm	92 cm
36G33QCT*	7F	D ●	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO™ 3 System PIU
CY1212CT	Catheter to Carto™ 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.



## WEBSTER™ HIS Catheter



Catheter handle

European Ordering No.	Diameter	Curve	Poles	Spacing	Insertion Length
D1085413	5F	HIS-F	4	5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to Carto™ 3 System PIU
CY1212CT	Catheter to CARTO™ 3 System PIU
D128709	Auto ID Catheters up to 10 poles to the EP recording system pinbox

## WEBSTER™ Quadrapolar Catheter Fixed



5 mm



Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
37C53RCT	5F	F ●	4	5 mm	110 cm
37C03RCT	5F	A ●	4	5 mm	110 cm
F6QF005CT	6F	F ●	4	5 mm	115 cm
F6QA005CT	6F	A ●	4	5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to Carto™ 3 System PIU
CY1212CT	Catheter to CARTO™ 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

## WEBSTER™ Decapolar Catheter Fixed



2-8-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
37S94RCT	5F	P	10	2-8-2 mm	60 cm

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to CARTO™ 3 System PIU
CY1212CT	Catheter to CARTO™ 3 System PIU
D128709	Auto ID Catheter up to 10 poles to the EP recording system pinbox

### HALO™ XP Catheter



2-8-2 mm



Catheter handle

European Ordering No.	Diameter	Number of Electrodes	Curve	Spacing	Insertion Length
36J13RCT*	7F	20	Tricuspid	2-8-2 mm	110 cm

Associated Cable:

Catalog No.:	Connects
CB3434CT	Catheter to CARTO™ 3 System PIU
D128622	Auto ID Catheter up to 20 poles to the EP recording system pinbox

### ISMUS™ Catheter



2-12-2 mm



Catheter handle


European Ordering No.	Diameter	Number of Electrodes	Curve	Spacing	Insertion Length
36Y32RCT*	7F	20	ISMUS™ Cath	2-12-2 mm	110 cm

Associated Cable:


Catalog No.:	Connects
CB3434CT	Catheter to CARTO™ 3 System PIU
D128622	Auto ID Catheter up to 20 poles to the EP recording system pinbox

\* Please contact your Biosense Webster representatives for further details on the product.

## LASO™ 2515 Variable Mapping Catheter



End View



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35O2515R*	7F	D ●	15-25 mm	10	8 mm	115 cm
35T2515R*	7F	D ●	15-25 mm	20	2-6-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>1</sup> or 39F54R <sup>1</sup>	ECG Recording System CARTO™ 3 System
D128624	Catheter to CARTO™ 3 System

## LASO™ Mapping Catheter

European Ordering No.	Diameter	Curve	Loop Diameter	Number of Electrodes	Spacing	Insertion Length
35O16R*	7F	D ●	12 mm	10	4 mm	115 cm
35O26R*	7F	D ●	15 mm	10	4.5 mm	115 cm
35O36R*	7F	D ●	20 mm	10	6 mm	115 cm
35O46R*	7F	D ●	25 mm	10	8 mm	115 cm
35O66R*	7F	D ●	35 mm	15	7 mm	115 cm
35T26R*	7F	D ●	15 mm	20	4.5 mm	115 cm
35T36R*	7F	D ●	20 mm	20	6 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>1</sup> or 39F54R <sup>1</sup>	ECG Recording System CARTO™ 3 System
D128624	Catheter to CARTO™ 3 System

<sup>1</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.

\*Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

### WEBSTER™ CS Bi-Directional Catheter



2-8-2 mm



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36D35R	7F	D-F ●●	2 mm	10	2-8-2 mm	115 cm
36D5JR	7F	F-J ●●	2 mm	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>1</sup>	ECG Recording System / CARTO™ XP System
CB3410CT	CARTO™ 3 System

### WEBSTER™ Octopolar Catheter Deflectable



Irregular Spacing due to anchor placement



Pairs at 10 mm (2-8-2 mm)



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35Y13Q*	5F	D ●	1 mm	8	2-8-2 mm	92 cm
36S31Q	6F	D ●	1 mm	8	2 mm	92 cm
36F30Q*	7F	D ●	1 mm	8	2.5-3.75-2.5-2.5-2.5-2.5-2.5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F44R	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

<sup>1</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

## WEBSTER™ Hexapolar Catheter Deflectable



Irregular Spacing due to anchor placement



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36E30Q	6F	D ●	1 mm	6	2.5-4-2.5-2.5-2.5 mm	92 cm

Associated Cable:

Catalog No.:	Connects
39F43R	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

## WEBSTER™ Quadrapolar Catheter Deflectable



2-5-2 mm



5 mm



10 mm



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36B38R	6F	D ●	1 mm	4	2-5-2 mm	92 cm
36B33R	6F	D ●	1 mm	4	5 mm	92 cm
36B35R	6F	D ●	1 mm	4	10 mm	92 cm
36C38R	7F	D ●	1 mm	4	2-5-2 mm	115 cm
36C33R	7F	D ●	1 mm	4	5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F40R <sup>1</sup>	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

<sup>1</sup> The interface cable 39F40R is 0.9 m long. It is also available in 1.8 m as 39F55R.

### WEBSTER™ Decapolar Catheter Deflectable




European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
35X38Q	6F	D ●	1 mm	10	2-5-2 mm	92 cm
35X33Q	6F	D ●	1 mm	10	2-8-2 mm	92 cm
35X58Q	6F	F ●	1 mm	10	2-5-2 mm	92 cm
36G33Q	7F	D ●	1 mm	10	2-8-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>1</sup>	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

<sup>1</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.

## WEBSTER™ Decapolar Catheter Fixed




European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
37S94R	5F	P-CS	1 mm	10	2-8-2 mm	60 cm
37S95R	5F	P-CS	-	10	2-8-2 mm	110 cm
37G58R	6F	F-Cournand Type <sup>1</sup> ●	1 mm	10	2-5-2 mm	115 cm
37G68R	6F	G-curve ●	1 mm	10	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F41R <sup>2</sup>	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

## WEBSTER™ Hexapolar Catheter Fixed



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
37U00R	6F	Josephson Type <sup>1</sup> ●	1 mm	6	2.5-2.5-2.5 2.5-247.5 mm	115 cm


Associated Cable:

Catalog No.:	Connects
39F43R	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

<sup>1</sup> The third party trademarks listed herein are trademarks of their respective owners.

<sup>2</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.

### WEBSTER™ Quadrapolar Catheter Fixed



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
37C08R	5F	A-Josephson Type <sup>1</sup> ●	1 mm	4	2-5-2 mm	110 cm
37C03R	5F	A-Josephson Type <sup>1</sup> ●	1 mm	4	5 mm	110 cm
37C53R	5F	F-Cournand Type <sup>1</sup> ●	1 mm	4	5 mm	110 cm
37T03R	5/6F <sup>2</sup>	A-Josephson Type <sup>1</sup> ●	1 mm	4	5 mm	115 cm
37D00R	6F	A-Josephson Type <sup>1</sup> ●	1 mm	4	2 mm	115 cm
37D08R	6F	A-Josephson Type <sup>1</sup> ●	1 mm	4	2-5-2 mm	115 cm
37D03R	6F	A-Josephson Type <sup>1</sup> ●	1 mm	4	5 mm	115 cm
37D05R	6F	A-Josephson Type <sup>1</sup> ●	1 mm	4	10 mm	115 cm
37D38R	6F	D-Damato Type <sup>1</sup> ●	1 mm	4	2-5-2 mm	115 cm
37D33R	6F	D-Damato Type <sup>1</sup> ●	1 mm	4	5 mm	115 cm
37D35R	6F	D-Damato Type <sup>1</sup> ●	1 mm	4	10 mm	115 cm
37D53R	6F	F-Cournand Type <sup>1</sup> ●	1 mm	4	5 mm	115 cm
37D55R	6F	F-Cournand Type <sup>1</sup> ●	1 mm	4	10 mm	115 cm
37D58R	6F	F-Cournand Type <sup>1</sup> ●	1 mm	4	2-5-2 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F40R <sup>2</sup>	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

<sup>1</sup> The third party trademarks listed herein are trademarks of their respective owners.

<sup>2</sup> The interface cable 39F40R is 0.9 m long. It is also available in 1.8 m as 39F55R.



## WEBSTER™ Duo-Decapolar Catheter



Catheter handle

European Ordering No.	Diameter	Loop diameter	Tip Electrodes	Spacing	Insertion Length
D728260RT	7F	60 mm	2 mm	2-8-2...60...2-8-2 mm	110 cm

Associated Cable:

Catalog No.:	Connects
2 x 39F54R	ECG Recording System CARTO™ 3 System
D128624	Catheter to CARTO™ 3 System

## HALO™ XP Tricuspid Mapping Catheter



HALO™ XP 2-8-2 mm (w/20 mm tail)



HALO™ XP 2-13-2 mm



HALO™ XP 2-18-2-8-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36J13R	7F	Tricuspid	20	(20)-2-8-2 mm	110 cm
36J16R	7F	Tricuspid	20	2-13-2 mm	110 cm
36J14R	7F	Tricuspid	20	2-18-2-8-2 mm	110 cm


Associated Cable:

Catalog No.:	Connects
2 x 39F41R <sup>1</sup>	ECG Recording System CARTO™ 3 System
D128624	Catheter to CARTO™ 3 System

## ISMUS™ Catheter



2-12-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Number of Electrodes	Spacing	Insertion Length
36Y32R	7F	ISMUS™	20	2-12-2 mm	110 cm

Associated Cable:

Catalog No.:	Connects
2 x 39F41R <sup>1</sup>	ECG Recording System CARTO™ 3 System
D128624	Catheter to CARTO™ 3 System

<sup>1</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.

### AVAIL™ Bipolar Fixed Shape Catheter




European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
37P0RR*	6F	A-Josephson Type <sup>1</sup>	na	2	10 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F40R <sup>2</sup>	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

### AVAIL™ Quadrapolar Fixed Shape Catheter



European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
37C0CR	5F	A-Josephson Type	●	4	2-5-2 mm	115 cm
37C3FR	5F	D-Damato Type <sup>1</sup>	●	4	10 mm	115 cm
37C5RR	5F	F-Cournand Type <sup>1</sup>	●	4	2-5-2 mm	115 cm
37C5CR	5F	F-Cournand Type <sup>1</sup>	●	4	5 mm	115 cm
37D0CR	6F	A-Josephson Type <sup>1</sup>	●	4	2-5-2 mm	115 cm
37D5CR	6F	F-Cournand Type <sup>1</sup>	●	4	2-5-2 mm	115 cm
37D0RR	6F	A-Josephson Type <sup>1</sup>	●	4	5 mm	115 cm
37D3FR	6F	D-Damato Type <sup>1</sup>	●	4	10 mm	115 cm
37D5RR	6F	F-Cournand Type <sup>1</sup>	●	4	5 mm	115 cm

Associated Cable:

Catalog No.:	Connects
39F40R <sup>2</sup>	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

\*Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

<sup>1</sup> The third party trademarks listed herein are trademarks of their respective owners.

<sup>2</sup> The interface cable 39F40R is 0.9 m long. It is also available in 1.8 m as 39F55R.

## CRISTACATH™ Catheter



1-3-1 mm



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36Y39R	7F	D-Type/Crista	3 mm	20	1-3-1 mm	115 cm

Associated Cable:

Catalog No.:	Connects
2 x 39F41R <sup>1</sup>	ECG Recording System CARTO™ 3
D128624	Catheter to CARTO™ 3 System

## PARAHISIAN™ Catheter



2-10-2-2-2-5-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length
36L00Q*	7F	D ●	1 mm	8	2-10-2-2-2-2-2 mm	90 cm

Associated Cable:

Catalog No.:	Connects
39F44R	ECG Recording System / CARTO™ XP

## Pediatric Catheter



2-10-2-2-2-5-2 mm



Catheter handle

European Ordering No.	Diameter	Curve	Tip Electrode	Number of Electrodes	Spacing	Insertion Length	Color Code	Compatibility
36S30Q*	6F	D	1 mm	8	2-10-2-2-2-5-2 mm	92 cm	blue	CARTO™ XP System/ CARTO™ 3 System/ ECG Recording System

Associated Cable:

Catalog No.:	Connects
39F44R	ECG Recording System / CARTO™ XP System
CY1210CT	CARTO™ 3 System

<sup>1</sup> The interface cable 39F41R is 0.9 m long. It is also available in 1.8 m as 39F54R.

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# Access



## Guiding Sheaths - Fixed Guiding Sheaths

PREFACE™ Fixed Guiding Sheaths for EP .....	56
PREFACE™ Fixed Guiding Sheaths .....	56
CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath .....	57

## Guiding Sheaths - Steerable Guiding Sheaths


CARTO VIZIGO™ Sheaths .....	58
MOBICATH™ Sheaths .....	58

## Transseptal Needles


HEARTSPAN™ Transseptal Needles .....	59
--------------------------------------	----



## PREFACE™ Fixed Guiding Sheaths for EP



European Ordering No.	Size	Internal Diameter	Length	Curve
301803MS	8F <sup>1</sup>	0.11 inch	62 cm	Multipurpose Short



European Ordering No.	Size	Internal Diameter	Length	Curve
301803M	8F <sup>1</sup>	0.11 inch	62 cm	Multipurpose
301805M	8F <sup>1</sup>	0.11 inch	77 cm	Multipurpose

## PREFACE™ Fixed Guiding Sheaths



European Ordering No.	Size	Internal Diameter	Length	Curve
301803A	8F <sup>1</sup>	0.11 inch	62 cm	Anterior



European Ordering No.	Size	Internal Diameter	Length	Curve
301803P	8F <sup>1</sup>	0.11 inch	62 cm	Posterior

<sup>1</sup> Maximum Acceptable Lead Outer Diameter.

### CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath



European Ordering No.	Curve Type	Inner Diameter	Length
FCL20100	MLO	8.5F	63 cm
FCL20101	ML1	8.5F	63 cm
FCL20102	MLO	8.5F	81 cm
FCL20103	ML1	8.5F	81 cm

# GUIDING SHEATHS - STEERABLE GUIDING SHEATHS

## CARTO VIZIGO™ Sheaths



European Ordering No.	Usable Length	Inner Diameter	Outer Diameter	Curve size	Description
D138501	71 cm	8.5F	11.5F	17 mm	Small Curve
D138502	71 cm	8.5F	11.5F	22 mm	Medium Curve
D138503	71 cm	8.5F	11.5F	50 mm	Large Curve

Associated Cable:

Catalog No.:	Connects
CB3412CT	Catheter to Carto3™ System PIU
CY1212CT	Catheter to Carto3™ System PIU


## MOBICATH™ Sheaths



European Ordering No.	Usable Length	Inner Diameter	Outer Diameter	Curve size	Description
D140010*	72 cm	8.5F	11.5F	17 mm	Sheath Alone Small Curve
D140011*	72 cm	8.5F	11.5F	30 mm	Sheath Alone Large Curve

\*Biosense Webster will be discontinuing this product code end of 2022.  
Please contact your Biosense Webster representative to identify an alternative product.

## HEARTSPAN™ Transseptal Needles



European Ordering No.	Description	Curve Angle	Needle ID
FND01901	18 gauge X 71 cm	50°	.015 inch Minimum
FND01902	18 gauge X 71 cm	86°	.015 inch Minimum
FND01903	18 gauge X 89 cm	50°	.015 inch Minimum
FND01905	18 gauge X 98 cm	50°	.015 inch Minimum
FND01906	18 gauge X 98 cm	86°	.015 inch Minimum



# Implantable Devices



## The COHEREX WAVECREST™ LAA Occlusion System

Device Description.....	61
Implant Characteristics .....	61
COHEREX WAVECREST™ LAA Occlusion System.....	62
Sheath ordering information.....	62



## The COHEREX WAVECREST™ LAA Implant


### Device Description

The COHEREX WAVECREST™ LAA Implant includes a self-expanding Nitinol frame and an ePTFE covered occluder. The COHEREX WAVECREST™ LAA Implant comes pre-attached to the delivery system which consists of a delivery catheter, loading device and a proximal control handle. The delivery system is designed to allow a single operator to position and deploy the Implant through The COHEREX WAVECREST™ LAA Occlusion System Delivery Sheath.

Delivery Sheath Usable Length: 70 cm



### Implant Characteristics

The laser cut Nitinol frame is designed to allow the device to be deployed independent of the anchoring system. Once position is optimized, 20 interlocking anchors are actuated to engage the surrounding tissue and complete the implant deployment. The implant can be fully recaptured at any time prior to release. The implant frame struts and markers delineate the implant for correct alignment and visualization under fluoroscopic imaging. A rim of polyurethane foam acts as a substrate to stimulate tissue in-growth.



# THE COHEREX WAVECREST™ LAA IMPLANT

## COHEREX WAVECREST™ LAA Implant




WCR1503	22 mm
---------	-------

WCR1513	27 mm
---------	-------

WCR1523	32 mm
---------	-------

## Sheath ordering information




WCR1530	60 curve
---------	----------

WCR1540	75 curve
---------	----------

WCR1541	90 curve
---------	----------

WCR1551	90S curve
---------	-----------

# UltraSound Catheters



## Ultrasound Catheters

SOUNDSTAR™ 3D eco Catheter.....	64
SOUNDSTAR™ eco 8F Catheter .....	64
SOUNDSTAR™ 3D Ultrasound Catheter .....	65
ACUNAV™ 10F Ultrasound Catheter .....	65
ACUNAV™ 8F Ultrasound Catheter .....	66

## System Compatibility

CARTOSOUND® Module System Compatibility.....	67
ACUNAV™ Catheter System Compatibility* .....	67

## SOUNDSTAR™ 3D eco Catheter



European Ordering No.	Outer Diameter	Insertion Length	System Compatibility
10438577	10F	90 cm	Siemens Sequoia™, SC 2000™ Cypress™, Acuson X300™
10439072	10F	90 cm	GE-Vivid i™ Vivid q™

Associated Cable:

Catalog No.:	Connects
CG2025CT	Catheter to CARTO™ 3 System PIU
EM5050011	SOUNDSTAR™ 3D eco Cable compatible with CARTO™ 3 System

## SOUNDSTAR™ eco 8F Catheter



European Ordering No.	Outer Diameter	Insertion Length	System Compatibility
10439011	8F	90 cm	Siemens Sequoia™, SC 2000™ Cypress™, Acuson X300™
10439236	8F	90 cm	GE-Vivid i™ Vivid q™

Associated Cable:

Catalog No.:	Connects
CG2025CT	Catheter to CARTO™ 3 System PIU
EM5050011	SOUNDSTAR™ 3D eco Cable compatible with CARTO™ 3 System

## SOUNDSTAR™ 3D Ultrasound Catheter

NAV



European Ordering No.	Outer Diameter	Insertion Length	System Compatibility
SNDSTR10	10F	90 cm	Siemens Sequoia™, SC 2000™ Cypress™, Acuson X300™
SNDSTR10G	10F	90 cm	GE-Vivid i™ Vivid q™

Associated Cable:

Catalog No.:	Connects
CG2025CT	Catheter to CARTO™ 3 System PIU

## ACUNAV™ Ultrasound Catheter\*\*



European Ordering No.	Outer Diameter	Length	System Compatibility*
08255790	10F	90 cm	Siemens Sequoia™, SC 2000™, Cypress™, Acuson X300™, Acuson X700™, S Series™
10043342	10F	90 cm	GE-Vivid i™, GE-Vivid q™

\* Please contact your Biosense Webster representatives for further details on exact compatible versions.  
The third party trademarks used herein are trademarks of their respective owners.

\*\* Trademark of Siemens Medical Solutions USA, Inc.

## ACUNAV™ Ultrasound Catheter\*\*



European Ordering No.	Outer Diameter	Length	System Compatibility*
10135936	8F	90 cm	Siemens Sequoia™, SC 2000™, Cypress™, Acuson X300™, Acuson X700™, S Series™
10135910	8F	90 cm	GE-Vivid i™, GE-Vivid q™

\* Please contact your Biosense Webster representatives for further details on exact compatible versions.  
The third party trademarks used herein are trademarks of their respective owners.

\*\* Trademark of Siemens Medical Solutions USA, Inc.

## SYSTEM COMPATIBILITY

### CARTOSOUND® Module System Compatibility\*

	Sequoia™	Cypress™	Acuson X300™	Vivid i™	Vivid q™
AcuNav™ Catheter Option	Required	Required	Required	Required	Required
BWI Comm Kit	BWI COMM KIT Prespective Router and Software 4-6 weeks leadtime	No Kit Required	BWI COMM KIT Prespective Router and Software 4-6 weeks leadtime	VGA SPLITTER KIT VGA Splitter Box Power Cord 2 DVI cables Software 1 week leadtime	Required VGA Splitter Box Power Cord 2 DVI cables Software 1 week leadtime
Ultrasound Software Level	10, 11, and 12.221	14 and 20	5, 6 and 7	BT10, BT11, BT12 with CARTOSOUND® Module Interface Software	BT10, BT11, BT12 with CARTOSOUND® Module Interface Software
Carto™ System Software Required	Carto™ XP System or Carto™ 3 System	Carto™ XP System or Carto™ 3 System	Acuson X300™ V7 and below: Carto™ XP System and Carto™ 3 System Acuson X300™ V7: Carto™ 3 System only	Carto™ XP System or Carto™ 3 System	Carto™ XP System or Carto™ 3 System

### ACUNAV™ Catheter System Compatibility\*\*

Siemens	Sequoia™	Cypress™	Acuson X300™	SC2000™
GE	Vivid i™	Vivid q™		

### ACUNAV™ Catheter System Compatibility\*\*

Compatibility	Ultrasound Catheters	Product Codes	Label Color	Handle Color
Siemens	ACUNAV™ 10Fr 90 cm	08255790	Blue	White Handle, Blue Knobs
Siemens	ACUNAV™ 8Fr 90 cm	10135936	Blue	Blue Handle, White Knobs
Siemens	SOUNDSTAR™ eco 10Fr 90 cm	SNDSTR10	Green	Black Handle, Blue Knobs
Siemens	SOUNDSTAR™ 10Fr 90 cm	10438577	Yellow	White Handle, Blue Knobs
Siemens	SOUNDSTAR™ eco 8Fr 90 cm	10439011		Blue Handle, White Knobs
GE	ACUNAV™ 10Fr 90 cm	10043342	Orange	White Handle, Blue Knobs
GE	ACUNAV™ 8Fr 90 cm	10135910	Orange	Blue Handle, White Knobs
GE	SOUNDSTAR™ 10Fr 90 cm	SNDSTR10G	Purple	Black Handle, Blue Knobs
GE	SOUNDSTAR™ eco 10Fr 90 cm	10439072	Grey	White Handle, Blue Knobs
GE	SOUNDSTAR™ eco 8Fr 90 cm	10439236		Blue Handle, White Knobs

\* Please contact your Biosense Webster representatives for further details on exact compatible versions.

\*\* ACUNAV™ is a trademark of Siemens Medical Solutions USA, Inc.





## RF Generators and Pump Systems



### **nGEN™ GENERATOR PLATFORM**

nGEN™ Generator Platform ..... 69

### **SMARTABLATE™ System**

SMARTABLATE™ System ..... 70

SMARTABLATE™ System Diagram ..... 71

SMARTABLATE™ System Service Contracts ..... 72

### **Stockert EP Shuttle Generator**

Indifferent Electrodes and Corresponding Connection Cables ..... 73

Spare Parts and Others ..... 74

Stockert EP Shuttle Accessories ..... 75

Service Contracts ..... 75

### **COOLFLOW™ Pump**

COOLFLOW™ Pump Tubing<sup>1</sup> ..... 76

COOLFLOW™ Pump Accessories ..... 76

Diagram - Fully Integrated COOLFLOW™ Pump System ..... 77

<sup>1</sup> Tubing to be used with COOLFLOW™ Pump.

# nGEN™ GENERATOR PLATFORM

## nGEN™ Generator Platform



Core System Kit		
<b>D138401</b>	nGEN™ Generator Kit, including Console, Power Supply, Monitor, and Accessory Kit	
<b>D138801</b>	nGEN™ Monitor	
<b>D139701</b>	nGEN™ Pump	
Installation Accessories		
<b>EM540200</b>	CARTO™ 3 System Cart	o
<b>MA5402300</b>	CARTO™ 3 System Cart Upgrade Kit for nGEN™	o
<b>M583196SP</b>	nGEN™ Monitor to Bed Rail Clamp	
<b>M549104SP</b>	nGEN™ Pump to IV Pole Clamp	+
<b>M581802SP</b>	nGEN™ Foot Pedal, 4.9 m	+
Catheter Interface Cables		
<b>D140102</b>	nGEN™ Generator TX eco Cable (Dongle for QDOT MICRO™)	o
<b>D135703</b>	nGEN™ Generator TX eco Extension Cable (for QDOT MICRO™)	o
<b>D133701</b>	HELIOSTAR™ Catheter to nGEN™ Generator Interface Cable	o
<b>D133703</b>	CELSIUS™ Catheter to nGEN™ Generator Interface Cable	o
System Interface Cables		
<b>M581101SP</b>	ECG Stimulator, 2 m (6.5 ft) Cable	+
<b>M583199SP</b>	nGEN™ Console to Recording System Cable, 3 m (10 ft)	
<b>M5831100SP</b>	nGEN™ Console to Recording System Cable, 7.6 m (25 ft)	
<b>M5831101SP</b>	nGEN™ Console to Recording System Cable, 13.7 m (45 ft)	
<b>M583197SP</b>	nGEN™ Monitor to Recording System Cable, 3 m (10 ft)	+
<b>M583198SP</b>	nGEN™ Monitor to Recording System Cable, 7.6 m (25 ft)	
<b>M589704SP</b>	Ethernet cable, 10 m (33 ft)	+(Qty.2)
<b>M589705SP</b>	Ethernet cable, 23 m (75 ft)	+
<b>M600203SP</b>	nGEN™ Console to nGEN™ Power Supply Unit Cable, 1.1 m (3.6 ft)	+
<b>M600202SP</b>	nGEN™ Console to nGEN™ Power Supply Unit Cable, 3 m (10 ft)	
<b>M600201SP</b>	nGEN™ Console to nGEN™ Power Supply Unit Cable, 5 m (16.4 ft)	
<b>M600301SP</b>	nGEN™ Console to nGEN™ Pump Cable, 5 m (16.4 ft)	+
<b>M581007SP</b>	nGEN™ Console to CARTO™ PIU Cable TX, 3 m (10 ft)	+
<b>M581006SP</b>	nGEN™ Console to CARTO™ PIU Cable TX, 3 m (10 ft) - for use with HELIOSTAR™	+
System Power Adapters		
<b>M621410SP</b>	nGEN™ Power Supply Unit AC Power Cord, 3 m (10 ft) - EU Standard	+
<b>M621402SP</b>	nGEN™ Power Supply Unit AC Power Cord, 3 m (10 ft) - United Kingdom	
<b>M621406SP</b>	nGEN™ Power Supply Unit AC Power Cord, 3 m (10 ft) - Israel	
<b>M621407SP</b>	nGEN™ Power Supply Unit AC Power Cord, 3 m (10 ft) - Switzerland	
<b>M529723SP</b>	nGEN™ Pump AC Adapter	+
<b>M649010SP</b>	nGEN™ Pump AC Power Cord, 3 m (10 ft) - EU Standard	+
<b>M649002SP</b>	nGEN™ Pump AC Power Cord, 3 m (10 ft) - United Kingdom	
<b>M649006SP</b>	nGEN™ Pump AC Power Cord, 3 m (10 ft) - Israel	
<b>M649007SP</b>	nGEN™ Pump AC Power Cord, 3 m (10 ft) - Switzerland	
<b>M583191SP</b>	nGEN™ Monitor Power Supply, 1.7 m (5.5 ft)	+(Qty.2)
Irrigation Tubing Set		
<b>SAT001</b>	SMARTABLATE™ Irrigation Tubing Set	o

+ = included in Core System Kit

o = not included and recommended

## SMARTABLATE™ System




European Ordering No. System	Product Description	Included in the SMARTABLATE™
M4900106*	SMARTABLATE™ System: SMARTABLATE™ Generator, SMARTABLATE™ Irrigation Pump, SMARTABLATE™ Remote control, Power Cord Cable, foot Pedal	X
SAT001	SMARTABLATE™ System Irrigation tubing set	
M490005	Foot pedal	X
M490012	SMARTABLATE™ Generator-Rem Cable (25 m)	X
M490010	SMARTABLATE™ Generator-Rem Cable (10 m)	
M490029	SMARTABLATE™ Generator-Rem Cable (15 m)	
M490013	SMARTABLATE™ Generator-Rem Cable (30 m)	
M490026	SMARTABLATE™ System Foot pedal ext. cable (7 m)	
M490027	SMARTABLATE™ Generator-Pump Cable (1 m)	
M490016	SMARTABLATE™ Generator-Pump Cable (3 m)	
M490028	SMARTABLATE™ Generator-Pump Cable (5 m)	X
M490030	SMARTABLATE™ Generator-Pump Cable (10 m)	
M490021	SMARTABLATE™-System-SCom Recording/CARTO™ System Cable (3 m)	
M490023	SMARTABLATE™-System-SCom Recording/CARTO™ System Cable (15 m)	X
M490025	SMARTABLATE™-System-SCom Recording/CARTO™ System Cable (5 m)	X
M490014	SMARTABLATE™ System Ground Cable (3 m)	X
M490017	Connection cable from Stockert Testbox III to an Indifferent Electrode	
M490018	SMARTABLATE™ System ECG Cable (3 m)	X
M490037	IV Pole Clamp	X
M4900602**	SMARTABLATE™ System Power Cord-UK	
M4900603**	SMARTABLATE™ System Power Cord-Russia	
M4900605**	SMARTABLATE™ System Power Cord-Switzerland	
M4900606**	SMARTABLATE™ System Power Cord-Israel	
MA5402100	SMARTABLATE™ System cart upgrade kit (needed for Carto™ 3 System Cart installed before Feb. 2017)	
M490034	SMARTABLATE™ System Erbe Indifferent Electrode Cable (3 m)	
M490035	SMARTABLATE™ System 3M Indifferent Electrode Cable (3 m)	

\* If intended to be connected to CARTO™ 3 System, please order D129704 or D129707.

For related cables, please see page 112.

\*\* Please order 3 times.

## SMARTABLATE™ System Diagram



## SMARTABLATE™ System Service Contracts

European Ordering No.	Product Description
SMAPRIV1	Privilege 1 Year SMARTABLATE™ System Complete Maintenance
SMAPRIV2	Privilege 2 Year SMARTABLATE™ System Complete Maintenance
SMAPRIV3	Privilege 3 Year SMARTABLATE™ System Complete Maintenance

### • Privilege

The **Privilege** Service Contract provides the following services:

- On-failure response time commitment
- Spare parts sent the same day (order before 1PM CET for loaners and 4PM CET for peripherals)
- Phone / email support
- Repair Labor
- Loaner for the repair time
- Faulty Parts including peripherals and damaged parts (excluding fire or water debris)
- One yearly Preventive Maintenance

The **First Year Warranty** includes the following:

- All terms of a **Privilege** Service Contract excluding the peripherals and damaged parts

# STOCKERT EP SHUTTLE GENERATOR

## Stockert EP Shuttle RF Accessories










## Indifferent Electrodes and Corresponding Connection Cables

European Ordering No.	Product Description
39D12X	Connection Cable for "3M" indifferent patch electrode (Redel-2 to "3M"-clip connector)
39D12X5	Connection Cable for "3M" indifferent patch electrode (Redel-2 to "3M"-clip connector, 5m)
39D16X	Connection cable for "Erbe" indifferent electrode (Redel-2 to "Erbe" - Plug, 3 m)
39D16X5*	Connection cable for "Erbe" indifferent electrode (Redel-2 to "Erbe" - Plug, 5 m)
39D82X	Connection cable for Valleylab™ indifferent electrode (Redel-2 to Valleylab™ - Plug, 3 m)
39F27R*	Multipurpose connection cable indifferent electrode (4 mm Banana - Plug)

\* For availability, please contact your local Biosense Webster representative.



## Spare Parts and Others

European Ordering No.	Product Description	Package 1	Package 2
39D03X	Foot pedal 	X	X
39D06X*	Serial interface cable Workstation to Global Port 		X
39D14X*	Grounding cable 	X	X
39D21X	Connection cable Stockert EP Shuttle RF Generator - Global Port (SubD+2xLWL to SubD+2xLWL), 3 m 		X
39D83X	ECG output cable 	X	X
39D91X	Connection cable EP- Shuttle to EP/IO, 2.9 m		
39D12X	Connection cable for "3M" indifferent patch electrode 	X	X
39E43R	Cable for CELSIUS™ Ablation Catheter (TC/THR) 	X	X

\* For availability, please contact your local Biosense Webster representative.

# STOCKERT EP SHUTTLE GENERATOR

## Stockert EP Shuttle Accessories

European Ordering No.	Product Description
39F28X10 <sup>1</sup>	Connection cable remote control to Stockert EP Shuttle RF Generator, 10 meters
39F28X15	Connection cable remote control to Stockert EP Shuttle RF Generator, 15 meters
39F28X20 <sup>1</sup>	Connection cable remote control to Stockert EP Shuttle RF Generator, 20 meters
39F28X30	Connection cable remote control to Stockert EP Shuttle RF Generator, 30 meters
39D60X	COOLFLOW™ Pump Control Interface Cable (6 m)
39D60X10	COOLFLOW™ Pump Control Interface Cable (10 m)

## Service Contracts

European Ordering No.	Product Description
SHUTBC1	Basic 1 Year EP Shuttle Preventive Maintenance
SHUTPV1 <sup>2</sup>	Privilege 1 Year EP Shuttle Complete Maintenance
SHUTPVP1 <sup>3</sup>	Privilege Plus 1 Year EP Shuttle Complete Maintenance with Peripherals

<sup>1</sup> For availability, please contact your local Biosense Webster representative.

<sup>2</sup> Preventive Maintenance, Corrective Maintenance, Spare Parts, and Loaners.

<sup>3</sup> Preventive Maintenance, Corrective Maintenance, Spare Parts, Loaners, and Peripherals.




## COOLFLOW™ Pump Tubing<sup>1</sup>

European  
Ordering No.

Product Description

CFT003

COOLFLOW™ Pump Tubing (1 per pack)




## COOLFLOW™ Pump Accessories

European  
Ordering No.

Product Description


39D03X

Foot pedal



39D14X\*

Grounding cable



39D60X

COOLFLOW™ Pump Control Interface Cable

CFP003\*


Poleclamp

CFP006\*

SMARTABLATE™/COOLFLOW™ Pump Flow Rate Verification Kit

\* For availability, please contact your local Biosense Webster representative.  
1 Tubing to be used with COOLFLOW™ Pump.

## Diagram - Fully Integrated COOLFLOW™ Pump System





## Systems



### CARTO™ 3 System

CARTO™ 3 System .....	79
CARTO™ 3 System Accessories.....	79
CARTO™ 3 System Software Module.....	81
CARTO™ 3 Systems software Module OBL (available in Middle East only) .....	82
CARTO™ 3 System Service Contracts.....	83-84
CARTO™ 3 System External Reference Patch .....	85

## CARTO™ 3 System



Cart Configuration

European Ordering No.	Product Description	Included in new CARTO™ 3 System
FG540000	CARTO™ 3 System (PIU, Location Pad, Workstation, Accessories) <sup>3</sup>	X
FG540000U	CARTO™ 3 System (PIU, Location Pad, Workstation, Accessories) <sup>3</sup> OBL*	X
C3RMT**	CARTO™ System RMT Module <sup>1</sup>	
EM540200	CARTO™ 3 System Cart <sup>1</sup>	
EM5402020	CARTO™ 3 Workstation Cart <sup>1</sup>	
MA5402100***	SMARTABLATE™ System cart upgrade kit included in EM540200	
MA540013	Patient Interface Unit (PIU) Table Arm <sup>1</sup>	
CP242101F	24" Computer Monitor <sup>2-3</sup>	XX
KT5400101	GE Recording System-Cables kit (4 cables, KIT needed for customers with GE System) <sup>3</sup>	
KT5400102	Generic Recording System-Cables kit (2 cables, ICOUT & BSECG Cable) <sup>3</sup>	

## CARTO™ 3 System Accessories

European System Ordering No.	Product Description	Included in new CARTO™ 3
EM540003	Pin Box	XX
EM540004	Yellow Pin Box	
KT5400103	DVI-to-DVI Fiber-Optic 98Ft (30m) Shield	X
CW417831F	Fiber-Optica Duplex SC-to-MTRJ 62.5/125W	X
CA540006	Ablation Adaptor Cable	X
CA540011	24V Power Cable	X
CA5400161	Two Pairs STIM Cable - 3m (part of GE Recording System Cable Kit)	

\* OBL refers to Own Brand Labeling. For availability, please contact your local Biosense Webster representative.  
Contact your customer service for lead time information.

\*\* This is make-to-order so expect lead times. Contact your customer service for lead time information.

\*\*\* For availability, please contact your local Biosense Webster representative.

1 This item is optional. For further details, refer to your customer regarding lab set up preferences.

2 Please note that two are needed.

3 Please note that these items are mandatory and needed for a full functional CARTO™ 3 System.

## CARTO™ 3 System Accessories - continued

European Ordering No.	Product Description	Included in CARTO™ 3 System
CA540019	BW-GE BS-ECG Cable (part of GE Recording System)	
CA540021	BW-GE ICEG Cable (part of GE Recording System Kit)	
CA540017	STIM SPLITTER CARTO™ - CARDIOLAB (part GE Recording System Kit)	
EM540015	Foot Pedal Assy	X
EM540020	ECG OUT CONNECTION UNIT (part of Generic Recording System Kit)	
EM5400201	Dual ECG Out Connection Unit	
CA540012	IC Out Cable (part of Generic Recording System Kit)	
CA540151	LOCATION PAD EXTENTION CABLE	X
KT0001429*	Patch Unit Holder Upgrade Kit*	
CA404103	JB PACING AUX.	X
CA404104	PACING EXTENTION CABLE	
CA400811	GROUNDING CABLE	X
CA401214*	POWER CORD (EUR)	X
KT0001427*	Cable Isolator Clamps Kit	
MA5400501	550 mm Location Pad Holder (molded)	X
MA5400801*	750 mm Location Pad Holder (molded, optional)*	
CW-3902-30F	BNC-to-DVI-I Cable 75ohm 30m W/Shrink-Cover <sup>1</sup>	
CW-3903-30F	T-Type15(F) to DVI-I Video Cable 75ohm 30m Molded <sup>1</sup>	
KT0001438*	LPH Mounting Kit (for 90 mm height)	
KT0001439*	LPH Mounting Kit (for 65 mm height)	
KT0001440*	LPH Mounting Kit (for 160 mm width)	
CA5400081	RMT Map Adaptor Cable (part of the CR RMT)	
CW540915N	BS-ECG 5-Lead Set: 5xVi 1.5m	X
KT5401110	ECG card upgrade for PIU (pacing through EP record system)	

\* For availability, please contact your local Biosense Webster representative.

<sup>1</sup> These items are optional and are part of the C3 Sound.

## CARTO™ 3 System Software Module

European Ordering No.	Product Description	Included in new CARTO™ 3 System
KT5400124	CARTOUNIVU™ Software Module	
KT5400130	VISITAG™ Software Module	X
KT5400234	VISITAG SURPOINT™ Software Module for V7 Windows 10	
C3SOUND	CARTOSOUND® Image Integration Software Module	
C3MERGE	CARTOMERGE™ Image Integration Module	
C3CFAE	CFAE Software Module	
C3TOUCH	CARTO SMARTTOUCH™ Module	X
C3PASO	PASO™ Software Module	
C3REPLAY	CARTOREPLAY™ MODULE	
C3CONFD	CARTO™ 3 CONFIDENSE™ Mapping Module	
C3CNFDRPLY	CARTO™ 3 CONFIDENSE™ Mapping Module & CARTOREPLAY™ Module	X
KT5400144	Dual Monitor arm for CARTO™ 3 System Cart (only available for CARTOREPLAY™ Module users)	
MA5402201	Dual Monitor arm for CARTO™ 3 Workstation Cart (only available for CARTOREPLAY™ Module users)	
KT540050	Ablation Index software CD with HP Z4 G4	
KT5400490	Ablation Index software CD	
KT5400155	Carto™ 3 System V6. upgrade kit	X
KT5400170	Ripple Mapping Module	
KT5400160	CARTOSEG™ CT Module	
KT5400165	CARTO™ FAM DX	
KT5400221	Base Upgrade Kit	
KT5400221U	Carto™ 3 System V7 Base Upgrade Kit OBL*	
KT5400233	HELIOSTAR™ Upgrade Kit for V7 (Windows 10)	
KT5400250	CARTO™ 3 System V7 Windows 10 Base Upgrade Kit for OCTARAY™, HELIOSTAR™ and QDOT MICRO™	
KT5400230	QDOT MICRO™ Module Installation Kit for V7	
KT5400235	QDOT MICRO™ Module Installation Kit for V7 Windows 10	
KT5400222X	CARTO PRIME™ Module including CARTO PRIME™ Workstation and V7 Base	
KT5400222U	CARTO PRIME™ Module including CARTO PRIME™ Workstation and V7 Base OBL*	
KT0001420	Green Patch Sensor cable	
KT0001421	Yellow Patch Sensor cable	
CNETSTANDARD	CARTONET™ Integrated CARTO™ Data Management Cloud Solution (annual subscription)	
C3NETRENT	CARTONET™ Integrated CARTO™ Data Management Cloud Solution (monthly fee)	
C3NETGEN	CARTONET™ Integrated CARTO™ Data Management Cloud Solution + nGEN™ bundle	

\* OBL refers to Own Brand Labeling.

## CARTO™ 3 Systems Software Module OBL (available in Middle East only)

Ordering No.	Product Description
FG540000U	Carto™ 3 System (PIU, Location Pad, Workstation, Accessories)
EM5050060U	Catheter to Carto™ 3 System PIU 20 pole eco Cable OBL
EM5050011U	SOUNDSTAR™ 3D eco Cable compatible with Carto™ 3 System OBL
KT5400101U	GE Recording System-Cables kit (4 cables, KIT needed for customers with GE System) OBL
KT5400102U	Generic Recording System-Cables kit (2 cables, ICOUT & BSECG Cable) OBL
KT5400155U	Carto™ 3 System V6 Upgrade Kit OBL
KT5400221U	Carto™ 3 System V7 Base Upgrade Kit OBL
KT5400136U	Carto™ 3 CONFIDENSE™ Mapping Module Upgrade Kit OBL
KT5400145U	Carto™ 3 CONFIDENSE™ Mapping Module & CARTOREPLAY™ Module OBL
KT5400222U	CARTO PRIME™ Module including CARTO PRIME™ Workstation and V7 Base OBL
KT540050U	Ablation Index software CD with HP Z4 G4 OBL
KT5400170U	Ripple Mapping Module OBL
KT5400160U	CARTOSEG™ CT Module OBL
KT5400130U	VISITAG™ Software Module OBL
KT5400124U	CARTOUNIVU™ Software Module OBL
KT5400165U	CARTO™ FAM DX OBL
KT5400111U	CARTOSOUND® Image Integration Software Module OBL
KT5400104U	CARTO SMARTTOUCH™ Module OBL
KT5400262U	VISITAG SURPOINT™ with V7 WS OBL
KT5400162U	VISITAG SURPOINT™ with V6 WS OBL
KT5400490U	Ablation Index Software CD OBL

*Note: Please contact your local Biosense Webster representatives or the technical team to order these codes.*

\* OBL refers to Own Brand Labeling.

## CARTO™ 3 System Service Contracts

European Ordering No.	Product Description
<b>C3PREMY1</b>	CARTO™ 3 System PREMIUM Service Contract 1 Year
<b>C3PREMUPY1</b>	Upgrade to CARTO™ 3 System PREMIUM Service Contract 1 Year - under warranty period
<b>C3PREMY2</b>	CARTO™ 3 System PREMIUM Service Contract 2 Years
<b>C3PREMUPY2</b>	Upgrade to CARTO™ 3 System PREMIUM Service Contract 2 Years
<b>C3PREMY3</b>	CARTO™ 3 System PREMIUM Service Contract 3 Years
<b>C3PREMUPY3</b>	Upgrade to CARTO™ 3 System PREMIUM Service Contract 3 Years
<b>C3EXTY1</b>	Extended Warranty Policy for CARTO™ 3 System 1 year
<b>C3EXTY2</b>	Extended Warranty Policy for CARTO™ 3 System 2 years
<b>C3EXTY3</b>	Extended Warranty Policy for CARTO™ 3 System 3 years
<b>C3LSBASE</b>	Lab Solution Base (Service contract on Carto™ 3, generator & pump as well as existing and future software technology)
<b>C3LSADV</b>	Lab Solution Advanced (Service contract on Carto™ 3, generator & pump and existing as well as future software technology incl. CARTONET license)
<b>C3LSPLUS</b>	Lab Solution Plus (Service contract on Carto™ 3, generator & pump as well as existing and future software technology incl. CARTONET™ license and generator upgrade)

*Note: Please contact your local Biosense Webster representatives or the technical team to order these codes.*



## CARTO™ 3 System Service Contracts

---

Biosense Webster delivers high quality products and we also do our best to exceed your expectations with our high level of service provided by our professionally trained Field Service Engineers. The Technical Service Team has built a reputation of maintaining a high level of system performance through fast and efficient service.

### • Premium Service Contract

---

The **Premium Service Contract** ensures that your CARTO™ 3 System is protected against any accidental damage and available for use at all times. In addition to the standard warranty terms, the upgrade to the **Premium Service** provides the following enhanced coverage:

- Replacement / Repair for any accidental damage (including any system parts)
- Replacement of any kind of system cable
- 2 business days of service response time guaranteed
- Multiple system relocations
- Priority in repair and software updates over customers without *Premium service contract*

After the warranty period, the **Premium Service Contract** will be available to cover - besides the above mentioned services -all the manufacturer defects and access to loaner unit.

### • Warranty Extension

---

The **Warranty Extension** includes the following services:

- All repairs for system breakdown (including labor, parts and travel)
- 12-Hour Call Center telephone support
- Software updates\*
- One relocation
- One breakage of each defined system cable and accessories per year
- One accidental hardware (Patient Interface Unit, Location Pad, Workstation) breakage per year


The **Warranty Extension** is available for 1 year, 2 years or 3 years.

Please check prices with your Biosense Webster sales representative.

---

\* Software update as part of Biosense Webster continuous dedication to enhance system software development. This does not cover software update including new software functionalities/ module such as CARTO™ 3 SMARTTOUCH™ Software Module, etc.

CARTO™ 3 System External Reference Patch



Set of CARTO™ 3 System External Reference Patches

European  
Ordering No.

Product Description

CREFP6

CARTO™ 3 System External Reference Patches - Set of 6



## Cables

### Charts

CONNECTION OF CHOICE™ Chart.....	87-88
----------------------------------	-------

### Cables

CARTO™ 3 System Cables.....	89
Interface Cable EP Shuttle RF Generator .....	90
Interface SMARTABLATE™ System Cables.....	90

### Non-Navigation


Non-Navigation Diagnostic Cables.....	91
Non-Navigation Ablation Cables.....	91

### Reverse Cables


Reverse Cables .....	92
----------------------	----

# CARTO™ 3 CONNECTION OF CHOICE™ Chart

## System



# CONNECTION OF CHOICE™ Chart



## CY1212CT

WEBSTER™ Quadrapolar Catheters with Auto ID - Fixed and Deflectable

WEBSTER™ Octapolar Catheters with Auto ID - Deflectable

LASSO™ Decapolar Catheters with Auto ID - Deflectable

CARTO VIZIGO™ Bi-Directional Guiding Sheath

## CY1210CT

Existing WEBSTER™ Quadrapolar Catheters - Fixed and Deflectable

Existing WEBSTER™ Diagnostic Catheters with up to 10 poles

Existing WEBSTER™ Octapolar Catheters - Deflectable

Existing LASSO™ Decapolar Catheters - Deflectable



20 POLE A/B interface cable from SPU to PIU.

For more OCTARAY™ connections refer to the OCTARAY™ SPU Connectology Deck



Connect Carto™ Signal Processing Unit to front of PIU with 20 Pole A/B interface cable



## CARTO™ 3 System Cables

### Diagnostic CARTO™ 3 System Cables

European Ordering No.	Length	Product Description	Connector	
			at catheter	at PIU
CB3434CT	3 m	Cable for LASSO™ 2515 NAV Catheter, HALO™, and ISMUS™ Catheter with Auto ID, LASSO™ Auto ID, DECANAV™ Catheter	34 pin dark blue	34 pin light blue
CB3412CT	3 m	Cable for WEBSTER™ Decapolar Catheter (Fixed and Deflectable) and WEBSTER™ CS Bi-Directional Catheter Auto ID	12 pin dark blue	34 pin light blue
CB3410CT	3 m	Cable for WEBSTER™ CS Decapolar Catheter	10 pin black	34 pin light blue
CY1212CT	3 m	Cable for WEBSTER™ Octapolar with Auto ID, Webster Quadrapolar with Auto ID, Webster CS, Webster HIS and LASSO™ Decapolar with Auto ID Catheter	12 pin dark blue	12 pin yellow
CY1210CT	3 m	Cable for WEBSTER™ Catheter with up to 10 poles	10 pin black	12 pin yellow

European Ordering No.	Length	Product Description	Connector	
			at catheter	at CARTO™ 3 System
D128624	3 m	Cable for 20-pole catheters, non Auto ID	2 x 10 pin black	2 x 10 pin black

### Diagnostic CARTO™ 3 eco System Cables



European Ordering No.	Length	Product Description	Connector	
			at catheter	at PIU
D134401	2.7 m	Cable for LASSO™ NAV eco Circular Mapping Catheter, LASSO™ 2515 NAV eco Circular Mapping Catheter and the PENTARAY™ NAV eco High Density Mapping Catheter		
D134402		OCTARAY™ Mapping Catheter connection cable		
EM5050060		34 pin eco Adaptor		

### UltraSound CARTO™ 3 System Cables

European Ordering No.	Length	Product Description	Connector	
			at catheter	at PIU
CG2025CT	3 m	Cable for SOUNDSTAR™ Catheter	25 pin black	20 pin green



### UltraSound eco CARTO™ 3 System Cables

European Ordering No.	Length	Product Description	Connector	
			at catheter	at PIU
EM5050011	6 m	Permanent Cable	25 pin black	20 pin green

### Therapeutic CARTO™ 3 System Cables

European Ordering No.	Length	Product Description	Connector	
			at catheter	at PIU
CR3425CT	3 m	Cable for existing NAVISTAR™ Catheter	25 pin black	34 pin red
CR3434CT	3 m	Cable for THERMOCOOL SMARTTOUCH™ Catheter	34 pin red	34 pin red

## Interface SMARTABLATE™ System Cables

European Ordering No.	Product Description	Sterile
<b>M490034*</b>	SMARTABLATE™ Erbe Indifferent Electrode Cable (3 m)	
<b>M490035*</b>	SMARTABLATE™ 3M Indifferent Electrode Cable (3 m)	
<b>VALLEYLAB</b>	REM Polyhesive II E7507 (Return Patch)	
<b>D129704**</b>	CARTO™ 3 System Interface Cable / Connecting Generator to the CARTO™ System, 10ft (305 cm)	
<b>D129707**</b>	CARTO™ 3 System Interface Cable / Connecting Generator to the CARTO™ 3 System, 4ft (122 cm)	
<b>D130302</b>	Ablation Cable from SMARTABLATE™ System to CELSIUS™ Catheter or CARTO™ XP System PIU (TC/THR)	X
<b>D130303</b>	Ablation Cable from SMARTABLATE™ System to CELSIUS™ Catheter or CARTO™ XP System PIU (DS)	X

\* This product is not regulatory approved and commercially available. Please contact your local BW representative for availability.

\*\* Please order one unit of either D129704 or D129707 per SMARTABLATE™ System that will connect to a CARTO™ 3 System.

## Non-Navigation Diagnostic Cables

European Ordering No.	Length	Product Description	Auto- clavable	Connector	
				at catheter	at generator
39F40R	0.9 m	For 2 and 4 pole catheters	yes	10-Pin	4 Protected Tip Pins
39F55R	1.8 m	For 2 and 4 pole catheters	yes	10-Pin	4 Protected Tip Pins
39F43R	0.9 m	For 2, 4 and 6 pole catheters	yes	10-Pin	6 Protected Tip Pins
39F44R	0.9 m	For 2, 4, 6 and 8 pole catheters	yes	10-Pin	8 Protected Tip Pins
39F41R	0.9 m	For 2, 4, 6, 8 and 10 pole catheters	yes	10-Pin	10 Protected Tip Pins
39F54R	1.8 m	For 2, 4, 6, 8 and 10 pole catheters	yes	10-Pin	10 Protected Tip Pins
39F47D	0.9 m	For use with Porterfield Catheter			

## Non-Navigation Ablation Cables

### For the Stockert EP Shuttle (with Biosense Webster Catheter)

European Ordering No.	Length	Product Description	Connector	
			at catheter	at generator
39E43R	3 m	Cable for CELSIUS™ Ablation Catheter (TC/THR)	10-Pin	10-Pin
39E68R	3 m	Cable for CELSIUS™ DS Ablation Catheter (TC)	10-Pin	10-Pin

### For the Stockert EP Shuttle (with foreign catheter)

European Ordering No.	Length	Product Description	Connector	
			at catheter	at generator
39E46R	3 m	Cable for Medtronic™ <sup>1</sup> Catheter, TC	10-Pin	10-Pin
39E53S	3 m	Cable for EPT Catheter, THR	9-Pin	10-Pin

### For connection to a foreign generator with Celsius Catheter

European Ordering No.	Length	Product Description	Connector	
			at catheter	at generator
39E56S	3 m	Cable for Osypka™ <sup>1</sup> RF Generator HAT 300 smart, THR	10-Pin	9-Pin
39F23R	3 m	Cable for CELSIUS™ DS to Osypka™ <sup>1</sup> HAT 300 smart	10-Pin	9-Pin
39E69R	3 m	Cable for EPT™ <sup>1</sup> (Box APM), THR	10-Pin	9-Pin

<sup>1</sup> The third party trademarks used herein are trademarks of their respective owners.



## Reverse Cables

European Ordering No.	Product Description	Connector	
		at catheter	at PIU
D128709	Interface cable for use with BW catheters and EP recording system or PIU pinbox	12-Pin	12-Pin
D128622	Interface cable for use with BW catheters and EP recording system or PIU pinbox	34-Pin	20-Pin

# Index



# PRODUCT NAME INDEX

ACUNAV™ Ultrasound Catheter .....	08255790 .....	65
ACUNAV™ Ultrasound Catheter .....	10135910 .....	66
ACUNAV™ Ultrasound Catheter .....	10043342 .....	65
ACUNAV™ Ultrasound Catheter .....	10135936 .....	66
AVAIL™ Bipolar Fixed Shape Catheter .....	37P0RR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37C0CR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37C3FR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37C5RR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37C5CR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37D0CR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37D5CR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37D0RR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37D3FR .....	53
AVAIL™ Quadrapolar Fixed Shape Catheter .....	37D5RR .....	53
CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath .....	FCL20100 .....	57
CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath .....	FCL20101 .....	57
CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath .....	FCL20102 .....	57
CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath .....	FCL20103 .....	57
CARTO™ 3 System .....	FG540000 .....	79
CARTO™ 3 System .....	FG540000U .....	79
CARTO™ 3 System .....	C3RMT .....	79
CARTO™ 3 System .....	EM540200 .....	79
CARTO™ 3 System .....	EM5402020 .....	79
CARTO™ 3 System .....	MA540013 .....	79
CARTO™ 3 System .....	CP242101F .....	79
CARTO™ 3 System .....	KT5400101 .....	79
CARTO™ 3 System .....	KT5400102 .....	79
CARTO™ 3 System Accessories .....	EM540003 .....	79
CARTO™ 3 System Accessories .....	EM540004 .....	79
CARTO™ 3 System Accessories .....	KT5400103 .....	79
CARTO™ 3 System Accessories .....	CW417831F .....	79
CARTO™ 3 System Accessories .....	CA540006 .....	79
CARTO™ 3 System Accessories .....	CA540011 .....	79
CARTO™ 3 System Accessories .....	CA5400161 .....	79
CARTO™ 3 System Accessories .....	CA540019 .....	80
CARTO™ 3 System Accessories .....	CA540021 .....	80
CARTO™ 3 System Accessories .....	CA540017 .....	80
CARTO™ 3 System Accessories .....	EM540015 .....	80
CARTO™ 3 System Accessories .....	EM540020 .....	80
CARTO™ 3 System Accessories .....	EM5400201 .....	80
CARTO™ 3 System Accessories .....	CA540012 .....	80
CARTO™ 3 System Accessories .....	CA540151 .....	80
CARTO™ 3 System Accessories .....	KT0001429 .....	80
CARTO™ 3 System Accessories .....	CA404103 .....	80
CARTO™ 3 System Accessories .....	CA404104 .....	80
CARTO™ 3 System Accessories .....	CA400811 .....	80
CARTO™ 3 System Accessories .....	CA401214 .....	80
CARTO™ 3 System Accessories .....	KT0001427 .....	80
CARTO™ 3 System Accessories .....	MA5400501 .....	80
CARTO™ 3 System Accessories .....	MA5400801 .....	80
CARTO™ 3 System Accessories .....	CW-3902-30F .....	80
CARTO™ 3 System Accessories .....	CW-3903-30F .....	80
CARTO™ 3 System Accessories .....	KT0001438 .....	80
CARTO™ 3 System Accessories .....	KT0001439 .....	80
CARTO™ 3 System Accessories .....	KT0001440 .....	80
CARTO™ 3 System Accessories .....	CA5400081 .....	80
CARTO™ 3 System Accessories .....	CW540915N .....	80
CARTO™ 3 System Accessories .....	CW540906N .....	80
CARTO™ 3 System Cables .....	CB3434CT .....	89
CARTO™ 3 System Cables .....	CB3412CT .....	89
CARTO™ 3 System Cables .....	CB3410CT .....	89
CARTO™ 3 System Cables .....	CY1212CT .....	89
CARTO™ 3 System Cables .....	CY1210CT .....	89
CARTO™ 3 System Cables .....	D128624 .....	89
CARTO™ 3 System Cables .....	D134401 .....	89
CARTO™ 3 System Cables .....	D13440 .....	89
CARTO™ 3 System Cables .....	EM5050060 .....	89
CARTO™ 3 System Cables .....	CG2025CT .....	89
CARTO™ 3 System Cables .....	EM5050011 .....	89
CARTO™ 3 System Cables .....	CR3425CT .....	89
CARTO™ 3 System Cables .....	CR3434CT .....	89
CARTO™ 3 System Software Module .....	KT5400124 .....	81
CARTO™ 3 System Software Module .....	KT5400127 .....	81

# PRODUCT NAME INDEX

CARTO™ 3 System Software Module	KT5400128	81
CARTO™ 3 System Software Module	KT5400129	81
CARTO™ 3 System Software Module	KT5400134	81
CARTO™ 3 System Software Module	KT5400134U	81
CARTO™ 3 System Software Module	KT5400130	81
CARTO™ 3 System Software Module	C3SOUND	81
CARTO™ 3 System Software Module	C3MERGE	81
CARTO™ 3 System Software Module	C3CFAE	81
CARTO™ 3 System Software Module	C3TOUCH	81
CARTO™ 3 System Software Module	C3PASO	81
CARTO™ 3 System V7 Base Upgrade Kit	KT5400250	81
CARTO™ 3 System V7 Base Upgrade Kit	KT5400221	81
CARTO™ 3 System V7 Base Upgrade Kit OBL	KT5400221U	81
CARTO™ 3 Software System Module	C3REPLAY	81
CARTO™ 3 Software System Module	C3CONFD	81
CARTO™ 3 Software System Module	C3CNFDRPLY	81
CARTO™ 3 Software System Module	C3V4	81
CARTO™ 3 System Service Contracts	C3PREMY1	83
CARTO™ 3 System Service Contracts	C3PREMUPY1	83
CARTO™ 3 System Service Contracts	C3PREMY2	83
CARTO™ 3 System Service Contracts	C3PREMUPY2	83
CARTO™ 3 System Service Contracts	C3PREMY3	83
CARTO™ 3 System Service Contracts	C3PREMUPY3	83
CARTO™ 3 System Service Contracts	C3EXTY1	83
CARTO™ 3 System Service Contracts	C3EXTY2	83
CARTO™ 3 System Service Contracts	C3EXTY3	83
CARTO™ 3 System External Reference Patch	CREFP6	85
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	FG540000U	96
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	EM5050060U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	EM5050011U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400101U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400102U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400155U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400221U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400136U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400145U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400222U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT540050U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400170U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400160U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400130U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400124U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400165U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400111U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400104U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400262U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400162U	82
CARTO™ 3 Systems Software Module OBL (available in Middle East only)	KT5400490U	82
CARTONET (ANNUAL SUBSCRIPTION)	CNETSTANDARD	81
CARTONET (MONTHLY FEE)	C3NETRENT	81
CARTONET + NGEN BUNDLE	C3NETGEN	81
CARTOSEG™ CT Module	KT5400160	81
CARTO PRIME™ Module including CARTO PRIME™ Workstation and V7 Base	KT5400222X	81
CARTO PRIME™ Module including CARTO PRIME™ Workstation and V7 Base OBL	KT5400222U	81
CARTO VIZIGO™ Sheaths	D138501	67
CARTO VIZIGO™ Sheaths	D138502	58
CARTO VIZIGO™ Sheaths	D138503	58
CARTO VIZIGO™ Sheaths	CB3412CT	58
CARTO VIZIGO™ Sheaths	CY1212CT	58
CELSIUS™ 4 mm Braided Tip Catheter	34D37R	33
CELSIUS™ 4 mm Braided Tip Catheter	34D47R	33
CELSIUS™ 4 mm Braided Tip Catheter	34D57R	33
CELSIUS™ 4 mm Braided Tip Catheter	D130302	33
CELSIUS™ 4 mm Catheter Thermistor	36O07R	33
CELSIUS™ 4 mm Catheter Thermistor	36O17R	33
CELSIUS™ 4 mm Catheter Thermistor	36O27R	33
CELSIUS™ 4 mm Catheter Thermistor	36O37R	33
CELSIUS™ 4 mm Catheter Thermistor	36O47R	33
CELSIUS™ 4 mm Catheter Thermistor	36O57R	33
CELSIUS™ 4 mm Catheter Thermistor	D130302	33
CELSIUS™ 4 mm Catheter Thermocouple	35G07R	32
CELSIUS™ 4 mm Catheter Thermocouple	35G17R	32
CELSIUS™ 4 mm Catheter Thermocouple	35G37R*	32

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# PRODUCT NAME INDEX

CELSIUS™ 4 mm Catheter Thermocouple.....	36N07R*	32
CELSIUS™ 4 mm Catheter Thermocouple.....	36N17R	32
CELSIUS™ 4 mm Catheter Thermocouple.....	36N27R	32
CELSIUS™ 4 mm Catheter Thermocouple.....	36N37R	32
CELSIUS™ 4 mm Catheter Thermocouple.....	36N47R*	32
CELSIUS™ 4 mm Catheter Thermocouple.....	36N57R	32
CELSIUS™ 4 mm Catheter Thermocouple.....	D130302	32
CELSIUS™ DS Catheter.....	35S17R	35
CELSIUS™ DS Catheter.....	35S37R	35
CELSIUS™ DS Catheter.....	35S57R	35
CELSIUS™ DS Catheter.....	D130303	35
CELSIUS FLTR™ 8 mm Uni-Directional Catheter.....	D135501*	34
CELSIUS FLTR™ 8 mm Uni-Directional Catheter.....	D130303	34
CELSIUS FLTR™ 8 mm Uni-Directional Catheter.....	D135502*	34
CELSIUS FLTR™ 8 mm Bi-Directional Catheter.....	D135901	34
CELSIUS FLTR™ 8 mm Bi-Directional Catheter.....	D130303	34
CELSIUS™ RMT Catheter.....	CR7TCS4RT	35
CELSIUS™ RMT Catheter.....	D130302	35
CELSIUS™ RMT THERMOCOOL™ Catheter.....	CR7TCSIRT	32
CELSIUS™ RMT THERMOCOOL™ Catheter.....	D130302	32
CELSIUS™ RMT THERMOCOOL™ Catheter.....	SAT001	32
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	35Q13R*	31
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	35Q33R	31
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	35Q53R	31
CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	35Q73R*	31
COHEREX WAVECREST™ LAA Implant.....	WCR1503	62
COHEREX WAVECREST™ LAA Implant.....	WCR1513	62
COHEREX WAVECREST™ LAA Implant.....	WCR1523	62
COHEREX WAVECREST™ LAA Implant.....	WCR1530	62
COHEREX WAVECREST™ LAA Implant.....	WCR1540	62
COHEREX WAVECREST™ LAA Implant.....	WCR1541	62
COHEREX WAVECREST™ LAA Implant.....	WCR1551	62
COOLFLOW™ Pump Accessories.....	39D03X	76
COOLFLOW™ Pump Accessories.....	39D14X	76
COOLFLOW™ Pump Accessories.....	39D60X	76
COOLFLOW™ Pump Accessories.....	39D60X10	75
COOLFLOW™ Pump Accessories.....	CFP003	76
COOLFLOW™ Pump Accessories.....	CFP006	76
COOLFLOW™ Pump Tubing.....	CFT003	76
CRISTACATH™ Catheter.....	36Y39R	54
DECANAV™ Catheter.....	R7D282CT	39
DECANAV™ Catheter.....	R7F282CT	39
EZ STEER™ 4 mm Bi-Directional Catheter.....	35I35R	30
EZ STEER™ 4 mm Bi-Directional Catheter.....	35I5JR	30
EZ STEER™ 4 mm Bi-Directional Catheter.....	35I33R	30
EZ STEER™ 4 mm Bi-Directional Catheter.....	35I55R	30
EZ STEER™ DS Bi-Directional Catheter.....	35U35R	31
EZ STEER™ DS Bi-Directional Catheter.....	35U5JR	31
EZ STEER™ DS Bi-Directional Catheter.....	35U33R	31
EZ STEER™ DS Bi-Directional Catheter.....	35U55R	31
EZ STEER™ DS Bi-Directional Catheter.....	D130302	31
EZ STEER™ DS Bi-Directional Catheter.....	D130303	31
EZ STEER™ DS Bi-Directional Catheter.....	SAT001	31
EZ STEER™ DS Bi-Directional Catheter.....	35UJJR	31
EZ STEER™ THERMOCOOL™ Catheter.....	36Q33M*	30
EZ STEER™ THERMOCOOL™ Catheter.....	36Q35M	30
EZ STEER™ THERMOCOOL™ Catheter.....	36Q55M	30
EZ STEER™ THERMOCOOL™ Catheter.....	36Q5JM	30
EZ STEER™ THERMOCOOL™ Catheter.....	D130302	30
EZ STEER™ THERMOCOOL™ Catheter.....	SAT001	30
FAM DX Module.....	KT5400165	40
FAM DX Module OBL.....	KT5400165U	40
GREEN PATCH SENSOR CABLE.....	KT0001420	81
GUIDESTAR™ Deflectable Sheath.....	D141101	36
GUIDESTAR™ Deflectable Sheath.....	D141102	36
GUIDESTAR™ Deflectable Sheath.....	D141103	36
HALO™ XP Catheter.....	36J13RCT	45
HALO™ XP Catheter.....	D128622	45
HALO™ XP Tricuspid Mapping Catheter.....	36J13R	52
HALO™ XP Tricuspid Mapping Catheter.....	36J16R	52
HALO™ XP Tricuspid Mapping Catheter.....	36J14R	52
HEARTSPAN™ Transseptal Needles.....	FND01901	59
HEARTSPAN™ Transseptal Needles.....	FND01902	59

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# PRODUCT NAME INDEX

HEARTSPAN™ Transseptal Needles .....	FND01903 .....	59
HEARTSPAN™ Transseptal Needles .....	FND01905 .....	59
HEARTSPAN™ Transseptal Needles .....	FND01906 .....	59
HELIOSTAR™ Balloon Ablation Catheter .....	D138904 .....	36
HELIOSTAR™ Upgrade Kit .....	KT5400233 .....	81
ISMUS™ Catheter .....	36Y32RCT .....	45
ISMUS™ Catheter .....	D128622 .....	45
ISMUS™ Catheter .....	36Y32R .....	52
Indifferent Electrodes and Corresponding Connection Cables .....	39D12X .....	73
Indifferent Electrodes and Corresponding Connection Cables .....	39D12X5 .....	73
Indifferent Electrodes and Corresponding Connection Cables .....	39D16X .....	73
Indifferent Electrodes and Corresponding Connection Cables .....	39D16X5 .....	73
Indifferent Electrodes and Corresponding Connection Cables .....	39D82X .....	73
Indifferent Electrodes and Corresponding Connection Cables .....	39F27R .....	73
Interface Cable EP Shuttle RF Generator .....	D129705 .....	90
Interface SMARTABLATE™ System Cables .....	M490034 .....	90
Interface SMARTABLATE™ System Cables .....	M490035 .....	90
Interface SMARTABLATE™ System Cables .....	M490043 .....	90
Interface SMARTABLATE™ System Cables .....	M490044 .....	90
Interface SMARTABLATE™ System Cables .....	D129704 .....	90
Interface SMARTABLATE™ System Cables .....	D129707 .....	90
Interface SMARTABLATE™ System Cables .....	D130302 .....	90
Interface SMARTABLATE™ System Cables .....	D130303 .....	90
LAB SOLUTION BASE .....	C3LSBASE .....	82
LAB SOLUTION ADVANCED .....	C3LSADV .....	82
LAB SOLUTION PLUS .....	C3LSPLUS .....	82
LASSO™ 2515 NAV eco Circular Mapping Catheter .....	D134302 .....	38
LASSO™ 2515 NAV eco Circular Mapping Catheter .....	D134301 .....	38
LASSO™ 2515 Variable Mapping Catheter .....	35O2515R* .....	46
LASSO™ 2515 Variable Mapping Catheter .....	35T2515R* .....	46
LASSO™ Decapolar Catheter .....	35O26RCT .....	41
LASSO™ Decapolar Catheter .....	35O36RCT .....	41
LASSO™ Duo-Decapolar Catheter .....	D7L2015CT* .....	41
LASSO™ Duo-Decapolar Catheter .....	D7L2020CT* .....	41
LASSO™ Mapping Catheter .....	35O16R* .....	46
LASSO™ Mapping Catheter .....	35O26R* .....	46
LASSO™ Mapping Catheter .....	35O36R* .....	46
LASSO™ Mapping Catheter .....	35O46R* .....	46
LASSO™ Mapping Catheter .....	35O66R* .....	46
LASSO™ Mapping Catheter .....	35T26R* .....	46
LASSO™ Mapping Catheter .....	35T36R* .....	46
LASSO™ NAV eco Circular Mapping Catheter .....	D134901 .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	D134903 .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	D134905 .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	D134902 .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	D134904 .....	38
LASSO™ NAV eco Circular Mapping Catheter .....	D134906 .....	38
LASSOSTAR™ Circular Catheter .....	D139001 .....	36
LASSOSTAR™ Circular Catheter .....	D139002 .....	36
LASSOSTAR™ Circular Catheter .....	D139003 .....	36
LASSOSTAR™ NAV Circular Catheter .....	D140401 .....	36
LASSOSTAR™ NAV Circular Catheter .....	D140402 .....	36
LASSOSTAR™ NAV Circular Catheter .....	D140403 .....	36
MOBICATH™ Sheaths .....	D140010* .....	58
MOBICATH™ Sheaths .....	D140011* .....	58
NAVISTAR™ 4 mm Catheter .....	34A15M .....	26
NAVISTAR™ 4 mm Catheter .....	34A25M* .....	26
NAVISTAR™ 4 mm Catheter .....	34A35M .....	26
NAVISTAR™ 4 mm Catheter .....	34A45M* .....	26
NAVISTAR™ 4 mm Catheter .....	34A55M .....	26
NAVISTAR™ 4 mm Catheter .....	34AJ5M* .....	26
NAVISTAR™ 4 mm Catheter .....	D129704 .....	26
NAVISTAR™ DS Catheter .....	34J27M .....	26
NAVISTAR™ DS Catheter .....	34J37M .....	26
NAVISTAR™ DS Catheter .....	34J57M .....	26
NAVISTAR™ DS Catheter .....	34JJ7M .....	26
NAVISTAR™ DS Catheter .....	D129704 .....	26
NAVISTAR™ RMT 4 mm Catheter .....	NR7TCS4Y .....	27
NAVISTAR™ RMT 4 mm Catheter .....	D129704 .....	27
NAVISTAR™ RMT THERMOCOOL™ Catheter .....	NR7TCSIY .....	27
NAVISTAR™ RMT THERMOCOOL™ Catheter .....	D129704 .....	27
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	36H33M .....	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	36H35M .....	25

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# PRODUCT NAME INDEX

NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	36H55M*	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	36H5JM	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	36HJMM	25
NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter .....	D129704	25
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	34H17M	25
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	34H27M	25
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	34H37M	25
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	34H57M	25
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	34HJ7M	25
NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter .....	D129704	25
nGEN™ Generator Kit .....	D138401	80
nGEN™ Monitor .....	D138801	69
nGEN™ Pump .....	D139701	69
nGEN™ Installation Accessories .....	EM540200	69
nGEN™ Installation Accessories .....	MA5402300	69
nGEN™ Clamp .....	M583196SP	69
nGEN™ Clamp .....	M549104SP	69
nGEN™ Foot Pedal .....	M581802SP	69
nGEN™ Generator Catheter Interface Cables .....	D140102	69
nGEN™ Generator Catheter Interface Cables .....	D135703	69
nGEN™ Generator Catheter Interface Cables .....	D133701	69
nGEN™ Generator Catheter Interface Cables .....	D133703	69
nGEN™ System Interface Cables .....	M581101SP	69
nGEN™ System Interface Cables .....	M583199SP	69
nGEN™ System Interface Cables .....	M5831100SP	69
nGEN™ System Interface Cables .....	M5831101SP	69
nGEN™ System Interface Cables .....	M583197SP	69
nGEN™ System Interface Cables .....	M583198SP	69
nGEN™ System Interface Cables .....	M589704SP	69
nGEN™ System Interface Cables .....	M589705SP	69
nGEN™ System Interface Cables .....	M600203SP	69
nGEN™ System Interface Cables .....	M600202SP	69
nGEN™ System Interface Cables .....	M600201SP	69
nGEN™ System Interface Cables .....	M600301SP	69
nGEN™ System Interface Cables .....	M581007SP	69
nGEN™ System Interface Cables .....	M581006SP	69
nGEN™ AC Power Cord .....	M621410SP	69
nGEN™ AC Power Cord .....	M621402SP	69
nGEN™ AC Power Cord .....	M621406SP	69
nGEN™ AC Power Cord .....	M621407SP	69
nGEN™ Pump AC Adapter .....	M529723SP	69
nGEN™ Pump AC Power Cord .....	M649010SP	69
nGEN™ Pump AC Power Cord .....	M649002SP	69
nGEN™ Pump AC Power Cord .....	M649006SP	69
nGEN™ Pump AC Power Cord .....	M649007SP	69
nGEN™ Monitor Power Cord .....	M583191SP	69
nGEN™ Irrigation Tubing Set .....	SAT001	69
Non-Navigation ABLATION CABLES .....	39E43R	91
Non-Navigation ABLATION CABLES .....	39E68R	91
Non-Navigation ABLATION CABLES .....	39E46R	91
Non-Navigation ABLATION CABLES .....	39E53S	91
Non-Navigation ABLATION CABLES .....	39E56S	91
Non-Navigation ABLATION CABLES .....	39F23R	91
Non-Navigation ABLATION CABLES .....	39E69R	91
Non-Navigation ABLATION CABLES .....	39E23R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F40R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F55R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F43R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F44R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F41R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F54R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F24R	91
Non-Navigation DIAGNOSTIC CABLES .....	39F47D	91
OCTARAY™ Catheter small .....	D160901	39
OCTARAY™ Catheter large .....	D160902	39
OCTARAY™ Catheter large .....	D160903	39
OCTARAY™ Catheter small .....	D160904	39
OCTARAY™ Catheter large .....	D160905	39
OCTARAY™ Catheter large .....	D160906	39
PARAHISIAN™ Catheter .....	36L00Q*	54
Pediatric Catheter .....	36S30Q*	54
PENTARAY™ NAV eco High Density Mapping Catheter .....	D128207	39
PENTARAY™ NAV eco High Density Mapping Catheter .....	D128208	39

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.



# PRODUCT NAME INDEX

PENTARAY™ NAV eco High Density Mapping Catheter .....	D128210 .....	39
PENTARAY™ NAV eco High Density Mapping Catheter .....	D128211 .....	39
PREFACE™ Fixed Guiding Sheaths .....	301803A .....	56
PREFACE™ Fixed Guiding Sheaths .....	301803P .....	56
PREFACE™ Fixed Guiding Sheaths for EP .....	301803MS .....	56
PREFACE™ Fixed Guiding Sheaths for EP .....	301803M .....	56
PREFACE™ Fixed Guiding Sheaths for EP .....	301805M .....	56
QDOT MICRO™ Uni-directional Catheter .....	D139401 .....	20
QDOT MICRO™ Uni-directional Catheter .....	D139402 .....	20
QDOT MICRO™ Uni-directional Catheter .....	D139403 .....	20
QDOT MICRO™ Bi-directional Catheter .....	D139501 .....	20
QDOT MICRO™ Bi-directional Catheter .....	D139502 .....	20
QDOT MICRO™ Bi-directional Catheter .....	D139503 .....	20
QDOT MICRO™ Bi-directional Catheter .....	D139504 .....	20
QDOT MICRO™ Bi-directional Catheter .....	D139505 .....	20
QDOT MICRO™ Module Installation Kit for V7 .....	KT5400230 .....	81
QDOT MICRO™ Module Installation Kit for V7 Windows 10 .....	KT5400235 .....	81
Ripple Mapping Module .....	KT5400170 .....	81
REVERSE CABLES .....	D128709 .....	92
REVERSE CABLES .....	D128622 .....	92
Service Contracts .....	SHUTBC1 .....	75
Service Contracts .....	SHUTPV1 .....	75
Service Contracts .....	SHUTPVP1 .....	75
Service Contracts .....	CFPPV1 .....	77
Service Contracts .....	CFPPVP1 .....	77
SMARTABLATE™ System .....	M4900106 .....	70
SMARTABLATE™ System .....	SAT001 .....	70
SMARTABLATE™ System .....	M490005 .....	70
SMARTABLATE™ System .....	M490012 .....	70
SMARTABLATE™ System .....	M490010 .....	70
SMARTABLATE™ System .....	M490029 .....	70
SMARTABLATE™ System .....	M490013 .....	70
SMARTABLATE™ System .....	M490026 .....	70
SMARTABLATE™ System .....	M490027 .....	70
SMARTABLATE™ System .....	M490016 .....	70
SMARTABLATE™ System .....	M490028 .....	70
SMARTABLATE™ System .....	M490030 .....	70
SMARTABLATE™ System .....	M490021 .....	70
SMARTABLATE™ System .....	M490023 .....	70
SMARTABLATE™ System .....	M490025 .....	70
SMARTABLATE™ System .....	M490014 .....	70
SMARTABLATE™ System .....	M490017 .....	70
SMARTABLATE™ System .....	M490018 .....	70
SMARTABLATE™ System .....	M490037 .....	70
SMARTABLATE™ System .....	M4900602 .....	70
SMARTABLATE™ System .....	M4900603 .....	70
SMARTABLATE™ System .....	M4900605 .....	70
SMARTABLATE™ System .....	M4900606 .....	70
SMARTABLATE™ System .....	MA5402100 .....	70
SMARTABLATE™ System .....	M490034 .....	70
SMARTABLATE™ System .....	M490035 .....	70
SOUNDSTAR™ 3D eco Catheter .....	10438577 .....	64
SOUNDSTAR™ 3D eco Catheter .....	10439072 .....	64
SOUNDSTAR™ 3D Ultrasound Catheter .....	SNDSTR10 .....	65
SOUNDSTAR™ 3D Ultrasound Catheter .....	SNDSTR10G .....	65
SOUNDSTAR™ eco 8F Catheter .....	10439011 .....	64
SOUNDSTAR™ eco 8F Catheter .....	10439236 .....	64
Spare Parts and Others .....	39D03X .....	74
Spare Parts and Others .....	39D06X .....	74
Spare Parts and Others .....	39D14X .....	74
Spare Parts and Others .....	39D21X .....	74
Spare Parts and Others .....	39D83X .....	74
Spare Parts and Others .....	39D12X .....	74
Spare Parts and Others .....	39E43R .....	74
Stockert EP Shuttle Accessories .....	39F28X10 .....	75
Stockert EP Shuttle Accessories .....	39F28X15 .....	75
Stockert EP Shuttle Accessories .....	39F28X20 .....	75
Stockert EP Shuttle Accessories .....	39F28X30 .....	75
Stockert EP Shuttle Accessories .....	39D60X .....	75
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	BNI35BDH* .....	23
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	BNI35BFH* .....	23
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	BNI35DFH .....	23
THERMOCOOL™ SF NAV Bi-Directional Catheter .....	BNI35FFH* .....	23

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.



# PRODUCT NAME INDEX

THERMOCOOL™ SF NAV Bi-Directional Catheter.....	BNI35FJH*	23
THERMOCOOL™ SF NAV Bi-Directional Catheter.....	D129704	23
THERMOCOOL™ SF NAV Uni-Directional Catheter.....	D129704	22
THERMOCOOL™ SF NAV Uni-Directional Catheter.....	D131802*	22
THERMOCOOL™ SF NAV Uni-Directional Catheter.....	D131803	22
THERMOCOOL™ SF NAV Uni-Directional Catheter.....	D131804*	22
THERMOCOOL™ SF NAV Bi-Directional Catheter with curve visualization.....	BNI35BDCT*	24
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization.....	39E43R	24
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization.....	CFT003	24
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization.....	D131504*	24
THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization.....	SAT001	24
THERMOCOOL™ SF Bi-Directional Catheter.....	BDI35BDRT*	29
THERMOCOOL™ SF Bi-Directional Catheter.....	BDI35DDRT*	29
THERMOCOOL™ SF Bi-Directional Catheter.....	BDI35DFRT	29
THERMOCOOL™ SF Bi-Directional Catheter.....	BDI35FJRT*	29
THERMOCOOL™ SF Uni-Directional Catheter.....	D131602*	29
THERMOCOOL™ SF Uni-Directional Catheter.....	D131603	29
THERMOCOOL™ SF Uni-Directional Catheter.....	D131604*	29
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter.....	D134801	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter.....	D134802	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter.....	D134803	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter.....	D134804	21
THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter.....	D134805	21
THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter.....	D134701	21
THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter.....	D134702	21
THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter.....	D134703	21
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	D132701	22
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	D132702	22
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	D132703	22
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	D132704	22
THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter.....	D132705	22
THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter.....	D133601	21
THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter.....	D133602	21
THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter.....	D133603	21
VISITAG SURPOINT™ Software Module.....	KT5400234	81
WEBSTER™ CS Bi-Directional Catheter.....	36G35QCT	42
WEBSTER™ CS Bi-Directional Catheter.....	36G5JMCT	42
WEBSTER™ CS Bi-Directional Catheter.....	36D35R	47
WEBSTER™ CS Bi-Directional Catheter.....	CY1212CT	47
WEBSTER™ CS Bi-Directional Catheter.....	D128709	42
WEBSTER™ CS Bi-Directional Catheter.....	36D5JR	42
WEBSTER™ CS Uni-Directional Catheter.....	D135303	42
WEBSTER™ CS Uni-Directional Catheter.....	D135304	42
WEBSTER™ Decapolar Catheter Deflectable.....	35X33QCT*	43
WEBSTER™ Decapolar Catheter Deflectable.....	36G33QCT*	43
WEBSTER™ Decapolar Catheter Deflectable.....	35X38Q	49
WEBSTER™ Decapolar Catheter Deflectable.....	35X33Q	49
WEBSTER™ Decapolar Catheter Deflectable.....	35X58Q	49
WEBSTER™ Decapolar Catheter Deflectable.....	36G33Q	49
WEBSTER™ Decapolar Catheter Deflectable.....	CY1212CT	43
WEBSTER™ Decapolar Catheter Deflectable.....	D128709	43
WEBSTER™ Decapolar Catheter Fixed.....	37S94RCT	44
WEBSTER™ Decapolar Catheter Fixed.....	37S94R	50
WEBSTER™ Decapolar Catheter Fixed.....	37S95R	50
WEBSTER™ Decapolar Catheter Fixed.....	37G58R	50
WEBSTER™ Decapolar Catheter Fixed.....	37G68R	50
WEBSTER™ Duo-Decapolar Catheter.....	D728260RT	52
WEBSTER™ Hexapolar Catheter Deflectable.....	36E30Q	48
WEBSTER™ Hexapolar Catheter Fixed.....	37U00R	50
WEBSTER™ HIS Catheter.....	CB3412CT	44
WEBSTER™ HIS Catheter.....	D1085413	44
WEBSTER™ Octopolar Catheter Deflectable.....	35Y13Q*	47
WEBSTER™ Octopolar Catheter Deflectable.....	36S31Q	47
WEBSTER™ Octopolar Catheter Deflectable.....	36F30Q*	47
WEBSTER™ Quadrapolar Catheter Deflectable.....	36B38RCT	43
WEBSTER™ Quadrapolar Catheter Deflectable.....	36B33RCT	43
WEBSTER™ Quadrapolar Catheter Deflectable.....	36C38RCT	43
WEBSTER™ Quadrapolar Catheter Deflectable.....	36B38R	48
WEBSTER™ Quadrapolar Catheter Deflectable.....	36B33R	48
WEBSTER™ Quadrapolar Catheter Deflectable.....	36B35R	48
WEBSTER™ Quadrapolar Catheter Deflectable.....	36C38R	48
WEBSTER™ Quadrapolar Catheter Deflectable.....	36C33R	48
WEBSTER™ Quadrapolar Catheter Deflectable.....	CB3412CT	43

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# CODE INDEX

---

WEBSTER™ Quadrapolar Catheter Fixed.....	37C53RCT.....	44
WEBSTER™ Quadrapolar Catheter Fixed.....	37C03RCT.....	44
WEBSTER™ Quadrapolar Catheter Fixed.....	F6QF005CT.....	44
WEBSTER™ Quadrapolar Catheter Fixed.....	F6QA005CT.....	44
WEBSTER™ Quadrapolar Catheter Fixed.....	CB3412CT.....	44
WEBSTER™ Quadrapolar Catheter Fixed.....	D128709.....	44
WEBSTER™ Quadrapolar Catheter Fixed.....	37C08R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37C03R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37C53R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37T03R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D00R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D08R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D03R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D05R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D38R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D33R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D35R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D53R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D55R.....	51
WEBSTER™ Quadrapolar Catheter Fixed.....	37D58R.....	51
YELLOW PATCH SENSOR CABLE .....	KT0001421 .....	81

08255790.....	ACUNAV™ Ultrasound Catheter.....	65
10043342.....	ACUNAV™ Ultrasound Catheter.....	65
10135910.....	ACUNAV™ Ultrasound Catheter.....	66
10135936.....	ACUNAV™ Ultrasound Catheter.....	66
10438577.....	SOUNDSTAR™ 3D eco Catheter.....	67
10439011.....	SOUNDSTAR™ eco 8F Catheter.....	64
10439072.....	SOUNDSTAR™ 3D eco Catheter.....	67
10439236.....	SOUNDSTAR™ eco 8F Catheter.....	64
301803A.....	PREFACE™ Fixed Guiding Sheaths.....	56
301803M.....	PREFACE™ Fixed Guiding Sheaths for EP.....	56
301803MS.....	PREFACE™ Fixed Guiding Sheaths for EP.....	56
301803P.....	PREFACE™ Fixed Guiding Sheaths.....	56
301805M.....	PREFACE™ Fixed Guiding Sheaths for EP.....	56
34A15M.....	NAVISTAR™ 4 mm Catheter.....	26
34A25M*.....	NAVISTAR™ 4 mm Catheter.....	26
34A35M.....	NAVISTAR™ 4 mm Catheter.....	26
34A45M*.....	NAVISTAR™ 4 mm Catheter.....	26
34A55M.....	NAVISTAR™ 4 mm Catheter.....	26
34AJ5M*.....	NAVISTAR™ 4 mm Catheter.....	26
34D37R.....	CELSIUS™ 4 mm Braided Tip Catheter.....	33
34D47R.....	CELSIUS™ 4 mm Braided Tip Catheter.....	33
34D57R.....	CELSIUS™ 4 mm Braided Tip Catheter.....	33
34H17M.....	NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
34H27M.....	NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
34H37M.....	NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
34H57M.....	NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
34HJ7M.....	NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter.....	25
34J27M.....	NAVISTAR™ DS Catheter.....	26
34J37M.....	NAVISTAR™ DS Catheter.....	26
34J57M.....	NAVISTAR™ DS Catheter.....	26
34JJ7M.....	NAVISTAR™ DS Catheter.....	26
35G07R.....	CELSIUS™ 4 mm Catheter Thermocouple.....	32
35G17R.....	CELSIUS™ 4 mm Catheter Thermocouple.....	32
35G37R*.....	CELSIUS™ 4 mm Catheter Thermocouple.....	32
35I33R.....	EZ STEER™ 4 mm Bi-Directional Catheter.....	30
35I35R.....	EZ STEER™ 4 mm Bi-Directional Catheter.....	30
35I55R.....	EZ STEER™ 4 mm Bi-Directional Catheter.....	30
35I5JR.....	EZ STEER™ 4 mm Bi-Directional Catheter.....	30
35O2515R*.....	LASSO™ 2515 Variable Mapping Catheter.....	46
35O16R*.....	LASSO™ Mapping Catheter.....	46
35O26R*.....	LASSO™ Mapping Catheter.....	46
35O26RCT.....	LASSO™ Decapolar Catheter.....	41
35O36R*.....	LASSO™ Mapping Catheter.....	46
35O36RCT.....	LASSO™ Decapolar Catheter.....	41
35O46R*.....	LASSO™ Mapping Catheter.....	46
35O66R*.....	LASSO™ Mapping Catheter.....	46
35Q13R*.....	CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	31
35Q33R.....	CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	31
35Q53R.....	CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	31
35Q73R*.....	CELSIUS™ THERMOCOOL™ Uni-Directional Catheter.....	31
35S17R.....	CELSIUS™ DS Catheter.....	35
35S37R.....	CELSIUS™ DS Catheter.....	35
35S57R.....	CELSIUS™ DS Catheter.....	35
35T2515R*.....	LASSO™ 2515 Variable Mapping Catheter.....	46
35T26R*.....	LASSO™ Mapping Catheter.....	46
35T36R*.....	LASSO™ Mapping Catheter.....	46
35U33R.....	EZ STEER™ DS Bi-Directional Catheter.....	31
35U35R.....	EZ STEER™ DS Bi-Directional Catheter.....	31
35U55R.....	EZ STEER™ DS Bi-Directional Catheter.....	31
35U5JR.....	EZ STEER™ DS Bi-Directional Catheter.....	31
35UJR.....	EZ STEER™ DS Bi-Directional Catheter.....	31
35X33Q.....	WEBSTER™ Decapolar Catheter Deflectable.....	49
35X33QCT*.....	WEBSTER™ Decapolar Catheter Deflectable.....	43
35X38Q.....	WEBSTER™ Decapolar Catheter Deflectable.....	49
35X58Q.....	WEBSTER™ Decapolar Catheter Deflectable.....	49
35Y13Q*.....	WEBSTER™ Octopolar Catheter Deflectable.....	54
36B33R.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	48
36B33RCT.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	43
36B35R.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	48
36B38R.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	48
36B38RCT.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	43
36C33R.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	48
36C38R.....	WEBSTER™ Quadrapolar Catheter Deflectable.....	48

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# CODE INDEX

36C38RCT	WEBSTER™ Quadrapolar Catheter Deflectable	43
36D35R	WEBSTER™ CS Bi-Directional Catheter	47
36D5JR	WEBSTER™ CS Bi-Directional Catheter	47
36E30Q	WEBSTER™ Hexapolar Catheter Deflectable	48
36F30Q*	WEBSTER™ Octopolar Catheter Deflectable	47
36G33Q	WEBSTER™ Decapolar Catheter Deflectable	49
36G33QCT*	WEBSTER™ Decapolar Catheter Deflectable	43
36G35QCT	WEBSTER™ CS Bi-Directional Catheter	42
36G5JMCT	WEBSTER™ CS Bi-Directional Catheter	42
36H33M	NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter	25
36H35M	NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter	25
36H55M*	NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter	25
36H5JM	NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter	25
36HJM	NAVISTAR™ THERMOCOOL™ Bi-Directional Catheter	25
36J13R	HALO™ XP Tricuspid Mapping Catheter	52
36J13RCT	HALO™ XP Catheter	45
36J14R	HALO™ XP Tricuspid Mapping Catheter	52
36J16R	HALO™ XP Tricuspid Mapping Catheter	52
36L00Q*	PARAHISIAN™ Catheter	54
36N07R*	CELSIUS™ 4 mm Catheter Thermocouple	32
36N17R	CELSIUS™ 4 mm Catheter Thermocouple	32
36N27R	CELSIUS™ 4 mm Catheter Thermocouple	32
36N37R	CELSIUS™ 4 mm Catheter Thermocouple	32
36N47R*	CELSIUS™ 4 mm Catheter Thermocouple	32
36N57R	CELSIUS™ 4 mm Catheter Thermocouple	32
36O07R	CELSIUS™ 4 mm Catheter Thermistor	33
36O17R	CELSIUS™ 4 mm Catheter Thermistor	33
36O27R	CELSIUS™ 4 mm Catheter Thermistor	33
36O37R	CELSIUS™ 4 mm Catheter Thermistor	33
36O47R	CELSIUS™ 4 mm Catheter Thermistor	33
36O57R	CELSIUS™ 4 mm Catheter Thermistor	33
36Q33M*	EZ STEER™ THERMOCOOL™ Catheter	30
36Q35M	EZ STEER™ THERMOCOOL™ Catheter	30
36Q55M	EZ STEER™ THERMOCOOL™ Catheter	30
36Q5JM	EZ STEER™ THERMOCOOL™ Catheter	30
36S30Q*	Pediatric Catheter	54
36S31Q	WEBSTER™ Octopolar Catheter Deflectable	47
36Y32R	ISMUS™ Catheter	52
36Y32RCT	ISMUS™ Catheter	45
36Y39R	CRISTACATH™ Catheter	54
37C0CR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37C3FR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37C5CR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37C5RR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37C03R	WEBSTER™ Quadrapolar Catheter Fixed	51
37C03RCT	WEBSTER™ Quadrapolar Catheter Fixed	44
37C08R	WEBSTER™ Quadrapolar Catheter Fixed	51
37C53R	WEBSTER™ Quadrapolar Catheter Fixed	51
37C53RCT	WEBSTER™ Quadrapolar Catheter Fixed	44
37D00R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D03R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D05R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D08R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D33R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D35R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D38R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D53R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D55R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D58R	WEBSTER™ Quadrapolar Catheter Fixed	51
37D0CR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37D3FR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37D5CR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37D0RR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37D5RR	AVAIL™ Quadrapolar Fixed Shape Catheter	53
37G58R	WEBSTER™ Decapolar Catheter Fixed	50
37G68R	WEBSTER™ Decapolar Catheter Fixed	50
37PORR	AVAIL™ Bipolar Fixed Shape Catheter	53
37S94R	WEBSTER™ Decapolar Catheter Fixed	50
37S94RCT	WEBSTER™ Decapolar Catheter Fixed	44
37S95R	WEBSTER™ Decapolar Catheter Fixed	50
37T03R	WEBSTER™ Quadrapolar Catheter Fixed	51
37U00R	WEBSTER™ Hexapolar Catheter Fixed	50
39D03X	Spare Parts and Others	74

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

39D03X	COOLFLOW™ Pump Accessories	76
39D06X	Spare Parts and Others	74
39D12X	Spare Parts and Others	74
39D12X5	Indifferent Electrodes and Corresponding Connection Cables	73
39D14X	COOLFLOW™ Pump Accessories	76
39D14X	Spare Parts and Others	74
39D16X	Indifferent Electrodes and Corresponding Connection Cables	73
39D16X5	Indifferent Electrodes and Corresponding Connection Cables	73
39D21X	Spare Parts and Others	74
39D60X	Stockert EP Shuttle Accessories	75
39D60X10	Stockert EP Shuttle Accessories	75
39D60X	COOLFLOW™ Pump Accessories	76
39D82X	Indifferent Electrodes and Corresponding Connection Cables	73
39D83X	Spare Parts and Others	74
39E43R	Spare Parts and Others	74
39E43R	Non-Navigation ABLATION CABLES	91
39E46R	Non-Navigation ABLATION CABLES	91
39E53S	Non-Navigation ABLATION CABLES	91
39E56S	Non-Navigation ABLATION CABLES	91
39E68R	Non-Navigation ABLATION CABLES	91
39E69R	Non-Navigation ABLATION CABLES	91
39E23R	Non-Navigation ABLATION CABLES	91
39F23R	Non-Navigation ABLATION CABLES	91
39F24R	Non-Navigation DIAGNOSTIC CABLES	91
39F27R	Indifferent Electrodes and Corresponding Connection Cables	73
39F28X10	Stockert EP Shuttle Accessories	75
39F28X15	Stockert EP Shuttle Accessories	75
39F28X20	Stockert EP Shuttle Accessories	75
39F28X30	Stockert EP Shuttle Accessories	75
39F40R	Non-Navigation DIAGNOSTIC CABLES	91
39F41R	Non-Navigation DIAGNOSTIC CABLES	91
39F43R	Non-Navigation DIAGNOSTIC CABLES	91
39F44R	Non-Navigation DIAGNOSTIC CABLES	91
39F47D	Non-Navigation DIAGNOSTIC CABLES	91
39F54R	Non-Navigation DIAGNOSTIC CABLES	91
39F55R	Non-Navigation DIAGNOSTIC CABLES	91
BDI35BDRT*	THERMOCOOL™ SF Bi-Directional Catheter	29
BDI35DDRT*	THERMOCOOL™ SF Bi-Directional Catheter	29
BDI35DFRT	THERMOCOOL™ SF Bi-Directional Catheter	29
BDI35FJRT*	THERMOCOOL™ SF Bi-Directional Catheter	29
BNI35BDCT*	THERMOCOOL™ SF NAV Bi-Directional Catheter with curve visualization	24
BNI35BDH*	THERMOCOOL™ SF NAV Bi-Directional Catheter	23
BNI35BFH*	THERMOCOOL™ SF NAV Bi-Directional Catheter	23
BNI35DFH	THERMOCOOL™ SF NAV Bi-Directional Catheter	23
BNI35FFH*	THERMOCOOL™ SF NAV Bi-Directional Catheter	23
BNI35FJH*	THERMOCOOL™ SF NAV Bi-Directional Catheter	23
C3CFAE	CARTO™ 3 System Software Module	81
C3EXTY1	CARTO™ 3 System Service Contracts	83
C3EXTY2	CARTO™ 3 System Service Contracts	83
C3EXTY3	CARTO™ 3 System Service Contracts	83
C3MERGE	CARTO™ 3 System Software Module	83
C3PASO	CARTO™ 3 System Software Module	83
C3PREMY1	CARTO™ 3 System Service Contracts	83
C3PREMUPY1	CARTO™ 3 System Service Contracts	83
C3PREMY2	CARTO™ 3 System Service Contracts	83
C3PREMUPY2	CARTO™ 3 System Service Contracts	83
C3PREMY3	CARTO™ 3 System Service Contracts	83
C3PREMUPY3	CARTO™ 3 System Service Contracts	83
C3RMT	CARTO™ 3 System	79
C3SOUND	CARTO™ 3 System Software Module	81
C3TOUCH	CARTO™ 3 System Software Module	81
C3CNFDRPLY	CARTO™ 3 Software System Module	81
C3CONFD	CARTO™ 3 Software System Module	81
C3REPLAY	CARTO™ 3 Software System Module	81
C3V4	CARTO™ 3 Software System Module	81
CA400811	CARTO™ 3 System Accessories	80
CA401214	CARTO™ 3 System Accessories	80
CA404103	CARTO™ 3 System Accessories	80
CA404104	CARTO™ 3 System Accessories	80
CA540006	CARTO™ 3 System Accessories	79
CA5400081	CARTO™ 3 System Accessories	80
CA540011	CARTO™ 3 System Accessories	79
CA540012	CARTO™ 3 System Accessories	80

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# CODE INDEX

CA5400161	CARTO™ 3 System Accessories	79
CA540017	CARTO™ 3 System Accessories	80
CA540019	CARTO™ 3 System Accessories	80
CA540021	CARTO™ 3 System Accessories	80
CA540151	CARTO™ 3 System Accessories	80
CB3410CT	CARTO™ 3 System Cables	89
CB3412CT	CARTO™ 3 System Cables	89
CB3412CT	WEBSTER™ HIS Catheter	44
CB3412CT	WEBSTER™ Quadropolar Catheter Deflectable	89
CB3412CT	WEBSTER™ Quadropolar Catheter Fixed	44
CB3412CT	CARTO VIZIGO™ Sheaths	58
CB3434CT	CARTO™ 3 System Cables	89
CFP003	COOLFLOW™ Pump Accessories	76
CFP006	COOLFLOW™ Pump Accessories	76
CFPPV1	Service Contracts	77
CFPPVP1	Service Contracts	77
CFT003	THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization	24
CFT003	COOLFLOW™ Pump Tubing	76
CG2025CT	CARTO™ 3 System Cables	89
CP242101F	CARTO™ 3 System	79
CR3425CT	CARTO™ 3 System Cables	87
CR3434CT	CARTO™ 3 System Cables	87
CR7TCS4RT	CELSIUS™ RMT Catheter	35
CR7TCSIRT	CELSIUS™ RMT THERMOCOOL™ Catheter	32
CREFP6	CARTO™ 3 System External Reference Patch	98
CW-3902-30F	CARTO™ 3 System Accessories	80
CW-3903-30F	CARTO™ 3 System Accessories	80
CW417831F	CARTO™ 3 System Accessories	79
KT5400103	CARTO™ 3 System Accessories	79
CW540906N	CARTO™ 3 System Accessories	80
CW540915N	CARTO™ 3 System Accessories	80
CY1212CT	CARTO™ 3 System Cables	89
CY1210CT	CARTO VIZIGO™ Sheaths	58
CY1212CT	CARTO™ 3 System Cables	89
CY1212CT	WEBSTER™ Bi-Directional Catheter	42
CY1212CT	WEBSTER™ Decapolar Catheter Deflectable	43
D1085413	WEBSTER™ HIS Catheter	44
D128207	PENTARAY™ NAV eco High Density Mapping Catheter	39
D128208	PENTARAY™ NAV eco High Density Mapping Catheter	39
D128210	PENTARAY™ NAV eco High Density Mapping Catheter	39
D128211	PENTARAY™ NAV eco High Density Mapping Catheter	39
D128622	REVERSE CABLES	92
D128622	Halo™ XP Catheter	45
D128622	ISMUS™ Catheter	45
D128624	CARTO™ 3 System Cables	89
D128709	WEBSTER™ CS Bi-Directional Catheter	42
D128709	REVERSE CABLES	90
D128709	WEBSTER™ Decapolar Catheter Deflectable	43
D128709	WEBSTER™ Quadropolar Catheter Fixed	44
D129704	Interface SMARTABLATE™ System Cables	90
D129704	THERMOCOOL™ SF NAV Uni-Directional Catheter	22
D129704	NAVISTAR™ THERMOCOOL™ Uni-Directional Catheter	25
D129704	NAVISTAR™ 4 MM CATHETER	26
D129704	NAVISTAR™ 4 DS CATHETER	26
D129707	Interface SMARTABLATE™ System Cables	90
D129709	WEBSTER™ Quadropolar Catheter Deflectable	43
D130302	EZ STEER™ THERMOCOOL™ Catheter	30
D130302	CELSIUS™ 4 mm Catheter Thermocouple	32
D130302	CELSIUS™ RMT THERMOCOOL™ Catheter	32
D130302	EZ STEER™ THERMOCOOL™ Catheter	32
D130302	EZ STEER™ DS Bi-Directional Catheter	31
D130302	CELSIUS™ 4 mm Catheter Thermistor	33
D130302	CELSIUS™ RMT Catheter	35
D130302	CELSIUS™ 4 mm Braided Tip Catheter	33
D130303	CELSIUS™ DS Catheter	35
D130303	EZ STEER™ DS Bi-Directional Catheter	31
D130303	Interface SMARTABLATE™ System Cables	90
D130303	CELSIUS FLTR™ 8 mm Uni-Directional Catheter	34
D130303	CELSIUS FLTR™ 8 mm Bi-Directional Catheter	34
D131504*	THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization	24
D131602*	THERMOCOOL™ SF Uni-Directional Catheter	29
D131603	THERMOCOOL™ SF Uni-Directional Catheter	29
D131604*	THERMOCOOL™ SF Uni-Directional Catheter	29

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.



D131802*	THERMOCOOL™ SF NAV Uni-Directional Catheter	22
D131803	THERMOCOOL™ SF NAV Uni-Directional Catheter	22
D131804*	THERMOCOOL™ SF NAV Uni-Directional Catheter	22
D132701	THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter	22
D132702	THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter	22
D132703	THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter	22
D132704	THERMOCOOL SMARTTOUCH™ Bi-Directional Catheter	22
D133601	THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter	21
D133602	THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter	21
D133603	THERMOCOOL SMARTTOUCH™ Uni-Directional Catheter	21
D133701	nGEN™ Generator Catheter Interface Cables	69
D133703	nGEN™ Generator Catheter Interface Cables	69
D134301	LASO™ 2515 NAV eco Circular Mapping Catheter	38
D134302	LASO™ 2515 NAV eco Circular Mapping Catheter	38
D134401	CARTO™ 3 System Cables	89
D134402	CARTO™ 3 System Cables	39
D134801	THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter	21
D134802	THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter	21
D134803	THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter	21
D134804	THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter	21
D134805	THERMOCOOL SMARTTOUCH™ SF Bi-Directional Catheter	21
D134701	THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter	21
D134702	THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter	21
D134703	THERMOCOOL SMARTTOUCH™ SF Uni-Directional Catheter	21
D134901	LASO™ NAV eco Circular Mapping Catheter	38
D134902	LASO™ NAV eco Circular Mapping Catheter	38
D134903	LASO™ NAV eco Circular Mapping Catheter	38
D134904	LASO™ NAV eco Circular Mapping Catheter	38
D134905	LASO™ NAV eco Circular Mapping Catheter	38
D134906	LASO™ NAV eco Circular Mapping Catheter	38
D135303	WEBSTER™ CS Uni-Directional Catheter	42
D135304	WEBSTER™ CS Uni-Directional Catheter	42
D135501*	CELSIUS FLTR™ 8 mm Uni-Directional Catheter	34
D135502*	CELSIUS FLTR™ 8 mm Uni-Directional Catheter	34
D135703	nGEN™ Generator Catheter Interface Cables	69
D135901	CELSIUS FLTR™ 8 mm Bi-Directional Catheter	34
D138401	nGEN™ Generator Kit	69
D138501	CARTO VIZIGO™ Sheath	58
D138502	CARTO VIZIGO™ Sheath	58
D138503	CARTO VIZIGO™ Sheath	58
D138801	nGEN™ Monitor	69
D138904	HELIOSTAR™ Balloon Ablation Catheter	36
D139001	LASSOSTAR™ Circular Catheter	36
D139002	LASSOSTAR™ Circular Catheter	36
D139003	LASSOSTAR™ Circular Catheter	36
D139401	QDOT MICRO™ Uni-directional Catheter	20
D139402	QDOT MICRO™ Uni-directional Catheter	20
D139403	QDOT MICRO™ Uni-directional Catheter	20
D139501	QDOT MICRO™ Bi-directional Catheter	20
D139502	QDOT MICRO™ Bi-directional Catheter	20
D139503	QDOT MICRO™ Bi-directional Catheter	20
D139504	QDOT MICRO™ Bi-directional Catheter	20
D139505	QDOT MICRO™ Bi-directional Catheter	20
D139701	nGEN™ Pump	69
D140010*	MOBICATH™ Sheath	58
D140011*	MOBICATH™ Sheath	58
D140102	nGEN™ Generator Catheter Interface Cables	69
D140401	LASSOSTAR™ NAV Circular Catheter	36
D140402	LASSOSTAR™ NAV Circular Catheter	36
D140403	LASSOSTAR™ NAV Circular Catheter	36
D141101	GUIDESTAR™ Deflectable Sheath	36
D141102	GUIDESTAR™ Deflectable Sheath	36
D141103	GUIDESTAR™ Deflectable Sheath	36
D129704	NAVISTAR™ RMT THERMOCOOL™ Catheter	27
D129704	NAVISTAR™ RMT 4 mm Catheter	26
D160901	OCTARAY™ Catheter small	39
D160902	OCTARAY™ Catheter large	39
D160903	OCTARAY™ Catheter large	39
D160904	OCTARAY™ Catheter small	39
D160905	OCTARAY™ Catheter large	39
D160906	OCTARAY™ Catheter large	39
D728260RT	WEBSTER™ Duo-Decapolar Catheter	52
D7L2015CT*	LASSO™ Duo-Decapolar Catheter	41
D7L2020CT*	LASSO™ Duo-Decapolar Catheter	41

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

# CODE INDEX

EM5050011	CARTO™ 3 System Cables	89
EM5050060	CARTO™ 3 System Cables	89
EM540003	CARTO™ 3 System Accessories	79
EM540004	CARTO™ 3 System Accessories	79
EM540015	CARTO™ 3 System Accessories	80
EM540020	CARTO™ 3 System Accessories	80
EM5400201	CARTO™ 3 System Accessories	80
EM540200	CARTO™ 3 System	79
EM5050060U	Systems Software Module OBL (available in Middle East only)	82
EM5050011U	Systems Software Module OBL (available in Middle East only)	82
EM5402020	CARTO™ 3 System	79
FCL20100	CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath	57
FCL20101	CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath	57
FCL20102	CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath	57
FCL20103	CARDIAGUIDE™ Fixed Curved Braided Transseptal Sheath	57
FG540000	CARTO™ 3 System	79
FG540000U	CARTO™ 3 System	79
F6QF005CT	WEBSTER™ Quadrapolar Catheter Fixed	44
F6QA005CT	WEBSTER™ Quadrapolar Catheter Fixed	44
FND01901	HEARTSPAN™ Transseptal Needles	59
FND01902	HEARTSPAN™ Transseptal Needles	59
FND01903	HEARTSPAN™ Transseptal Needles	59
FND01905	HEARTSPAN™ Transseptal Needles	59
FND01906	HEARTSPAN™ Transseptal Needles	59
KT0001420	Green Patch Sensor cable	81
KT0001421	Yellow Patch Sensor cable	81
KT0001427	CARTO™ 3 System Accessories	80
KT0001429	CARTO™ 3 System Accessories	80
KT0001438	CARTO™ 3 System Accessories	80
KT0001439	CARTO™ 3 System Accessories	80
KT0001440	CARTO™ 3 System Accessories	80
KT5400101	CARTO™ 3 System	79
KT5400101U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400102	CARTO™ 3 System	79
KT5400102U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400155U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400221U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400136U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400145U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400222U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT540050U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400170U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400160U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400130U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400124U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400165U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400111U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400104U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400262U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400162U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400162U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400165	FAM DX Module	40
KT5400165U	FAM DX Module OBL	40
KT5400490U	CARTO™ 3 Systems Software Module OBL (available in Middle East only)	82
KT5400124	CARTO™ 3 System Software Module	81
KT5400127	CARTO™ 3 System Software Module	81
KT5400128	CARTO™ 3 System Software Module	81
KT5400129	CARTO™ 3 System Software Module	81
KT5400130	CARTO™ 3 System Software Module	81
KT5400134	CARTO™ 3 System Software Module	81
KT5400170	Ripple Mapping Module	81
KT5400160	CARTOSEG™ CT Module	81
KT5400230	QDOT MICRO™ Module Installation Kit for V7	81
KT5400233	HELIOSTAR™ Upgrade Kit	81
KT5400234	VISITAG SURPOINT™ Software Module	81
KT5400235	QDOT MICRO™ Module Installation Kit for V7 Windows 10	81
KT5400250	CARTO™ 3 System V7 Windows 10 Base Upgrade Kit	81
MA540013	CARTO™ 3 System	79
MA5400501	CARTO™ 3 System Accessories	80
MA5400801	CARTO™ 3 System Accessories	80
MA5402100	SMARTABLATE™ System	79
MA5402300	CARTO™ 3 System Accessories	69
M490005	SMARTABLATE™ System	70



M490010	SMARTABLATE™ System	70
M4900106	SMARTABLATE™ System	70
M490012	SMARTABLATE™ System	70
M490013	SMARTABLATE™ System	70
M490014	SMARTABLATE™ System	70
M490016	SMARTABLATE™ System	70
M490017	SMARTABLATE™ System	70
M490018	SMARTABLATE™ System	70
M490021	SMARTABLATE™ System	70
M490023	SMARTABLATE™ System	70
M490025	SMARTABLATE™ System	70
M490026	SMARTABLATE™ System	70
M490027	SMARTABLATE™ System	70
M490028	SMARTABLATE™ System	70
M490029	SMARTABLATE™ System	70
M490030	SMARTABLATE™ System	70
M490034	Interface SMARTABLATE™ System Cables	90
M490034	SMARTABLATE™ System	70
M490035	Interface SMARTABLATE™ System Cables	90
M490035	SMARTABLATE™ System	70
M490037	SMARTABLATE™ System	70
M4900602	SMARTABLATE™ System	70
M4900603	SMARTABLATE™ System	70
M4900605	SMARTABLATE™ System	70
M4900606	SMARTABLATE™ System	70
M529723SP	nGEN™ Pump AC Adapter	69
M549104SP	nGEN™ Clamp	69
M581006SP	nGEN™ System Interface Cables	69
M581007SP	nGEN™ System Interface Cables	69
M581101SP	nGEN™ System Interface Cables	69
M5831100SP	nGEN™ System Interface Cables	69
M5831101SP	nGEN™ System Interface Cables	69
M583191SP	nGEN™ Monitor Power Cord	69
M583196SP	nGEN™ Clamp	69
M583197SP	nGEN™ System Interface Cables	69
M583198SP	nGEN™ System Interface Cables	69
M583199SP	nGEN™ System Interface Cables	69
M589704SP	nGEN™ System Interface Cables	69
M589705SP	nGEN™ System Interface Cables	69
M600201SP	nGEN™ System Interface Cables	69
M600202SP	nGEN™ System Interface Cables	69
M600203SP	nGEN™ System Interface Cables	69
M600301SP	nGEN™ System Interface Cables	69
M621402SP	nGEN™ AC Power Cord	69
M621406SP	nGEN™ AC Power Cord	69
M621407SP	nGEN™ AC Power Cord	69
M621410SP	nGEN™ AC Power Cord	69
M649002SP	nGEN™ Pump AC Power Cord	69
M649006SP	nGEN™ Pump AC Power Cord	69
M649007SP	nGEN™ Pump AC Power Cord	69
M649010SP	nGEN™ Pump AC Power Cord	69
M649002SP	nGEN™ System Interface Cables	69
NR7TCSIY	NAVISTAR™ RMT THERMOCOOL™ Catheter	27
NR7TCS4Y	NAVISTAR™ RMT 4 mm Catheter	27
R7D282CT	DECANAV™ Catheter	39
R7F282CT	DECANAV™ Catheter	39
SAT001	THERMOCOOL™ SF NAV Uni-Directional Catheter with curve visualization	24
SAT001	SMARTABLATE™ System	70
SAT001	EZ STEER™ THERMOCOOL™ Catheter	30
SAT001	SMARTABLATE™ System Irrigation Tubing Set	31
SAT001	CELSIUS™ 4 mm Catheter Thermocouple	32
SAT001	nGEN™ Irrigation Tubing Set	69
SHUTBC1	Service Contracts	75
SHUTPV1	Service Contracts	75
SHUTPVP1	Service Contracts	75
SNDSTR10	SOUNDSTAR™ 3D Ultrasound Catheter	65
SNDSTR10G	SOUNDSTAR™ 3D Ultrasound Catheter	65
WCR1503	COHEREX WAVECREST™ LAA Implant	62
WCR1513	COHEREX WAVECREST™ LAA Implant	62
WCR1523	COHEREX WAVECREST™ LAA Implant	62
WCR1530	COHEREX WAVECREST™ LAA Implant	62
WCR1540	COHEREX WAVECREST™ LAA Implant	62
WCR1541	COHEREX WAVECREST™ LAA Implant	62
WCR1551	COHEREX WAVECREST™ LAA Implant	62

\* Biosense Webster will be discontinuing this product code end of 2022. Please contact your Biosense Webster representative to identify an alternative product.

The third-party trademarks used herein are trademarks of their respective owners.  
Capitalized product names are trademarks of Biosense Webster.

Important Information: For product details such as indications, contraindications, warnings and precautions please consult the IFU.  
Please refer to the instructions for use accompanying each medical device herein mentioned. This promotional document is intended for healthcare professionals only.  
As part of Biosense Webster policy of continuous development, we reserve the right to change product specifications without prior notification. This publication is not intended for distribution outside of the EMEA region.  
©Johnson & Johnson Medical NV/SA 2022.

**EC Representative / Biosense Webster**  
**A Division of Johnson & Johnson Medical NV/SA**  
Leonardo da Vincilaan 15 / 1831 Diegem, Belgium  
Tel: +32-2-7463-401 / Fax: +32-2-7463-403

**Biosense Webster, Inc.**  
31 Technology Drive, Suite 200  
Irvine, CA 92618 USA  
Tel: +1-909-839-8500  
Tel: +1-800-729-9010