

---

# Getinge Assured Wash Monitor

## Product Specification

**GETINGE**  
GETINGE GROUP

---

The Getinge Assured Wash Monitor is a single use synthetic soil test for routine monitoring of automated wash processes.

The test is inserted into a Wash Monitor Holder and placed into the washer disinfectant.



---

### Performance Characteristics

The synthetic test soil (red diamond) simulates human blood and tissue. When inserted into the Wash Monitor Holder the test soil is half obscured, thus simulating the joint of a surgical instrument. The removal of all test soil indicates effective cleaning.

---

### Physical Characteristics

The Wash Monitor consists of a metal strip with synthetic test soil.  
Dimensions: 89 mm x 25 mm (3.5" X 1")  
Substrate Material: Metal

---

### Environmental Limitations

Storage: 10° - 38°C (50° - 100°F) and 10-60% Relative Humidity. Keep away from disinfectants and chemical fumes

---

### Shelf Life

24 months.

---

### Standards Classification

Consistent with ISO 15883, routine monitoring of washer disinfectors.

---

Customer

Reference

---

## Ordering Information

<b>Getinge Assured</b>	<b>Order Number</b>
Wash Monitor 25 pcs	503864500
Wash Monitor 50 pcs	503864600
Wash Monitor 100 pcs	503864700
Wash Monitor Holder	503864800

---

---

**GETINGE**  
GETINGE GROUP

Getinge Infection Control AB  
P O Box 69, SE-305 05 Getinge, Sweden  
Phone: +46 10 335 00 00  
info@getinge.com  
  
www.getinge.com

---

## GETINGE GROUP

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, Getinge and Maquet. ArjoHuntleigh focuses on patient mobility and wound management solutions. Getinge provides solutions for infection control within healthcare and contamination prevention within life sciences. Maquet specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.

---