



DESILET CATHETER INTRODUCERS

1129 Desilet - 13 cm

1137 Desilet - 15 cm

10/150

Introducer (PE) for silicone or cardiovascular catheters.

Device placed in the vascular system by the Seldinger technique for introducing cardiovascular or silicone catheters.

Presented in a rigid tubular pack containing:
- 1 introducer needle Ø 1.06 x 1.26 mm
- 1 straight guidewire Ø 0.97 mm, 39 cm long
- 1 dilator + sheath assembly.



Desilet

Code		Colour of dilator	For Catheter		Needle Length mm	Sheath length cm	
Code 1129	Code 1137		Fr	Ø mm		1129	1137
1129.04*	-	Purple	04	1.3	38	5	-
1129.05	1137.05	Grey	05	1.7	70	13	15
1129.06	1137.06	Green	06	2.0	70	13	15
1129.07	1137.07	Pink	07	2.3	70	13	15
1129.08	1137.08	Blue	08	2.7	70	13	15
1129.09	1137.09	Yellow	09	3.0	70	13	15
1129.10	-	White	10	3.3	70	13	15

* Guidewire Length 30 cm, Ø 0.5 mm - Needle Ø 0.66 x 0.9 mm

1149 Desivalve


10/40

Catheter introducer with haemostatic valve.

PE-XRO catheter introducer with haemostatic valve for introducing cardiovascular catheters such as Swan-Ganz or temporary pacing catheters.

- XRO polyethylene sheath provided with colour coded stress-relief collar and fixation wing
- Side arm for injecting medication or venous pressure monitoring
- Nontouch catheter sleeve (AseptCath) incorporates a bayonet lock mechanism for secure and rapid attachment to the body of the Desivalve
- Distal compression hub minimises the risk of catheter migration
- Proximal compression hub acts as a barrier to bacterial contamination.

Tray presentation containing:
- 1 dilator and sheath
- 1 non-touch handling sleeve (AseptCath)
- 1 introducer needle
- 1 "J" guidewire.



Desivalve

Code	For catheter*			Sheath (PE) Length cm	Needle			Guidewire	
	Colour of dilator	Fr	Ø mm		Ø mm	G	Length mm	Ø mm	Length cm
1149.05	Grey	05	1.7	10	1.4	17	70	0.96	46
1149.06	Green	06	2.0	10	1.4	17	70	0.96	46
1149.07	Pink	07	2.3	10	1.4	17	70	0.96	46
1149.08	Blue	08	2.7	10	1.4	17	70	0.96	46
1149.09	Yellow	09	3.0	10	1.4	17	70	0.96	46

* For catheters with balloon, use a Desivalve the next Fr size larger than the catheter being placed.

DESILET CATHETER INTRODUCERS



1146 Peel-apart Desilet

5/45


Catheter introducer with peel-away sheath.

Tray presentation containing:
- 1 introducer needle Ø 1.06 x 1.26 mm - Length 70 mm
- 1 "J" guidewire - Length 53 cm, Ø 0.97 mm
- 1 dilator
- 1 Desilet with peel-apart sheath
- 1 10 ml syringe.

Code	For catheter **		Sheath	
	Fr	Ø mm	Length cm	Int. Ø - Ext. Ø mm
1146.04*	04	1.3	7	1.4 - 1.95
1146.05	05	1.7	14	1.8 - 2.3
1146.06	06	2.0	14	2.1 - 2.6
1146.07	07	2.3	14	2.4 - 3.1
1146.08	08	2.7	14	2.8 - 3.4
1146.09	09	3.0	14	3.1 - 3.7
1146.10	10	3.3	14	3.4 - 4.1
1146.12	12	4.0	14	4.1 - 4.7
1146.14	14	4.6	14	4.8 - 5.4

* Needle Ø 0.6 x 0.9 mm length 38 mm - Guidewire Ø 0.5 mm length 300 mm, FEP sheath and additional dilator. Peel-pack presentation.

** For catheters with balloon, use a Desilet the next Fr size larger than the catheter being placed.



Peel-apart Desilet

Prie 19 - desileto rinkinys: 4 Fr (tinkamas kateterio išoriniam diametrui iki 1,3 mm)

ACCESSORIES FOR DESILET

1135.08 Haemostatic valve with side outlet

10/400

Allows heparinisation and change of intra-cardiac* catheters without Desilet withdrawal and without loss of blood.



Code 1135

* Intracardiac catheters 8 Fr without balloon or 7 Fr with balloon.