

Shimadzu Europa GmbH · Postfach 29 02 60 · D-47262 Duisburg F.R.G.

To WHOM IT MAY CONCERN

Shimadzu AOC-6000 plus Specifications

Dear Sir/Madam,

We would like to inform you that Shimadzu is constantly working to improve specifications and performance of mass spectrometer instrument line.

In that scope, we would like to confirm the following supplementary information for Shimadzu AOC-6000 plus systems:

- Possible injection volume: **0.01-250 µl**
- Repeatability of liquid injection **<0.3% RSD***.

*using FID as GC detector

Sincerely Yours, *i. A. W. Weber*

**Dr. Waldemar Weber | Product Manager GCMS,
Marketing Europe | Analytical Business Unit
Shimadzu Europa GmbH**

Address: Albert-Hahn-Straße 6-10, D-47269 Duisburg
Tel: +49 (0)203/7687 - 463 | Mobile: +49(0) 172 2562102

Follow us   