

Vascular Instruments – Gefäßinstrumente

Vascular Clips and Bulldog Clamps

Gefäßclips und Bulldogklemmen




JOHNS-HOPKINS

BH018R

Forceps for applying and removing bulldog clamps BH011R-BH031R and FB420R-FB427R

Zange zum Anlegen und Abnehmen der Bulldogklemmen BH011R-BH031R und FB420R-FB427R

240 mm, 9 3/4"



BH020R

BH021R

BH026R

BH027R

BH022R

BH023R

DIEFFENBACH

	Total length:	Length of jaw serration:
	Gesamtlänge:	Länge der Zahnung:
BH020R	38 mm, 1 1/2"	12 mm, 1/2"
BH022R	48 mm, 2"	16 mm, 5/8"
BH021R	38 mm, 1 1/2"	12 mm, 1/2"
BH023R	48 mm, 2"	16 mm, 5/8"
BH026R	58 mm, 2 1/8"	16 mm, 5/8"
BH027R	58 mm, 2 1/8"	16 mm, 5/8"



BH011R-BH017R

JOHNS-HOPKINS

	Total length:	Length of jaw serration:
	Gesamtlänge:	Länge der Zahnung:
BH011R	38 mm, 1 1/2"	12 mm, 1/2"
BH012R	48 mm, 2"	17 mm, 5/8"
BH013R	56 mm, 2 1/8"	25 mm, 1"
BH014R	63 mm, 2 1/2"	30 mm, 1 3/16"
BH015R	70 mm, 2 3/4"	33 mm, 1 1/4"
BH016R	80 mm, 3 1/8"	42 mm, 1 3/4"
BH017R	91 mm, 3 1/2"	52 mm, 2"