

Straight Connector P. d. 1.1

Meets all connection requirements of the extracorporeal circulation system

Unique design and precision moulding techniques allow a smooth transition at tubing and connector junctions reducing blood turbulence

Produced from crystal clear polycarbonate allowing maximum visibility during de-airing procedures



Ref	Size	
	Luer	Size
750 0111 1	751 0111 1	3/16" x 3/16"
750 0121 1	751 0121 1	3/16" x 1/4"
750 0221 1	751 0221 1	1/4" x 1/4"
750 0321 1	751 0321 1	1/4" x 3/8"
750 0331 1	751 0331 1	3/8" x 3/8"
750 0421 1	751 0421 1	1/4" x 1/2"
750 0431 1	751 0431 1	3/8" x 1/2"
750 0441 1	751 0441 1	1/2" x 1/2"

"Y" Connector P. d. 1.2

Meets all connection requirements of the extracorporeal circulation system

Unique design and precision moulding techniques allow a smooth transition at tubing and connector junctions reducing blood turbulence

Produced from crystal clear polycarbonate allowing maximum visibility during de-airing procedures



Ref	Size		
	Luer	①	② ③
754 2221 1	755 2221 1	1/4"	1/4" 1/4"
754 3221 1	755 3221 1	3/8"	1/4" 1/4"
754 3321 1	755 3321 1	3/8"	1/4" 3/8"
754 3331 1	755 3331 1	3/8"	3/8" 3/8"
754 4321 1	755 4321 1	1/2"	3/8" 1/4"
754 4331 1	755 4331 1	1/2"	3/8" 3/8"
754 4421 1	755 4421 1	1/2"	1/2" 1/4"
754 4431 1	755 4431 1	1/2"	1/2" 3/8"
754 4441 1	755 4441 1	1/2"	1/2" 1/2"



