



85378

2 p.d Rosen tipo viola

[🏠](#) / [PRODUCTS](#) / [ROSEN WIRE GUIDE](#)

Rosen Wire Guide

[SPECIFICATIONS](#) [IMAGES](#)

This product line is serviced by the following clinical division(s):

[Aortic Intervention »](#)

[Peripheral Intervention »](#)

**LEARN HOW TO
ORDER »**

Set up a customer account
and order this or other
Cook products

**ASK A SUPPORT
QUESTION »**

Ask Customer Support
about our products or
services

**FIND PRODUCT
TRAINING »**

Connect with our Vista®
Education & Training
program

**VIEW YOUR ACCOUNT
INFORMATION »**

Register your customer
number on our Customer
Portal

 **Specifications**

[Give Us Feedback](#)


COLLAPSED EXPANDED

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Diameter inch	Length cm	Taper Length cm	More Info
PTFE-Coated Stainless Steel							
G01261	THSCF-35-145-1.5-ROSEN		N/A	.035	145	4.5	Expand »
G01264	THSCF-35-180-1.5-ROSEN		N/A	.035	180	4.5	Expand »
G01253	THSCF-35-260-1.5-ROSEN		N/A	.035	260	4.5	Expand »

ADDITIONAL SPECS

Description	-
Floppy Tip Length cm	1.5
Tip Curve Radius mm	1.5

 **Image Library**



Not all products shown on this website may be approved or available in all jurisdictions.
Consult with your local Cook representative, distribution company or customer support centre for details.

© Cook 2023

 **+353 61 334 440**

[Products](#)
[Patient Resources](#)
[Support](#)
[About](#)
[Divisions](#)
[Contact](#)
[Education & Training](#)
[Customer Portal](#)

[Privacy Statement](#)
[Human Rights](#)
[Responsibility Statement](#)
[Terms of Use](#)
[Cookie Compliance](#)
[Sitemap](#)

**Sign Up. Get Email about News, Products,
and Events from Cook Medical.**

[SIGN UP](#)