

Thermo Scientific™

Finnpipette™ F2 Good Laboratory Pipetting (GLP) Kits

Catalog number: 4700870

Related applications: [Pipettes & Pipette Tips](#)



Product Overview

Videos

Documents

Catalog Number

4700870

Unit Size

Each

Description

Finnpipette F2 GLP Kit 1: 1-10µL micro, 10-100µL and 100-1000µL

Price (USD)

Price: 977.00

Your Price: 779.00

Availability

In stock



Full specifications ▼

Quantity

Save to list

Catalog Number [4700880](#)
Unit Size Each
Description Finnpipette F2 GLP Kit 2: 0.2-2µL micro, 2-20µL, 20-200µL and 100-1000µL
Price (USD) **Price:** 1,199.00
Your Price: 959.00
Availability **In stock**



Full specifications ▾

Quantity

Add to cart

Save to list

Catalog Number [4700885](#)
Unit Size Each
Description Finnpipette F2 GLP Kit 3: 10-100µL, 100-1000µL and 1-10mL
Price (USD) **Price:** 997.00
Your Price: 795.00
Availability Estimated availability date
Pending 



Full specifications ▾

Quantity

Add to cart

Save to list

Description

Finnpipette F2 GLP Kit 4: 2-20µL, 20-200µL and 100-1000µL

Price (USD)**Price:** 977.00**Your Price:** 779.00**Availability****In stock**

Full specifications ▾

Quantity

[Add to cart](#)[Save to list](#)

Showing 4 of 4

For decades, Thermo Scientific Finn timer pipetting systems have delivered unmatched productivity and ergonomics across thousands of labs and applications worldwide.

Finn pipette F2: Proven reliable for long-term intensive use**Features:**

- Fully autoclavable for protection
- Tough PVDF components that stand up to harsh chemicals and the damaging effects of UV light
- Easy maintenance
- Color coded together with Finntips for easier recognition of compatible tips

Kit 1 includes:

- Three Thermo Scientific Finn timer pipette F2 Pipettes from 1 to 1000 µL:
 - 1-10 µL micro
 - 10-100 µL
 - 100-1000 µL
- Flex tips: 10 micro: 1 x 96; 200: 1 x 96; 1000: 1 x 96
- F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

Kit 2 includes:

- Four Thermo Scientific Finn timer pipette F2 Pipettes from 0.2 to 1000 µL:
 - 0.2-2 µL
 - 2-20 µL univ
 - 20-200 µL
 - 100-1000 µL
- Flex tips: 10 micro: 1 x 96; 200: 2 x 96; 1000: 1 x 96

Kit 3 includes:

- Three Thermo Scientific Finnpiquette F2 Pipettes from 10 to 10,000 μL :
 - 10-100 μL
 - 100-1000 μL
 - 1000-10,000 μL
- Flex tips: 200: 1 x 96; 1000: 1 x 96; FT 10 mL: 1 x 24
- F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

Kit 4 includes:

- Three Thermo Scientific Finnpiquette F2 Pipettes from 20 to 1000 μL :
 - 2-20 μL
 - 20-200 μL
 - 100-1000 μL
- Flex tips: 200: 2 x 96; 1000: 1 x 96
- F-stand
- Good Laboratory Pipetting Guide
- Reagent reservoir demo pack

Ordering Information:

The kit contains samples of Finntip Flex, the most advanced traditional tip offering very low attachment and ejection forces. The Good Laboratory Pipetting Guide contains helpful tips for pipetting including pipetting techniques, ergonomic considerations and how to select tips for best performance.

Warranty: 5 years with web registration.

Videos

Certificates

Search by lot number or partial lot number



Search

[Request a Certificate](#)

Scientific Resources

Application Notes

-  [Technical Bulletin: Ergonomics and usability study of Finnpiquette F1 and F2 Pipet](#)
-  [Application Note: Pipetting performance of Finnpiquette F2 with Finntip Flex](#)
-  [Improve pipetting performance with pipette systems: Finnpiquette F2 pipettes and Finntip Flex pipette tips](#)
-  [Technical Bulletin: Effect of Liquid Properties in Pipetting Liquid Handling Note No. 1](#)
-  [Application Note: Dispense Liquids Containing Proteins More Reliably with Reverse Pipetting](#)
-  [Application Note: Finnpiquette F2-durabilty you can rely on](#)

Technical Notes

-  [Tech Note: Emphazing accuracy and precision of pipette Systems for human identification applications \(EMEA version\)](#)
-  [Tech Note: Emphazing accuracy and precision of pipette Systems for human identification applications](#)

Brochures

-  [Finnpiquette Manual Pipetting Systems Brochure](#)

 [Finnpipette F2移液器](#)

Product Information

Specification Sheets —

-  [Specification Sheet: ClipTip Pipette Tip Specification Sheet](#)
-  [E1-ClipTip Electronic Pipettes Specifications Sheet](#)

Manuals —

-  [Instructions for Safe Use of IVD/CE, Finnpipette and Finntip, Finnpipette F1, F2, F3 and Digital with Finntip](#)
-  [Calibrating your Finnpipettes](#)
-  [Finnpipette Service Request - Decontamination Declaration](#)
-  [Finnpipette F2 Single Channel Variable and Fixed Volume Multichannel Instructions for Use](#)

Instructions —

-  [Introduction:Finnpipette F2 Single Channel Variable and Fixed Volume Multichannel Instructions for Use](#)
-  [Instructions for Safe Use of IVD/CE marked \(Europe\) Finnpipette and Finntip](#)