



[🏠](#) / [PRODUCTS](#) / [EXTRA LARGE CHECK-FLO® INTRODUCERS...](#)

Extra Large Check-Flo® Introducers

[SPECIFICATIONS](#) [IMAGES](#)

Introducers are intended for introduction of balloons, closed and non-tapered end catheters or other diagnostic and interventional devices.

This product line is serviced by the following clinical division(s):
[Aortic Intervention »](#)

[LEARN HOW TO ORDER »](#)

Set up a customer account and order this or other Cook products

[ASK A SUPPORT QUESTION »](#)

Ask Customer Support about our products or services

[FIND PRODUCT TRAINING »](#)

Connect with our Vista® Education & Training program

[VIEW YOUR ACCOUNT INFORMATION »](#)

Register your customer number on our Customer Portal

Specifications

[Give Us Feedback](#)

COLLAPSED EXPANDED

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Introducer Fr*	Accepts Wire Guide Diameter Inch	Sheath Length cm	More Info
G28999	XVCFW-20.0-35-25		N/A	20	.035	25	Expand »

ADDITIONAL SPECS

Description	-
Dilator Length cm	41

G12830	XVCFW-20.0-35-40		N/A	20	.035	40	Expand »
--------	------------------	---	-----	----	------	----	--------------------------

ADDITIONAL SPECS

Description	-
Dilator Length cm	56

G12831	XVCFW-22.0-35-25		N/A	22	.035	25	Expand »
--------	------------------	---	-----	----	------	----	--------------------------

ADDITIONAL SPECS

Description	-
--------------------	---

Order Number	Reference Part Number	Instructions for Use (IFU)	MR Status	Introducer Fr*	Accepts Wire Guide Diameter Inch	Sheath Length cm	More Info
	Dilator Length cm	41					
G12896	XVCFW-24.0-35-25		N/A	24	.035	25	Expand »

ADDITIONAL SPECS

Description	-
Dilator Length cm	41

*Determine the maximum diameter of the instrument or catheter to be introduced to ensure its passage through the sheath. Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Image Library



Not all products shown on this website may be approved or available in all jurisdictions.
Consult with your local Cook representative, distribution company or customer support centre for details.

© Cook 2023

 **+353 61 334 440**

[Products](#)
[Patient Resources](#)
[Support](#)
[About](#)
[Divisions](#)
[Contact](#)
[Education & Training](#)
[Customer Portal](#)

[Privacy Statement](#)
[Terms of Use](#)
[Human Rights](#)
[Responsibility Statement](#)
[Cookie Compliance](#)
[Sitemap](#)

**Sign Up. Get Email about News, Products,
and Events from Cook Medical.**

[SIGN UP](#)