

SECTION 3 - PRODUCT OVERVIEW

Unpacking and Inspection

- Unpack the equipment.
- Check the packing list to insure all accessories are there.
- Examine the packing material carefully for separately packed items.
- Check the incubator for external surface damage.
- Check that the display is not broken or damaged.

Accessories

Standard Accessories

24 volts main adapter (PSU)
Order code: K51122



Mains Cable

Order code: 1212807



Optional Equipment

Optional equipment consists of a range of warming blocks. Please contact your local CooperSurgical distributor or CooperSurgical, Inc. for details about the full range of warming blocks.

Warming block B01

Holds seven tubes Ø17 mm
Falcon 352001
Order code: 26001



Warming block B06

Holds twelve tubes of
Ø12mm
Order code: 26006



Warming block B03

Holds twelve tubes of Ø12mm
Falcon 352003
Order code: 26003



Warming block B08

Holds ten tubes of
Ø17mm
Falcon 352001
Order code: 26008



Warming block B04

Holds four tubes of
25ml – Ø23,5mm
Order code: 26004



Warming block B40

Holds one IVF
MediaFlask of 125ml
Order code: 26055

