


Polyp trap – 2 drawers



Description

- Recovery of polyps,
- 2 removal drawers,
- Magnifying glass.



Primary packaging : Individual box
Secondary packaging: cardboard (with vegetal ink) of 24
Minimum Delivery Quantity (MDQ) in units: 24

References	Technical specifications
VTL/P-TRAP	2 removable drawers polyp trap with retrieval stick

Indication

The single-use polyp trap is used to retrieve polyps removed by endoscopy.

Technical informations

Materials :

Suction tube : Silicone
Trap body : ABS
Filter : TPE

Contains :

Latex :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No	Bisphenol A :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No
Product of animal origin :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No	Nanoparticules :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No
PVC :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No	Lead :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No
Phtalate and DEHP :	<input type="checkbox"/>	Yes	<input checked="" type="checkbox"/>	No					

Sterility :

Sterile : Yes No
Sterilization mode : N/A
Number of use : Single use

Storage :

Storage conditions : Keep away from the light.
Use by date : 5 years

Other data :

GMDN code: 34858
Barcode : No GS1-128 Standard HIBC Standard
Si oui, avec les informations : Reference Lot number Use by date

Use :

Correct and safe use : Cf. Instructions for use
Format : Electronical Paper
User advices : Instructions for use, precaution and contraindications – Cf. Instructions for use

Manufacturer:

Manufacturer : Pergo Medikal
Country : Turkey

CE Mark :

Notified body : Autocertification
Standard: European Regulation (UE) 2017/745
Class : Class I medical device

IMPORTER

VYTIL

161 avenue Gallieni – 93170 Bagnolet, France

Phone : +33 1 49 88 01 01

www.vytil.com / info@vytil.com

Order :

E-mail : sales@vytil.com

Vigilance :

E-mail : quality@vytil.com

