



HI-TORQUE INFILTRAC™

Guide Wire



**Information contained herein for DISTRIBUTION outside of the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

STATE-OF-THE-ART TIP TECHNOLOGY FOR PENETRATING RESISTANT CTO CAPS

TAPERED TIP


for high penetration powers

PRE-SHAPED MICRO-J TIP

for enhanced steering

UNIQUE MICRO-TEXTURED TIP


for CTO cap engagement



Data on file at Abbott.

**Information contained herein for DISTRIBUTION outside of the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

©2023 Abbott. All rights reserved. MAT-2302347 v2.0



ORDERING INFORMATION



HI-TORQUE INFILTRAC™ Guide Wire

NON-CE* PART NUMBER	CE PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	PENETRATION POWER (kg/in ²)	TIP DIAMETER (in)	TIP LOAD (g)
1030001J	1030005J	190	Micro-J	167	0.009	11
1030002J	1030006J	300	Micro-J	167	0.009	11

HI-TORQUE INFILTRAC™ PLUS Guide Wire

NON-CE* PART NUMBER	CE PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	PENETRATION POWER (kg/in ²)	TIP DIAMETER (in)	TIP LOAD (g)
1030003J	1030007J	190	Micro-J	224	0.009	14
1030004J	1030008J	300	Micro-J	224	0.009	14

*Same as U.S. FDA part numbers applicable to select OUS countries.
Data on file at Abbott.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

**Information contained herein for DISTRIBUTION outside of the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Abbott International BVBA
Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies

www.cardiovascular.abbott

©2023 Abbott. All rights reserved. MAT-2302347 v2.0

