

HI-TORQUE PILOT™

Guide Wires

THE RIGHT TOOL FOR A SOFT APPROACH

A CTO INDICATED GUIDE WIRE

CORE-TO-TIP DESIGN

for excellent torque response¹

FULL POLYMER COVER

for smooth vessel tracking

POLYMER COVER AND HYDROPHILIC COATING

for outstanding lubricity¹

3 OPTIONS OF INCREASING TIP LOADS

for CTO lesion navigation

Data on file at Abbott.

1. Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645–652.

**Information contained herein for DISTRIBUTION outside the U.S. ONLY.
Check the regulatory status of the device in areas where CE marking is not the regulation in force.**

© 2021 Abbott. All rights reserved. MAT-2102003 v1.0



HI-TORQUE PILOT™ GUIDE WIRE ORDERING INFORMATION



HI-TORQUE PILOT™ 50


PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE PILOT™ 50	1010480-H	190	Straight	3	1.3
HI-TORQUE PILOT™ 50	1010480-HJ	190	J	3	1.3
HI-TORQUE PILOT™ 50	1010483-H	300	Straight	3	1.3
HI-TORQUE PILOT™ 50	1010483-HJ	300	J	3	1.3

HI-TORQUE PILOT™ 150

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE PILOT™ 150	1010481-H	190	Straight	3	2.5
HI-TORQUE PILOT™ 150	1010481-HJ	190	J	3	2.5
HI-TORQUE PILOT™ 150	1010484-H	300	Straight	3	2.5
HI-TORQUE PILOT™ 150	1010484-HJ	300	J	3	2.5

HI-TORQUE PILOT™ 200

PRODUCT DESCRIPTION	PART NUMBER	WIRE LENGTH (cm)	TIP SHAPE	RADIOPAQUE TIP LENGTH (cm)	TIP LOAD* (g)
HI-TORQUE PILOT™ 200	1010482-H	190	Straight	3	3.9
HI-TORQUE PILOT™ 200	1010482-HJ	190	J	3	3.9
HI-TORQUE PILOT™ 200	1010485-H	300	Straight	3	3.9
HI-TORQUE PILOT™ 200	1010485-HJ	300	J	3	3.9



*Tests performed by and data on file at Abbott.

1. Tóth GG, Yamane M, Heyndrickx GR. Heart 2015;101:645–652.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at vascular.eifu.abbott or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. This material is intended for use with healthcare professionals only.

Information contained herein for DISTRIBUTION outside the U.S. ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

Data on file at Abbott.

Abbott International BVBA

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

www.cardiovascular.abbott

© 2021 Abbott. All rights reserved. MAT-2102003 v1.0