


MagicTouch PTA

SIROLIMUS COATED PTA BALLOON CATHETER



TECHNOLOGY

An innovative proprietary Nanolute technology providing better bioavailability of the drug.



DRUG

A proven drug with a safe profile



COATING

A uniform coating in the inflated position of the balloon catheter.



CARRIER

A phospholipid carrier which is biocompatible in nature

MagicTouch PTA

SIROLIMUS COATED PTA BALLOON CATHETER

0.014" OTW

Shaft Length: 90 cm,
120 cm & 150 cm

Baliono diametrij intervalas turi būti ne mažesnis nei nuo 1,5 mm ikimm

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.014" OTW)			NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft Length							
		90 cm	120 cm	150 cm					
1.5	20	CPTW14025	CPT14150020	CPTW14061	8	16	0.014"	5F	6F
2	40	CPTW14026	CPT14200040	CPTW14062	8	14	0.014"	5F	6F
2	60	CPTW14027	CPT14200060	CPTW14063	8	14	0.014"	5F	6F
2	80	CPTW14028	CPT14200080	CPTW14064	8	14	0.014"	5F	6F
2	100	CPTW14029	CPT14200100	CPTW14065	6	14	0.014"	5F	6F
2	120	CPTW14030	CPT14200120	CPTW14066	6	14	0.014"	6F	7F
2	150	CPTW14031	CPT14200150	CPTW14067	6	14	0.014"	6F	7F
2	200	CPTW14032	CPT14200200	CPTW14068	6	14	0.014"	6F	7F
2.5	40	CPTW14033	CPT14250040	CPTW14069	8	14	0.014"	5F	6F
2.5	60	CPTW14034	CPT14250060	CPTW14070	8	14	0.014"	5F	6F
2.5	80	CPTW14035	CPT14250080	CPTW14071	8	14	0.014"	5F	6F
2.5	100	CPTW14036	CPT14250100	CPTW14072	6	12	0.014"	5F	6F
2.5	120	CPTW14037	CPT14250120	CPTW14073	6	12	0.014"	6F	7F
2.5	150	CPTW14038	CPT14250150	CPTW14074	6	12	0.014"	6F	7F
2.5	200	CPTW14039	CPT14250200	CPTW14075	6	12	0.014"	6F	7F
3	40	CPTW14040	CPT14300040	CPTW14076	8	14	0.014"	5F	6F
3	60	CPTW14041	CPT14300060	CPTW14077	8	14	0.014"	5F	6F
3	80	CPTW14042	CPT14300080	CPTW14078	8	14	0.014"	5F	6F
3	100	CPTW14043	CPT14300100	CPTW14079	6	12	0.014"	5F	6F
3	120	CPTW14044	CPT14300120	CPTW14080	6	12	0.014"	6F	7F
3	150	CPTW14045	CPT14300150	CPTW14081	6	12	0.014"	6F	7F
3	200	CPTW14046	CPT14300200	CPTW14082	6	12	0.014"	6F	7F
3.5	40	CPTW14047	CPT14350040	CPTW14083	8	14	0.014"	5F	6F
3.5	60	CPTW14048	CPT14350060	CPTW14084	8	14	0.014"	5F	6F
3.5	80	CPTW14049	CPT14350080	CPTW14085	8	14	0.014"	5F	6F
3.5	100	CPTW14050	CPT14350100	CPTW14086	6	12	0.014"	5F	6F
3.5	120	CPTW14051	CPT14350120	CPTW14087	6	12	0.014"	6F	7F
3.5	150	CPTW14052	CPT14350150	CPTW14088	6	12	0.014"	6F	7F
3.5	200	CPTW14053	CPT14350200	CPTW14089	6	12	0.014"	6F	7F
4	40	CPTW14054	CPT14400040	CPTW14090	8	14	0.014"	6F	7F
4	60	CPTW14055	CPT14400060	CPTW14091	8	14	0.014"	6F	7F
4	80	CPTW14056	CPT14400080	CPTW14092	8	14	0.014"	6F	7F

MagicTouch PTA

SIROLIMUS COATED PTA BALLOON CATHETER

0.014" OTW

Shaft Length: 90 cm,
120 cm & 150 cm

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.014" OTW)			NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft Length							
		90 cm	120 cm	150 cm					
4	100	CPTW14057	CPT14400100	CPTW14093	6	12	0.014"	6F	7F
4	120	CPTW14058	CPT14400120	CPTW14094	6	12	0.014"	6F	7F
4	150	CPTW14059	CPT14400150	CPTW14095	6	12	0.014"	6F	7F
4	200	CPTW14060	CPT14400200	CPTW14096	6	12	0.014"	6F	7F

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.018" OTW)			NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft Length							
		90 cm	120 cm	150 cm					
2	40	CPTW18099	CPT18200040	CPTW18156	8	14	0.018"	5F	6F
2	60	CPTW18100	CPT18200060	CPTW18157	8	14	0.018"	5F	6F
2	80	CPTW18101	CPT18200080	CPTW18158	8	14	0.018"	5F	6F
2	100	CPTW18102	CPT18200100	CPTW18159	6	12	0.018"	5F	6F
2	120	CPTW18103	CPT18200120	CPTW18160	6	12	0.018"	6F	7F
2	150	CPTW18104	CPT18200150	CPTW18161	6	12	0.018"	6F	7F
2	200	CPTW18105	CPT18200200	CPTW18162	6	12	0.018"	6F	7F
2.5	40	CPTW18106	CPT18250040	CPTW18163	8	14	0.018"	5F	6F
2.5	60	CPTW18107	CPT18250060	CPTW18164	8	14	0.018"	5F	6F
2.5	80	CPTW18108	CPT18250080	CPTW18165	8	14	0.018"	5F	6F
2.5	100	CPTW18109	CPT18250100	CPTW18166	6	12	0.018"	5F	6F
2.5	120	CPTW18110	CPT18250120	CPTW18167	6	12	0.018"	6F	7F
2.5	150	CPTW18111	CPT18250150	CPTW18168	6	12	0.018"	6F	7F
2.5	200	CPTW18112	CPT18250200	CPTW18169	6	12	0.018"	6F	7F
3	40	CPTW18113	CPT18300040	CPTW18170	8	14	0.018"	5F	6F
3	60	CPTW18114	CPT18300060	CPTW18171	8	14	0.018"	5F	6F
3	80	CPTW18115	CPT18300080	CPTW18172	8	14	0.018"	5F	6F
3	100	CPTW18116	CPT18300100	CPTW18173	6	12	0.018"	5F	6F
3	120	CPTW18117	CPT18300120	CPTW18174	6	12	0.018"	6F	7F
3	150	CPTW18118	CPT18300150	CPTW18175	6	12	0.018"	6F	7F
3	200	CPTW18119	CPT18300200	CPTW18176	6	12	0.018"	6F	7F
3.5	40	CPTW18120	CPT18350040	CPTW18177	8	14	0.018"	6F	7F
3.5	60	CPTW18121	CPT18350060	CPTW18178	8	14	0.018"	6F	7F
3.5	80	CPTW18122	CPT18350080	CPTW18179	8	14	0.018"	6F	7F
3.5	100	CPTW18123	CPT18350100	CPTW18180	6	12	0.018"	6F	7F
3.5	120	CPTW18124	CPT18350120	CPTW18181	6	12	0.018"	6F	7F
3.5	150	CPTW18125	CPT18350150	CPTW18182	6	12	0.018"	6F	7F
3.5	200	CPTW18126	CPT18350200	CPTW18183	6	12	0.018"	6F	7F
4	40	CPTW18127	CPT18400040	CPTW18184	8	14	0.018"	6F	7F
4	60	CPTW18128	CPT18400060	CPTW18185	8	14	0.018"	6F	7F
4	80	CPTW18129	CPT18400080	CPTW18186	8	14	0.018"	6F	7F

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.018" OTW)			NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft Length							
		90 cm	120 cm	150 cm					
4	100	CPTW18130	CPT18400100	CPTW18187	6	12	0.018"	6F	7F
4	120	CPTW18131	CPT18400120	CPTW18188	6	12	0.018"	6F	7F
4	150	CPTW18132	CPT18400150	CPTW18189	6	12	0.018"	6F	7F
4	200	CPTW18133	CPT18400200	CPTW18190	6	12	0.018"	6F	7F
5	40	CPTW18134	CPT18500040	CPTW18191	8	14	0.018"	7F	8F
5	60	CPTW18135	CPT18500060	CPTW18192	8	14	0.018"	7F	8F
5	80	CPTW18136	CPT18500080	CPTW18193	8	14	0.018"	7F	8F
5	100	CPTW18137	CPT18500100	CPTW18194	6	12	0.018"	7F	8F
5	120	CPTW18138	CPT18500120	CPTW18195	6	12	0.018"	7F	8F
5	150	CPTW18139	CPT18500150	CPTW18196	6	12	0.018"	7F	8F
5	200	CPTW18140	CPT18500200	CPTW18197	6	12	0.018"	7F	8F
5.5	40	CPTW18141	CPT18550040	CPTW18198	8	14	0.018"	7F	8F
5.5	60	CPTW18142	CPT18550060	CPTW18199	8	14	0.018"	7F	8F
5.5	80	CPTW18143	CPT18550080	CPTW18200	8	14	0.018"	7F	8F
5.5	100	CPTW18144	CPT18550100	CPTW18201	6	12	0.018"	7F	8F
5.5	120	CPTW18145	CPT18550120	CPTW18202	6	12	0.018"	7F	8F
5.5	150	CPTW18146	CPT18550150	CPTW18203	6	12	0.018"	7F	8F
5.5	200	CPTW18147	CPT18550200	CPTW18204	6	12	0.018"	7F	8F
6	40	CPTW18148	CPT18600040	CPTW18205	8	14	0.018"	7F	8F
6	60	CPTW18149	CPT18600060	CPTW18206	8	14	0.018"	7F	8F
6	80	CPTW18150	CPT18600080	CPTW18207	8	14	0.018"	7F	8F
6	100	CPTW18151	CPT18600100	CPTW18208	6	12	0.018"	7F	8F
6	120	CPTW18152	CPT18600120	CPTW18209	6	12	0.018"	7F	8F
6	150	CPTW18153	CPT18600150	CPTW18210	6	12	0.018"	7F	8F
6	200	CPTW18154	CPT18600200	CPTW18211	6	12	0.018"	7F	8F
7	40	CPTW18155	CPT18700040	CPTW18212	8	14	0.018"	7F	8F

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.035" OTW)	NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft length: 45 cm					
3	20	CPTW35113	8	14	0.035"	6F	7F
3	40	CPTW35114	8	14	0.035"	6F	7F
3	60	CPTW35115	8	14	0.035"	6F	7F
3	80	CPTW35116	8	14	0.035"	6F	7F
4	20	CPTW35117	8	14	0.035"	6F	7F
4	40	CPTW35118	8	14	0.035"	6F	7F
4	60	CPTW35119	8	14	0.035"	6F	7F
4	80	CPTW35120	8	14	0.035"	6F	7F
4	100	CPTW35121	6	12	0.035"	6F	7F
4	120	CPTW35122	6	12	0.035"	6F	7F
5	20	CPTW35123	8	14	0.035"	7F	8F
5	40	CPTW35124	8	14	0.035"	7F	8F
5	60	CPTW35125	8	14	0.035"	7F	8F
5	80	CPTW35126	8	14	0.035"	7F	8F
5	100	CPTW35127	6	12	0.035"	7F	8F
5	120	CPTW35128	6	12	0.035"	7F	8F
6	20	CPTW35129	8	14	0.035"	7F	8F
6	40	CPTW35130	8	14	0.035"	7F	8F
6	60	CPTW35131	8	14	0.035"	7F	8F
6	80	CPTW35132	8	14	0.035"	7F	8F
6	100	CPTW35133	6	12	0.035"	7F	8F
6	120	CPTW35134	6	12	0.035"	7F	8F
7	20	CPTW35135	8	13	0.035"	7F	8F
7	40	CPTW35136	8	13	0.035"	7F	8F
7	60	CPTW35137	8	13	0.035"	7F	8F
7	80	CPTW35138	8	13	0.035"	7F	8F
7	100	CPTW35139	6	10	0.035"	7F	8F
8	20	CPTW35140	8	12	0.035"	8F	9F
8	40	CPTW35141	8	12	0.035"	8F	9F
8	60	CPTW35142	8	12	0.035"	8F	9F

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.035" OTW)	NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft length: 45 cm					
8	80	CPTW35143	8	12	0.035"	8F	9F
8	100	CPTW35144	5	10	0.035"	8F	9F
9	20	CPTW35145	8	12	0.035"	8F	9F
9	40	CPTW35146	8	12	0.035"	8F	9F
9	60	CPTW35147	8	12	0.035"	8F	9F
9	80	CPTW35148	8	12	0.035"	8F	9F
10	20	CPTW35149	6	10	0.035"	8F	9F
10	40	CPTW35150	6	10	0.035"	8F	9F
10	60	CPTW35151	6	10	0.035"	8F	9F
10	80	CPTW35152	6	10	0.035"	8F	9F

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.035" OTW)				NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft Length								
		80 cm	90 cm	130 cm	150 cm					
3	20	CPTW35001	CPTW35057	CPTW35153	CPTW35209	8	14	0.035"	6F	7F
3	40	CPTW35002	CPTW35058	CPTW35154	CPTW35210	8	14	0.035"	6F	7F
3	60	CPTW35003	CPTW35059	CPTW35155	CPTW35211	8	14	0.035"	6F	7F
3	80	CPTW35004	CPTW35060	CPTW35156	CPTW35212	8	14	0.035"	6F	7F
4	20	CPTW35005	CPTW35061	CPTW35157	CPTW35213	8	14	0.035"	6F	7F
4	40	CPTW35006	CPTW35062	CPTW35158	CPTW35214	8	14	0.035"	6F	7F
4	60	CPTW35007	CPTW35063	CPTW35159	CPTW35215	8	14	0.035"	6F	7F
4	80	CPTW35008	CPTW35064	CPTW35160	CPTW35216	8	14	0.035"	6F	7F
4	100	CPTW35009	CPTW35065	CPTW35161	CPTW35217	6	12	0.035"	6F	7F
4	120	CPTW35010	CPTW35066	CPTW35162	CPTW35218	6	12	0.035"	6F	7F
4	150	CPTW35011	CPTW35067	CPTW35163	CPTW35219	6	12	0.035"	6F	7F
4	200	CPTW35012	CPTW35068	CPTW35164	CPTW35220	6	12	0.035"	6F	7F
5	20	CPTW35013	CPTW35069	CPTW35165	CPTW35221	8	14	0.035"	7F	8F
5	40	CPTW35014	CPTW35070	CPTW35166	CPTW35222	8	14	0.035"	7F	8F
5	60	CPTW35015	CPTW35071	CPTW35167	CPTW35223	8	14	0.035"	7F	8F
5	80	CPTW35016	CPTW35072	CPTW35168	CPTW35224	8	14	0.035"	7F	8F
5	100	CPTW35017	CPTW35073	CPTW35169	CPTW35225	6	12	0.035"	7F	8F
5	120	CPTW35018	CPTW35074	CPTW35170	CPTW35226	6	12	0.035"	7F	8F
5	150	CPTW35019	CPTW35075	CPTW35171	CPTW35227	6	12	0.035"	7F	8F
5	200	CPTW35020	CPTW35076	CPTW35172	CPTW35228	6	12	0.035"	7F	8F
6	20	CPTW35021	CPTW35077	CPTW35173	CPTW35229	8	14	0.035"	7F	8F
6	40	CPTW35022	CPTW35078	CPTW35174	CPTW35230	8	14	0.035"	7F	8F
6	60	CPTW35023	CPTW35079	CPTW35175	CPTW35231	8	14	0.035"	7F	8F
6	80	CPTW35024	CPTW35080	CPTW35176	CPTW35232	8	14	0.035"	7F	8F
6	100	CPTW35025	CPTW35081	CPTW35177	CPTW35233	6	12	0.035"	7F	8F
6	120	CPTW35026	CPTW35082	CPTW35178	CPTW35234	6	12	0.035"	7F	8F
6	150	CPTW35027	CPTW35083	CPTW35179	CPTW35235	6	12	0.035"	7F	8F
6	200	CPTW35028	CPTW35084	CPTW35180	CPTW35236	6	12	0.035"	7F	8F
7	20	CPTW35029	CPTW35085	CPTW35181	CPTW35237	8	13	0.035"	7F	8F
7	40	CPTW35030	CPTW35086	CPTW35182	CPTW35238	8	13	0.035"	7F	8F
7	60	CPTW35031	CPTW35087	CPTW35183	CPTW35239	8	13	0.035"	7F	8F
7	80	CPTW35032	CPTW35088	CPTW35184	CPTW35240	8	13	0.035"	7F	8F
7	100	CPTW35033	CPTW35089	CPTW35185	CPTW35241	6	10	0.035"	7F	8F
7	120	CPTW35034	CPTW35090	CPTW35186	CPTW35242	6	10	0.035"	7F	8F
7	150	CPTW35035	CPTW35091	CPTW35187	CPTW35243	6	10	0.035"	7F	8F
7	200	CPTW35036	CPTW35092	CPTW35188	CPTW35244	6	10	0.035"	7F	8F

Kateterio ilgiai turi apimti ne mažesnj intervalą nuo 80 cm iki 150 cm

Balloon Diameter (mm)	Balloon Length (mm)	Magic Touch PTA (0.035" OTW)				NP (atm)	RBP (atm)	Recommended Guiding wire	Recommended Sheath Compatibility	Guiding Catheter Compatibility
		Usable Shaft Length								
		80 cm	90 cm	130 cm	150 cm					
8	20	CPTW35037	CPTW35093	CPTW35189	CPTW35245	8	12	0.035"	8F	9F
8	40	CPTW35038	CPTW35094	CPTW35190	CPTW35246	8	12	0.035"	8F	9F
8	60	CPTW35039	CPTW35095	CPTW35191	CPTW35247	8	12	0.035"	8F	9F
8	80	CPTW35040	CPTW35096	CPTW35192	CPTW35248	8	12	0.035"	8F	9F
8	100	CPTW35041	CPTW35097	CPTW35193	CPTW35249	5	10	0.035"	8F	9F
8	120	CPTW35042	CPTW35098	CPTW35194	CPTW35250	5	10	0.035"	8F	9F
8	150	CPTW35043	CPTW35099	CPTW35195	CPTW35251	5	10	0.035"	8F	9F
8	200	CPTW35044	CPTW35100	CPTW35196	CPTW35252	5	10	0.035"	8F	9F
9	20	CPTW35045	CPTW35101	CPTW35197	CPTW35253	8	12	0.035"	8F	9F
9	40	CPTW35046	CPTW35102	CPTW35198	CPTW35254	8	12	0.035"	8F	9F
9	60	CPTW35047	CPTW35103	CPTW35199	CPTW35255	8	12	0.035"	8F	9F
9	80	CPTW35048	CPTW35104	CPTW35200	CPTW35256	8	12	0.035"	8F	9F
10	20	CPTW35049	CPTW35105	CPTW35201	CPTW35257	6	10	0.035"	8F	9F
10	40	CPTW35050	CPTW35106	CPTW35202	CPTW35258	6	10	0.035"	8F	9F
10	60	CPTW35051	CPTW35107	CPTW35203	CPTW35259	6	10	0.035"	8F	9F
10	80	CPTW35052	CPTW35108	CPTW35204	CPTW35260	6	10	0.035"	8F	9F
12	20	CPTW35053	CPTW35109	CPTW35205	CPTW35261	6	10	0.035"	8F	9F
12	40	CPTW35054	CPTW35110	CPTW35206	CPTW35262	6	10	0.035"	8F	9F
12	60	CPTW35055	CPTW35111	CPTW35207	CPTW35263	6	10	0.035"	8F	9F
12	80	CPTW35056	CPTW35112	CPTW35208	CPTW35264	6	10	0.035"	8F	9F

Baliono diametrų intervalas turi būti ne mažesnis nei nuo mm iki 12 mm