



4-Level Wash Cart

Article No 503284900

Article No 503288500 US Wash Cart

Especially designed to take advantage of the taller chamber in the 46-5 series, this wash cart has four shelves for loading larger goods (three removable shelves and spray arms).

Washer-Disinfector Model:

Getinge 46-5

Getinge WD46-5 Turbo

CAPACITY

8 SPRI II (340x250x70 mm / 13.3x9.8x2.8 in) or

8 DIN (480x250x50 mm / 18.9x9.8x1.97 in)

Free loading height from shelf to spray arm: 104 mm / 4.1 in

Maximum available shelf space for loading: 592x488 mm / 23.3x19.2 in



5-Level Wash Cart

Article No 503284100

This wash cart is designed to maximize productivity in processing standard instrument trays with 5 shelves for loading (four removable shelves and spray arms).

Washer-Disinfector Model:

Getinge 46-5

Getinge WD46-5 Turbo

CAPACITY

10 SPRI II (340x250x70 mm / 13.3x9.8x2.8 in) or

10 DIN (480x250x50 mm / 18.9x9.8x1.97 in))

Free loading height from shelf to spray arm: 72 mm / 2.8 in

Maximum available shelf space for loading: 592x488 mm / 23.3x19.2 in



Manifold Connections

Article No 503393200 Silicon cap

Article No 503300000 Hose nipple

The manifolds on all instrument wash carts are prepared for 6 hose nipple connections. These connections are used to clean and disinfect a lower number of tubular instruments on a regular wash cart without the need of a MIS-module. Always use the silicon cap when no instrument is connected to the hose nipple.