



PTA BALLOON CATHETER

Ordering Information

| Balloon Diameter (mm) | Balloon Working Length (mm) | | | | | New Sizes | 0.014" GW |
|-----------------------|-----------------------------|------------|------------|------------|------------|------------|------------|
| | 12 | 20 | 40 | 80 | 120 | | |
| 1.2 | 515-015-41 | 515-020-41 | 515-040-41 | 515-080-41 | 515-120-41 | 515-180-41 | 515-240-41 |
| 2.0 | 520-015-41 | 520-020-41 | 520-040-41 | 520-080-41 | 520-120-41 | 520-180-41 | 520-240-41 |
| 2.5 | - | 525-020-41 | 525-040-41 | 525-080-41 | 525-120-41 | 525-180-41 | 525-240-41 |
| 3.0 | 530-015-41 | 530-020-41 | 530-040-41 | 530-080-41 | 530-120-41 | 530-180-41 | 530-240-41 |
| 3.5* | - | 535-020-41 | 535-040-41 | 535-080-41 | 535-120-41 | 535-180-41 | 535-240-41 |
| 4.0 | 540-015-41 | 540-020-41 | 540-040-41 | 540-080-41 | 540-120-41 | 540-180-41 | 540-240-41 |
| 4.5* | - | 545-020-41 | 545-040-41 | 545-080-41 | 545-120-41 | - | - |
| 5.0 | 550-015-41 | 550-020-41 | 550-040-41 | 550-080-41 | 550-120-41 | 550-180-41 | 550-240-41 |
| 5.5* | - | 555-020-41 | 555-040-41 | 555-080-41 | 555-120-41 | - | - |
| 6.0 | 560-015-41 | 560-020-41 | 560-040-41 | 560-080-41 | 560-120-41 | 560-180-41 | 560-240-41 |


*Available only in selected countries. Please contact our local sales representative for further information.

Compliance Chart

| Pressure atm (kPa) | 1.5 | 2.0 | 2.5 | 3.0 | 3.5 | 4.0 | 4.5 | 5.0 | 5.5 | 6.0 |
|--------------------|------|------|------|------|------|------|------|------|------|------|
| 2 (203) | 1.41 | 1.88 | 2.37 | 2.82 | 3.27 | 3.70 | 4.22 | 4.67 | 5.09 | 5.60 |
| 4 (405) | 1.43 | 1.90 | 2.40 | 2.86 | 3.32 | 3.76 | 4.27 | 4.74 | 5.17 | 5.68 |
| 6 (608) | 1.45 | 1.93 | 2.42 | 2.89 | 3.36 | 3.82 | 4.33 | 4.80 | 5.26 | 5.76 |
| 8 (811) | 1.47 | 1.95 | 2.45 | 2.93 | 3.41 | 3.88 | 4.39 | 4.87 | 5.34 | 5.84 |
| 10 (1013) | 1.48 | 1.98 | 2.47 | 2.96 | 3.45 | 3.94 | 4.44 | 4.93 | 5.42 | 5.92 |
| NOM* 12 (1216) | 1.50 | 2.00 | 2.50 | 3.00 | 3.50 | 4.00 | 4.50 | 5.00 | 5.50 | 6.00 |
| 14 (1419) | 1.52 | 2.02 | 2.53 | 3.04 | 3.55 | 4.06 | 4.56 | 5.07 | 5.58 | 6.08 |
| 16 (1621) | 1.53 | 2.05 | 2.55 | 3.07 | 3.59 | 4.12 | 4.61 | 5.13 | 5.66 | 6.16 |
| 18 (1824) | 1.55 | 2.07 | 2.58 | 3.11 | 3.64 | 4.18 | 4.67 | 5.20 | 5.74 | 6.24 |
| 20 (2027) | 1.57 | 2.10 | 2.60 | 3.14 | 3.68 | 4.24 | 4.73 | 5.26 | 5.83 | 6.32 |
| RBP** 22 (2229) | 1.59 | 2.12 | 2.63 | 3.18 | 3.73 | 4.30 | 4.78 | 5.33 | 5.91 | 6.40 |
| 24 (2432) | 1.60 | 2.14 | 2.65 | 3.22 | 3.78 | 4.36 | 4.84 | 5.40 | - | - |
| 26 (2634) | 1.62 | 2.17 | 2.68 | 3.25 | 3.82 | 4.42 | - | - | - | - |

*Nominal Pressure. The nominal in vitro device specifications do not take into account any lesion resistance.

**Rated Burst Pressure - Do Not Exceed (1.5-4.0mm: 22atm, 4.5-5.0mm: 20atm, 5.5-6.0mm: 18atm)



For more information please visit our website at www.OrbusNeich.com or contact us:

EMEA HEADQUARTERS

Drs. W. van Royenstraat 5
3871 AN Hoevelaken
The Netherlands

Phone +31-33.254.1150

Toll Free Phone 00800.0254.1150*

* Only for Belgium, Denmark, France, Germany, Ireland, Netherlands, Norway, Sweden and UK

©2018 OrbusNeich Medical. All rights reserved.

The listed catalog numbers are not available for sale in the USA.

#G-70-0547 Rev02

LIFE and OrbusNeich are registered trademarks of OrbusNeich Medical, Inc. All data and photos on file. Illustrations are not to be considered as engineering drawings or photos.



OrbusNeich®
Pioneers in life-changing technologies



180 mm & 240 mm
balloon length
NOW
AVAILABLE!

Optimal Combination of
Robustness and Superb Deliverability