

s e n s e u s TM

The New Endo-Sensation

DENTSPLY
MAILLEFER



Sensory

Sensitivity

Success

Sensation

Sensuality

- Experience a new standard of working comfort and tactile feel
- Explore a new level of control and precision
- Enjoy a new level of clinical performance and security
- Experience a new World of Endo-Sensations


Start Enjoying Your Endo

s e n s e u s TM



The large, soft, ergonomic silicone handle transforms your working sensations and improves your clinical performance.

- Greater Working Comfort
- Improved Torque and Force Transmission, due to :
 - Reduction of effort achieved by large handle
 - Excellent grip when used with gloves
- Improved Tactile Feedback
- Reduction of Fatigue



X-ray visible calibration rings will improve your clinical performance during endodontic procedures, through :

- Improved Diagnostics
- Increased Accuracy
- Greater Visibility



FlexoFile® FlexoReamer®

The highly flexible instrument blades with guiding action preserve canal morphology and improve clinical security.

- More Reliable Instrumentation of Canals
- Reduced Risks of Transportation, Ledging and Perforation
- Reduced Risk of Fracture of Instruments



Triangular Cross-Section*


- High Flexibility
- High Cutting Efficiency



Batt® Tip

- Non-Cutting 55° Batt® Tip** with Rounded End
- Safe and Secure Guiding Action

* ISO 015 and up; square cross-section for smaller sizes



Hedstroem

- Improved Preservation of Canal Morphology and Reduced Transportation and Ledging
- Excellent Cutting Efficiency



High Quality Blade


- Proven Instrument Design
- High Quality Swedish Steel



Guiding Tip

- Non-Cutting 55° Tip**
- Guiding Action

** 75° tip for instrument sizes ISO 055 and up




ProFinder™

Root Canal Catheterization File

ProFinder™ features an exclusive minimal, variable taper, giving you unmatched root canal catheterization capabilities.

- Easier to Establish the Glide Path
- Faster and More Efficient Negotiation of Difficult or Calcified Canals
- Easier Penetration of Canal
- Better Tactile Feel



Variable & Minimal Taper Blade

- Square Cross-Section with Increased Taper in Tip Region
 - Increased buckling resistance
- Minimal Taper Behind Tip Region
 - Higher instrument flexibility

Guiding Tip



- Non-Cutting 65° Tip with Rounded End
- Guiding Action

*“ If you own the glide path,
you can consistently win the rotary
game of endodontics. ”*

Dr. Clifford J. Ruddle




	Ø ISO	006	008	010	015	020	025	030	035	040	045	050	055	060	070	080
Refills	6x															
Assortments	015-040															
	045-080															



	Ø ISO	008	010	015	020	025	030	035	040	045	050	055	060	070	080
Refills	6x														
Assortments	015-040														
	045-080														



	Ø ISO	008	010	015	020	025	030	035	040	045	050	055	060	070	080
Refills	6x														
Assortments	015-040														
	045-080														



		010	013	017
Refills	6x			
Assortments	010-017			

* Available Autumn 2005

