

DDMS™ Diamond Dusted Membrane Scrapers

The Diamond-dusted, soft, silicone tip has been designed to identify and lift the edge of the membrane quickly and easily for your forceps. Available in several styles, including a retractable version, the patented DDMS has become a staple in operating rooms around the world.

DDMS benefits:

- Aids in the atraumatic removal of immature membranes and pigmented cells from the retina
- Assists atraumatic separation of hyaloids from neurosensory retina
- Especially beneficial for relieving traction from the rim of macular holes or retinal breaks in an atraumatic fashion

Tano

2004 20G   = 5
 2004.23 23G
 2004.25S (Stiff) 25G

NEW!

2004.27 27G



Stiff Tano

2007 20G   = 5



Retractable Tano

2016 20G   = 3
 2016.23 23G   = 5
 2016.25 25G

98

97



Adjustable Tano with Pick

2018 20G   = 3



Disposable Infusion Cannulas

6mm Tip

30.32 20G   = 5
 • Ideal aphakic length



6mm Tip

30.32.23 23G   = 6
 30.32.25 25G
 • For severe vitreous hemorrhaging
 • For use with Synergetics® steel cannula sets (sold separately)



Bipolar Bipolar Cautery

53.00.23 23G   = 6
 53.00.25 25G

NEW!

53.00.27 27G



Iris Retractors Disposable

40.02-2   = 2
 40.02-6   = 6
 • 5 retractors/set in case one is lost or damaged



Visco Dissector Hubbard Visco Dissector

12.06.25 25G 28G Tip   = 6

NEW!

12.06.27 27G

- Injects viscoelastic for delaminating large membranes from the retina



VRL100 DK-Line® 5ml   = 1

VRL110 DK-Line® 7ml   = 1



VRL200 Okta-Line® 5ml   = 1



VRL500 Oxane® 1300, 10ml Syringe   = 1

VRL600 Oxane® 5700, 10ml Syringe   = 1



VRL700 Oxane® HD, 10ml Syringe   = 1



CX5710 Viscous Fluid Extraction Kits   = 6



CX5720 Viscous Fluid Injection Kits for Oxane range   = 6



3090 Silicone Oil collector chamber   = 6



NEW!

BL7620 High Flow VFI Cannula 20G   = 6



BL7623 High Flow VFI Cannula 23G   = 6



BL7625 High Flow VFI Cannula 25G   = 6

