

INSTRUCTIONS FOR USE

14. 230 V, 50 Hz elektros tinklas

Transformer specifications

| | |
|-------------------------|---|
| Table top standard size | 380 x 500 mm |
| Power supply voltage | 100V / 120V / 230V / 240V CA ±10% |
| Fuses: 5x20 mm | 100-120V CA --- 1 A 230-240V CA --- 0.5A |
| Mains frequency | 50-60Hz |
| Maximum power absorbed | 25 VA |

Caratteristiche Switching

| | | |
|------------------------|-----------------------------|-----------------------------|
| Table Switching | 50,8 x 76,2 mm | 50,8 x 76,2 mm |
| Power supply voltage | 100V/120V/230V/240V CA ±10% | 100V/120V/230V/240V CA ±10% |
| FUSES: 5x20 mm | 1,6 A | 1,6 A |
| Maximum power absorbed | 25 VA | 25 VA |

Other features

| | |
|---------------------|---------------------|
| Lamp size | 299x313x(644±15) mm |
| Lamp weight | 8,7 Kg |
| Digital lamp weight | 9,4 Kg |

6.5. STEREOSCOPIC MICROSCOPE

| Microscope Model | 2x | | | | 3x | | | 5x | | | | Zoom | |
|---------------------------|---|-----|-----|------|--|-----|-----|---------------------------|-----|-----|-----|------------|---------------------|
| Type | Greenough convergent binocular microscope | | | | 6° Galileian convergent binocular microscope | | | | | | | | Input change system |
| | Two-position rotating lens | | | | Three-position rotary drum | | | Five-position rotary drum | | | | Continuous | |
| Eyepieces | 10x st. | | 16x | | 12.5x | | | 12.5x | | | | 12.5x | |
| Ingrandimenti totali | 10x | 16x | 16x | 25x | 10x | 16x | 25x | 6x | 10x | 16x | 25x | 40x | 6x-33x |
| Actual visual field: (mm) | 18.5 | 12 | 16 | 10.5 | 24 | 14 | 8 | 37 | 24 | 14 | 8 | 5.2 | 32 - 6,2 |
| Interpupillary distance | 51.5-87 mm | | | | 48.5-80 mm | | | | | | | | |

3.
37 mm, kai didinimas 6x
24 mm, kai didinimas 10x
14 mm, kai didinimas 16x
8 mm, kai didinimas 25x
5.2 mm, kai didinimas 40x