



Aesculap Surgical Instruments

Aesculap chirurginiai instrumentai


General Catalog
Pagrindinis katalogas



Think Surgical Instruments. Think Aesculap.




LEXER
BC283R
 schmal
 narrow
 165 mm, 6 1/2"



MAYO-LEXER
BC284R
 165 mm, 6 1/2"

8.1




BC256R-BC257R
 115 mm, 4 1/2"

Scheren – Scissors
DUROTIP® Feine Scheren
 DuroTIP® Delicate Scissors



BC



DUROTIP®-Scheren
 mit eingesetzten Hartmetall-Schneiden


- 3 Jahre Garantie!
- Hervorragende Schnitthaltigkeit
- Lange Lebensdauer
- Vergoldete Ring

DUROTIP® Scissors
 with inserted carbide cutting edges

- 3 years guarantee!
- Excellent cutting accuracy
- Long working life
- Gold-plated rings

DUROTIP žirkės su karbido
 įdėklais kerpančiose dalyse

Puikus kirpimo tikslumas
 Ilgas eksploatacinis laikotarpis
 Paausuoti žiedai



BC210R

BC211R

BC210R-BC211R

Fadenscheren
 Ligature scissors
 110 mm, 4 3/8"




1 Schneide gezahnt
 1 blade serrated

GOLDMAN-FOX

DO250R

130 mm, 5 1/8"



DO208R

DO209R

DO208R-DO209R

120 mm, 4 3/4"



■ We can only guarantee reproducible Aesculap quality over decades thanks to the technical skill, experience and care of Aesculap's surgical instrument makers. Expertise that is responsibly and conscientiously shared with others.



Tungsten Carbide Instruments

For our instruments with tungsten carbide inserts we use only high quality nickel-bound tungsten carbide tips, well known for their high corrosion resistance and low wear.

Instrumentams su volframo karbido įdėklais mes naudojame tik aukštos kokybės nikeliu sutvirtintą volframo karbidą, kuris gerai žinomas dideliu atsparumu korozijai ir dėvėjimuisi



Hemostatic Forceps – Olive in the Jaw

The olive forms the transition from the working end or jaw to the "joint" of the instrument. This special design makes the jaw more elastic and the instrument less liable to break.




How do we ensure that every Aesculap instrument conforms exactly to the shape and dimensions of the master pattern? We inspect every single one!



| BC313R | BC314R | BC315R | BC316R | BC317R |
|--------|--------|--------|--------|--------|
| 130 mm | 145 mm | 150 mm | 165 mm | 175 mm |
| 5 1/8" | 5 3/4" | 6" | 6 1/2" | 7" |

BC313R-BC317R
Standard-Modelle
Standard patterns

8.4



| BC411R | BC413R | BC414R | BC415R | BC416R | BC417R |
|--------|--------|--------|--------|--------|--------|
| 115 mm | 130 mm | 145 mm | 150 mm | 165 mm | 175 mm |
| 4 1/2" | 5 1/8" | 5 3/4" | 6" | 6 1/2" | 7" |

COOPER
BC411R-BC417R
Standard-Modelle
Standard patterns

8.19 8.2 8.3



Scheren – Scissors
 Präparierscheren
 Dissecting Scissors
 Preparacines zirkles

BC



BABY-METZENBAUM
 BC600R-BC603R
 145 mm, 5 3/4"



METZENBAUM
 BC680R-BC683R




| | | |
|--------|--------|--------|
| BC601R | BC602R | BC604R |
| 145 mm | 180 mm | 200 mm |
| 5 3/4" | 7" | 8" |

| | | |
|--------|--------|--------|
| BC605R | BC606R | BC607R |
| 145 mm | 180 mm | 200 mm |
| 5 3/4" | 7" | 8" |

METZENBAUM
 BC601R-BC607R

8.5



| BC576R | BC578R |
|--------|--------|
| 140 mm | 170 mm |
| 5 1/2" | 6 3/4" |

| BC577R | BC579R |
|--------|--------|
| 140 mm | 170 mm |
| 5 1/2" | 6 3/4" |


MAYO-STILLE

BC576R-BC579R

mit abgechrägten Schneiden

with chamfered blades *su nusklembtomis briaunomis*

8.6



| BC562R | BC564R |
|--------|--------|
| 190 mm | 215 mm |
| 7 1/2" | 8 1/2" |

| BC563R | BC565R |
|--------|--------|
| 190 mm | 215 mm |
| 7 1/2" | 8 1/2" |

MAYO

BC562R-BC565R



mit abgechrägten Schneiden

with chamfered blades *su nusklembtomis briaunomis*




Scheren – Scissors
Nagelscheren
Nail Scissors


BC



HF400R
Nagelschere
Nail scissors
90 mm, 3 1/2"



HF409R
Nagelschere,
extra feines Modell
Nail scissors,
especially delicate pattern
90 mm, 3 1/2"



HF032R
Baby-Nagelschere
Baby nail scissors
85 mm, 3 1/8"



HF014R–HF021R
Nagelscheren
Nail scissors





HF028R–HF031R
Nagelscheren
Nail scissors

8.7

Scheren – Scissors
 Feine Scheren
 Delicate Scissors

BC



BC100R

BC101R

BC100R-BC101R

90 mm, 3 1/2"



BC056R

BC057R

BC056R-BC057R

90 mm, 3 3/4"

8.8




BC060R

BC061R

BC060R-BC061R

105 mm, 4 1/8"



BC064R

BC065R

BC064R-BC065R

120 mm, 4 3/4"



Pinzetten – Dissecting- and Tissue Forceps

Chirurgische Pinzetten

Tissue Forceps

Chirurginiai pincetai

BD



| BD552R | BD553R | BD555R | BD557R | BD559R | BD560R | BD561R | BD562R | BD563R | BD564R |
|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 105 mm | 115 mm | 130 mm | 145 mm | 160 mm | 180 mm | 200 mm | 250 mm | 300 mm | 350 mm |
| 4 1/8" | 4 1/2" | 5 1/8" | 5 3/4" | 6 1/4" | 7" | 8" | 10" | 12" | 14" |


8.9 8.10

BD552R–BD564R

Standard-Modelle

Standard patterns

Standartinis modelis



| BD533R | BD535R | BD537R | BD527R | BD539R | BD540R | BD541R | BD542R | BD543R |
|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 115 mm | 130 mm | 145 mm | 145 mm | 160 mm | 180 mm | 200 mm | 250 mm | 300 mm |
| 4 1/2" | 5 1/8" | 5 3/4" | 5 3/4" | 6 1/4" | 7" | 8" | 10" | 12" |

BD527R–BD543R

mittelbreit

medium

Vidutinio stambumo



BC313R-BC317R
Standard-Modelle
Standard patterns

8.11 8.12




COOPER
BC411R-BC417R
Standard-Modelle
Standard patterns



Scheren – Scissors
 Präparierscheren
 Dissecting Scissors
 Preparacines zirkles

BC



BABY-METZENBAUM
 BC600R-BC603R
 145 mm, 5 3/4"




METZENBAUM
 BC680R-BC683R



METZENBAUM
 BC601R-BC607R

8.13



8.14


Scheren – Scissors

Präparierscheren

Dissecting Scissors

Preparacīnēs žirkļs

BC



| BC584R | BC585R | 8.15 |
|--------|--------|------|
| 165 mm | 200 mm | |
| 6 1/2" | 8" | |



| BC587R | BC588R | BC589R |
|--------|--------|--------|
| 165 mm | 200 mm | 210 mm |
| 6 1/2" | 8" | 8 1/4" |

MAYO

BC584R-BC589R

Scheren – Scissors
 Chirurgische Scheren
 Surgical Scissors

BC



| BC341R | BC343R | BC344R | BC345R | BC346R |
|--------|--------|--------|--------|--------|
| 115 mm | 130 mm | 145 mm | 150 mm | 165 mm |
| 4 1/2" | 5 1/8" | 5 3/4" | 6" | 6 1/2" |

BC341R-BC346R

Standard-Modelle
 Standard patterns

8.18 8.17



| BC441R | BC443R | BC444R | BC445R |
|--------|--------|--------|--------|
| 115 mm | 130 mm | 145 mm | 150 mm |
| 4 1/2" | 5 1/8" | 5 3/4" | 6" |


BC441R-BC445R

Standard-Modelle
 Standard patterns



| BC313R | BC314R | BC315R | BC316R | BC317R |
|--------|--------|--------|--------|--------|
| 130 mm | 145 mm | 150 mm | 165 mm | 175 mm |
| 5 1/8" | 5 3/4" | 6" | 6 1/2" | 7" |

BC313R-BC317R
Standard-Modelle
Standard patterns



| BC411R | BC413R | BC414R | BC415R | BC416R | BC417R |
|--------|--------|--------|--------|--------|--------|
| 115 mm | 130 mm | 145 mm | 150 mm | 165 mm | 175 mm |
| 4 1/2" | 5 1/8" | 5 3/4" | 6" | 6 1/2" | 7" |


COOPER
BC411R-BC417R
Standard-Modelle
Standard patterns

8.19



Scheren – Scissors
 Präparierscheren
 Dissecting Scissors
 Preparacinės žirklys

BC



| BC576R | BC578R |
|--------|--------|
| 140 mm | 170 mm |
| 5 1/2" | 6 3/4" |

| BC577R | BC579R |
|--------|--------|
| 140 mm | 170 mm |
| 5 1/2" | 6 3/4" |

MAYO-STILLE


BC576R-BC579R

mit abgescrãgten Schneiden

with chamfered blades su nusklembtomis briaunomis

8.20

8.21



| BC562R | BC564R |
|--------|--------|
| 190 mm | 215 mm |
| 7 1/2" | 8 1/2" |

| BC563R | BC565R |
|--------|--------|
| 190 mm | 215 mm |
| 7 1/2" | 8 1/2" |

MAYO

BC562R-BC565R

mit abgescrãgten Schneiden

with chamfered blades su nusklembtomis briaunomis






Scheren – Scissors

Feine Scheren - große Ringform
Delicate Scissors - Large Ring Design

BC




BC004R-BC005R
115 mm, 4 1/2"



8.22



JOSEPH
BC006R-BC007R
150 mm, 6"



MILLIMETERS