

REGULAR INDICATING LAEBL SEAL PLUS HUB HOT SEAL

Packing: 100 pcs / box, 50 boxes / case

Ref.	Description
1A2004-B21	27G (0.40) x 21 mm
1A2004-B25	27G (0.40) x 25 mm
1A2004-B34	27G (0.40) x 34 mm
1A2005-B12	30G (0.30) x 12 mm
1A2005-B16	30G (0.30) x 16 mm
1A2005-B21	30G (0.30) x 21 mm
1A2005-B25	30G (0.30) x 25 mm
1A2005-B38	30G (0.30) x 38 mm

Other sizes are available upon requests.



HUB HOT SEAL ONLY 2.1.P.D. ir 2.2. P.D. ECONOMIC TYPE

Packing: 100 pcs / box, 50 boxes / case

Ref.	Description
1A2004-A21	27G (0.40) x 21 mm
1A2004-A25	27G (0.40) x 25 mm
1A2004-A34	27G (0.40) x 34 mm
1A2005-A12	30G (0.30) x 12 mm
1A2005-A16	30G (0.30) x 16 mm
1A2005-A21	30G (0.30) x 21 mm
1A2005-A25	30G (0.30) x 25 mm
1A2005-A38	30G (0.30) x 38 mm



ASPIRATING SYRINGE BRASS MADE

Packing: 1 pc / bag, 200 bags / case

Ref.	Description
1A2010	1.8ml, Metric thread & Imperial screw thread, cuspidal head
1A2011	1.8ml, Metric thread & Imperial screw thread, hook head

TEMPORARY CROWN & BRIDGE MATERIAL

Automix Cartridge System

Master-Dent® Automix Temporary Crown & Bridge Material is a quick and easy-to-use automix system that prevents cross contamination. Low brittleness, strong when bending/removing, so it will not break. With no oxygen inhibition layer, wiping residual layer is not needed. Less exothermic reaction eliminates heat sensitivity. No need to polish with this smooth finish!



REF		
60-203	Shade A1	1 x 50 ml Cartridge w/ Tips
60-200	Shade A2	1 x 50 ml Cartridge w/ Tips
60-201	Shade A3	1 x 50 ml Cartridge w/ Tips
60-202	Shade B1	1 x 50 ml Cartridge w/ Tips
20-515	HP Mixing Gun	
20-555	HP Mixing Tips (Blue/Clear)	(25/pk)

TOPICAL ANESTHETIC

Master-Dent® Topical Anesthetic is a rapid-acting, long-lasting gel with no systemic absorption. The benzocaine provides temporary relief of pain and discomfort during procedures, including local anesthetic injections, periodontal curettage, impression taking, scaling, intra-oral radiographs, root planing and prophylaxis. The gel comforts patients with ulcers, wounds or other minor mouth irritations. Master-Dent® Topical Anesthetic has no bitter aftertaste and is now available in four great-tasting flavors.



REF	
47-149	Bubble Gum, 1 Oz. (30ml)
47-150	Mint, 1 Oz. (30ml)
47-151	Cherry, 1 Oz. (30ml)
47-152	Strawberry, 1 Oz. (30ml)

ULTIMATE™ BOND

Light Cure Adhesive, One Component

Master-Dent® Ultimate™ Bond creates a chemical bond between dentin/enamel and all restorative materials. A single, light cured component that does not require mixing. Easy to use. Microleakage is minimized. Secondary caries are reduced. Excellent compatibility. Relieves hypersensitivity. Provides root surface desensitization. Combines mono and polymethacrylate resins with HEMA (hydroxyethylmethacrylate) to produce desired hydrophilic properties.



REF	
19-215	7ml Bottle

5.1., 5.2. P.D.

Absorbent Paper Points

.02/.04/.06/.08/ Protaper(F1~F5) Paper Points




Advantage

- Highly Absorbent
- Sterilized Individual Packages
- Variety Sizes and Package Available
- 5yrs Shelf life

Length Marked

Accurate insertion of the working length



Size & Specification

- .02 Taper: 15~140
- .04, .06, .08 Taper: 15~80
- .08 Taper: 15~40
- Non Standard
 - Non-Standard: XX-Fine, X-Fine, Fine, Medium, Coarse, X-Coarse
 - Protaper Size(F1~F5)

Specification

Size	Taper(points)				Packing Type
	02(200)	04(60)	06(60)	08(60)	
10					
15					
20					
25					
30					
35					
40					
45					
50					
55					
60					
70					
80					
90					
100					
110					
120					
130					
140					

※ Please check Page12

6.1., 6.2., 6.3. P.D.

Gutta Percha Points

.02/.04/.06/.08/ Protaper(F1~F5) Gutta Percha Points



Advantage

- Long Shelf Life 4yrs
- Stiff and Flexible
- Excellent Radiopacity
- Various Size, Taper, Package

Length Marked



Size & Specification

- 02 Taper: 10~140
- 04, 06 Taper: 15~80
- 08 Taper: 15~40
- Non Standard
 - Non-Standard: XF, FF, MF, F, FM, M, ML, L, XL
 - Japanese Type: JF-S, JF-M, JF-L, JS, JM, JL, JL-S, JL-M, JL-L
 - De Trey Type: XX-F, X-F, Fine, Medium
 - Protaper Size(F1~F5)

Specification

Size	Taper(points)				Packing Type
	02(120)	04(60)	06(60)	08(60)	
10					
15					
20					
25					
30					
35					
40					
45					
50					
55					
60					
70					
80					
90					
100					
110					
120					
130					
140					

※ Please check Page12

NiTi Protaper
(Hand use)



Package: 6pcs/box

SIZE	SX	S1	S2	F1	F2	F3	SX-F3
Ref.#	NSX19	NS121	NS221	NF121	NF221	NF321	NSXF321
SIZE		S1	S2	F1	F2	F3	SX-F3
Ref.#		NS125	NS225	NF125	NF225	NF325	NSXF325

NiTi Protaper
(Engine Use)



Rec.R.P.M.: 150-350
SX: 19MM

SIZE	SX	S1	S2	F1	F2	F3	SX-F3
Ref.#	ENSX19	ENS121	ENS221	ENF121	ENF221	ENF321	ENSXF321
SIZE		S1	S2	F1	F2	F3	SX-F3
Ref.#		ENS125	ENS225	ENF125	ENF225	ENF325	ENSXF325

Paste Carriers
(Stainless Steel)



Rec.R.P.M.: 250-500
Package: 4cs/box

- Made of imported Stainless steel(Medical grade)
- Precision processing technology, good flexibility, high resistance to broken.

SIZE	25	30	35	40	25-40
Ref.#	PC2521	PC3021	PC3521	PC4021	PC254021
SIZE	25	30	35	40	25-40
Ref.#	PC2525	PC3025	PC3525	PC4025	PC254025
SIZE	25	30	35	40	25-40
Ref.#	PC2529	PC3029	PC3529	PC4029	PC254029

Gates Drills

Rec.R.P.M.: 700-1100
Package: 6pcs/box



- High quality Stainless steel(Medical grade)
- For the provision of a root canal orifice, remove the root canal filling material and preparation.

SIZE	1	2	3	4	5	6	1-6
Ref.#	G0128	G0228	G0328	G0428	G0528	G0628	G1628
SIZE	1	2	3	4	5	6	1-6
Ref.#	G0132	G0232	G0332	G0432	G0532	G0632	G1632

Peeso Reamers

Rec.R.P.M.: 700-1100
Package: 6pcs/box




- High quality Stainless Steel(Medical grade)
- Used for crown one-third of root canal preparation, remove the root canal filling material and preparation.

SIZE	1	2	3	4	5	6	1-6
Ref.#	P0128	P0228	P0328	P0428	P0528	P0628	P1628
SIZE	1	2	3	4	5	6	1-6
Ref.#	P0132	P0232	P0332	P0432	G0532	P0632	P1632

7.1. P.D. ir 7.2. P.D.

Endo File Foam Inserts

Autoclavable



Ref.#	Description	Package
ED002-1	50*50*20(mm) White	48pcs/bag
ED002-2	100*69*13(mm) Blue	100pcs/bag

Ref.#	Description	Package
ED004-1	1big ring + 1small ring + 1foam foam only	1set/bag
ED004-2		50pcs/bag

Endo Ruler

autoclavable



Ref.#	Description	Package
ED003-1	Autoclavable, stainless steel	1pcs/bag

Ref.#	Description	Package
ED003-2	Aluminum, R ng	1pcs/box
ED003-3	Aluminum, Straight	1pcs/box

Endo File holder

Autoclavable



35-holes

Ref.#	Description	Package
ED003-4	aluminum, 118*83*50mm	1pc/box

Ref.#	Description	Package
ED003-5	aluminum, 105*47*55mm	1pc/box

7.2. P.D.

7.1. P.D.

8. P.D.

Root Canal Sealer / ADSEAL

• 15

ADSEAL™

Resin Based Sealer



Applications

- Permanent root canal sealing
- Thermal obturation

Composition

- Base-Epoxy oligomer resin, Ethylene glycol salicylate, Bismuth subcarbonate
- Catalyst-Poly butanediol aminobenzoate, Calcium Phosphate, Bismuth subcarbonate

Advantages


- Excellent biocompatibility
- Easy to mix paste-paste
- Hermetic sealing ability
- Nonstaining to the teeth
- Insoluble in tissue fluids
- Good radiopacity

Components


- 13.5g dual syringe
 - Base 9g (Epoxy resin, Calcium phosphate)
 - Catalyst 4.5g (Amines, Bismuth subcarbonate)
- Spatula
- Mixing Plate
- Dimension : 55x180x40(mm)
- Weight : 98g

Physical and Chemical Properties


• Flow




• Solubility



• Radiopacity



• Material Properties



- Working Time : 35 Minutes at 23°C (73°F)
- Setting Time: 45 Minutes at 37°C (99°F)
- Film Thickness: 3.3 μm

Shelf life

- Two years from the production date.

MD-ChelCream

Viscous Chelator / EDTA Gel



Applications

- Root canal preparation with files and reamers
- Presence of calcifications in the canal
- Narrow canals in posterior teeth
- Smear layer removal

Composition

- EDTA, PEG

Advantages

- Removal of inorganic substances by chelation (EDTA-Ca²⁺)
- Easier filing with lubricating action
- Produces the cleanest canal wall after smear layer removal

Components

- MD-ChelCream 7g(5.4ml)×2
- Disposable tips
- Dimension : 186×45×30(mm)
- Weight : 59g

Shelf life

- Two years from the production date.

MD-Cleanser™

Smear Layer Remover / EDTA Solution



Components

- MD-Cleanser : 100 ml
- 1 cc syringe
- Disposable tips
- Dimension : 106 x 50 x 53 mm
- Weight : 138g

Clinical Case

Applications

- Presence of calcification in the canal
- Narrow canals on posterior teeth
- Smear layer removal

Composition


- EDTA, Water, Ammonia water

Advantages


- Removal of inorganic substances by chelation (EDTA-Ca²⁺)
- Easier instrumentation by softening dentin canal wall
- Produces the cleanest dentinal walls after smear layer removal




Application Method



After root canal preparation



After MD-Cleanser application (x1000)



After MD-Cleanser application (x7000)

Shelf life

- Two years from the production date.

CORE BUILD-UP (PASTE-PASTE)

Master-Dent® Core Build-up is a self cure product designed for core build-ups and cementation of posts and pins. Creamy smooth formulation, radiopaque and non-slumping. Exclusive matched-viscosity pastes ensure easy, reliable mixing even if proportions are varied to adjust setting time. Sets hard in minutes so finishing can take place without delay. High compressive strength with easy application.



REF		
19-000-B	Blue	14g Base Paste / 14g Catalyst Paste
19-000-N	Natural	14g Base Paste / 14g Catalyst Paste
19-000-W	White	14g Base Paste / 14g Catalyst Paste

CALCIUM HYDROXIDE PASTE With Barium Sulfate, Root Canal Filling Material

Master-Dent® Calcium Hydroxide Paste with Barium Sulfate provides the full benefits of calcium hydroxide in root canal therapy. It is used as a root canal sealer for the temporary filling of root canals after conventional endodontic treatment between office visits to disinfect the canal. It is also indicated for direct pulp capping, pulpotomy, pulpectomy, cavity lining and indirect pulp caps to stimulate dentin bridging. Water soluble for easy removal, it has an excellent antibacterial effect. Barium Sulfate provides high radiopacity. This material has a high pH (>12.0) and its convenient premixed syringe with disposable tips makes delivery and application simple.



REF	
11-280	4 x 1g Syringes w/ Tips



CALCIUM HYDROXIDE POWDER

Master-Dent® Calcium Hydroxide Powder is ideal for those who prefer to formulate their own direct and indirect pulp capping paste or sub-bases for use under restorations.

REF	
11-276	1 Lb. (454g)




EDTA SOLUTION, 17%

Master-Dent® EDTA Solution irrigates root canals. Also assists in softening root canal walls for ease of instrumentation.

REF	
12-750	2 Oz. bottle (60ml)
12-800	8 Oz. bottle (240ml)

Endo File Holder



- Autoclavable
- Material: PP

Ref.#	Description	Package
ED004A-1	10X2.4X3.2cm,millimeter scale	1pcs/box

Endo File Holder with Sponge



- Autoclavable
- Material: Polyamide

Ref.#	Description	Package
ED004A-3	45mm	1pc/box

Gutta Percha Points



ISO.02

Ref.#	Description	Package
ED005-115	15#	120pcs/box
ED005-120	20#	120pcs/box
ED005-125	25#	120pcs/box
ED005-130	30#	120pcs/box
ED005-135	35#	120pcs/box
ED005-140	40#	120pcs/box
ED005-1	15-40#	120pcs/box
ED005-2	45-80#	120pcs/box
ED005-3	90-140#	60pcs/box

ISO.02

- Material:aluminium
- 6 Holes
- Size:115*80*50mm

GP Holder



13.1., 13.2. P.D.

Absorbent Paper Points




ISO.02


Ref.#	Description	Package
ED006-115	15#	200pcs/box
ED006-120	20#	200pcs/box
ED006-125	25#	200pcs/box
ED006-130	30#	200pcs/box
ED006-135	35#	200pcs/box
ED006-140	40#	200pcs/box
ED006-1	15-40#	200pcs/box
ED006-2	45-80#	200pcs/box
ED006-3	90-140#	100pcs/box




Notes: All dimension to be measured to an accuracy 0.05mm.A tolerance of ± 0.05 mm applies to D3,D16.Taper proportion is 0.02mm per 1mm at uniform taper.Overall length is not less than 28mm ± 0.01 mm unless otherwise specified.



SHAPE



Non-Standard Points



Ref.#	Description	Dia.D3 3mm From tip(mm)	Dia.D16 16mm From tip(mm)
ED008-1	XF	0.20	0.45
ED008-2	FF	0.24	0.56
ED008-3	MF	0.27	0.68
ED008-4	F	0.30	0.80
ED008-5	FM	0.35	0.88
ED008-6	M	0.40	1.10
ED008-7	ML	0.43	1.25
ED008-8	L	0.49	1.55
ED008-9	XL	0.52	1.60



ISO.02

Ref.#	Description	Package
ED006M-115	15#	200pcs/box
ED006M-120	20#	200pcs/box
ED006M-125	25#	200pcs/box
ED006M-130	30#	200pcs/box
ED006M-135	35#	200pcs/box
ED006M-140	40#	200pcs/box
ED006M-1	15-40#	200pcs/box
ED006M-2	45-80#	200pcs/box
ED006M-3	90-140#	100pcs/box



27.1. P.D.

Syringe & Dental Needle

Syringe & Dental Needle



Aspirating Syringe (Brass Made)

SN002

- Material: Brass with chromium coated
- Capacity: 1.8ml



Ref. #	Description	Package
SN002-1	USA/Metric thread,triangle head	1pcs/box
SN002-2	USA/Metric thread,hook head	1pcs/box

12CC Utility Syringes

SN003

- Designed for surgical site and post-surgical home irrigation, similar to Monoject #412 syringe
- Can be used with silicone and rubber-based impression materials
- Package:50pcs/bag,24bag/ctn



- Curved Straight



Disposable Dental Needle(Sterilized)

SN001

- Size: 27G X 15/25/38mm
30G X 15/25/38mm


Ref. #	Description	Package
SN001-1	UK & Metric General, 27G (0.40)	100pcs/box,50boxes/ctn
SN001-2	UK & Metric General, 30G (0.30)	100pcs/box,50boxes/ctn

Luer Lock & Luer Slip Irrigation Syringe

Ref. #	Description	Package
SN004-1	1.2cc	100pcs/bag,50bags/ctn
SN004-2	3cc	100pcs/bag,38bags/ctn
SN004-3	5cc	100pcs/bag,30bags/ctn
SN004-4	10cc	100pcs/bag,15bags/ctn
SN004-5	12cc	100pcs/bag,15bags/ctn



27.2. P.D.



● Endo Irrigation Needle of **MEDIC** is the very useful products which can irrigate perfectly root canals.

● Endo Irrigation Needle has closed end tips and side hole opening. The needle is able to clean Probe canal.

Products	description	size	color
Irri - 1	close end tips, one side hole opening	30G 0.30 x 25mm	yellow
Irri - 2	close end tips, two side hole opening	27G 0.40 x 25mm	grey
Irri - 3	opening tips, half cutting tips	27G 0.40 x 25mm	grey

- Sterilized by EO gas.
- Unit package: plastic container.



Polishing Strips **28. P.D.**



No 1.050

Strips for Gross Reduction
gross reduction / contouring
25 pcs.



No 1.051

Strips for Contouring
contouring / finishing
25 pcs.



No 1.052

Strips for Polishing
finishing / polishing
25 pcs.



No 1.055

Strips. Universal Kit
75 pcs.:
No 1.050 25 pcs.
No 1.051 25 pcs.
No 1.052 25 pcs.

Finishing & Polishing Kits



No 1.020

Finishing & Polishing Kit

stem polishing discs 48 pcs.:
No 1.631 12 pcs.
No 1.632 12 pcs.
No 1.633 12 pcs.
No 1.634 12 pcs.

polishing strips 75 pcs.:
No 1.050 25 pcs.
No 1.051 25 pcs.
No 1.052 25 pcs.

mandrel No 1.111 1 pc.



No 1.021

Finishing & Polishing Kit


polishing discs 48 pcs.:
No 1.731(14) 12 pcs.
No 1.732(14) 12 pcs.
No 1.733(14) 12 pcs.
No 1.734(14) 12 pcs.

polishing strips 75 pcs.:
No 1.050 25 pcs.
No 1.051 25 pcs.
No 1.052 25 pcs.


mandrel No 1.121 1 pc.

Three types:


- for gross reduction / contouring



- for contouring / finishing



- for finishing / polishing



Catalogue was prepared using materials of
K. Aleinikov (Moscow, Russia),
K. Asimakopolus (Athens, Greece),
V. Bezruchko (Dmitrov, Russia),
A. Volom (Budapest, Hungary),
V. Durov (Tver', Russia),
Yu. Nekipelov (Vologda, Russia),
A. Nikolaev (Smolensk, Russia),
D. Nikolaev (Bryansk, Russia),
A. Salova (St. Peterburg, Russia),
I. Makeeva (Moscow, Russia),
A. Sidorov (Moscow, Russia),
V. Tishchenko (Krasnodar, Russia)

**CUSTOM
MADE**

All the products can also be assorted in a different than standard way.

Complete kit of implantology wear



TA021-4A

- Each Kit contains:
- 1 surgical gown, SMS blue, 55g/m², impermeable with knit cuffs
 - 1 assistant gown, polypropylene, 50g/m², impermeable with cotton cuffs
 - 50g/m², impermeable with knit cuffs
 - 2 pair of non-slippery overshoes
 - 1 pair of patient overshoes
 - 2 hand towels
 - 3 tie on masks
 - 2 surgical hats
 - 3 bouffant caps

TA021-4B

- Each Kit contains:
- 2 drapes 75 x 75 cm
 - 1 instrument (bridge) table cover 90 x 120cm, 43g/m², impermeable
 - 1 instrument table cover, 100 x 180 cm, 43g/m², impermeable for work plan.
 - 1 non woven "U" drape with adhesive 100 x 150 cm
 - 1 non woven surgical drape with adhesive 70 x 100 cm
 - 3 plastic sheets 8 x 120 cm with adhesive
 - 2 adhesive films 5 x 80 cm
 - 2 elastic wrist protectors
 - 1 plastic bag 30 x 80 cm
 - 1 protection for panel control with adhesive
 - 1 elastic cap radio cover

Dental Surgical kit



TA021-4C

- Contains:
- 1 patient's drape 75x90cm
 - 1 transparent hose cover: 10x200cm
 - 1 transparent hose cover: 10x180cm
 - 1 reflector drape 15x115cm
 - 1 table drape DL 60x75cm
 - 1 bundle cover 50x50cm



TA021-4D

- Contains:
- 3 drapes 60*80cm
 - 2 instrument (bridge) table cover 120*120cm, 50g/m², impermeable
 - 1 non-woven fenestrated drape 160*170cm
 - 1 suture bag 30*16cm
 - 4 double side tape 2*20cm
 - 2 surgical hats
 - 2 tie-on face masks
 - 1 surgical gown, SMS blue, 50g/m², impermeable with cotton cuffs
 - 2 plastic sheets 8*120cm with adhesive

TA021-4E

- Contains:
- 2 Surgical Gown Standard-Large with in overwrap
 - 2 Absorbent Towel
 - 1 Drape (175*150 cm), circular fenestration (7) with adhesive (10*10cm)
 - 2 sleeve cover (7.5*120cm)
 - 1 Back table cover (80*100cm)

Propphy Cup

100pcs/box



TA023-1A

TA023-1B

TA023-1C

Propphy Angle Plus

144pcs/box



TA023-2

Propphy Paste Ring

TA023-1F

- Fit for all 2g propphy paste cups
- Package: 200pcs/bag

Propphy Brush

100pcs/bag



Prophylaxis Powder



TA025-1

- Mesh 200, Mint-taste
- 40g/bag

TA025-2

- Mesh: 180-200, 130g/bottle
- Package: 20bottle/box, 600x/ctn



TA025-2A Original
TA025-2B Orange
TA025-2A Mint

ZINC OXIDE POWDER USP GRADE

Master-Dent® USP Grade Zinc Oxide Powder can be used to formulate temporary filling materials, cements and root canal sealers. ZOE bases are compatible with hard and soft tissue and have a local anesthetic effect.



REF	
12-700	1 Lb. (454g) Can

BLUE GEL ETCH

Master-Dent® Blue Gel Etch offers 37% phosphoric acid gel in various packaging options. Vivid blue for visual confirmation of application and rinsing. Luer lock tips are angled for easy access, offering better control.

REF	
18-120	12g luer lock syringe w/ 20 black luer lock needle tips
18-125	4 x 1.2ml syringes w/ 20 blue luer lock needle tips
18-130	60g syringe, 5 empty 3ml syringes, 50 blue luer lock needle tips
ACCESSORIES	
11-258	Blue Luer lock needle tip (20 pkg)
11-259	Blue Luer lock needle tip (100 pkg)
18-106	Empty 3ml syringes (50 per pkg)



GREEN GEL ETCH

Master-Dent® Green Gel Etch offers a sophisticated 37% phosphoric acid formulation containing wetting agents, easily soluble thickeners and FDA certified dye. It is easy to apply and to control etched areas, as well as easy to rinse off. Dispensing system is precise and simple, thus preventing the possibility of contamination and waste. The economy size 12 gram package is priced to allow substantial savings and the 60 gram saves even more.



REF	
18-100	12g Syringe with 25 clear straight plastic tips
18-199	Straight plastic tips for 18-100 (50pkg)
18-200	Angled plastic tips for 18-100 (50pkg)
18-114	12g Syringe (luer-lock) with 20 blue luer-lock tips
18-105	60g Syringe, 5 empty 3ml syringes, 50 black luer lock needle tips

PRODUCT CATEGORIES



Dentistry

- Endodontics
- Prosthetics
- Surgery and periodontology
- Whitening of teeth
- Prophylaxis and cosmetics
- Dermatology
- Musculoskeletal system
- Cardiovascular system
- Vitamin D and generics
- Veterinary products

ALUSTIN



Composition

Aluminum chloride hexahydrate, mint essence, purified water, preservative (quinosol).
1 g of the product contains 0,2 g of aluminum chloride hexahydrate.

Application

Dressing for perigingival bleeding control that is a result of intraoral procedures.

Action description

Alustin® applied on the bleeding area creates superficial protective layer. Bleeding stops after several minutes, which enables single session treatment without need for temporary dressing.

Product handling

Inhibit the bleeding using a cotton pad moistened with Alustin®. After hemorrhage arrested, rinse harmed area with purified water. In case of bleeding repeat procedure.
Do not use the product in case of hypersensitivity against any of constituents.

Available packaging

Plastic bottle sealed with a cap with dropper containing 10 g of the product.

Storage conditions

Keep containers tightly closed. **Keep out of the reach of children.**
Store in temperature 5°C to 25°C.

Sectional Contoured Metal Matrices

All matrices are available of 4 types

- 0.050 mm, hard
- 0.035 mm, hard
- 0.050 mm, soft
- 0.035 mm, soft

Specify please the necessary type, e.g. № 1.098(0.050 mm, soft).

There are two variants of kiting of matrices of one shape - per 10 or per 50 pcs.

When ordering kit with matrices of one shape specify please the necessary kiting variant and type, e.g. № 1.0971(10) (0.035 mm, hard) or № 1.0976(50) (0.050 mm, soft).



Sectional Contoured Matrices. Small with Ledge

№ 1.0971(10) 10 pcs.
№ 1.0971(50) 50 pcs.



Sectional Contoured Matrices. Small

№ 1.0972(10) 10 pcs.
№ 1.0972(50) 50 pcs.



Sectional Contoured Matrices. Large

№ 1.0973(10) 10 pcs.
№ 1.0973(50) 50 pcs.



Sectional Contoured Matrices. Large with Ledge

№ 1.0974(10) 10 pcs.
№ 1.0974(50) 50 pcs.



Sectional Contoured Matrices. Medium with Ledge

№ 1.0975(10) 10 pcs.
№ 1.0975(50) 50 pcs.



Sectional Contoured Matrices. Medium

№ 1.0976(10) 10 pcs.
№ 1.0976(50) 50 pcs.



№ 1.1972 Kit of Small Matrices 4 types

30 pcs.:
 № 1.0972(0.050 mm, hard) 10 pcs.
 № 1.0972(0.035 mm, hard) 10 pcs.
 № 1.0972(0.050 mm, soft) 5 pcs.
 № 1.0972(0.035 mm, soft) 5 pcs.



№ 1.1976 Kit of Medium Matrices 4 types

30 pcs.:
 № 1.0976(0.050 mm, hard) 10 pcs.
 № 1.0976(0.035 mm, hard) 10 pcs.
 № 1.0976(0.050 mm, soft) 5 pcs.
 № 1.0976(0.035 mm, soft) 5 pcs.



№ 1.1973 Kit of Large Matrices 4 types

30 pcs.:
 № 1.0973(0.050 mm, hard) 10 pcs.
 № 1.0973(0.035 mm, hard) 10 pcs.
 № 1.0973(0.050 mm, soft) 5 pcs.
 № 1.0973(0.035 mm, soft) 5 pcs.

- Anatomically pre-shaped
- Provide tight anatomical contact
- Stainless steel

Matrices of 0.035 mm thickness are more effective for restoration of small decays.

Matrices of 0.05 mm thickness are intended for restoration of large decays in subgingival area. Such matrices retain the shape well.

Kits **№ 1.1972, 1.1973, 1.1976** are intended for advanced users of sectional systems or for dentists who desire to test the difference between hard and soft matrices.

Hard matrices are easy to install.

Soft matrices are easily adapted for contact point formation.

Clinical Case



Operating Positions



Sectional Contoured Metal Matrices

All matrices are available of 4 types

- 0.050 mm, hard
- 0.035 mm, hard
- 0.050 mm, soft
- 0.035 mm, soft

Specify please the necessary type, e.g. № 1.098(0.050 mm, soft).

There are two variants of kiting of matrices of one shape - per 10 or per 50 pcs.

When ordering kit with matrices of one shape specify please the necessary kiting variant and type, e.g. № 1.0971(10) (0.035 mm, hard) or № 1.0976(50) (0.050 mm, soft).



Sectional Contoured Matrices. Small with Ledge

№ 1.0971(10) 10 pcs.
№ 1.0971(50) 50 pcs.



Sectional Contoured Matrices. Small

№ 1.0972(10) 10 pcs.
№ 1.0972(50) 50 pcs.



Sectional Contoured Matrices. Large

№ 1.0973(10) 10 pcs.
№ 1.0973(50) 50 pcs.



Sectional Contoured Matrices. Large with Ledge

№ 1.0974(10) 10 pcs.
№ 1.0974(50) 50 pcs.



Sectional Contoured Matrices. Medium with Ledge

№ 1.0975(10) 10 pcs.
№ 1.0975(50) 50 pcs.



Sectional Contoured Matrices. Medium

№ 1.0976(10) 10 pcs.
№ 1.0976(50) 50 pcs.



№ 1.1972 Kit of Small Matrices 4 types

30 pcs.:
 № 1.0972(0.050 mm, hard) 10 pcs.
 № 1.0972(0.035 mm, hard) 10 pcs.
 № 1.0972(0.050 mm, soft) 5 pcs.
 № 1.0972(0.035 mm, soft) 5 pcs.



№ 1.1976 Kit of Medium Matrices 4 types

30 pcs.:
 № 1.0976(0.050 mm, hard) 10 pcs.
 № 1.0976(0.035 mm, hard) 10 pcs.
 № 1.0976(0.050 mm, soft) 5 pcs.
 № 1.0976(0.035 mm, soft) 5 pcs.



№ 1.1973 Kit of Large Matrices 4 types

30 pcs.:
 № 1.0973(0.050 mm, hard) 10 pcs.
 № 1.0973(0.035 mm, hard) 10 pcs.
 № 1.0973(0.050 mm, soft) 5 pcs.
 № 1.0973(0.035 mm, soft) 5 pcs.

- Anatomically pre-shaped
- Provide tight anatomical contact
- Stainless steel

Matrices of 0.035 mm thickness are more effective for restoration of small decays.

Matrices of 0.05 mm thickness are intended for restoration of large decays in subgingival area. Such matrices retain the shape well.

Kits **№ 1.1972, 1.1973, 1.1976** are intended for advanced users of sectional systems or for dentists who desire to test the difference between hard and soft matrices.

Hard matrices are easy to install.

Soft matrices are easily adapted for contact point formation.

Clinical Case



Operating Positions



Saddle Matrices

39.2. P.D.



Edges of saddle matrices are performed as cylinder-shape tubes where tines of ring or springclip are inserted upon matrix installation.

All matrices are available of three shapes and two thicknesses (0.050 mm and 0.035 mm). Specify please the necessary type, for example if you need small saddle matrices of shape 2 and thickness 0.035 mm order please - №1.301(2)(0.035 mm).

	Small Matrices 12 pcs. per box	Medium Matrices 12 pcs. per box	Large Matrices 12 pcs. per box
Standard shape	 № 1.301(1)	 № 1.302(1)	 № 1.303(1)
Shape with flexible central part	 № 1.301(2)	 № 1.302(2)	 № 1.303(2)
Shape with enlarged undergingiva ledge	 № 1.301(3)	 № 1.302(3)	 № 1.303(3)

- Easy-to-use matrix system
- Optimal approximal contact without difficulties
- Three lengths
- Stainless steel

Important!

To provide perfect contact choose appropriate matrix length



incorrect



correct



incorrect

Saddle Contoured Matrices

39.1. P.D.



Edges of saddle matrices are performed as cylinder-shape tubes where tines of ring or springclip are inserted upon matrix installation.

All matrices are available of three shapes and two thicknesses (0.050 mm and 0.035 mm). Specify please the necessary type, for example if you need medium saddle contoured matrices of shape 1 and thickness 0.050 mm order please - №1.312(1)(0.050 mm).

	Small Matrices 12 pcs. per box	Medium Matrices 12 pcs. per box	Large Matrices 12 pcs. per box
Standard shape	 № 1.311(1)	 № 1.312(1)	 № 1.313(1)
Shape with flexible central part	 № 1.311(2)	 № 1.312(2)	 № 1.313(2)
Shape with enlarged undergingiva ledge	 № 1.311(3)	 № 1.312(3)	 № 1.313(3)

Installation of saddle matrix with springclip № 1.003

1. Insert tines of springclip outside the patient mouth into the matrix side tubes



2. Press springclip holders with matrix and install matrix onto the tooth



Accessories for Installation of Saddle Matrices

39.3. P.D.



№ 1.033

Ring
1 pc.

№ 1.033c

Ring with Silicone Safety Cover
1 pc.
Autoclavable

№ 1.033A

Ring
1 pc.

Ring № 1.033A has the same construction of tines as standard ring № 1.033 for saddle matrices. Addition of extra loop on the ring makes it more flexible and allows to operate in wider range of openings without increase of residual distance between tines after use.

№ 1.003

Springclip small
1 pc.

№ 1.005

Springclip large
1 pc.
More tight than № 1.003



Tines of rings & springclips are narrowed

WARNING!!!

- When operating with ring:
- 1) make sure of integrity of the ring and absence of rust on it;
 - 2) avoid high-speed opening;
 - 3) open the ring for width not more than 10 - 12 mm;
 - 4) open the ring for necessary width outside the mouth before placing it on a tooth

Installation of saddle matrix with ring № 1.033

1. Insert tines of the ring into the matrix side tubes



2. Open the ring assembled with saddle matrix via forceps № 1.099-1 (page 11) or via standard rubber dam clamp forceps and install matrix onto the tooth



Matrices can be installed on distal and on mesial tooth surface

Wooden Wedges



No 1.181 **Wooden Wedges**
superthin, short (orange)
100 pcs.



No 1.182 **Wooden Wedges**
thin, short (white)
100 pcs.



No 1.183 **Wooden Wedges**
thin, short (green)
100 pcs.



No 1.184 **Wooden Wedges**
thin, long (yellow)
100 pcs.



No 1.185 **Wooden Wedges**
medium, short (blue)
100 pcs.



No 1.186 **Wooden Wedges**
medium, long (pink)
100 pcs.



No 1.187 **Wooden Wedges**
thick, long (violet)
100 pcs.



No 1.080 **Wooden Wedges.**
2 Types
100 pcs.
white (thin, short) 50 pcs.
blue (medium, long) 50 pcs.



No 1.083 **Wooden Wedges.**
4 Types
200 pcs.
white (thin, short) 50 pcs.
yellow (thin, long) 50 pcs.
blue (medium, short) 50 pcs.
pink (medium, long) 50 pcs.



No 1.085 **Wooden Wedges.**
6 Types
400 pcs.
orange (superthin, short) 100 pcs.
white (thin, short) 100 pcs.
yellow (thin, long) 50 pcs.
blue (medium, short) 50 pcs.
pink (medium, long) 50 pcs.
violet (thick, long) 50 pcs.



No 1.281(2) **Triangle Wooden Wedges**
size 13x1,7x2,7 mm, 100 pcs.

Sizes available upon request:
No 1.281(1) 16x2,4x2,5 mm, 100 pcs.
No 1.281(3) 15x 2 x3,4 mm, 100 pcs.
No 1.281(4) 15x 3 x 3 mm, 100 pcs.

Anatomically shaped



Made of maple wood

Operating positions

Wooden wedges



Elastic wedges



Delta ring (page 12) and add-on wedges (page15)



Transparent plastic wedge



Non-transparent plastic wedges



Triangular shaped



Made of maple wood

42.1. P.D.

Folded dental bibs (33cm x 45cm)

Pastel colors



- Excellent absorbency, wide range of colors, tear resistance, nice folding

- **Manufactured by 2 plies of Tissue + 1 ply of PE**
- **Available in 1 ply Tissue + 1 ply of PE**



4 x 125 pcs in plastic bag with the sticker



1 box / 500 pcs 1 EUR pallet / 80 boxes



20': 1600 boxes
40'HQ: 3600 boxes

Folded dental bibs (33cm x 45cm)

Premium colors



- Great absorbency, attractive premium colors, nice folding

- **Manufactured by 1 ply of Tissue + 1 ply of PE**
- **Dimension: 33cm x 45cm**



4 x 125 pcs in plastic bag with the sticker



1 EUR pallet / 80 boxes



1 box / 500 pcs



20': 1600 boxes
40'HQ: 3600 boxes

Dark blue	●	KA/500 2 DB
Green lime	●	KA/500 2 GI
Orange	●	KA/500 2 O
Bordeaux	●	KA/500 2 Bord
Black	●	KA/500 2 Black

Light blue	●	KA/500 3 B
Green mint	●	KA/500 3 GM
White	●	KA/500 3 W
Pink	●	KA/500 3 P
Apricot	●	KA/500 3 A
Yellow	●	KA/500 3 Y
Lavender	●	KA/500/3 L

Color:



Non Woven Bouffant cap

PB012

- Material: PP, 14g/㎡
- Latex-free
- Size: 18", 21"
- Package: 100pcs/bag, 20bag/ctn

regular vacuumed

Polygard Dental bib

2ply paper + 1ply poly

PB017

- 100% wooden paper, medical grade PE film.
- Paper: 17g/㎡, PE film: 13g/㎡
- Package: 125pcs/bag, 4bags/ctn



Color:

Clips for bib



PB023

- Material: Steel
- Length: 18"
- Package: 1pcs/bag



PB024

- Material: Plastic
- Length: 15-1/2"
- Package: 20pcs/bag



PB025

- Material: Plastic
- Length: 14-1/5"
- Package: 1pcs/bag



PB025A

- Material: Silicone tie, SS tip, autoclavable
- Color:
- Length: 14-1/5"
- Package: 1pcs/polybag



Bibs Holder



PB026-1

- Material: Acrylic
- Length(mm): L:180 X W:133 X D:155
- Package: 1pcs/box



PB026-2

- Material: Stainless Steel, autoclavable
- Length(mm): L:177 X W:130 X D:155
- Package: 1pcs/box

Disposable Isolation Gown

PB015


- Hand sewn ties at waist
- Material: PP, 30g/35g/㎡
- Size(cm): small 115X137 medium 120X140 large 130X150
- Package: 5pcs/bag, 10bags/ctn

regular vacuumed


Tie back




Knit cuff



Velcro



Elastic cuff



Lab Gown


PB016

- Material: PP, 40g/㎡
- Size(cm): small 98X118 medium 103X123 large 107X128 extra large 110X134
- Package: 5pcs/bag, 10bags/ctn


regular vacuumed



Shirt collar



Knit collar



Elastic cuff



Knit cuff



Pocket

Shoe Cover

PB013

- Material: Non-woven
- Size: regular 16X40cm
- Color: blue
- 40g/㎡ Skid-proof
- 30g/㎡ Non Skid-proof



PB014

- 2g/3g/㎡
- Size: Regular 15X39cm
- CPE skid-proof
- PE



PRODUCT CATEGORIES



Dentistry

Endodontics

Prosthetics

Surgery and periodontology

Whitening of teeth

Prophylaxis and cosmetics

Dermatology

Musculoskeletal system

Cardiovascular system

Vitamin D and generics

Veterinary products

CITRIC ACID Plus



Application

Preparation for root canal irrigation: cleans root canals, including smear layer removal, in order to expose dentinal tubules before root canal filling.

Action description

The product consists citric acid, which dissolves well inorganic debris of a smear layer. Low surface tension of Citric Acid Plus provides much better penetration inside narrow root canals.

Product handling


Apply Citric Acid Plus on depth 1,5-2 mm shorter than working length. Apply this preparation using syringe with endodontic needle (with lateral perforation). After reaming action, place 1 - 2 ml of preparation inside root canal for about 1 minute (2 minutes max.). Then rinse accurately, preferably with purified water. Do not rinse again and again, especially when rinsing repeatedly alternately with sodium hypochlorite. Do not use directly before/after sodium hypochlorite. This may cause chlorine liberation. Use purified water between irrigating by any of root canal agents.

Storage conditions

Keep container tightly closed, **out of the reach of children**.
Store in temperature below 25°C.
Keep away from sunlight.

Filling instruments

48.1., 48.2. P.D.



Cement Spatulas Heidemann Separating Spatulas

48.3., 48.4. P.D.



980



981

48.4. P.D.



982



983



984



985



986



987

48.3. P.D.




988
Fig. 1



989
Fig. 2




990
Fig. 3



48.6. P.D.

EXPLORERS ENDODONTIC INSTRUMENTS



48.7. P.D.



DENTAL Instruments


- [Extracting Forceps](#)
- [English Pattern](#)
- [Tooth Extracting Forceps For Children](#)
- [Tooth Extracting Forceps American Pattern](#)
- [Root Splinter Forceps](#)
- [Bone Rongeurs](#)
- [Root Elevators](#)
- [Pliers For Orthodontics](#)
- [Orthodontic Instruments](#)
- [Dental Scissors](#)
- [Haemostatic Forceps](#)
- [Sponge & Dressing Forceps](#)
- [Needle Holders](#)
- [Towel Clamps](#)
- [Skin Hook & Retractors](#)
- [Retractors](#)
- [Cheek Retractors](#)
- [Elevators](#)
- [Chisel & Gouges](#)
- [Rubber Dam Instruments](#)
- [Tissue Pliers](#)
- [Cotton & Dressing Pliers](#)
- [Mouth Mirrors & Handles](#)
- [Scalpel Handles](#)
- [Handles Reference Guide](#)
- [Endodontic Instruments](#)
- [Periodontal Pocket Probes](#)
- [Explorers](#)
- [Filling Instruments](#)
- [Plastic Filling Instruments](#)
- [Amalgam Carvers](#)
- [Cavity Preparation Inst](#)
- [Excavators](#)
- [Periodontia Instruments](#)
- [Bone Files](#)
- [Bone Curettes](#)
- [Periodontia Instruments](#)
- [Gracey Curettes](#)
- [Scalers](#)
- [Cement Spatulas](#)
- [Wax & Modelling Carvers](#)
- [Wax Knives](#)
- [Plaster Spatulas](#)
- [Amalgam Carriers](#)
- [Dental Cartridge Syringes](#)
- [Napkin & Film Holders](#)
- [Matrix Retainers](#)
- [Crown Instruments](#)
- [Impression Trays](#)
- [Instrument Trays](#)
- [Holloware](#)
- [Rubber Dam Clamps](#)

ZONA INDUSTRIES
WWW.ZONADENT.COM

Excavators English Pattern

Page :- [69](#) [70](#) [71](#)

71



ROOT ELEVATORS

50.1., 50.2. P.D.



239
Fig. 1B
Bein



240
Fig. 2B
Bein



241
Fig. 3B
Bein



242
Fig. 4B
Bein



243
Fig. 5B
Bein



244
Fig. 6B
Bein



245
Fig. 4SA
Schmeckebeer-Apexo



246
Fig. 5SA
Schmeckebeer-Apexo

50.1. P.D.



247
Fig. 34S
Seldin



248
Fig. 34
Seldin



249
Fig. 2S
Seldin



250
Fig. 3S
Seldin



251
Fig. 1F
Flohr




252
Fig. 2F
Flohr



253
Fig. 3F
Flohr

50.2. P.D.

50.3. P.D.



101
Fig.1
Upper centrals
& canines



102
Fig.2
Upper laterals
& bicuspids



103
Fig.3
Upper incisors
& canines



104
Fig.4
Lower incisors
& canines



105
Fig.7
Upper primolars




106
Fig.8
Lower primolars



107
Fig.13
Lower primolars

EXTRACTING FORCEPS
 ENGLISH PATTERN



109
 Fig.18
 Left upper molars

110
 Fig.18A
 Upper molars

111
 Fig.19
 Upper wisdoms

108
 Fig.17
 Right upper molars

112
 Fig.20
 Lower wisdoms

113
 Fig.21
 Lower molars

114
 Fig.22
 Lower molars
 Hawk's bills




115
Fig.23
Right lower molars




116
Fig.24
Left lower molars



117
Fig.29
Upper root & incisors



118
Fig.29N
Roots of upper incisors



119
Fig.30
Upper roots & incisors



120
Fig.31
Lower roots




121
Fig.32
Lower molars



122
Fig.33
Lower roots



123
Fig.33A
Lower roots




124
Fig.33C
Lower roots




125
Fig.33M
Lower roots with long beak



126
Fig.41
Upper roots



127
Fig.44
Upper roots



128
Fig.45
Lower roots




129
Fig.51
Upper roots




130
Fig.51A
Upper roots narrow beak



131
Fig.51C
Upper roots



132
Fig.52
Upper roots & canines
with wider beak



133
Fig.52A
Upper teeth & resection




134
Fig.55
Separating upper molars roots



135
Fig.56
Separating lower molars roots




136
Fig.67
Upper wisdoms




137
Fig.67A
Upper wisdoms


EXTRACTING FORCEPS
 ENGLISH PATTERN




138
 Fig.68
 Lower roots




139
 Fig.73
 Lower molars (hawk's bill)




140
 Fig.74
 Lower incisors & roots




141
 Fig.74M
 Lower incisors & roots




142
 Fig.74N
 Lower incisors & roots



143
 Fig.75
 Lower premolars & canines




144
 Fig.76
 Upper roots




145
 Fig.76N
 Upper roots
 Narrow beak


EXTRACTING FORCEPS
 ENGLISH PATTERN




146
 Fig.79
 Lower wisdoms



147
 Fig.79A
 Lower wisdoms




148
 Fig.86A
 Lower molars




149
 Fig.86B
 Lower molars



150
 Fig.86C
 Lower molars



151
 Fig.87
 Cowhorn beak




152
 Fig.89
 Right upper molars




153
 Fig.9
 Left upper molars

EXTRACTING FORCEPS


ENGLISH PATTERN




154
 Fig.94
 Right upper molars




155
 Fig.95
 Left upper molars




156
 Fig.99
 Lower molars
 (broken-down)




157
 Fig.99 1/4
 Lower molars
 (broken-down)




158
 Fig.121
 Upper wisdoms



159
 Fig.136
 Upper centrals
 & roots







160
 Fig.137
 Lower centrals
 & roots



161
 Fig.167
 Universal pincer
 for lower teeth
 molars & premolars


50.4. P.D.




Curved, saw edge



413
Fig. 1
La Grange
11.5cm



414
Fig. 1
Quinby
12.5cm




415
Fig. 2
Quinby
12.5cm



416
Fig. 3
Goldman-Fox (Special)
12.5cm




417
Fig. 1
Kelly
16cm



418
Fig. 2
Kelly
16cm




419
Fig. 1
Locklin
16.5 cm




420
Fig. 2
Locklin
16.5cm


50.5. P.D.




50.6. P.D.



Designed for Frontal & Palate Suture



50.7. P.D.



56. P.D.

Expanding Dental Floss

HILO DENTAL
FIL DENTAIRE



30m



Ref. 48

EAN: 8411732100487

INNER / ENVASE / PACK: 12
 OUTER / CAJA / BOITE: 144

DESIGNED TO EXPAND TO REMOVE MORE PLAQUE .GB
 MADE WITH POLYESTER WITH ALCO VERA .ES
 AROMA A MENTA
 WAX COATING

SE EXPANDE PARA ELIMINAR MÁS PLACA DENTAL .FR
 FABRICADO EN POLIÉSTER CON ALCO VERA
 AROMA A MENTA
 CONTIENE CERA

IL S'ÉLARGIT POUR UN NETTOYAGE EN PROFONDEUR .FR
 FABRIQUÉ EN POLYESTER A L'ALCOES VERA
 AROMA A MENTHE
 AVEC CIRE

Aloe Vera Dental Floss

HILO DENTAL
FIL DENTAIRE



50m



Ref. 41

EAN: 8411732104119

INNER / ENVASE / PACK: 12
 OUTER / CAJA / BOITE: 144

MADE WITH NYLON 6.6 .GB
 WITH ALCO VERA
 MINT AROMA
 WAX COATING

FABRICADO EN NYLON 6.6 .ES
 CON ALCO VERA
 AROMA A MENTA
 CONTIENE CERA

FABRIQUÉ AVEC NYLON 6.6 .FR
 A L'ALCOES VERA
 AROMA A MENTHE
 AVEC CIRE

Strawberry Dental Floss

HILO DENTAL
FIL DENTAIRE



50m



Ref. 43

EAN: 8411732000435

INNER / ENVASE / PACK: 12
 OUTER / CAJA / BOITE: 144

MADE WITH POLYESTER .GB
 STRAWBERRY AROMA
 WAX COATING

FABRICADO EN POLIÉSTER .ES
 AROMA A FRESA
 CONTIENE CERA

FABRIQUÉ EN POLYESTER .FR
 AROMA A FRAISE
 AVEC CIRE



Search our catalog



BASKET (0)



OUR COMPANY

OUR BRANDS

PRODUCTS

INTERNATIONAL

STRATEGIC PARTNERS

CAREER

NEWS

[Home](#) - [Our Products](#) - [Surgery](#) - [wound treatment](#) - [ribbon gauze \(tamponades\)](#) - [NOBADENT®](#)

NOBADENT®

Dental Cotton Rolls

Size PU REF: 149408 **57.1. P.D.**

PZN: 02477143

(PU = Packaging Unit)

Size PU REF: 149410 **57.2. P.D.**

PZN: 02477166

(PU = Packaging Unit)

The dental rolls are made of pure cotton, bleached without chlorine. Dental rolls are used for dental treatment. They serve for the absorption of saliva, secretion and blood. NOBADENT® dental rolls have a high capacity of absorption, at the same time having a padding function and holding off function. They are packed in a convenient folding box.

Palatinalkopfbürsten mit Kunststoffkern

cup brush - brosse champignon - cepillo de copa spazzole testa

Art.-Nr. rund round	Art.-Nr. zylindrisch cylindric	in mm in mm
---------------------------	--------------------------------------	----------------

Naturborste schwarz, natural bristle - soie nature - cerda naturaleza - setola naturale

140-30P		30
140-35P		35
140-40	140-40Z	40
140-50	140-50Z	50
140-60	140-60Z	60
140-70	140-70Z	70
140-80	140-80Z	80



Messingdraht gewellt 0,08, brass wire crimped - fil de laiton ondulé - láton ondulado - ottone ondo 0,08

143-135		35
143-140	143-140Z	40
143-150	143-150Z	50
143-170	143-170Z	70
143-180	143-180Z	80


Bohrerreinigungsbürsten,

burs cleaning brushes - brosses pou le nettoyage des fraises - limpia-fresas - spazzole per la pulitura di frese

Art.-Nr. art.-no.	Stärke thickness	Maß size
----------------------	---------------------	-------------


Messingdraht gewellt, brass wire crimped - fil de laiton ondulé - láton ondulado - ottone ondo

318	0,15	55 x 23
319	0,15	70 x 28



Micro Applicator

Disposable



Bendable Brush Applicators

TA026C

- Chemical resistant fibers hold up to the harshest dental materials
- Flexible and bendable to reach all areas
- Package: 144pcs/bag, 100bags/case

Composite Brush

TA026D

- Made of rigid plastic to withstand repeated autoclaving (up to 120°C)
- Tip can be replaced, save cost



Ref. #	Description	Package
TA026D-1	Handle+Tip	100pcs/bag, 40bag/case
TA026D-2	Tip (black, white)	100pcs/bag

Valustick Sticky Applicators

TA026A

- Non-absorbent
- Non-linting fibers
- Bendable neck
- Different color for different solutions
- Package: 100pcs/bottle, 4bottle/box, 60box/ctn

Ref. #	Size	Color
TA026A-1	Regular	Green, Purple, Orange, Blue
TA026A-2	Fine	Yellow, Mauve
TA026A-3	Super fine	White
TA026A-4	Extra fine	Black



Brush Applicator

TA026B

- Chemical resistant fibers hold up to the harshest dental materials
- Bends close to tip for improved control and easily reaching
- Usage: Whitening, Fluoride, Etchants, Bonding Agents, Sealants, Desensitizers, Adhesives / Resins

Ref. #	Description	Package
TA026B-1	Single head brush	100pcs/bundle 128bundles/ctn
TA026B-2	Double-head brush	100pcs/bag 4bags/box 32boxes/ctn

Fixing Plastic Wedges



TA028

- Anatomically designed

Ref. #	Description	Package
TA028A-1	L (16X2.6X2.2mm)	100pcs/box
TA028A-2	M (14X2.2X2mm)	100pcs/box
TA028A-3	S (12X1.9X1.9mm)	100pcs/box
TA028A-4	XS (10X1.6X1.5mm)	100pcs/box

Fixing Plastic Wedges with hole



TA028B

- Anatomically designed

Ref. #	Description	Package
TA028B-1	L (17.8X2.8mm)	100pcs/box
TA028B-2	M (17.8X2.45mm)	100pcs/box
TA028B-3	S (17.8X2.1mm)	100pcs/box
TA028B-4	XS (17.8X1.7mm)	100pcs/box




LB002-3

- Material: PP
- Handle type: 185*33mm
- Package: 40pcs/box, 20bags/ctn



LB002-4


- Material: PP
- Slim type: 110X4mm, white
- Package: 25pcs/box, 16boxes/ctn



Ortho Study Model Molds

LB003A

- Material: TPR
- Package: 6pcs/box, 30box/ctn



Ortho Study Model Molds

LB003B

- Material: TPE Medical grade
- Package: 6pcs/box, 30box/ctn



Alcohol Lamp

Ref.#	Description	Package	Ref.#	Description	Package
LB015	Blow lamp 100ml	1pc/box	LB016	Glass 125ml	1pc/box



Neaten Case

LB017

- Material: PP
- Slim type: 152X134X205mm
- Packing: 50pcs/ctn

Poly-coated Mixing Pad

- Used in conjunction with the spatula
- Coated paper 128g/m²
- 50pcs/book



Ref.#	Size	Package
LB022-1	2"*2"	10books/box, 60boxes/ctn
LB022-2	3"*3"	10books/box, 24boxes/ctn
LB022-3	3"*5"	5books/box, 32boxes/ctn
LB022-4	4"*5.5"	5books/box, 24boxes/ctn
LB022-5	4"*7"	5books/box, 20boxes/ctn
LB022-6	5.5"*7"	5books/box, 18boxes/ctn

Soft/Hard Plastic Sheet

LB023

- Mainly used in: Dental medical retainer and other products blister.
- Size: 125X125mm



Ref.#	Size (mm)	Package
LB023-1	Soft, Thickness: 1.0	20pcs/bag
LB023-2	Soft, Thickness: 1.5	15pcs/bag
LB023-3	Soft, Thickness: 2.0	10pcs/bag
LB023-4	Soft, Thickness: 3.0	10pcs/bag
LB023-5	Hard, Thickness: 1.0	20pcs/bag
LB023-6	Hard, Thickness: 1.5	15pcs/bag
LB023-7	Hard, Thickness: 2.0	10pcs/bag



Disposable Articulator




Ref.#	Description	Package
LB024-1B	Black	100pcs/bag
LB024-1C	Color	20bags/ctn

69. P.D.

Transparent Strips

Transparent strips are 0.05 mm thick



№ 1.040 **Transparent Strips**
8 mm wide, 100 mm long
100 pcs.

№ 1.041 **Transparent Strips**
10 mm wide, 100 mm long
100 pcs.

№ 1.240 **Transparent Striproll**
8 mm wide, 10 m in length
(round box Ø 32 mm)

№ 1.241 **Transparent Striproll**
10 mm wide, 10 m in length
(round box Ø 32 mm)

№ 1.242 **Transparent Striproll**
8 mm wide, 10 m in length
(squared box 50x50 mm)

№ 1.243 **Transparent Striproll**
10 mm wide, 10 m in length
(squared box 50x50 mm)

№ 1.440 **Transparent Stopstrips**
8 mm wide, 100 mm long
12 pcs.

№ 1.441 **Transparent Stopstrips**
10 mm wide, 100 mm long
12 pcs.


Operating Positions



Striprolls of alternative lengths and widths can be supplied upon special request.




Stopstrip offers the possibility to operate with one hand only



Transparent Cervical Matrix System



№ 1.141 **Transparent Cervical Matrix Bands**
matrix bands 10x37mm 100 pcs.
cervical former 1 pc.
stopper 1 pc.



The system operates as follows

1. Insert matrix into the cervical former (fig. 1).

2. Cut out matrix in accordance with the tooth shape (fig. 2). Regulate the matrix position according to the tooth size.

3. Fix matrix in the cervical former by the stopper (fig. 3).

4. Perform the tooth restoration (figs. 4 - 6).



70.2. P.D.

https://www.swedent.se/products/endodontics/

Search

Endodontics



RUBBER DAM

In 3 colours, dark, light and green. Withstand for autoclaving upto 118°C (224°F). Vacuum packed for prolonged shelf life.

IN ROLLS:

- No. 3080 – Thin dark, roll of 450 x 15 cm
- No. 3081 – Thin light, roll of 450 x 15 cm
- No. 3060 – Thin green, roll of 450 x 15 cm
- No. 3082 – Medium dark, roll of 450 x 15 cm
- No. 3083 – Medium light, roll of 450 x 15 cm
- No. 3061 – Medium green, roll of 450 x 15 cm
- No. 3084 – Heavy dark, roll of 225 x 15 cm



RUBBER DAM

In 3 colours, dark, light and green. Withstand for autoclaving upto 118°C (224°F). Vacuum packed for prolonged shelf life.

IN CUT PIECES:

- No. 3085 – Thin dark, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3086 – Thin light, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3070 – Thin green, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3087 – Medium dark, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3088 – Medium light, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3071 – Medium green, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3089 – Heavy dark, 5" x 5", 12,5 x 12,5 cm, 52 pcs
- No. 3090 – Thin dark, 6" x 6", 15 x 15 cm, 36 pcs
- No. 3091 – Thin light, 6" x 6", 15 x 15 cm, 36 pcs
- No. 3075 – Thin green, 6" x 6", 15 x 15 cm, 36 pcs
- No. 3092 – Medium dark, 6" x 6", 15 x 15 cm, 36 pcs
- No. 3093 – Medium light, 6" x 6", 15 x 15 cm, 36 pcs
- No. 3076 – Medium green, 6" x 6", 15 x 15 cm, 36 pcs
- No. 3094 – Heavy dark, 6" x 6", 15 x 15 cm, 36 pcs

Teeth Brushing Model With tongue

A0002



A0002



A0003
Times Brushing Model

Removable Standard Model



A0001
● 28pcs, hard gum



A0011
● 28pcs, soft gum

Ortho Metal & Ceramic



A0012



A0013



A0014

Teeth Model

Standard Model



A0005
● 28pcs, hard gum, wax fix



A0006
● 28pcs, hard gum

Standard Children Model



A0009
● 24pcs, soft gum

Restorion With Implant



A0010

● Clear, Blue, Orange, Pink, Yellow can be choose



A0015

Pathology Model with half implant



A0016

Clear primary teeth model



A0017



A0018

Transparent Standard Model



A0019




70.1. P.D.


Calmp forceps



596
Ivory light weight
17cm



597
Brewer
17.5cm



598
Stokes
17.5cm


Punch forceps



599
Ainsworth
17cm



600
Ivory
17cm



Rubber dam frame

601

Molar Clamps Winged



1156
Fig. 30



1157
Fig. 31



1158
Fig. 51



1159
Fig. 138



1160
Fig. 139



1161
Fig. 15



1162
Fig. 18



1163
Fig. 24



1164
Fig. 25



1165
Fig. 26



1166
Fig. 28

Bicuspid Clamps Winged



1167
Fig. 206



1168
Fig. 207



1169
Fig. 208



1170
Fig. 209



1171
Fig. 22



1172
Fig. 27



1173
Fig. 29

Molar Clamps Winged



1174
Fig. 26A



1175
Fig. 28A



1176
Fig. 30A



1177
Fig. 31A



1178
Fig. 51A



1179
Fig. 27



1180
Fig. 29



1181
Fig. 60

Bicuspid Clamps Winged



1182
Fig. 200



1183
Fig. 201



1184
Fig. 202



1185
Fig. 203



1186
Fig. 204



1187
Fig. 205

Molar Clamps Winged



1188
Fig. 210




1189
Fig. 211



1190
Fig. 212

Labial Clamps Winged

76. P.D.



Dokumentą elektroniniu
parašu pasirašė JULIJA
LIČKŪNIENĖ

Data: 2020-06-22 14:37:39

Paskirtis: 492219

Vieta: Kaunas

Kontaktinė informacija:

Viešųjų pirkimų specialistė

Julija Ličkūnienė UAB

"SKIRGESA" Energetikų g. 8,

LT-52461 Kaunas Tel.: +370

37457746 El.p.:

konkursai@skirgesa.lt