

68 pozicija

Pressure Connecting Tube

Prailginimo linija pritaikyta naudoti su Ulrich Medical CT motion XD8000 automatiniu švirkštu



P/N: 601150

Length: 250cm

Ilgis: 250cm

Pressure: 350PSI

Slėgis: 350 psi = 24 Bar

Description: CT/MR Straight Tube with 2 safety dual check valves Su dviem apsauginiais vožtuvais

Used for Ulrich CT motion XD8000 injection system

Pritaikyta naudoti su Ulrich Medical CT motion XD8000 automatiniu švirkštu

Packing: 200PCS/CASE

filling volume: 5ml/m Sistemos užpildymo tūris 5 ml/1 m

Latex-free, Sterile, Disposable, DEHP free, Transparent

Be latekso

Sterili

Vienkartinė

Be DEHP

Permatoma



antmed

Partner For Life



Shenzhen Ant Medical Devices Co., Ltd.

18 Jinhui Ave., Pingshan New District, Shenzhen 518122, China

Tel: +86 755 8606 0992

Fax: +86 755 8606 0991

E-mail: info@szat.com.cn

http://www.szat.com.cn

STATEMENT OF COMPATIBILITY

As a manufacturer Shenzhen Ant Medical Devices Co., Ltd., 18 Jinhui AVE., Pingshan New District, 518122 Shenzhen, China, we declare that product **P/N 601150** is fully compatible with **Ulrich Medical system**.

Shenzhen Ant Medical Devices Co., Ltd.

Pingshan New District, 518122 Shenzhen, China

Tel: +86-755 86060992



Lareta

2022.3.10